

BLACKline

Die BLACKline-Instrumente erhalten mit einer DLC – Diamond-like Carbon – Beschichtung ein High-End Finish.

Durch das matte Finish werden störende Reflexionen eliminiert. Dies ist ein wichtiger Aspekt für die Arbeit unter dem Mikroskop, die in der Endodontie immer mehr an Bedeutung gewinnt.

Die BLACKline verbindet Arbeitskomfort und Ergonomie und bietet in allen Bereichen der Diagnostik überlegende Vorteile.

The BLACKline instruments receive with a DLC (Diamond-like carbon) coating a high-end finish.

The matte finish eliminates disturbing reflections. This is an important aspect for working under the microscope, which is becoming increasingly important in endodontics.

The BLACKline increases working comfort while adding ergonomic elements and offers superior advantages for every purpose.

BLACKline

ULTRA FS
MEGA FS Rhodium


- ▶ DLC-Beschichtung (Diamond-like carbon) für ein High-End-Finish
- ▶ keine Reflektion durch extra matte, tief-schwarze Schicht
- ▶ augenschonendes Arbeiten bei starkem Licht und unter Vergrößerung
- ▶ DLC coating (Diamond-Like Carbon) for a high-end finish
- ▶ antireflective, dark black coating
- ▶ avoids visual fatigue by working under light and magnification



BLACKline

ULTRAmicro
MICROflex ultra
MEGAmicro
MICROflex rhodium

- ▶ Microspiegel mit biegbarem Stiel für eine punktgenaue Ausrichtung des Spiegels
- ▶ auch für die Parodontologie, insbesondere für Aufklappungen
- ▶ micro mirrors with bendable stem for a precise alignment of the mirror
- ▶ can be also used for Periodontics, especially for flap openings



BLACKline


Handles + Explorers

- ▶ keine Reflektion durch extra matte, tief-schwarze DLC-Schicht
- ▶ optimiertes Griff-Design für eine angenehme Haptik
- ▶ Sonde mit hoher Bruchfestigkeit und Elastizität

Die Sonden sind permanent montiert.

- ▶ antireflective, dark black DLC coating
- ▶ Handle design with maximum grip
- ▶ Explorer tip with high bending strength flexibility







The Explorers are permanently fixed.




BLACKline ULTRA FS 3-Pack			SS-thread	CS-thread
BLACK ULTRA	4-5 plane	Ø 22-24 mm	B6160-B6161	B6165-B6166
BLACK ULTRAduo	4-5 plane	Ø 22-24 mm	B6640-B6641	B6645-B6646
BLACK ULTRAretract, bifurcated stem	4-5 plane	Ø 22-24 mm	B6690-B6691	B6695-B6696

BLACKline MEGA FS Rhodium 3-Pack			SS-thread	CS-thread
BLACK MEGA	4-5 plane	Ø 22-24 mm	B6060-B6061	B6070-B6071
BLACK MEGAduo	4-5 plane	Ø 22-24 mm	B6100-B6101	B6170-B6171
BLACK MEGAretract, bifurcated stem	4-5 plane	Ø 22-24 mm	B6680-B6681	B6685-B6686

BLACKline each	BLACK ULTRA FS		BLACK Rhodium FS	
	ULTRAmicro	MICROflex ultra	MEGAmicro	MICROflex rhodium
Ø 3 mm	B6181	B5181	B6200	B5200
Ø 5 mm	B6182	B5182	B6220	B5220
3x6 mm oval	B6183	B5183	B6232	B5232
3x6 mm drop-shaped	B6184	B5184	B6240	B5240

BLACKline ERGOtouch Explorers each	Form A 010 Fig. 23 REF	Form C 030 Fig. 8 REF	Form D 040 Fig. 9 REF	Form F 060 Fig. 6 REF	Form G 070 Fig. 17 REF
					
	B600-A	B600-C	B600-D	B600-F	B600-G

BLACKline ERGOtouch Handles each	SS-thread	CS-thread
	B442	B443