


premium plus international



ISO13485: 2016



# dental products 2020



innovators, designers, manufacturers,  
partners and much, much more...

premium  
plus™

## P.1-2 INTRODUCTION

Who are we  
Why choose us  
What do we make  
Where are we

## P.3 DIAMOND BURS FG

Diamond Burs FG Group *I*  
Diamond Burs FG Group *II*  
Diamond Burs FG Group *III*  
Diamond Burs FG Group *IV*  
Diamond Polishing Paste

## P.4-10 SMALL EQUIPMENT

Apex Locators  
Curing Lights  
Curing Light Meter  
High and Low Speed Handpieces

## P.11-17 PROPHY

Prophy Angles  
Prophy Brushes  
Prophy Cups  
Prophy Handpieces  
Prophy Points

## P.18 SYRINGE TIPS

Disposable Syringe Tips  
Disposable Syringe Tip Shields  
Metal Syringe Tip

## P.19 APPLICATORS

Brush Applicators  
Brush Tips and Handle  
Composite Gun and Tips  
Micro Applicators  
Micro Applicator Dispensers  
Micro Applicator Tips

## P.20 EVACUATION

HVE Suction Tubes  
Oral Evacuators  
Saliva Ejectors  
Surgical Aspirators

## P.21-22 IMPRESSION

Alginate and Plaster Spatulas  
Autoclavable Impression Trays  
Autoclavable Mixing Bowls  
Bite Registration Material  
Disposable 3-in-1 Trays  
Disposable Impression Syringes  
Disposable Impression Trays  
Disposable Mixing Tips  
Intraoral Syringe Tips

## P.23-24 BARRIERS

Barrier Film  
Bite Block Sleeves  
Chair Sleeves  
Curing Light Guide Sleeves  
Curing Light Sleeves  
Flowable Composite Sleeves  
Handpiece & Pen Sleeves  
Headrest Covers Non Woven  
Headrest Sleeves  
Impression Gun Sleeves  
Instrument Tray Sleeves

## P.23-24 BARRIERS

Intraoral Camera Sleeves  
Keyboard Sleeves  
Lab Transport Bags  
Mouth Mirror Sleeves  
PC Mouse Sleeves  
PSP Barrier Envelopes  
Suction Tube Sleeves  
Syringe Sleeves  
T-Style Light Handle Sleeves  
Tubing Sleeves  
X-Ray Head Sleeves  
X-Ray Sensor Sleeves

## P.25 DISPOSABLE

Disposable Bibs  
Disposable Masks  
Disposable Napkin Holders  
Disposable Tray Cover Papers

## P.26 PRACTICE BUILDERS

Disposable Bibs  
Disposable Cups  
Disposable Masks  
Interdental Brushes

## P.27-30 ENDO INSTRUMENTS

Aluminium Endo Holder  
Endo Brush Applicators  
Endo Cushions & Sponges  
Endo Dispensers and Holders  
Endo Instruments  
Endo Measure Round  
Endo Measuring Blocks  
Endo Rulers  
Endo Stands  
Endo Stops  
Rubber Dam Frames  
Rubber Dam Punches and Forceps

## P.31-34 SUNDRIES

Acrylic Mixing Cups  
Aluminium Bur Blocks  
Articulating Mounting Plates  
Aspirator Brushes  
Autoclavable Cheek Retractors  
Autoclavable Mouth Lippers  
Autoclavable Napkin Holders  
Autoclavable Prophy Rings  
Bite Tabs and Loops  
Brass Bur Brush  
Cotton Roll Dispenser  
Dental Appliance Boxes  
Dental Bib Holder and Clip  
Dental Mandrels  
Disposable Cheek Retractors  
Disposable Composite Tips  
Disposable Curved Tip Syringes  
Disposable Dappen Dishes  
Disposable Evacuation Traps  
Disposable Eye Shields  
Disposable Fluoride Trays  
Disposable Mixing Pads  
Disposable Mixing Wells  
Disposable Mouth Lippers  
Disposable Mouth Mirrors  
Disposable Plastic Wedges  
Disposable Prophy Rings  
Disposable Sticky Applicators

## P.31-34 SUNDRIES

Face Shields  
Fracture Detector  
Glass Dappen Dishes  
Glass Mixing Slab  
Implant Scaler  
Jumping Teeth  
Magnetic Bur Stands  
Mouth Props  
Patient Chair Cushions  
Patient Tongue Cleaners  
Pen/Toothbrush Holder  
Plastic Bur Dispensers  
Prophy Ring Stainless Steel  
Protective Eyewears  
Retainer/Mouth Boxes/Denture Baths  
Sterilization Pouches  
Suction Tip Sleeves and Adaptors  
X-Ray Film Holders

## P.35-36 ACRYLIC PRODUCTS

Barrier Film Dispenser  
Bib and C-Fold Towel Dispensers  
Chain Dispenser  
Composite Dispensers and Organizers  
Composite Syringe Holders  
Cup Dispensers  
Glove and Tissue Dispensers  
Grinding Box  
Impression Gun Holders  
Light Protection Covers  
Light Protection Hand Shield  
Mask Dispenser  
Micro Applicator Dispensers  
Pouch Dispensers

## P.37-42 ORTHODONTIC

Brackets  
Buccal Tubes  
Orthodontic Accessories  
Orthodontic Pliers  
Orthodontic Wires

## P.43-46 PIEZO SCALERS

LED Scalers  
Scaler Handpieces  
Scaler Tips  
Scaler Wrenches and Accessories

## P.47-52 HAND INSTRUMENTS

Amalgam Instruments  
Filling Instruments  
Impression Trays  
Mouth Mirrors and Handle  
Operation Instruments  
Probes  
Scalers and Curettes  
Sundry Instruments  
Surgical Instruments

## P.53 INSTRUMENTS

Autoclavable Cabinet Trays  
Autoclavable Bur/Endo Organizers  
Autoclavable Instrument Trays  
Instrument Tray Racks  
Needle Station  
Silicone Instrument Mat  
Silicone Instrument Rings

## Who are we ?



At **Premium Plus**, manufacturing dental products is what we do. Only for the past decade we've not just made dental products... we've made better dental products.

Through innovation, research, investment in manufacturing, technical know-how and over 40 years experience in the dental trade, we've learned how to control and improve every step of the manufacturing process... and that's a huge advantage to our customers.


With a global view of the dental industry and knowledge of the different geographical markets, we have learned to understand what you expect and are committed to far exceeding your requirements.



## Why choose us ?

Whether its under our own **Premium Plus** brand, private labels or OEM projects, our emphasis is on quality for not only the products, but also its presentation. We not only manufacture in-house – We innovate and develop products, research raw material selection, manufacture moulds, construct production lines, create and design artwork and packaging for our partners, warehouse goods and turn our customers' ideas and requests into reality... providing the complete package.

Our continual investment in purpose-built machinery, technical training and market research along with our flexible approach keeps us at the forefront of dental manufacturers.



# What do we make ?

Over the past decade our range has increased to cover small equipment and the majority of plastic and paper disposable dental products. This gives us, and more importantly you, the ability to offer a huge range that includes LED curing lights, apex locators, handpieces, disposable prophy cups and angles, micro applicators, syringe tips, barrier sleeves, saliva ejectors and suction tubes etc.



disposable oral evacuators

disposable saliva ejectors



LED curing lights



disposable syringe tips



disposable trays for implants



disposable 3-in-1 trays



disposable prophy angles



disposable prophy cups



disposable prophy brushes



disposable micro applicators



disposable PSP barrier envelopes

# Where are we ?

**Premium Plus** now supplies to over 90 countries and regions throughout the world and has regional offices in the USA, Japan, Germany and the United Kingdom. Our corporate head office in Hong Kong oversees the production of our range in our purpose-built manufacturing site. It is also home to our dedicated regulatory affairs team that are constantly working on not only maintaining but exceeding standards set by ISO, CE and the FDA. We are proud of the fact that we hold **ISO13485:2016** certification for the **Quality Management System of Medical Devices** and have all of our class IIa products certified with the **CE<sub>0123</sub>** mark by **TUV SUD**. We know that with our experience, our integrated range of production processes and our ever-expanding range we have the ability to cover all our customers' requirements.



ISO13485:2016



# Diamond Burs FG and Polishing Paste

## DIAMOND POLISHING PASTE EXTRA FINE

This diamond standard polishing paste, 2.0 g extra fine diamond grit, 1-2 microns produces an excellent high gloss and natural brightness on **all restoration materials, composites, ceramics, gold etc.** and **tooth enamel.**

The grey coloured (flavourless) paste is specially formulated for intra and extra oral use and is hydro-soluble for easy cleaning.



PASD002  
Diamond Polishing  
Paste 2.0 gm  
Extra Fine

### DISPOSABLE PROPHY BRUSHES RA REGULAR



Disposable  
Prophy Brushes  
100 pcs/box

### RECIPROCATING CONTRA ANGLE HANDPIECES STANDARD OR MINI HEAD



L20-4PY ( R )  
60° Turn Reciprocating  
Push Button Standard Head  
CA Handpiece  
4:1 Reduction




L20-4PY ( RM )  
60° Turn Reciprocating  
Push Button Mini Head  
CA Handpiece  
4:1 Reduction

**THIS PASTE CONTAINS EXTRA FINE DIAMOND GRIT FOR HIGH SHINE POLISHING. FOR BEST RESULTS, POLISH WITH A DISPOSABLE PROPHY BRUSH (REGULAR HARDNESS) TOGETHER WITH RECIPROCATING HANDPIECE (NO SPLATTER, NO HEAT).**

## DIAMOND BURS FG 25 PIECES BULK PACK

### ROUND



ISO No.	801	801	801	801	801	801	801	801	801	801	801L	805	805	805	806	806	806	806	807	830	830	830	
Size Ø 1/10 mm	008	009	010	012	014	016	018	021	023	025	031	016	012	014	018	009	010	012	014	014	010	012	014
Length (mm)	0.8	0.9	1.0	1.2	1.4	1.6	1.8	2.1	2.3	2.5	3.1	1.6	1.2	1.6	1.6	2.0	2.0	3.0	3.0	4.0	2.7	2.7	2.7
Medium Order No. Without Ring 115µm	008M	009M	010M	012M	014M	016M	018M	021M	023M	025M	031M	016M	012M	014M	018M	009M	010M	012M	014M	014M	010M	012M	014M
Coarse Order No. Green Ring 150µm				012C	014C	016C	018C	021C	023C	025C	031C	016C	012C	014C	018C				014C	014C	010C	012C	014C
Fine Order No. Red Ring 40µm							018F	021F	023F	025F	031F												
SuperFine Order No. Yellow Ring 20µm							018XF	021XF	023XF	025XF	031XF												

### INVERTED CONE



### INVERTED CONE WITH COLLAR




### PEAR




actual size

### CYLINDER



### TAPER FLAT END



### TAPER ROUND END (CHAMFER)



### CYLINDER ROUND END




### WHEEL ROUND EDGE



ISO No.	835	835	835	835L	835	835	836	836	836	847	847	848	848	848	856	856	850	850	850	880	880	881	909
Size Ø 1/10 mm	008	009	010	010	012	014	010	012	014	014	016	014	016	018	014	016	014	016	018	010	012	012	041
Length (mm)	3.0	3.0	3.0	4.0	4.0	4.0	6.0	6.0	6.0	8.0	8.0	10.0	10.0	10.0	8.0	10.0	10.0	10.0	10.0	8.0	8.0	10.0	1.0
Medium Order No. Without Ring 115µm	008M	009M	010M	010M	012M	014M	010M	012M	014M	014M	016M	014M	016M	018M	014M	016M	014M	016M	018M	010M	012M	012M	041M
Coarse Order No. Green Ring 150µm		009C	010C	010C	012C	014C	010C	012C	014C	014C	016C	014C	016C	018C	014C	016C	014C	016C	018C	010C	012C	012C	041C
Fine Order No. Red Ring 40µm															014F	016F				010F	012F	012F	

### POINTED CONE




### POINTED BUD



### EGG



### PEAR LONG



Diamond Burs FG  
25 pcs/pack Bulk  
Group I 008-014 Medium  
Group II 016-041 Medium  
Group III 008-014 C / F / SF  
Group IV 016-041 C / F / SF

ISO No.	858	858	859	859	859	859	859L	859L	862	862	862	863	863	368	368	379	379	830L	830L	830L
Size Ø 1/10 mm	010	012	010	012	014	016	012	016	010	012	014	012	014	018	023	018	023	014	016	022
Length (mm)	8.0	8.0	10.0	10.0	10.0	10.0	11.5	11.5	8.0	8.0	8.0	10.0	10.0	4.5	5.0	3.5	4.5	5.0	5.0	5.0
Medium Order No. Without Ring 115µm	010M	012M	010M	012M	014M	016M	012M	016M	010M	012M	014M	012M	014M	018M	023M	018M	023M	014M	016M	022M
Coarse Order No. Green Ring 150µm						016C	012C	016C	010C	012C	014C	012C	014C	018C	023C	018C	023C	014C	016C	022C
Fine Order No. Red Ring 40µm	010F	012F	010F	012F	014F	016F	012F	016F	010F	012F	014F	012F	014F	018F	023F	018F	023F	014F	016F	022F
SuperFine Order No. Yellow Ring 20µm	010XF	012XF	010XF	012XF	014XF	016XF	012XF	016XF	010XF	012XF	014XF	012XF	014XF	018XF	023XF	018XF	023XF	014XF	016XF	022XF

## C01-M / C02-M LED CURING LIGHTS MINI (with USB Charger)



- 360° Rotatable tips
- Automatic compensation circuit for constant power output
- Fully covered membrane panel for easy cleaning and infection control
- Light weight: approx. 85 gms
- Dual modes: FULL / SOFT
- 12 Operation times
- Fully charged battery provides over 100 x 10 second curing cycles

C01-M / C02-M SPECIFICATIONS		
Model	C01-M (F.O. Light Guide)	C02-M (Direct Light Source)
LIGHT TIP	8 mm ( dia. )	11 mm ( dia. )
LIGHT WAVELENGTH	440~480 nm	440~480 nm
MODE	FULL	1,200 mW/cm <sup>2</sup> for 5, 10, 15, 20, 25, 30 seconds±10%
	SOFT	1 <sup>st</sup> five seconds from 0-900 mW/cm <sup>2</sup> next 5, 10, 15, 20, 25 seconds 900 mW/cm <sup>2</sup> ±10%
WEIGHT	86 gms	84 gms
LENGTH	230 mm	220 mm


C01-M / C02-M CURING LIGHT DEPTH OF CURE							
MODE	MODEL	C01-M (F.O. Light Guide) Curing Time	Depth Measurements	ISO 4049 Depth Measurements (-50%)	C02-M (Direct Light Source) Curing Time	Depth Measurements	ISO 4049 Depth Measurements (-50%)
FULL		10 seconds	4.3 mm	2.15 mm	10 seconds	5.1 mm	2.55 mm
SOFT		15 seconds	4.6 mm	2.30 mm	15 seconds	4.8 mm	2.40 mm

**REMARKS:** Above data is based on test results using an A2 shade composite. The light tip was placed above and pointed perpendicular to the composite. Tests were done as per ISO 4049 Standard. The samples were measured and the data was divided by two ( 50% of the total composite depth cured ). ISO 4049 Standard requires that the results for class II restorative materials should be at least 1.5 mm for non-opaque shades and 1.0 mm for opaque shades. All resins and composites have different curing times. It is recommended that you should test your resins and composites performance to determine the curing times required.



C01-M Curing Light ( F.O. Light Guide )  
 C02-M Curing Light ( Direct Light Head )  
 C01-M-1 F.O. Light Guide 8 mm  
 C02-M-1 Direct Light Head 11 mm  
 C01/2-M-1 Replacement Battery  
 228 Light Sleeves Mini 500 pcs/box

**24 MONTH WARRANTY**



## C01 LED CURING LIGHTS with fibre optic light guides

C01 SPECIFICATIONS			
Model	C01-C (classic)	C01-S (super power)	C01-D (dual range)
F.O. LIGHT GUIDE	8 mm ( dia. )	11 mm ( dia. )	
LIGHT WAVELENGTH	440~480 nm	440~480 nm	390~430 nm 440~480 nm
APPLICATION	General	Ortho / General	General / Ortho
MODE	TACK	N/A	900 mW/cm <sup>2</sup> for 1 second±10%
	TURBO	1,500 mW/cm <sup>2</sup> for 5 seconds ±10%	2,000 mW/cm <sup>2</sup> for 4 seconds±10%
	3' + 3'	N/A	2,000 mW/cm <sup>2</sup> for 3 + 3 seconds±10% (with 1 second blink in between)
	6' + 6'	N/A	1,500 mW/cm <sup>2</sup> for 6 + 6 seconds±10% (with 1 second blink in between)
	FULL	1,200 mW/cm <sup>2</sup> for 5, 10, 15, 20 seconds±10%	
	RAMP	1 <sup>st</sup> five seconds from 0-1,200 mW/cm <sup>2</sup>	next 5, 10, 15, 20 seconds 1,200 mW/cm <sup>2</sup> ±10%
	SOFT CURE	1 <sup>st</sup> five seconds from 0-900 mW/cm <sup>2</sup>	next 5, 10, 15, 20, 25 seconds 900 mW/cm <sup>2</sup> ±10%
WEIGHT	205 gms	200 gms	200 gms
LENGTH	260 mm	265 mm	265 mm

CURING LIGHT MODE APPLICATIONS	
TACK	Tack cure, a very short initial curing for temporary processing before final cure. For example, veneer bonding or temporary cement for in-position and easy removal of excess material.
TURBO	Super fast curing for resin and composite under 2 mm thickness. Also recommended for ceramic brackets, veneer bonding etc.
3' + 3'	3 Seconds Mesial + 3 Seconds Distal for metal brackets resin bonding. (with 1 second blink in between)
6' + 6'	6 Seconds Mesial + 6 Seconds Distal for metal brackets with bonding resin requiring more curing time, and buccal tubes bonding. (with 1 second blink in between)
FULL	For general applications.
RAMP	With first 0-5 seconds ramping for general application.
SOFT CURE	With first 0-5 seconds ramping for more gentle processing, designed for the patients sensitive to the heat.

AUTOCLAVABLE FIBRE OPTIC LIGHT GUIDE CAN BE ROTATED 360°

AUTOMATIC COMPENSATION CIRCUIT FOR CONSTANT POWER OUTPUT

FULLY COVERED MEMBRANE PANEL FOR EASY CLEANING AND INFECTION CONTROL

SIZE: L265 x Dia.27 (mm)

REPLACEABLE BATTERY OPENING

A FULL CHARGE OF BATTERY CAN BE USED MORE THAN 100 TIMES WHEN SET AT 10 SECOND TIME MODE

ON/OFF SWITCH

DISPLAY SCREEN

TIME SWITCH

MODE INDICATORS

MODE SWITCH


C01 Complete Set with Built-In Light Intensity Indicator

063 Protective Hand Light Shield

C01-C Curing Light ( Classic ) with 8 mm Light Guide  
 C01-S Curing Light ( Super Power ) with 11 mm Light Guide  
 C01-D Curing Light ( Dual Range ) with 11 mm Light Guide  
 C01-C-1 8 mm Fibre Optic Light Guide  
 C01-S/D-1 11 mm Fibre Optic Light Guide  
 C01/2-1 Replacement Battery  
 123S Light Sleeves Small ( for C01-C ) 500 pcs/box  
 123L Light Sleeves Large ( for C01-S/C01-D ) 500 pcs/box



**24 MONTH WARRANTY**



C01-C / C01-S / C01-D CURING LIGHT DEPTH OF CURE									
MODEL	C01-C (classic) Curing Time	Depth Measurements	ISO 4049 Depth Measurements (-50%)	C01-S (super power) Curing Time	Depth Measurements	ISO 4049 Depth Measurements (-50%)	C01-D (dual range) Curing Time	Depth Measurements	ISO 4049 Depth Measurements (-50%)
TURBO	5 seconds	4.2 mm	2.10 mm	4 seconds	4.6 mm	2.30 mm	4 seconds	4.3 mm	2.15 mm
FULL	10 seconds	4.3 mm	2.15 mm	10 seconds	5.1 mm	2.55 mm	10 seconds	4.4 mm	2.20 mm
RAMP	10 seconds	4.1 mm	2.05 mm	10 seconds	5.0 mm	2.50 mm	10 seconds	4.2 mm	2.10 mm
SOFT CURE	15 seconds	4.6 mm	2.30 mm	15 seconds	4.8 mm	2.40 mm	15 seconds	4.6 mm	2.30 mm

**REMARKS:** Above data is based on test results using an A2 shade composite. The light guide tip was placed above and pointed perpendicular to the composite. Tests were done as per ISO 4049 Standard. The samples were measured and the data was divided by two (50% of the total composite depth cured). ISO 4049 Standard requires that the results for class II restorative materials should be at least 1.5 mm for non-opaque shades and 1.0 mm for opaque shades. All resins and composites have different curing times. It is strongly recommended that you should test your composite and resin performance to determine the curing times required.

Data is only for reference, always follow the instructions from your resin and composite manufacturers.

- There are many variables (distance and angle between the light tip and composite, area, shade, thickness and type etc.) affecting composite curing times. It is recommended to place and cure a maximum of 2 mm for each composite placement.
- Using a barrier sleeve on the curing light will reduce 5-10% of light intensity.
- The light guide tip should be placed approximately 1-2 mm above the composite or resin and pointed perpendicular to the tooth surface.
- High intensity curing lights produce more heat. This is a direct indication of the power of the light being emitted. When using the device for more than 10 seconds, do not keep the light tip in one position. It is always recommended to use a dental syringe to blow a cold air flow on the area during long periods of curing. This will maintain a lower temperature for the composite and surrounding area, particularly when using rubber dam so as to avoid any rubber dam breakage due to heat built-up.
- Do not expose soft tissue to the light being emitted for more than 2 seconds or burns may occur.
- Avoid the light guide tip coming into contact with any soft tissue in the mouth.

10 colours available:






## C02 LED CURING LIGHTS 90° right angle direct light source head

C02 SPECIFICATIONS			
Model	C02-C (classic)	C02-S (super power)	C02-D (dual range)
DIRECT LIGHT HEAD	10 mm ( dia. )		
LIGHT WAVELENGTH	440~480 nm	440~480 nm	390~430 nm 440~480 nm
APPLICATION	General	Ortho / General	General / Ortho
MODE	TACK	N/A	900 mW/cm <sup>2</sup> for 1 second±10%
	TURBO	1,500 mW/cm <sup>2</sup> for 5 seconds ±10%	2,500 mW/cm <sup>2</sup> for 3 seconds±10%
	3' + 3'	N/A	2,000 mW/cm <sup>2</sup> for 3 + 3 seconds±10% (with 1 second blink in between)
	6' + 6'	N/A	1,500 mW/cm <sup>2</sup> for 6 + 6 seconds±10% (with 1 second blink in between)
	FULL	1,200 mW/cm <sup>2</sup> for 5, 10, 15, 20 seconds±10%	
	RAMP	1 <sup>st</sup> five seconds from 0-1,200 mW/cm <sup>2</sup>	next 5, 10, 15, 20 seconds 1,200 mW/cm <sup>2</sup> ±10%
	SOFT CURE	1 <sup>st</sup> five seconds from 0-900 mW/cm <sup>2</sup>	next 5, 10, 15, 20, 25 seconds 900 mW/cm <sup>2</sup> ±10%
WEIGHT	185 gms	185 gms	185 gms
LENGTH	255 mm	255 mm	255 mm

CURING LIGHT MODE APPLICATIONS	
TACK	Tack cure, a very short initial curing for temporary processing before final cure. For example, veneer bonding or temporary cement for in-position and easy removal of excess material.
TURBO	Super fast curing for resin and composite under 2 mm thickness. Also recommended for ceramic brackets, veneer bonding etc.
3' + 3'	3 Seconds Mesial + 3 Seconds Distal for metal brackets resin bonding. (with 1 second blink in between)
6' + 6'	6 Seconds Mesial + 6 Seconds Distal for metal brackets with bonding resin requiring more curing time, and buccal tubes bonding. (with 1 second blink in between)
FULL	For general applications.
RAMP	With first 0-5 seconds ramping for general application.
SOFT CURE	With first 0-5 seconds ramping for more gentle processing, designed for the patients sensitive to the heat.




C02 Complete Set with Built-In Light Intensity Indicator



- C02-C ( Classic ) Curing Light with 10 mm Light Head
- C02-S ( Super Power ) Curing Light with 10 mm Light Head
- C02-D ( Dual Range ) Curing Light with 10 mm Light Head
- C02-C-1 ( Classic ) Direct Light Source Head
- C02-S-1 ( Super Power ) Direct Light Source Head
- C02-D-1 ( Dual Range ) Direct Light Source Head
- C01/2-1 Replacement Battery
- 125 Light Sleeves ( for C02-C/S/D ) 500 pcs/box
- C02-3 Heat Protect Rings, 6 pcs/pack



**24 MONTH WARRANTY**



C02-C / C02-S / C02-D CURING LIGHT DEPTH OF CURE									
MODEL	C02-C (classic) Curing Time	Depth Measurements	ISO 4049 Depth Measurements (-50%)	C02-S (super power) Curing Time	Depth Measurements	ISO 4049 Depth Measurements (-50%)	C02-D (dual range) Curing Time	Depth Measurements	ISO 4049 Depth Measurements (-50%)
TURBO	5 seconds	4.4 mm	2.20 mm	3 seconds	4.6 mm	2.30 mm	3 seconds	4.5 mm	2.25 mm
FULL	10 seconds	5.1 mm	2.55 mm	10 seconds	5.3 mm	2.65 mm	10 seconds	5.2 mm	2.60 mm
RAMP	10 seconds	4.8 mm	2.40 mm	10 seconds	5.3 mm	2.65 mm	10 seconds	4.9 mm	2.45 mm
SOFT CURE	15 seconds	5.0 mm	2.50 mm	15 seconds	5.4 mm	2.70 mm	15 seconds	5.1 mm	2.55 mm

**REMARKS:** Above data is based on test results using an A2 shade composite. The light tip was placed above and pointed perpendicular to the composite. Tests were done as per ISO 4049 Standard. The samples were measured and the data was divided by two (50% of the total composite depth cured). ISO 4049 Standard requires that the results for class II restorative materials should be at least 1.5 mm for non-opaque shades and 1.0 mm for opaque shades. All resins and composites have different curing times. It is strongly recommended that you should test your composite and resin performance to determine the curing times required.

Data is only for reference, always follow the instructions from your resin and composite manufacturers.

- There are many variables (distance and angle between the light tip and composite, area, shade, thickness and type etc.) affecting composite curing times. It is recommended to place and cure a maximum of 2 mm for each composite placement.
- Using a barrier sleeve on the curing light will reduce 5-10% of light intensity.
- The light tip should be placed approximately 1-2 mm above the composite or resin and pointed perpendicular to the tooth surface.
- High intensity curing lights produce more heat. This is a direct indication of the power of the light being emitted. When using the device for more than 10 seconds, do not keep the light tip in one position. It is always recommended to use a dental syringe to blow a cold air flow on the area during long periods of curing. This will maintain a lower temperature for the composite and surrounding area, particularly when using rubber dam so as to avoid any rubber dam breakage due to heat built-up.
- Do not expose soft tissue to the light being emitted for more than 2 seconds or burns may occur.
- Avoid the light tip coming into contact with any soft tissue in the mouth. As an option, you may cover the entire curing light with a disposable sleeve and then place the Heat Protective Ring (included inside the packing box) onto the light tip covered with sleeve as shown below. Ensure the Heat Protective Ring is securely positioned on the curing light before starting the cure.

10 colours available:





**L10-1**  
Straight  
Handpiece  
1:1 Standard

**L20-1L**  
Latch  
CA Handpiece  
1:1 Standard

**L20-1P**  
Push Button  
CA Handpiece  
1:1 Standard

**L20-1PM**  
Push Button Mini Head  
CA Handpiece  
1:1 Standard

## LOW SPEED HANDPIECES AND AIR MOTOR



**L20-4P**  
Push Button  
CA Handpiece  
4:1 Reduction

**L20-4PM**  
Push Button Mini Head  
CA Handpiece  
4:1 Reduction

**L30**  
Air Motor ( 2 or 4 holes )  
20,000 r.p.m.

**12 MONTH WARRANTY**




**L20-1PY ( SS )**  
Screw-In / Snap-On  
CA Handpiece  
1:1 Standard

**L20-4PY ( SS )**  
Screw-In / Snap-On  
CA Handpiece  
4:1 Reduction

**L20-4PY ( R )**  
60° Turn Reciprocating  
Push Button Standard Head  
CA Handpiece  
4:1 Reduction

**L20-4PY ( RM )**  
60° Turn Reciprocating  
Push Button Mini Head  
CA Handpiece  
4:1 Reduction

## PROPHY LOW SPEED HANDPIECES AND AIR MOTOR



**L10-4PY**  
Prophy Angle  
Straight Handpiece  
4:1 Reduction

**1002-4 Prophy Angle Air Motor**  
Complete with 3 Nose Cones  
4:1 Reduction 3,000 r.p.m.  
( 4 holes )

**12 MONTH WARRANTY**

## C10 LED & HALOGEN CURING LIGHT METER



- Simple, compact design for easy use
- Measures intensity up to 2500 mw/cm<sup>2</sup>
- Size: 110 mm x 68 mm x 33 mm

C10 Light Meter

**24 MONTH WARRANTY**



**A08 Apex Locator**  
A08-1 Battery  
A08-2 Spare Probe  
A08-3 Spare Hook  
A08-4 Spare Cord  
A08-5 Spare Clip

## A08 DENTAL APEX LOCATOR

- Bright LCD screen with the trajectory of the file being clearly indicated by different colours
- Foldable design with adjustable view angle, saves space and protects the screen
- The unit is designed with a smooth seamless body and membrane buttons to enhance cold sterilization and help avoid cross contamination
- Advanced multiple frequency technology with automatic calibrating ensures accurate measurement
- Measurement results are not affected by tooth type, size of root canal or the internal environment of root canal
- Different beep tones indicate the position of the file in the root canal
- Rechargeable lithium battery
- Accessories autoclavable up to 134°C
- Screen size: 110 mm x 70 mm
- Overall size: L122 mm x W118 mm x H25 mm
- Weight: 300 gms ( including battery )

**24 MONTH WARRANTY**

12 colours available:



## HIGH SPEED TURBINE WRENCH STANDARD CONNECTIONS



**T80 WRENCH**


Torque Power: 18 Watts  
Speed: 400,000 r.p.m. ±10%

- Ceramic Bearings
- Stainless steel head
- Ergonomic design
- Quiet operation noise level, less than 65 dB
- Autoclavable up to 134°C
- Midwest 4 holes or Borden 2 holes

*Turbine Wrench Standard Connection*  
T80W Turbine ( 2 or 4 holes )  
T80W-A Turbine Cartridge

**12 MONTH WARRANTY**

## HIGH SPEED TURBINE PUSH BUTTON STANDARD CONNECTIONS



**T90P HIGH TORQUE PUSH BUTTON**

Torque Power: 22 Watts  
Speed: 350,000 r.p.m. ±10%

- Anti-retraction
- Ceramic Bearings
- Stainless steel head
- Ergonomic design
- Quiet operation noise level, less than 65 dB
- Autoclavable up to 134°C
- Midwest 4 holes or Borden 2 holes

*Turbine Push Button Standard Connection*  
T90P High Torque Push Button ( 2 or 4 holes )  
T90P-A High Torque PB Cartridge

**12 MONTH WARRANTY**

## HIGH SPEED TURBINE QUICK DISCONNECT



**T90+QK HIGH TORQUE PUSH BUTTON**

Torque Power: 22 Watts  
Speed: 350,000 r.p.m. ±10%

**NON LED QUICK COUPLING**



QK KaVo<sup>1</sup> Type  
with Spray Adjustment

- Anti-retraction
- Ceramic Bearings
- Stainless steel head
- Quiet operation noise level, less than 65 dB
- Ergonomic design
- Autoclavable up to 134°C
- KaVo<sup>1</sup> type fitting

*NON LED High Torque Turbine Quick Disconnect*  
T90+QK KaVo<sup>1</sup> Type ( without Coupling )  
T90+-A High Torque PB Cartridge

*NON LED Quick Coupling*  
QK KaVo<sup>1</sup> Type ( non LED ) w Spray Adjustment

**24 MONTH WARRANTY**

1. KAVO is a registered trademark of Kavo Dental GmbH Germany

## LED HIGH SPEED TURBINE QUICK DISCONNECT



**T90+LK HIGH TORQUE PUSH BUTTON**

Torque Power: 22 Watts  
Speed: 350,000 r.p.m. ±10%

**LED QUICK COUPLING**



LK KaVo<sup>1</sup> Type  
with Spray Adjustment


- LED Fibre Optic
- Anti-retraction
- Ceramic Bearings
- Stainless steel head
- Quiet operation noise level, less than 65 dB
- Ergonomic design
- Autoclavable up to 134°C
- KaVo<sup>1</sup> type fitting

*LED High Torque Turbine Quick Disconnect*  
T90+LK KaVo<sup>1</sup> Type ( without Coupling )  
T90+-A High Torque PB Cartridge

*LED Quick Coupling*  
LK KaVo<sup>1</sup> Type ( LED ) w Spray Adjustment

**24 MONTH WARRANTY**

1. KAVO is a registered trademark of Kavo Dental GmbH Germany



## DISPOSABLE 60° TURN RECIPROCATING LATEX FREE PROPHY ANGLES

**NO SPLATTER** - The centrifugal force generated by a normal 360° turn prophy angle results in splattering and mess. Our **NEW 60° turn reciprocating** prophy angles eliminates potential splattering of prophy paste, blood and saliva, keeping the prophy paste in the cup where it belongs.

**NO HEAT** - Our **NEW 60° turn reciprocating** prophy angles build up much less heat during use, assuring a smooth, cool treatment for your patient.

**EFFECTIVE** - Due to less heat generated during use, if necessary more pressure may be applied to remove stubborn stains.

**FASTER** - Thanks to our innovatively designed **NEW 60° turn reciprocating** prophy angles, the prophy cup or brush can be used continuously on the tooth surface, saving valuable chair time.

**SAFER** - The rotational motion of traditional prophy cups and brushes can cause laceration and bleeding of the gums. Our **NEW 60° turn reciprocating** prophy angles provide a much safer and pain-free prophylaxis treatment.

**WARNING: USE OF THESE 60° TURN RECIPROCATING PROPHY ANGLES IS EXTREMELY ADDICTIVE**  
*Once you've tried them, you will never want to use any other prophy angles again*

## DISPOSABLE 60° TURN RECIPROCATING PROPHY ANGLES STANDARD CUPS LATEX FREE TURBINE BLADES

**REGULAR STANDARD CUPS TURBINE BLADES**

1,500 - 3,000 r.p.m.

<b>Regular Cup Colours</b>	White	Light Blue	Lavender	Orange
<b>Regular Single Colour Order No.</b>	CT3399-60WT	CT3399-60LB	CT3399-60LV	CT3399-60OR
<b>Regular Assorted Colours Order No.</b>	CT3399-60RA Assorted: White, L.Blue, Lavender, Orange			

**SOFT STANDARD CUPS TURBINE BLADES**

1,500 - 3,000 r.p.m.

<b>Soft Cup Colours</b>	Pink	Lime	Green	Dark Blue
<b>Soft Single Colour Order No.</b>	CT3399-60PK	CT3399-60LM	CT3399-60GR	CT3399-60DB
<b>Soft Assorted Colours Order No.</b>	CT3399-60SA Assorted: Pink, Lime, Green, D.Blue			

**X-SOFT STANDARD CUPS TURBINE BLADES**

1,500 - 3,000 r.p.m.

<b>X-Soft Cup Colours</b>	Red	Yellow
<b>X-Soft Single Colour Order No.</b>	CT3399-60RD	CT3399-60YL
<b>X-Soft Assorted Colours Order No.</b>	CT3399-60XA Assorted: Red, Yellow	

## DISPOSABLE 60° TURN RECIPROCATING PROPHY ANGLES SHORT CUPS LATEX FREE TURBINE BLADES

**REGULAR SHORT CUPS TURBINE BLADES**

1,500 - 3,000 r.p.m.

<b>Regular Cup Colours</b>	White	Light Blue	Lavender	Orange
<b>Regular Single Colour Order No.</b>	CTS3399-60WT	CTS3399-60LB	CTS3399-60LV	CTS3399-60OR
<b>Regular Assorted Colours Order No.</b>	CTS3399-60RA Assorted: White, L.Blue, Lavender, Orange			

**SOFT SHORT CUPS TURBINE BLADES**

1,500 - 3,000 r.p.m.

<b>Soft Cup Colours</b>	Pink	Lime	Green	Dark Blue
<b>Soft Single Colour Order No.</b>	CTS3399-60PK	CTS3399-60LM	CTS3399-60GR	CTS3399-60DB
<b>Soft Assorted Colours Order No.</b>	CTS3399-60SA Assorted: Pink, Lime, Green, D.Blue			

**X-SOFT SHORT CUPS TURBINE BLADES**

1,500 - 3,000 r.p.m.

<b>X-Soft Cup Colours</b>	Red	Yellow
<b>X-Soft Single Colour Order No.</b>	CTS3399-60RD	CTS3399-60YL
<b>X-Soft Assorted Colours Order No.</b>	CTS3399-60XA Assorted: Red, Yellow	

## DISPOSABLE 60° TURN RECIPROCATING PROPHY ANGLES LATEX FREE BRUSHES REGULAR / SOFT / X-SOFT

**REGULAR WHITE (nylon brush dia 0.2 mm)**

1,500 - 3,000 r.p.m.

<b>Regular Type</b>	Flat	Tapered	Mini Cup
<b>Regular White Colour Order No.</b>	C3399-60-R261F	C3399-60-R261T	C3399-60-R269
<b>Regular Assorted Colours Order No.</b>	C3399-60-R261FA C3399-60-R261TA C3399-60-R269A Assorted: White, L.Blue, Lavender, Orange		

**SOFT CLEAR (nylon brush dia 0.15 mm)**

1,500 - 3,000 r.p.m.


<b>Soft Type</b>	Flat	Tapered	Mini Cup
<b>Soft Clear Colour Order No.</b>	C3399-60-S261F	C3399-60-S261T	C3399-60-S269
<b>Soft Assorted Colours Order No.</b>	C3399-60-S261FA C3399-60-S261TA C3399-60-S269A Assorted: Pink, Lime, Green, D.Blue		

**X-SOFT YELLOW (nylon brush dia 0.10 mm)**

1,500 - 3,000 r.p.m.


<b>X-Soft Type</b>	Flat	Tapered	Mini Cup
<b>X-Soft Yellow Colour Order No.</b>	C3399-60-X261F	C3399-60-X261T	C3399-60-X269
<b>X-Soft Assorted Colours Order No.</b>	C3399-60-X261FA C3399-60-X261TA C3399-60-X269A Assorted: Yellow, Red		

## PROPHY ANGLE STRAIGHT HANDPIECE 4:1




L10-4PY Prophy Angle Handpiece Short 4:1 Reduction

## STANDARD STRAIGHT HANDPIECE 1:1



L10-1 Standard Straight Handpiece 1:1

## DISPOSABLE PROPHY ANGLES STANDARD CUPS LATEX FREE RIBBED AND WEBBED




Type	90° Right Angle		105° Contra Angle	
Cup Colours	Regular Lavender	Soft Green	Regular Lavender	Soft Green
Order No.	3399-090RL	3399-090SG	3399-105RL	3399-105SG

## DISPOSABLE PROPHY ANGLES BRUSHES LATEX FREE FLAT AND TAPERED




Type	90° Right Angle		105° Contra Angle	
Brushes	Regular Flat	Regular Tapered	Regular Flat	Regular Tapered
Order No.	3399-090F	3399-090T	3399-105F	3399-105T

## DISPOSABLE CLEAR PROPHY ANGLES STANDARD CUPS LATEX FREE RIBBED AND WEBBED




Type	90° Right Angle		90° Right Angle		Type	105° Contra Angle		105° Contra Angle	
Cup Colours	Regular Light Blue	Regular Lavender	Soft Green	Soft Pink	Cup Colours	Regular Light Blue	Regular Lavender	Soft Green	Soft Pink
Order No.	C3399-090RB	C3399-090RL	C3399-090SG	C3399-090SP	Order No.	C3399-105RB	C3399-105RL	C3399-105SG	C3399-105SP

## DISPOSABLE CLEAR PROPHY ANGLES SHORT CUPS LATEX FREE RIBBED AND WEBBED



Type	90° Right Angle		90° Right Angle		Type	105° Contra Angle		105° Contra Angle	
Cup Colours	Regular Light Blue	Regular Lavender	Soft Green	Soft Pink	Cup Colours	Regular Light Blue	Regular Lavender	Soft Green	Soft Pink
Order No.	CS3399-090RB	CS3399-090RL	CS3399-090SG	CS3399-090SP	Order No.	CS3399-105RB	CS3399-105RL	CS3399-105SG	CS3399-105SP

## DISPOSABLE CLEAR PROPHY ANGLES STANDARD CUPS LATEX FREE TURBINE BLADES



Type	90° Right Angle		90° Right Angle		Type	105° Contra Angle		105° Contra Angle	
Cup Colours	Regular Light Blue	Regular Lavender	Soft Green	Soft Pink	Cup Colours	Regular Light Blue	Regular Lavender	Soft Green	Soft Pink
Order No.	CT3399-090RB	CT3399-090RL	CT3399-090SG	CT3399-090SP	Order No.	CT3399-105RB	CT3399-105RL	CT3399-105SG	CT3399-105SP

# Prophy Cups Standard Shanks



## 60° TURN RECIPROCATING PROPHY HANDPIECE 4:1 REDUCTION PUSH BUTTON STANDARD HEAD

**NO SPLATTER** - The centrifugal force generated by a normal 360° turn prophy handpiece results in splattering and mess. Our **NEW 60° turn reciprocating** prophy handpiece eliminates potential splattering of prophy paste, blood and saliva, keeping the prophy paste in the cup where it belongs.

**NO HEAT** - Our **NEW 60° turn reciprocating** prophy handpieces build up much less heat during use, assuring a smooth, cool treatment for your patient.

**EFFECTIVE** - Due to less heat generated during use, if necessary more pressure may be applied to remove stubborn stains.

**FASTER** - Thanks to our innovatively designed **NEW 60° turn reciprocating** prophy handpiece, the prophy cup or brush can be used continuously on the tooth surface, saving valuable chair time.

**SAFER** - The rotational motion of traditional prophy cups and brushes can cause laceration and bleeding of the gums. Our **NEW 60° turn reciprocating** prophy handpieces provide a much safer and pain-free prophylaxis treatment.

### METAL SHANKS CUPS (8.5 mm)



Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	M222WT	M222LB	M222LV	M222OR
Regular Assorted Colours Order No.	M222RA Assorted: White, L.Blue, Lavender, Orange			

### METAL SHANKS CUPS (8.5 mm)



Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	M222PK	M222LM	M222GR	M222DB
Soft Assorted Colours Order No.	M222SA Assorted: Pink, Lime, Green, D.Blue			

### METAL SHANKS CUPS (8.5 mm)



X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	M222RD	M222YL
X-Soft Assorted Colours Order No.	M222XA Assorted: Red, Yellow	

**DISPOSABLE LATEX FREE STANDARD PROPHY CUPS CA METAL SHANKS RIBBED AND WEBBED**

### PLASTIC SHANKS CUPS (8.5 mm)



Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	P222WT	P222LB	P222LV	P222OR
Regular Assorted Colours Order No.	P222RA Assorted: White, L.Blue, Lavender, Orange			

### PLASTIC SHANKS CUPS (8.5 mm)



Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	P222PK	P222LM	P222GR	P222DB
Soft Assorted Colours Order No.	P222SA Assorted: Pink, Lime, Green, D.Blue			

### PLASTIC SHANKS CUPS (8.5 mm)



X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	P222RD	P222YL
X-Soft Assorted Colours Order No.	P222XA Assorted: Red, Yellow	

**DISPOSABLE LATEX FREE STANDARD PROPHY CUPS CA PLASTIC SHANKS RIBBED AND WEBBED**

### METAL SHANKS SHORT CUPS (7.0 mm)



Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	S222WT	S222LB	S222LV	S222OR
Regular Assorted Colours Order No.	S222RA Assorted: White, L.Blue, Lavender, Orange			

### METAL SHANKS SHORT CUPS (7.0 mm)



Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	S222PK	S222LM	S222GR	S222DB
Soft Assorted Colours Order No.	S222SA Assorted: Pink, Lime, Green, D.Blue			

### METAL SHANKS SHORT CUPS (7.0 mm)



X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	S222RD	S222YL
X-Soft Assorted Colours Order No.	S222XA Assorted: Red, Yellow	

**DISPOSABLE LATEX FREE SHORT PROPHY CUPS CA METAL SHANKS RIBBED AND WEBBED**

## DISPOSABLE LATEX FREE PROPHY BRUSHES CA REGULAR / SOFT / X-SOFT



Regular Type	Flat	Tapered	Micro	Junior Cup	Mini Cup
<b>Regular White Colour Order No.</b>	R-262F	R-262T	R-265	R-M268N	R-269N
<b>Regular Assorted Colours Order No.</b>	Assorted: White, L.Blue, Lavender, Orange				



Soft Type	Flat	Tapered	Micro	Junior Cup	Mini Cup
<b>Soft Clear Colour Order No.</b>	S-262F	S-262T	S-265	S-268	S-269
<b>Soft Assorted Colours Order No.</b>	Assorted: Pink, Lime, Green, D.Blue				



X-Soft Type	Flat	Tapered	Micro	Junior Cup	Mini Cup
<b>X-Soft Yellow Colour Order No.</b>	X-262F	X-262T	X-265	X-268	X-269
<b>X-Soft Assorted Colours Order No.</b>	Assorted: Red, Yellow				



Regular Cup Colours	White	Light Blue	Lavender	Orange
<b>Regular Single Colour Order No.</b>	225WT	225LB	225LV	225OR
<b>Regular Assorted Colours Order No.</b>	225RA Assorted: White, L.Blue, Lavender, Orange			



Soft Cup Colours	Pink	Lime	Green	Dark Blue
<b>Soft Single Colour Order No.</b>	225PK	225LM	225GR	225DB
<b>Soft Assorted Colours Order No.</b>	225SA Assorted: Pink, Lime, Green, D.Blue			



X-Soft Cup Colours	Red	Yellow
<b>X-Soft Single Colour Order No.</b>	225RD	225YL
<b>X-Soft Assorted Colours Order No.</b>	225XA Assorted: Red, Yellow	

**DISPOSABLE LATEX FREE PROPHY POINTS CA METAL SHANKS**  
[ for lingual and interproximal ]



Regular Cup Colours	White	Light Blue	Lavender	Orange
<b>Regular Single Colour Order No.</b>	T-M222 WT	T-M222 LB	T-M222 LV	T-M222 OR
<b>Regular Assorted Colours Order No.</b>	T-M222RA Assorted: White, L.Blue, Lavender, Orange			



Soft Cup Colours	Pink	Lime	Green	Dark Blue
<b>Soft Single Colour Order No.</b>	T-M222 PK	T-M222 LM	T-M222 GR	T-M222 DB
<b>Soft Assorted Colours Order No.</b>	T-M222SA Assorted: Pink, Lime, Green, D.Blue			



X-Soft Cup Colours	Red	Yellow
<b>X-Soft Single Colour Order No.</b>	T-M222 RD	T-M222 YL
<b>X-Soft Assorted Colours Order No.</b>	T-M222XA Assorted: Red, Yellow	

**DISPOSABLE LATEX FREE STANDARD PROPHY CUPS CA METAL SHANKS TURBINE BLADES**



Regular Cup Colours	White	Light Blue	Lavender	Orange
<b>Regular Single Colour Order No.</b>	T-P222 WT	T-P222 LB	T-P222 LV	T-P222 OR
<b>Regular Assorted Colours Order No.</b>	T-P222RA Assorted: White, L.Blue, Lavender, Orange			



Soft Cup Colours	Pink	Lime	Green	Dark Blue
<b>Soft Single Colour Order No.</b>	T-P222 PK	T-P222 LM	T-P222 GR	T-P222 DB
<b>Soft Assorted Colours Order No.</b>	T-P222SA Assorted: Pink, Lime, Green, D.Blue			



X-Soft Cup Colours	Red	Yellow
<b>X-Soft Single Colour Order No.</b>	T-P222 RD	T-P222 YL
<b>X-Soft Assorted Colours Order No.</b>	T-P222XA Assorted: Red, Yellow	

**DISPOSABLE LATEX FREE STANDARD PROPHY CUPS CA PLASTIC SHANKS TURBINE BLADES**



Regular Cup Colours	White	Light Blue	Lavender	Orange
<b>Regular Single Colour Order No.</b>	T-S222 WT	T-S222 LB	T-S222 LV	T-S222 OR
<b>Regular Assorted Colours Order No.</b>	T-S222RA Assorted: White, L.Blue, Lavender, Orange			



Soft Cup Colours	Pink	Lime	Green	Dark Blue
<b>Soft Single Colour Order No.</b>	T-S222 PK	T-S222 LM	T-S222 GR	T-S222 DB
<b>Soft Assorted Colours Order No.</b>	T-S222SA Assorted: Pink, Lime, Green, D.Blue			



X-Soft Cup Colours	Red	Yellow
<b>X-Soft Single Colour Order No.</b>	T-S222 RD	T-S222 YL
<b>X-Soft Assorted Colours Order No.</b>	T-S222XA Assorted: Red, Yellow	

**DISPOSABLE LATEX FREE SHORT PROPHY CUPS CA METAL SHANKS TURBINE BLADES**



## PROPHY HANDPIECE SNAP-ON / SCREW-IN 4:1 REDUCTION

# SNAP-ON SCREW-IN



L20-4PY (SS)  
Prophy Handpiece CA  
Snap On / Screw In  
4:1 Reduction

### DISPOSABLE LATEX FREE PROPHY POINTS SCREW-IN [ for lingual and interproximal ]

1,500 - 3,000 r.p.m.

**REGULAR POINTS SCREW-IN**

Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	226WT	226LB	226LV	226OR
Regular Assorted Colours Order No.	226RA Assorted: White, L.Blue, Lavender, Orange			

1,500 - 3,000 r.p.m.

**SOFT POINTS SCREW-IN**

Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	226PK	226LM	226GR	226DB
Soft Assorted Colours Order No.	226SA Assorted: Pink, Lime, Green, D.Blue			

1,500 - 3,000 r.p.m.

**X-SOFT POINTS SCREW-IN**

X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	226RD	226YL
X-Soft Assorted Colours Order No.	226XA Assorted: Red, Yellow	

### DISPOSABLE LATEX FREE PROPHY CUPS SNAP-ON RIBBED AND WEBBED

1,500 - 3,000 r.p.m.

**REGULAR CUPS RIBBED AND WEBBED**

SNAP-ON CUPS (8.0 mm)

Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	220WT	220LB	220LV	220OR
Regular Assorted Colours Order No.	220RA Assorted: White, L.Blue, Lavender, Orange			

1,500 - 3,000 r.p.m.

**SOFT CUPS RIBBED AND WEBBED**

SNAP-ON CUPS (8.0 mm)

Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	220PK	220LM	220GR	220DB
Soft Assorted Colours Order No.	220SA Assorted: Pink, Lime, Green, D.Blue			

1,500 - 3,000 r.p.m.

**X-SOFT CUPS RIBBED AND WEBBED**

SNAP-ON CUPS (8.0 mm)

X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	220RD	220YL
X-Soft Assorted Colours Order No.	220XA Assorted: Red, Yellow	

### DISPOSABLE LATEX FREE STANDARD PROPHY CUPS SCREW-IN RIBBED AND WEBBED

1,500 - 3,000 r.p.m.

**REGULAR STANDARD CUPS RIBBED AND WEBBED**

SCREW-IN CUPS (8.5 mm)

Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	221WT	221LB	221LV	221OR
Regular Assorted Colours Order No.	221RA Assorted: White, L.Blue, Lavender, Orange			

1,500 - 3,000 r.p.m.

**SOFT STANDARD CUPS RIBBED AND WEBBED**

SCREW-IN CUPS (8.5 mm)

Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	221PK	221LM	221GR	221DB
Soft Assorted Colours Order No.	221SA Assorted: Pink, Lime, Green, D.Blue			

1,500 - 3,000 r.p.m.

**X-SOFT STANDARD CUPS RIBBED AND WEBBED**

SCREW-IN CUPS (8.5 mm)

X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	221RD	221YL
X-Soft Assorted Colours Order No.	221XA Assorted: Red, Yellow	

### DISPOSABLE LATEX FREE SHORT PROPHY CUPS SCREW-IN RIBBED AND WEBBED

1,500 - 3,000 r.p.m.

**REGULAR SHORT CUPS RIBBED AND WEBBED**

SCREW-IN SHORT CUPS (7.0 mm)

Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	S221WT	S221LB	S221LV	S221OR
Regular Assorted Colours Order No.	S221RA Assorted: White, L.Blue, Lavender, Orange			

1,500 - 3,000 r.p.m.

**SOFT SHORT CUPS RIBBED AND WEBBED**

SCREW-IN SHORT CUPS (7.0 mm)

Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	S221PK	S221LM	S221GR	S221DB
Soft Assorted Colours Order No.	S221SA Assorted: Pink, Lime, Green, D.Blue			

1,500 - 3,000 r.p.m.

**X-SOFT SHORT CUPS RIBBED AND WEBBED**

SCREW-IN SHORT CUPS (7.0 mm)

X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	S221RD	S221YL
X-Soft Assorted Colours Order No.	S221XA Assorted: Red, Yellow	

# Prophy Cups and Brushes Snap-On / Screw-In



**REGULAR WHITE**  
(nylon brush dia. 0.20 mm)

Regular Type	Flat	Tapered
Regular White Brushes Order No.	R-260F	R-260T



**SOFT CLEAR**  
(nylon brush dia. 0.15 mm)

Soft Type	Flat	Tapered
Soft Clear Brushes Order No.	S-260F	S-260T



**X-SOFT YELLOW**  
(nylon brush dia. 0.10 mm)

X-Soft Type	Flat	Tapered
X-Soft Yellow Brushes Order No.	X-260F	X-260T

## DISPOSABLE LATEX FREE PROPHY BRUSHES SNAP-ON REGULAR / SOFT / X-SOFT

## DISPOSABLE LATEX FREE PROPHY BRUSHES SCREW-IN REGULAR / SOFT / X-SOFT



**REGULAR WHITE** (nylon brush dia. 0.20 mm)

Regular Type	Flat	Tapered	Micro	Junior Cup	Mini Cup
Regular White Colour Order No.	R-261F	R-261T	R-265-1	R-268-1	R-269-1
Regular Assorted Colours Order No.	Assorted: White, L.Blue, Lavender, Orange				



**SOFT CLEAR** (nylon brush dia. 0.15 mm)

Soft Type	Flat	Tapered	Micro	Junior Cup	Mini Cup
Soft Clear Colour Order No.	S-261F	S-261T	S-265-1	S-268-1	S-269-1
Soft Assorted Colours Order No.	Assorted: Pink, Lime, Green, D.Blue				



**X-SOFT YELLOW** (nylon brush dia. 0.10 mm)

X-Soft Type	Flat	Tapered	Micro	Junior Cup	Mini Cup
X-Soft Yellow Colour Order No.	X-261F	X-261T	X-265-1	X-268-1	X-269-1
X-Soft Assorted Colours Order No.	Assorted: Red, Yellow				



**SNAP-ON CUPS (8.0 mm)**



**REGULAR CUPS TURBINE BLADES**

Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	T-220WT	T-220LB	T-220LV	T-220OR
Regular Assorted Colours Order No.	T-220RA Assorted: White, L.Blue, Lavender, Orange			



**SNAP-ON CUPS (8.0 mm)**



**SOFT CUPS TURBINE BLADES**

Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	T-220PK	T-220LM	T-220GR	T-220DB
Soft Assorted Colours Order No.	T-220SA Assorted: Pink, Lime, Green, D.Blue			



**SNAP-ON CUPS (8.0 mm)**



**X-SOFT CUPS TURBINE BLADES**

X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	T-220RD	T-220YL
X-Soft Assorted Colours Order No.	T-220XA Assorted: Red, Yellow	

## DISPOSABLE LATEX FREE PROPHY CUPS SNAP-ON TURBINE BLADES



**SCREW-IN CUPS (8.5 mm)**



**REGULAR STANDARD CUPS TURBINE BLADES**

Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	T-221WT	T-221LB	T-221LV	T-221OR
Regular Assorted Colours Order No.	T-221RA Assorted: White, L.Blue, Lavender, Orange			



**SCREW-IN CUPS (8.5 mm)**



**SOFT STANDARD CUPS TURBINE BLADES**

Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	T-221PK	T-221LM	T-221GR	T-221DB
Soft Assorted Colours Order No.	T-221SA Assorted: Pink, Lime, Green, D.Blue			



**SCREW-IN CUPS (8.5 mm)**



**X-SOFT STANDARD CUPS TURBINE BLADES**

X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	T-221RD	T-221YL
X-Soft Assorted Colours Order No.	T-221XA Assorted: Red, Yellow	

## DISPOSABLE LATEX FREE STANDARD PROPHY CUPS SCREW-IN TURBINE BLADES



**SCREW-IN SHORT CUPS (7.0 mm)**



**REGULAR SHORT CUPS TURBINE BLADES**

Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	T-S221WT	T-S221LB	T-S221LV	T-S221OR
Regular Assorted Colours Order No.	T-S221RA Assorted: White, L.Blue, Lavender, Orange			



**SCREW-IN SHORT CUPS (7.0 mm)**



**SOFT SHORT CUPS TURBINE BLADES**

Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	T-S221PK	T-S221LM	T-S221GR	T-S221DB
Soft Assorted Colours Order No.	T-S221SA Assorted: Pink, Lime, Green, D.Blue			



**SCREW-IN SHORT CUPS (7.0 mm)**



**X-SOFT SHORT CUPS TURBINE BLADES**

X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	T-S221RD	T-S221YL
X-Soft Assorted Colours Order No.	T-S221XA Assorted: Red, Yellow	

## DISPOSABLE LATEX FREE SHORT PROPHY CUPS SCREW-IN TURBINE BLADES

# Prophy Cups and Brushes Short Shanks



## 60° TURN RECIPROCATING PROPHY HANDPIECE 4:1 REDUCTION PUSH BUTTON MINI HEAD

**NO SPLATTER** - The centrifugal force generated by a normal 360° turn prophy handpiece results in splattering and mess. Our NEW 60° turn reciprocating prophy handpiece eliminates potential splattering of prophy paste, blood and saliva, keeping the prophy paste in the cup where it belongs.

**NO HEAT** - Our NEW 60° turn reciprocating prophy handpieces build up much less heat during use, assuring a smooth, cool treatment for your patient.

**EFFECTIVE** - Due to less heat generated during use, if necessary more pressure may be applied to remove stubborn stains.

**FASTER** - Thanks to our innovatively designed NEW 60° turn reciprocating prophy handpiece, the prophy cup or brush can be used continuously on the tooth surface, saving valuable chair time.

**SAFER** - The rotational motion of traditional prophy cups and brushes can cause laceration and bleeding of the gums. Our NEW 60° turn reciprocating prophy handpieces provide a much safer and pain-free prophylaxis treatment.

L20-4PY (RM)  
60° Turn Reciprocating  
Push Button  
Mini Head  
Prophy Handpiece  
4:1 Reduction

**RIBBED AND WEBBED**  
1,500 - 3,000 r.p.m.

SHORT CUPS (7.0 mm)

REGULAR SHORT CUPS  
SHORT SHANKS  
RIBBED AND WEBBED

Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	SS-222WT	SS-222LB	SS-222LV	SS-222OR
Regular Assorted Colours Order No.	SS-222RA Assorted: White, L.Blue, Lavender, Orange			

**RIBBED AND WEBBED**  
1,500 - 3,000 r.p.m.

SHORT CUPS (7.0 mm)

SOFT SHORT CUPS  
SHORT SHANKS  
RIBBED AND WEBBED

Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	SS-222PK	SS-222LM	SS-222GR	SS-222DB
Soft Assorted Colours Order No.	SS-222SA Assorted: Pink, Lime, Green, D.Blue			

**RIBBED AND WEBBED**  
1,500 - 3,000 r.p.m.

SHORT CUPS (7.0 mm)

X-SOFT SHORT CUPS  
SHORT SHANKS  
RIBBED AND WEBBED

X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	SS-222RD	SS-222YL
X-Soft Assorted Colours Order No.	SS-222XA Assorted: Red, Yellow	

**DISPOSABLE LATEX FREE PROPHY CUPS CA METAL SHANKS SHORT SHANKS RIBBED AND WEBBED**  
[ fit mini head CA handpieces only ]

**TURBINE BLADES**  
1,500 - 3,000 r.p.m.

SHORT CUPS (7.0 mm)

REGULAR SHORT CUPS  
SHORT SHANKS  
TURBINE BLADES

Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	SS-T222WT	SS-T222LB	SS-T222LV	SS-T222OR
Regular Assorted Colours Order No.	SS-T222RA Assorted: White, L.Blue, Lavender, Orange			

**TURBINE BLADES**  
1,500 - 3,000 r.p.m.

SHORT CUPS (7.0 mm)

SOFT SHORT CUPS  
SHORT SHANKS  
TURBINE BLADES

Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	SS-T222PK	SS-T222LM	SS-T222GR	SS-T222DB
Soft Assorted Colours Order No.	SS-T222SA Assorted: Pink, Lime, Green, D.Blue			

**TURBINE BLADES**  
1,500 - 3,000 r.p.m.

SHORT CUPS (7.0 mm)

X-SOFT SHORT CUPS  
SHORT SHANKS  
TURBINE BLADES

X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	SS-T222RD	SS-T222YL
X-Soft Assorted Colours Order No.	SS-T222XA Assorted: Red, Yellow	

**DISPOSABLE LATEX FREE PROPHY CUPS CA METAL SHANKS SHORT SHANKS TURBINE BLADES**  
[ fit mini head CA handpieces only ]

## DISPOSABLE LATEX FREE SHORT SHANKS PROPHY BRUSHES CA [ fit mini head CA handpieces only ]

1,500 - 3,000 r.p.m.

REGULAR WHITE ( nylon brush dia. 0.20 mm )

Regular Type	Flat	Tapered	Micro	Junior Cup	Mini Cup
Regular White Colour Order No.	SS-R262F	SS-R262T	SS-R265	SS-R268	SS-R269
Regular Assorted Colours Order No.	SS-R262FA SS-R262TA SS-R265A SS-R268A SS-R269A Assorted: White, L.Blue, Lavender, Orange				

1,500 - 3,000 r.p.m.

SOFT CLEAR ( nylon brush dia. 0.15 mm )

Soft Type	Flat	Tapered	Micro	Junior Cup	Mini Cup
Soft Clear Colour Order No.	SS-S262F	SS-S262T	SS-S265	SS-S268	SS-S269
Soft Assorted Colours Order No.	SS-S262FA SS-S262TA SS-S265A SS-S268A SS-S269A Assorted: Pink, Lime, Green, D.Blue				

1,500 - 3,000 r.p.m.

X-SOFT YELLOW ( nylon brush dia. 0.10 mm )


X-Soft Type	Flat	Tapered	Micro	Junior Cup	Mini Cup
X-Soft Yellow Colour Order No.	SS-X262F	SS-X262T	SS-X265	SS-X268	SS-X269
X-Soft Assorted Colours Order No.	SS-X261FA SS-X261TA SS-X265A SS-X268A SS-X269A Assorted: Red, Yellow				

## DISPOSABLE AIR-WATER SYRINGE TIPS (SHORT) EURO TYPE



- Directly interchangeable with American type short syringe tips (syringes which have been converted to use American type short syringe tips)
- No sharp edges for patient comfort
- Unique quad air core design to ensure dry air

Syringe Tips Short, 250 pcs/pack  
028S Syringe Tips Short Spectrum  
028W Syringe Tips Short White



## DISPOSABLE AIR-WATER 'e' SYRINGE TIPS WITHOUT CORE

ready to use, no converter required for American style syringe

- Directly interchangeable with most American style metal air/water syringe tips
- Economically designed without separate water core, in exchange for a split second of moisture when switching from water to air
- Rounded tip end provides maximum patient comfort

Disposable 'e' Air-Water Syringe Tips without core 250 pcs/pack  
0120S Opaque Colour Tips  
012W Opaque White Tips




## DISPOSABLE AIR-WATER 'a' SYRINGE TIPS

ready to use, no converter required for American style syringe

- Directly interchangeable with most American style metal air/water syringe tips
- No sharp edges for patient comfort
- Unique quad air core design to ensure dry air

Disposable 'a' Air-Water Syringe Tips, 250 pcs/pack  
018S Spectrum Colour Tips/White Cores  
018C Clear Tips/Colour Cores  
018OS Opaque Colour Tips/White Cores  
018W Opaque White Tips/White Cores



## DISPOSABLE TRANSPARENT SYRINGE TIPS

- Directly interchangeable with \*Sani-Tip® type transparent tip (syringes which have been converted to use transparent tip)
- Costs less than sterilizing metal tips
- Size: Long( American type ) Short( European type )
- Colour: Clear Transparent


Clear Tips, 250 pcs/pack  
011A Standard Length ( American Type Syringe )  
011E Short Length ( European Type Syringe )

\*Sani-Tip is a registered trademark of Dentsply Sirona

## AUTOCLAVABLE METAL AIR-WATER SYRINGE TIP



- Designed with a positive locking feature between the inner and outer tubes for increased safety
- Improved spray pattern
- Special alloy is employed to ensure a virtually flawless finish
- Compatible with all American type standard quick change connectors



## DISPOSABLE SYRINGE TIP SHIELDS

- Prevent spattering of saliva and blood
- Increases the effectiveness of all rinsing and dry
- Fully disposable

015 Syringe Tip Shields  
100 pcs/pack



Regular  
Blue and  
Light Blue



Fine  
Green and  
Light Green



Ultrafine  
Lavender and  
Light Lavender



Cylinder  
Black and  
White

## DISPOSABLE MICRO APPLICATORS

- Economical and disposable, saving up to 80% by the reduction of excess materials used
- Bends into any angle for the placement of small controlled amounts of material in limited access areas




**TUBE PACK**

Available in :-  
Standard Handles  
100 mm x dia. 3 mm  
Mini Handles  
85 mm x dia. 2.5 mm

Standard Handles  
Tube Pack 400 pieces  
901T Regular  
902T Fine  
903T Ultrafine  
904T Cylinder

Mini Handles  
Tube Pack 400 pieces  
9901T Regular  
9902T Fine  
9903T Ultrafine



**BULK PACK**

Available in :-  
Standard Handles  
100 mm x dia. 3 mm  
Mini Handles  
85 mm x dia. 2.5 mm

Standard Handles  
Bulk Pack 1000 pieces  
901B Regular  
902B Fine  
903B Ultrafine  
904B Cylinder


Mini Handles  
Bulk Pack 1000 pieces  
9901B Regular  
9902B Fine  
9903B Ultrafine



**DISPENSER PACK**

Available in :-  
Standard Handles  
100 mm x dia. 3 mm

Standard Handles  
Dispenser Pack 400 pieces  
901D Regular  
902D Fine  
903D Ultrafine  
904D Cylinder



18-900  
STANDARD


18-9900  
MINI

## DISPOSABLE MICRO APPLICATOR DISPENSERS

- Small, sturdy, durable construction
- Provides aseptic barrier against operatory aerosols and accidental splatter

Applicator Dispenser  
18-900 Standard Handle  
18-9900 Mini Handle

## DISPOSABLE MICRO APPLICATOR TIPS



REGULAR

FINE

921  
Regular


922  
Fine

- Economical and disposable saving up to 80% by the reduction of excess materials used
- For the placement of small controlled amounts of material in limited access areas
- Can be used for etchants, bonding agents, cavity lining, fluoride varnishing, hemostatic and sealers etc

Micro Applicator Tips, 100 pcs/box  
(Blue, Green, Lavender or Yellow)  
921 Regular  
922 Fine

Brush Handle Autoclavable  
7502W - White  
7502T - Transparent

## DISPOSABLE BRUSH APPLICATORS



REGULAR

BENDABLE

MINI

- Perfect for placement of etchants, bonding agents, sealers and cements etc
- Bendable to reach around corners and the tightest spaces
- Handles are colour coded to eliminate confusion between materials
- Fully disposable, eliminates cross-contamination
- Available colours: Blue, Green, Lavender or Yellow

Brush Applicator Single Colour, 100 pcs/pack  
7500 Regular  
7510 Mini

## DISPOSABLE BRUSH TIPS AND HANDLES



REGULAR

MINI

- Perfect for placement of etchants, bonding agents, sealers and cements etc
- Available Colours: Blue, Green, Lavender, or Yellow

7501 Brush Tips Regular Standard Pack, 100 pcs/box  
7511 Brush Tips Mini Standard Pack, 100 pcs/box  
7502W Brush Handle Autoclavable White  
7502T Brush Handle Autoclavable Clear

## COMPOSITE TIP GUN AND TIPS



7600


High-viscosity

7605

7605 Syringe Gun  
7600 Syringe Tips, 100pcs/pack

- For placement of composites, cements etc
- Reusable handle with disposable tubes and plugs
- Autoclavable plastic handle with pistol grip design
- Tips are photoblock

## AUTOCLAVABLE HVE SUCTION TUBES




- Fits European style 16 mm or standard 11 mm suction tubing
- Autoclavable up to 135°C
- Special high quality plastic, sturdy enough for cheek retraction
- Fully polished internal surface and edge for quiet operation and patient comfort
- Assortment of bright, lively colours: Lavender, Light Blue, Light Green, Lime and Yellow

Assorted Colours, 10 pcs/pack  
 335L Adult Large (16 mm)  
 335M Adult Medium (16 mm)  
 336P Pedo (16 mm)  
 005L Adult Large (16 mm)  
 005M Adult Medium (16 mm)  
 006P Pedo (16 mm)  
 205M Adult Medium (11 mm)  
 206P Pedo (11 mm)

### VENTED

### SOFT TOUCH™



ST001 WHITE

## DISPOSABLE SOFT TOUCH™ ORAL EVACUATORS (SLOTTED) VENTED LATEX FREE

- Soft Touch edge eliminates the risk of lacerating the oral tissue
- Rigid tubing strong enough for cheek retraction
- Vented to reduce unnecessary tissue grab and increase evacuation cover area
- Thoroughly polished surface for smooth and quiet operation
- Fits standard 11 mm suction holders
- Colour: White (Blue or Lavender Tip)

PATENT PENDING

ST001 Soft Touch™ Oral Evacuators, 50 pcs/pack

## VENTED (SLOTTED) ORAL EVACUATORS



- Vented to reduce unnecessary tissue grab
- Soft rigid tubing, allows for greater patient comfort yet is strong enough for cheek retraction
- Thoroughly polished surface for smooth and quiet operation
- Fits standard 11 mm suction holders

Oral Evacuators, 100 pcs/pack  
 001 White Disposable  
 001A Single Colour Disposable  
 003 Grey Autoclavable



## VENTED/NON-VENTED ORAL EVACUATORS 'S' CUT

- Vented and non-vented combination
- Soft rigid tubing, allows for greater patient comfort yet is strong enough for cheek retraction
- Thoroughly polished surface for smooth and quiet operation
- Fits standard 11 mm suction holders

Oral Evacuators, 100 pcs/pack  
 301 White Disposable  
 301A Single Colour Disposable  
 303 Grey Autoclavable



## VENTED ORAL EVACUATORS 'S' CUT

- Vented to reduce unnecessary tissue grab
- Soft rigid tubing, allows for greater patient comfort yet is strong enough for cheek retraction
- Thoroughly polished surface for smooth and quiet operation
- Fits standard 11 mm suction holders
- Colours: Blue, Yellow, Green, Lime or Lavender

Oral Evacuators, 100 pcs/pack  
 308 White Disposable  
 308A Single Colour Disposable  
 309 Grey Autoclavable



## AUTOCLAVABLE SURGICAL ASPIRATOR TIPS

- Non-sterile
- Autoclavable up to 135°C
- Colours: Standard - White Large orifice - Blue



Autoclavable Surgical Tips, 25 pcs/pack  
 007 Mini Tip  
 009 Standard Non-Vented  
 009V Standard Vented  
 008 Large Orifice Non-Vented  
 008V Large Orifice Vented

## DISPOSABLE SALIVA EJECTORS



- Tip design that eliminates suction of tissue
- Smooth fixed tip for patient safety and comfort

Saliva Ejectors, 100 pcs/pack  
 300W White Tip/Clear Tube  
 300B Blue Tip/Clear Tube  
 300A White Tip/Coloured Tube


## DISPOSABLE SALIVA EJECTORS SOFT TIPS




- Flexible soft cushion tip for patient comfort
- Latex free tip design that eliminates suction of soft tissue
- Fixed tip for safety

Saliva Ejectors, 100 pcs/pack  
 333W White Tip/Coloured Tube  
 333B Blue Tip/Clear Tube  
 333A White Tip/Coloured Tube

## DISPOSABLE IMPRESSION TRAYS WITH RIM LOCK for Implant Post and General Impressions




base sections can be easily removed with a bar for implant posts



rim lock reduced thickness trim lines on labial wall for easy removal in inclination post cases

Disposable Impression Trays, 10 pcs/pack  
410 Intro (one size each)  
411-420 Single pack

- Save time and expensive laboratory charges
- Strong, rigid distortion free plastic construction
- Anatomically designed
- Smooth edges for patient comfort
- Excellent retention for impression materials whilst providing sufficient pressure for all details
- Colour: Blue



## AUTOCLAVABLE IMPRESSION TRAYS WITH RIM LOCK

- Autoclavable up to 135°C (275°F)
- Strong, rigid distortion free plastic construction
- Anatomically designed
- Smooth edges for patient comfort
- Excellent retention for impression materials whilst providing sufficient pressure for all details
- Colour: Light Grey

Autoclavable Impression Trays  
430 Intro 10 pcs/pack (one size each)  
431-440 Single 5 pcs/pack



## DISPOSABLE 3-IN-1 TRAYS

- Wide occlusal plane
- Strong and tough plastic
- Specially designed walls for extra retention
- Tear resistant quality mesh
- Colour: Blue

3-in-1 Trays  
18-N401 Posterior 36 pcs/box  
18-N402 Quadrant 36 pcs/box  
18-N403 Sideless 36 pcs/box  
18-N404 Anterior 36 pcs/pack  
18-N405 Full Arch 24 pcs/pack



## DISPOSABLE 3-IN-1 TRAYS Stainless Steel Wire Connectors

- Stainless steel wire connector for greater access in limited space
- Wide occlusal plane
- Special rigid walls designed for extra retention
- Tear resistant quality mesh
- Colour: Blue

3-in-1 Trays  
18-S401 Posterior 36 pcs/box  
18-S402 Quadrant 36 pcs/box  
18-S403 Sideless 36 pcs/box



## DISPOSABLE 3-IN-1 TRAYS Rigid Bar Connectors

- Rigid Aluminium Bar Connector
- Strong and tough
- Wide occlusal plane
- Special rigid walls designed for extra retention
- Tear resistant quality mesh
- Colour: Blue

3-in-1 Trays  
18-T401 Posterior 36 pcs/box  
18-T402 Quadrant 36 pcs/box  
18-T403 Sideless 36 pcs/box

## DISPOSABLE PLASTIC COATED METAL PLIABLE IMPRESSION TRAYS



- Fully disposable
- Pliable metal


451-460 Pliable Metal Tray



## AUTOCLAVABLE ALGINATE/PLASTER CURVED SPATULA

- For alginate or plaster
- Light weight, 58 gms only
- Fully autoclavable up to 135°C
- Stainless steel blade
- Handle colour: Blue

MS455 Curved Spatula



## DISPOSABLE IMPRESSION SYRINGES

- Pre-bent for easy application
- Fully disposable
- 5 ml capacity


089 Impression Syringes 50 pcs/pack

## SOLEBITE VINYL POLYSILOXANE BITE REGISTRATION MATERIAL

- Easy recognition of occlusal contacts
- Extremely low resistance to bite
- Long working (20 sec), short setting time in mouth (30 sec)
- Right hardness for positioning and not easy to break
- 0.1% or lower shrinkage

### Technical Data Hardening process image

Hardness



0 for reference only

Colour	Light Blue
Base/Catalyst mixing ratio	1:1
Working time	20 sec.
Setting time in mouth	30 sec.
Dimensional variation after 24hr	-0.1% or less
Hardness	90 Shore A



PJ8110  
Solebite Vinyl Polysiloxane  
Bite Registration Material  
2 x 50 ml Cartridges +  
3 x Mixing Tips Green



### 5 special features of the products

- **Just the right colour to see the occlusal contact**  
Blue, but not too light and not too dark either...
- **Just the right setting time ! 30 seconds**  
Working time is 20 seconds and the setting time is 30 seconds.  
Easy to apply for any technique, occlusal or buccal, due to minimal resistance when extruding.
- **Just enough shrinkage ! Less than 0.1%**  
Keeps the correct shape for long periods, ensure fine prosthesis work and minimizes time and mistakes.
- **90 Shore A is just the right hardness**  
Not too soft ensuring correct positioning of upper and lower models. Not too hard ensuring the strength of the impression. The right hardness also helps technicians with easy trimming.
- **Incredibly minimal resistance for extrusion**  
Easy extrusion of the material makes placement simple and precise.




### DISPOSABLE 3-IN-1 TRAYS

- Wide occlusal plane
- Strong and tough plastic
- Specially designed walls for extra retention
- Tear resistant quality mesh
- Individually wrapped

#### 3-in-1 Trays

N401 Posterior 36 pcs/box  
N402 Quadrant 36 pcs/box  
N403 Sideless 36 pcs/box  
N404 Anterior 36 pcs/pack  
N405 Full Arch 24 pcs/pack

RIGID  
Stainless  
Steel Wire  
Connector



S401  
Posterior  
(Light Blue)

RIGID  
Stainless  
Steel Wire  
Connector



S402  
Quadrant  
(Dark Green)

RIGID  
Stainless  
Steel Wire  
Connector



S403  
Sideless  
(Peach)

### DISPOSABLE 3-IN-1 TRAYS Wire Connectors Stainless Steel

- Stainless steel wire connector for greater access in limited space
- Wide occlusal plane
- Special rigid walls designed for extra retention
- Tear resistant quality mesh

#### 3-in-1 Trays

S401 Posterior 36 pcs/box  
S402 Quadrant 36 pcs/box  
S403 Sideless 36 pcs/box

RIGID  
Aluminium  
Bar  
Connector



T401  
Posterior  
(Lavender)

RIGID  
Aluminium  
Bar  
Connector



T402  
Quadrant  
(Light Blue)

RIGID  
Aluminium  
Bar  
Connector



T403  
Sideless  
(Yellow)

### DISPOSABLE 3-IN-1 TRAYS Rigid Bar Connectors

- Rigid Aluminium Bar Connector
- Strong and tough
- Wide occlusal plane
- Special rigid walls designed for extra retention
- Tear resistant quality mesh

#### 3-in-1 Trays

T401 Posterior 36 pcs/box  
T402 Quadrant 36 pcs/box  
T403 Sideless 36 pcs/box


### FLEXIBLE MIXING BOWLS



- For effective mixing of alginates and plasters
- Sizes: 7460 Small dia. 11 cm x H8 cm  
7470 Medium dia. 12.5 cm x H10 cm
- Available in Transparent Red, Green, Blue, Purple, Yellow and solid Blue

7460 Flexible Mixing Bowl Small  
7470 Flexible Mixing Bowl Medium


### AUTOCLAVABLE ALGINATE AND PLASTER SPATULAS



- Fully autoclavable up to 134°C
- Fit all flexible mixing bowls
- Will not scratch the surface of bowls and clean up easily
- Flexible control in mixing with easy grip handle
- Size: Length: 190 mm

450R Regular Blue  
450XS Extra Stiff Grey

### DISPOSABLE INTRAORAL SYRINGE TIPS



- Universal fitting
- Available in Yellow or White

Intraoral Tips, 100 pcs/pack  
085 Yellow for Yellow mixing tips  
086 Clear for Green and Pink mixing tips

### DISPOSABLE DYNAMIC MIXING TIPS



- Fits Pentamix<sup>®</sup> type mixing machine
- Economically priced
- Colours: Red or White

7055 Dynamic Mixing Tips, 50 pcs/pack

<sup>1</sup> Pentamix is a registered trademark of 3M or 3M ESPE AG.



## DISPOSABLE PSP BARRIER ENVELOPES OXY-BIODEGRADABLE



- Quick tear-off opening
- Fits most phosphor plates
- Fully disposable envelope
- Soft vinyl with rounded corners

### DURR® TYPE

#### Mid-Opening Barrier Envelopes

188-0 No. 0 (I.D. 39 mm x 28 mm) 100 pcs/box  
 188-1 No. 1 (I.D. 44 mm x 29 mm) 100 pcs/box  
 188-2 No. 2 (I.D. 45 mm x 35 mm) 100 pcs/box  
 188-3 No. 3 (I.D. 58 mm x 30 mm) 50 pcs/box  
 188-4 No. 4 (I.D. 82 mm x 60 mm) 50 pcs/box

#### Mid-Opening Barrier Envelopes Bulk

188-2B No. 2 (I.D. 45 mm x 35 mm) 500 pcs/box

#### Top-Opening Barrier Envelopes

198-0 No. 0 (I.D. 27 mm x 40 mm) 100 pcs/box  
 198-1 No. 1 (I.D. 28 mm x 45 mm) 100 pcs/box  
 198-2 No. 2 (I.D. 34 mm x 46 mm) 100 pcs/box

#### Top-Opening Barrier Envelopes Bulk

198-2B No. 2 (I.D. 33 mm x 44 mm) 500 pcs/box

### DIGORA® (SOREDEX®) TYPE

#### Top-Opening Barrier Envelopes

189-0 No. 0 (I.D. 24 mm x 34 mm) 100 pcs/box  
 189-1 No. 1 (I.D. 26 mm x 43 mm) 100 pcs/box  
 189-2 No. 2 (I.D. 33 mm x 44 mm) 100 pcs/box

### KODAK® (CARESTEAM®) TYPE

#### Top-Opening Barrier Envelopes

199-0 No. 0 (I.D. 31 mm x 36 mm) 100 pcs/box  
 199-1 No. 1 (I.D. 33 mm x 44 mm) 100 pcs/box  
 199-2 No. 2 (I.D. 40 mm x 46 mm) 100 pcs/box

#### Protective Paper Covers

178-0 (2 x 23 mm x 36 mm) 100 pcs/box  
 178-1 (2 x 24 mm x 42 mm) 100 pcs/box  
 178-2 (2 x 30 mm x 42 mm) 100 pcs/box



## DISPOSABLE PC MOUSE SLEEVES OXY-BIODEGRADABLE

- Universal size, fit both cordless and corded optical mouse
- Soft, pliable polypropylene sleeves
- Easy to apply

PC Mouse Sleeves, 500 pcs/box  
 133 15 cm x 8.5 cm (5 7/8" x 3 3/8")

## DISPOSABLE X-RAY SLEEVES OXY-BIODEGRADABLE



- Full x-ray unit coverage
- Disposable yet economical
- Soft, pliable polypropylene sleeve holds position until removed
- Easy application and removal


X-Ray Head Sleeves, 250 pcs/box  
 108 Small 20 cm + 17 cm x 60 cm  
 (7 7/8" + 6 5/8" x 23 1/2")  
 109 Large 35 cm + 24 cm x 75 cm  
 (13 3/4" + 9 1/2" x 29 1/2")



## DISPOSABLE HVE TUBE SLEEVES OXY-BIODEGRADABLE

- Pre-cut opening for easier application
- Made from very supple poly plastic for greater tactile feel

HVE Suction Tube Sleeves, 500 pcs/box  
 135S Small (11 mm) 35 cm x 5.5 cm (13 3/4" x 2 1/8")  
 135L Large (16 mm) 38 cm x 6 cm (15" x 2 3/8")



### A. STANDARD

140XS  
 140S  
 140M



### B. TWO STEP

183-1  
 183-2

## DISPOSABLE X-RAY SENSOR SLEEVES OXY-BIODEGRADABLE

- Made from very supple poly plastic for greater tactile feel
- Easy to apply

A. Standard Sensor Sleeve, 500 pcs/box  
 140XS X-Small 3.8 cm x 20 cm (1 1/2" x 7 7/8")  
 140S Small 4.2 cm x 20 cm (1 5/8" x 7 7/8")  
 140M Medium 4.5 cm x 20 cm (1 3/4" x 7 7/8")  
 140L Large 6 cm x 20 cm (2 3/8" x 10")

B. Two Step Sensor Sleeves, 500 pcs/box  
 183-1 Size 1 3.4 cm / 5 cm x 22 cm (1 3/8" / 2" x 8 5/8")  
 183-2 Size 2 3.8 cm / 5.5 cm x 22 cm (1 1/2" / 2 1/8" x 8 5/8")



## DISPOSABLE INSTRUMENT TRAY SLEEVES OXY-BIODEGRADABLE

- Instruments are placed directly on top of the disposable plastic sleeve
- Reduces the possibility of cross-contamination, saves cleaning time

Tray Sleeves, 500 pcs/box  
 101 Standard 'B' 27 cm x 36 cm (10 1/2" x 14")  
 102 Mini 'F' 19 cm x 28 cm (7 1/2" x 11")



## DISPOSABLE ISO INSTRUMENT TRAY SLEEVES OXY-BIODEGRADABLE

- Instruments are placed directly on top of the disposable plastic sleeve
- Reduces the possibility of cross-contamination, saves cleaning time

ISO Metal Tray Sleeves, 500 pcs/box  
 190 Standard 31 cm x 21 cm  
 (12 1/8" x 8 1/4")



## DISPOSABLE IMPRESSION GUN SLEEVES OXY-BIODEGRADABLE

- Plastic sleeve covers cartridge and handle
- Provides protection in case of cartridge failure or leakage

Impression Gun Sleeves, 500 pcs/pack  
 107 Universal 25 cm x 21 cm (9 7/8" x 8 1/4")




## DISPOSABLE MOUTH MIRROR SLEEVES OXY-BIODEGRADABLE

- Made from very supple poly plastic for greater tactile feel

Mouth Mirror Sleeves, 500 pcs/box  
 140XS X-Small 3.8 cm x 20 cm (1 1/2" x 7 7/8")

## DISPOSABLE BARRIER FILM



- With 5 mm non-adhesive edge for simple placement and removal
- Residue free
- Colours: Blue and Clear

Barrier Film, 1200 pcs/roll  
 8008-C Clear 10 cm x 15 cm (3 7/8" x 6")  
 8008-B Blue 10 cm x 15 cm (3 7/8" x 6")

## DISPOSABLE T-STYLE LIGHT HANDLE SLEEVES OXY-BIODEGRADABLE



- Universal sizing to fit most T-style handles
- Full handle coverage
- Disposable yet economical
- Soft, pliable polypropylene sleeve holds position until removed
- Easy application and removal

T-Shape Light Handle Sleeves, 500 pcs/box  
 106 Universal 10 cm / 15 cm x 8.5 cm (3 7/8" / 6" x 3 3/8")



## DISPOSABLE BITE BLOCK SLEEVES OXY-BIODEGRADABLE

- Polypropylene sleeves
- Easy application and removal
- Fully disposable

Bite Block Sleeves, 5 x 200 pcs/box  
 145 Universal 5.5 cm x 2.6 cm (2 1/8" x 1")

## DISPOSABLE PEN TYPE CURING LIGHT SLEEVES



- Made from very supple poly plastic for greater tactile feel
- Easy to apply

**Pen Type Curing Light Sleeves, 500 pcs/box**  
 228 Mini Short 25 cm x 3.8 cm (9" x 1 1/2")  
 123SS Small Short 28 cm x 5 cm (11" x 2")  
 123S Small 32 cm x 5 cm (12 1/2" x 2")  
 123L Large 35 cm x 6.5 cm (13 3/4" x 2 1/2")  
 123XL X-Large 38 cm x 7.2 cm (15" x 2 7/8")

## DISPOSABLE INTRAORAL CAMERA /RA CURING LIGHT SLEEVES



- Made from very supple poly plastic for greater tactile feel
- Easy to apply

**Intraoral Camera Sleeves / RA Curing Light Sleeves**  
 228 Mini Short 25 cm x 3.8 cm (9" x 1 1/2"), 500 pcs/box  
 125 Standard 33 cm x 4.8 cm (13" x 1 7/8"), 500 pcs/box  
 125SS Short 28 cm x 4.8 cm (11" x 1 7/8"), 500 pcs/box

## DISPOSABLE PEN AND HANDPIECE SLEEVES




### OXY-BIODEGRADABLE

- Disposable sleeves for pen, handpieces, and ultrasonic scalers
- Made from soft pliable poly plastic for greater tactile feel

**500 pcs/box**  
 103 Pen Shape Instrument 22 cm x 2.4 cm (8 1/2" x 1")  
 104S Low Speed/Universal Slim 22 cm x 3.8 cm (8 1/2" x 1 1/2")  
 104L Low Speed/Universal Large 22 cm x 4 cm (8 1/2" x 1 5/8")

## DISPOSABLE LAB TRANSPORT BAGS OXY-BIODEGRADABLE



- Eliminates cross-contamination of packing materials and shipping containers from exposure to potentially infectious materials
- Permanent adhesive closure
- Printed with easy to understand instructions and the required Biohazard symbol for easy identification and compliance
- Convenient two-pouch system keeps requisition slips separate from lab work
- Meets requirements and recommendations for the transport of dental cases
- White front provides for easy writing of pertinent information

**Lab Transport Bags, 100 pcs/box**  
 121 Universal 17.5 cm x 24 cm (6 7/8" x 9 1/2")

## DISPOSABLE PEEL-N-STICK LIGHT GUIDE SLEEVES



- Soft, pliable polypropylene Sleeve
- Easy to apply
- Economical

**Light Guide Sleeves, 200 pcs/box**  
 128 Small (fits 7-8 mm Guide)  
 129 Large (fits 9-11 mm Guide)



## DISPOSABLE SYRINGE/HVE SLEEVES OXY-BIODEGRADABLE

- Pre-cut opening for easier application
- Ideal for 3-way syringes, saliva ejectors and HVE valves etc

**Syringe Sleeves, 500 pcs/box**  
 100C Clear 26 cm x 6.5 cm (10" x 2 1/2")  
 100B Blue 26 cm x 6.5 cm (10" x 2 1/2")

## DISPOSABLE FLOWABLE COMPOSITE SYRINGE SLEEVES OXY-BIODEGRADABLE



- Made from very supple poly plastic for greater tactile feel
- Easy to apply
- Economical

**Flowable Composite Syringe Sleeves, 200 pcs/box**  
 138 Universal 1.8 cm / 5 cm x 12.5 cm (3/8" / 2" x 5")



## DISPOSABLE KEYBOARD SLEEVES OXY-BIODEGRADABLE

- Soft, pliable polypropylene sleeves

**Keyboard Sleeves, 500 pcs/box**  
 158 Standard 16 cm x 55 cm (6 1/4" x 21 5/8")  
 108 Large 20 cm + 17 cm x 60 cm (7 7/8" + 6 5/8" x 23 1/2")




## DISPOSABLE CHAIR SLEEVES OXY-BIODEGRADABLE

- Eliminates the need for separate headrest covers
- Avoids the need for harmful prolonging the life of dental chairs/stools

**Chair Sleeves, 125 pcs/box**  
 111 Full Chair 180 cm x 72 cm (70" x 28 3/8")  
 112 Half Chair 100 cm x 65 cm (39 3/8" x 25 5/8")

## DISPOSABLE HEADREST SLEEVES OXY-BIODEGRADABLE



- Polypropylene sleeves
- Reduces possible cross-contamination from infectious aerosol spray

**Headrest Sleeves, 500 pcs/box**  
 113 Small 26 cm x 27.5 cm (10 1/4" x 10 7/8")  
 114 Large 35 cm x 28 cm (13 3/4" x 11")



## DISPOSABLE NON-WOVEN HEADREST COVERS

- Breathable non-woven fabric
- Soft, provides maximum patient comfort
- Colours: Light Blue and White

**Non Woven Headrest Covers, 500 pcs/box**  
 780S Small 26 cm x 28 cm (10 1/4" x 11")  
 780M Medium 35 cm x 28 cm (13 3/4" x 11")

## LEVEL 3 DISPOSABLE 3-PLY MASKS EAR-LOOP WITH SHIELDS FLUID RESISTANT ULTRA SENSITIVE LATEX FREE



FS2288 Level 3 Ultra Sensitive Masks with Shields Ear-Loop, 25 pcs/box

- LEVEL 3 the most appropriate masks for dental treatment ( based on US Standard ASTM F2101-7, F2299 and European Standard 14683: 2014 Annex B)
- Water proof, excellent filtration (99.9% at 3.0 µm)
- Anti-fog shields
- Comfortable, soft and non-irritating
- Permits air to pass through mask to prevent shield from fogging
- Colours: Blue, Green, Pink or Lavender




## LEVEL 3 DISPOSABLE FLUID RESISTANT 3-PLY MASKS ULTRA SENSITIVE

Level 3 Ultra-Sensitive Dental Masks, 50 pcs/box  
2288 Ear-Loop  
2299 Tie-On


- LEVEL 3 the most appropriate masks for dental treatment ( based on US Standard ASTM F2101-7, F2299 and European Standard 14683: 2014 Annex B)
- Water proof, excellent filtration (99.9% at 3.0 µm)
- Comfortable, soft and non-irritating
- Permits air to pass through mask to prevent glasses from fogging
- Colours: Blue, Green, Pink, Lavender or White

## LEVEL 3 DISPOSABLE DENTAL MASKS EAR-LOOP WITH PATTERNS FLUID RESISTANT ULTRA SENSITIVE LATEX FREE



- LEVEL 3 the most appropriate masks for dental treatment ( based on US Standard ASTM F2101-7, F2299 and European Standard 14683: 2014 Annex B)
  - Waterproof, excellent filtration (99.9% at 3.0 µm)
  - Comfortable, soft and non-irritating
  - Permits air to pass through mask to prevent glasses from fogging
- 2288 Level 3 Dental Masks Ear-Loop with Patterns, 50 pcs/box

## DISPOSABLE FLUID RESISTANT DENTAL BIBS WITH PATTERNS



- Two plies of facial quality tissue plus a layer of polyethylene film offers quick absorbency with a positive moisture barrier
- Paper specially treated for wet strength
- Embossed horizontal pattern, water contact disperses horizontally instead of "channelling" down onto patients clothing
- Size: 45 cm x 33 cm

Dental Bibs with Patterns  
2200FA Flowers Assorted, 500 pcs/box  
2200LA Leaves Assorted, 500 pcs/box  
2200KID Kids, 500 pcs/box  
2200MON Monsters, 500 pcs/box


## DISPOSABLE PROPHY BIBS EXTRA LONG FLUID RESISTANT



- Ideal for prophy, where extra protection is needed
- Features contour neck design for additional coverage to shoulders and arm
- Two plies of facial quality tissue plus a layer of polyethylene
- Size: W45 cm x L60 cm
- Colour: Light Blue

2208 Prophy Bibs Extra Long  
100 pcs/box + Bib Holder XS


## DISPOSABLE PROPHY BIBS FLUID RESISTANT TRAY COVER PAPERS



- Two plies of facial quality tissue plus a layer of polyethylene film for positive barrier
- Colours: Blue and white barrier

2220 "B" Plastic Tray Covers 200'S (29.5 cm X 21.5 cm)  
2221 "F" Plastic Tray Covers 200'S (20 cm X 13.5 cm)  
2223 "ISO" Aluminium Tray Covers 200'S (27.5 cm X 17.5 cm)

## DISPOSABLE NAPKIN HOLDERS




- Fully disposable
- Easy to apply
- Width: 15 mm
- Length: 44 cm extendable up to 60 cm

2209 Napkin Holders, 100 pcs/pack



## ez INTERDENTAL BRUSHES



- ISO size (colour coded)
- 100% recycled plastic handle
- Plastic free packaging
- BPA free
- Vegan friendly

**G800 Interdenal Brushes**  
**ISO Size No.1-6**  
**Single or Assorted Sizes**  
**Standard Pack, 6 pcs/box**  
**Bulk Pack, 10 x 6 pcs/pack**



G800-A No.1-6 Assorted



## DISPOSABLE BAMBOO CUPS

- Leakproof construction, non-slip design
- Dispenses easily from cup dispenser due to their rolled lip

**6800BB Bamboo Cups 170 ml (6 oz) Blank**  
**6800BH Bamboo Cups 170 ml (6 oz) Healthy Teeth**  
**[ four assorted designs ]**



vegan friendly



sustainably sourced



## DISPOSABLE HEALTHY TEETH DENTAL MASKS, BIBS AND CUPS

### HEALTHY TEETH LEVEL 3 DENTAL MASKS EAR-LOOP

- Ultra sensitive, non-irritating
- LEVEL 3 the most appropriate masks for dental treatment (based on US Standard ASTM F2101-7, F2299 and European Standard 14683: 2014 Annex B)
- Water proof, excellent filtration (99.9% at 3.0 µm)
- Comfortable, soft and non-irritating
- Permits air to pass through mask to prevent glasses from fogging

**2288HT Healthy Teeth Level 3 Dental Masks, 50 pcs/box**  
**2200HT Healthy Teeth Dental Bibs, 500 pcs/pack**  
**6800HT Healthy Teeth Paper Cups 170 ml (6 oz), 1,000 pcs/box**

### HEALTHY TEETH DENTAL BIBS

- Fun way to promote dental health to patients
- Two plies of facial quality tissue plus a layer of polyethylene film offers quick absorbency with a positive moisture barrier
- Paper specially treated for wet strength
- Size: **33 cm x 45 cm**

### HEALTHY TEETH PAPER CUPS

- Leakproof construction, non-slip design
- Dispenses easily from cup dispenser due to their rolled lip
- Fun way to promote dental health to patients



## DISPOSABLE DENTAL MASKS, BIBS AND CUPS

### LEVEL 3 DENTAL MASKS EAR-LOOP

- Ultra sensitive, non-irritating
- LEVEL 3 the most appropriate masks for dental treatment (based on US Standard ASTM F2101-7, F2299 and European Standard 14683: 2014 Annex B)
- Water proof, excellent filtration (99.9% at 3.0 µm)
- Comfortable, soft and non-irritating
- Permits air to pass through mask to prevent glasses from fogging
- Available colours: **Light Blue, Green, Lavender Pink** or White

### DENTAL BIBS

- Two plies of facial quality tissue plus a layer of polyethylene film offers quick absorbency with a positive moisture barrier
- Paper specially treated for wet strength
- Size: **33 cm x 45 cm**
- Available colours: **Light Blue, Green, Lavender, Dusty Rose** or White

### PAPER CUPS

- Leakproof construction, non-slip design
- Dispenses easily from cup dispenser due to their rolled lip
- Available colours: **Light Blue, Green, Lavender, Rose** or White

**2288 Level 3 Dental Masks, 50 pcs/box**  
**2200 Dental Bibs, 500 pcs/pack**  
**6800 Paper Cups 170 ml (6 oz), 1,000 pcs/box**

# Endo Instruments


## AUTOCLAVABLE ENDO STAND



- Fully autoclavable plastic stand with disposable sponge insert
- Cleans and holds endo instruments
- Silicone ring base resists sliding
- Fully autoclavable up to 134°C
- Colour: Blue

4833 Endo Stand Large (dia. 45 mm)  
812 Endo Sponges, 50 pcs/pack

## DISPOSABLE ENDO SPONGES




811  
for round endo stand  
size: 2" dia. x 3/8" thick  
(50mm dia. x 8mm)

812  
for oval and round endo stand  
size: 2 1/2" dia. x 3/8" thick  
(65mm dia. x 8mm)

- Fully disposable blue polyurethane sponges
- Colour: Blue

Endo Sponges, 50 pcs/pack  
811 for Round Endo Stand  
812 for Oval and Large Endo Stand

## DISPOSABLE PEEL-AND-STICK ENDO CUSHIONS



802 Peel-and-Stick  
Endo Cushions  
50 pcs/box

- Hold, arrange and clean small instruments during endodontic or implant procedures
- Ideal for files, screws, drills and pluggers etc
- Peel-and-stick back secures to patient napkins, tray covers or anywhere handy
- Pre-arranging saves set up time
- Dense blue polyurethane foam 75 mm long

## AUTOCLAVABLE ENDO CUSHION HOLDER



- Holder for used files with disposable foam cushions
- Built in endo file ruler
- Stainless steel construction
- Fully autoclavable up to 134°C
- Sizes: 50 mm x 50 mm x 25 mm

1180 Endo Cushion Holder  
(without Cushion)  
801 Square Endo Cushions  
50 pcs/pack

## AUTOCLAVABLE TRIANGULAR ENDO RING AND CUSHIONS



- Fully disposable cushions
- Colour: Blue
- Autoclavable up to 134°C

088 Autoclavable Cushion Holders  
with Large and Small Rings  
803 Disposable Cushions 50 pcs/box

## AUTOCLAVABLE ROUND ENDO MEASURE STAINLESS STEEL



032  
Endo  
Measure

806  
Cylindrical  
Endo Cushions  
50 pcs/box

- Holds 8 instruments at preset measurement
- Solid stainless construction, fully autoclavable
- With Endo Cushion Holder
- Size: dia. 33 mm



## AUTOCLAVABLE ENDO HOLDER ALUMINIUM

2800 Endo Holder

- Specially designed for ultrasonic cleaning and autoclaving of root canal instruments
- Holds 18 hand or engine endo instrument up to size 100 all lengths
- Millimeter rule for rubber stop measurements for both right or left hand use
- Adjustable cover with silicone stop usage counter for sterilization
- Size: W65 mm x D30 mm x H55 mm
- Fully autoclavable up to 135°C
- Colour: Blue




## AUTOCLAVABLE ENDO DISPENSER PLASTIC

055 Endo Dispenser

- Specially designed for cleaning and autoclaving of endo instruments
- Holds 18 hand or engine instruments up to size 100 (all lengths)
- Medical grade plastic, autoclavable up to 134°C
- Size: 85 mm x 52 mm x 14 mm
- Colour: Blue

## AUTOCLAVABLE ENDO MEASURING BLOCK PLASTIC



4888  
Endo Measuring Block

- Accurate settings
- Durable construction
- Fully autoclavable up to 134°C
- Size: L100 mm x D24 mm x H32 mm
- Colour: Blue

## AUTOCLAVABLE ENDO RULER STAINLESS STEEL



- Durable Stainless steel
- Accurate setting
- Size: L80 mm x D12.5 mm x H3 mm

4899 Endo Ruler


## AUTOCLAVABLE ENDO RULER ALUMINIUM



- For both hand and engine instruments
- Gauges from 1 to 30 mm for right or left hand

4898 Endo Ruler

## AUTOCLAVABLE ENDO MEASURING BLOCK ALUMINIUM



10005  
Endo Measuring Block

- Durable aluminium construction
- Accurate settings
- Fully autoclavable up to 134°C
- Size: L150 mm x W24 mm x H33 mm
- Colour: Blue



## ENDO BRUSH APPLICATORS

- For cleaning and bonding material application in root canals
- Economical, saving up to 80% by the reduction of excess materials used
- Fully disposable

3223 Endo Brush Applicators  
100 pcs/pack

## AUTOCLAVABLE RUBBER DAM FRAME



- Size: 12.5 cm
- Stainless steel

58-165-002  
Rubber Dam Frame

## AUTOCLAVABLE RUBBER DAM FRAME



- Proven design
- Better grip
- Enhance patient comfort

4866  
Rubber  
Dam Frame

## ENDO STOPS



- Round with tail
- Assorted Colours: White, Yellow, Red, Blue, Green and Black

329  
Endo Stops Assorted Colours  
100 pcs/pack



58-110-170  
AINSWORTH  
160 mm



58-110-160  
IVORY  
170 mm




58-164-170  
BREWER  
175 mm



58-152-170  
IVORY  
170 mm

## RUBBER DAM PUNCHES AND FORCEPS

### CROSS-SECTION



## E01 SUPER FILES GOLD NiTi heat activation

- Smart varied taper design
- High cutting efficiency
- Good flexibility and fracture resistance
- Torque 2.0 Ncm - 2.5 Ncm
- Recommended speed 150-300 r.p.m.

### E01 SUPER FILES GOLD NiTi heat activation


Length	Size	Taper %	Tip (mm)
19 mm	SX	3.5-9	0.19
21 mm or 25 mm	S1	2-11	0.17
	S2	4-11.5	0.20
	F1	7	0.20
	F2	8	0.25
	F3	9	0.30
	F4	6	0.40
	F5	5	0.50

E01 Super Files Gold  
Single or Assorted Sizes  
6 pcs/pack

### E01 RECOMMENDED OPERATION PROCESS

PROBE	#10	
	#15	
OPENING	SX	
SHAPING	S1	
	S2	
FINISHING	F1	
	F2	
	F3	
ENLARGING	F4	
	F5	

### CROSS-SECTION



## E02 SUPER FILES BLUE NiTi heat activation

- Smart varied taper design
- High flexibility and pre-curved memory
- Safe tip, reduces the possibility of root canal deviation
- Torque 2.0 Ncm - 2.5 Ncm
- Recommended speed 150-300 r.p.m.

### E02 SUPER FILES GOLD NiTi heat activation

Length	Size	Taper %	Tip (mm)
19 mm	SX	3.5-9	0.19
21 mm or 25 mm	S1	2-11	0.17
	S2	4-11.5	0.20
	F1	7	0.20
	F2	8	0.25
	F3	9	0.30
	F4	6	0.40
	F5	5	0.50

E02 Super Files Blue  
Single or Assorted Sizes  
6 pcs/pack

### E02 RECOMMENDED OPERATION PROCESS

PROBE	#10	
	#15	
OPENING	SX	
SHAPING	S1	
	S2	
FINISHING	F1	
	F2	
	F3	
ENLARGING	F4	
	F5	

CROSS-SECTION



## E03 PRO-FLEXI FILES NiTi heat activation function


- Simple effective system
- Unique flexibility
- Torque **2.0 Ncm - 2.5 Ncm**
- Recommended speed **150-300 r.p.m.**

*E03 Pro-Flexi Files*  
Single or Assorted Sizes  
6 pcs/pack

E03 PRO-FLEXI FILES heat activation function			
Length	Size	Taper %	Tip (mm)
19 mm	#01/09 ○	9	0.25
21 mm or 25 mm	#20/04 ●	4	0.20
	#25/04 ●	4	0.25
	#30/04 ●	4	0.30
	#35/04 ●	4	0.35
	#25/06 ●	6	0.25

E03 RECOMMENDED OPERATION PROCESS		
OPENING	#01	
NARROW CANAL	#20 (-4%)	
	#25 (-4%)	
MEDIUM CANAL	#30 (-4%)	
WIDE CANAL	#35 (-4%)	
	#25 (-6%)	

CROSS-SECTION



## E05 ONE FILES NiTi heat activation

- Reciprocating movement
- Use only one instrument for root canal preparation, efficient and time saving
- 3 sizes for different size root canals
- Torque **2.0 Ncm - 2.5 Ncm**
- Recommended speed **150-300 r.p.m.**

*E05 One Files*  
Single or Assorted Sizes  
6 pcs/pack

E05 ONE FILES NiTi heat activation			
Length	Size	Taper %	Tip (mm)
21 mm or 25 mm	#R25 ●	8	0.25
	#R40 ●	6	0.40
	#R50 ●	5	0.50

E05 ONE FILES	
#R25	
#R40	
#R50	



## E04 ULTRASONIC FILES NiTi alloy

- For use with ultrasonic scaler
- Excellent tenacity
- Length: 33 mm

E04 ULTRASONIC FILES NiTi alloy	
#15 ○	#30 ●
#20 ●	#35 ●
#25 ●	#40 ●

*E04 Ultrasonic Files*  
Single Size  
6 pcs/pack



## E11 SPIRAL FILIERS stainless steel

E11 SPIRAL FILLERS stainless steel	
Length	Size
21 mm or 25 mm	#25 ●
	#30 ●
	#35 ●
	#40 ●

*E11 Spiral Fillers*  
Single or Assorted Sizes  
4 pcs/pack



## E12 GATES DRILLS stainless steel

- High cutting efficiency
- Recommended speed **800-1,200 r.p.m.**

E12 GATES DRILLS stainless steel		
Length	Size	Diameter (mm)
28 mm or 32 mm	#1	0.5
	#2	0.7
	#3	0.9
	#4	1.1
	#5	1.3
	#6	1.5

*E12 Gates Drills*  
Single or Assorted Sizes  
6 pcs/pack



## E15 PEESO REAMERS stainless steel


- High cutting efficiency
- Recommended speed **800-1,200 r.p.m.**

E15 PEESO REAMERS stainless steel		
Length	Size	Diameter (mm)
28 mm or 32 mm	#1	0.7
	#2	0.9
	#3	1.1
	#4	1.3
	#5	1.5
	#6	1.7









*E15 Peeso Reamers*  
Single or Assorted Sizes  
6 pcs/pack



E38 BARBED BROACHES stainless steel	
Length	Size
21 mm	XXXXF (0)
	XXXF (1)
	XXF (2)
	XF (3)
	F (4)
	M (5)
	C (6)



E48 C+ FILES stainless steel	
Length	Size
21 mm or 25 mm	#06
	#08
	#10
	#12.5
	#15







## DISPOSABLE EYE SHIELDS

- Fully disposable
- Available in tinted or clear lenses
- Light weight, optical plastic shields
- Assorted colours frames: Blue, Light Blue, Lavender, Green and White

**Disposable Eye Shield**  
20 Lenses + 20 Frames (Assorted Colours)  
720T, Tinted  
720C, Clear

**Replacement Eye Shield**  
Lens Only, 50 pcs/pack  
720T-R, Tinted  
720C-R, Clear



## FACE SHIELDS

- Light weight
- Adjustable design to suit every user
- High quality replaceable shields
- Colours: Blue, White and Lavender

1145 Face Shield Kit complete with 3 shields  
1146 Replacement Shields, 3's



## PROTECTIVE EYEWEARS

- Adjustable side arms
- Contemporary design
- Economical price

500  
500 Protective Eyewear Clear  
505 Protective Eyewear Tinted




## DISPOSABLE MOUTH MIRRORS

- Fully disposable plastic mouth mirrors for patients
- Anti-fog front surface fiber glass mirror
- Mirror size: dia. 23 mm

800 Disposable Mouth Mirrors, 10 pcs/box

## DISPOSABLE PLASTIC WEDGES



actual size  
7650X (17 mm)  
7650M (13 mm)  
7650L (15 mm)  
7650S (11 mm)

- Anatomically designed
  - Soft tough plastic
- Disposable Plastic Wedges, 100 pcs/box  
7650X X-Large (17 mm) White  
7650L Large (15 mm) Yellow  
7650M Medium (13 mm) Green  
7650S Small (11 mm) Blue



## TONGUE CLEANERS

- Effective in the treatment of halitosis
- Multiple cleaning ridges
- Assorted colours

9000 Tongue Cleaners  
20 pcs/pack

## RETAINER BOXES/MOUTH GUARD BOXES/DENTURE BATHS



600 Retainer Box W75 mm x D61 mm x H28 mm  
601 Mouth Guard Box W77 mm x D64 mm x H40 mm  
602 All Guard Box W77 mm x D77 mm x H38 mm  
603 Denture Bath W83 mm x D68 mm x H52 mm

Colours: Blue, Green, Yellow, White, Lime, Light Blue, Lavender or Assorted  
600 Retainer Boxes, 10 pcs/box  
601 Mouth Guard Boxes, 10 pcs/box  
602 All Guard Boxes, 10 pcs/box  
603 Denture Baths, 10 pcs/box

## DENTAL APPLIANCE BOXES WITH MIRRORS




PASTEL COLOURS	
Vanilla	Sea Green
Sand	Lilac
Light Mauve	Baby Blue
White	Assorted

- Sturdy strong plastic
- Size: W80 mm x D72 mm x H32 mm

612  
Dental Appliance Boxes  
10 pcs/box

## INTRAORAL CHEEK RETRACTORS



- Holds mouth open, depresses tongue movement
- Anatomically designed
- Available in both autoclavable and disposable

**Autoclavable Intraoral Cheek Retractor**  
2060A Adult Large, 1 pc/box  
2060B Adult Small, 1 pc/box  
2060C Child, 1 pc/box



**Disposable Intraoral Cheek Retractors**  
2080A Adult Large, 10 pcs/box  
2080B Adult Small, 10 pcs/box  
2080C Child, 10 pcs/box

## INTRAORAL MOUTH LIPPERS



- Expands the mouth vertically rather than horizontally
- Ideal for working in the gingival areas for anteriors, ortho procedures, bleaching etc...
- Available in both autoclavable and disposable

**Autoclavable Intraoral Lippers, 2 pcs/box**  
2040A Adult, Semi-Transparent  
2040C Child, Semi-Transparent  
**Disposable Intraoral Lippers, 10 pcs/box**  
2030A Adult, White  
2030C Child, White




## AUTOCLAVABLE CHEEK RETRACTORS



- Autoclavable up to 134°C

**Mouth Retractor, 2 pcs/pack**  
2050C-2 Child  
2050A-2 Adult  
2055C-2 Child  
2055A-2 Adult



## AUTOCLAVABLE MOUTH PROPS LATEX FREE

2090A Mouth Prop Adult  
2090P Mouth Prop Pado

- Made of chemical resistant, autoclavable silicone
- Durable
- Available in: - Adult (30 mm x 40 mm x 20 mm)  
Pado (30 mm x 25 mm x 20 mm)

### DISPOSABLE MULTI-PURPOSE DAPPEN DISHES



004B Standard Pack  
Assorted Colours  
1000 pcs/box

- Fully disposable to save time and prevent cross-contamination
- For prophy, etchant and bonding etc
- Colours: Lavender, White, Green, Light blue and Yellow

### DISPOSABLE PROPHY RINGS



- Fully disposable to save time and prevent cross-contamination
- For prophy, endo etc
- Colour: Blue

010 Disposable Propy Rings, 200 pcs/box

### AUTOCLAVABLE PROPHY RINGS



- Autoclavable up to 134°C
- For prophy, endo etc
- Colour: Grey

A010 Autoclavable Propy Rings, 10 pcs/pack

### GLASS DAPPEN DISH



- Colours: Light Blue, Green, Clear, Brown or Pink

040 Glass Dappen Dish

### DOUBLE ARCH



LARGE blue



MEDIUM yellow



SMALL green

### SINGLE ARCH



LARGE blue



MEDIUM yellow



SMALL green

### DISPOSABLE FLUORIDE TRAYS

- Extremely comfortable, anatomically designed trays
- Economical price
- Colour coded

Double Arch, 50 pcs/box  
7800S Small (green)  
7800M Medium (yellow)  
7800L Large (blue)

Single Arch, 100 pcs/box  
7801S Small (green)  
7801M Medium (yellow)  
7801L Large (blue)

### DISPOSABLE MIXING PADS



3220

3225

Mixing Pads, 3X100 pcs/pack  
3220 Small 35 mm X 45 mm  
3225 Large 70 mm X 80 mm

### DISPOSABLE MIXING WELLS



906

907

- Fully disposable
- Sturdy design

906 Mixing Wells 4 Wells 200 pcs/box  
907 Mixing Wells 2 Wells 500 pcs/box

### DISPOSABLE FLOWABLE COMPOSITE AND ETCHING GEL TIPS



actual size

081 ETCH-TIP

080 FLOW-TIP

082 LINER-TIP

- Universal Luer-lock design
- Bendable
- Gauges: 080 G.20 (0.9 mm)  
081 G.23 (0.6 mm)  
082 G.18 (1.2 mm)

100 pcs/pack  
080 Flow-Tips Black G.20  
081 Etch-Tips Light Blue G.23  
082 Liner-Tips Pink G.18

### DISPOSABLE VALUSTICK STICKY APPLICATORS



TWIST for removal

- Bends into any angle for placement of Veneers, Inlays, Onlays, Crowns etc
- Strong adhesive for multiple applications
- Fully disposable
- Length: 110 mm
- Colour: White

909 Valustick Sticky Applicators  
50 pcs/box

### AUTOCLAVABLE ASPIRATOR BRUSHES



020

021

022

023

- Autoclavable up to 134°C
- Rigid bristles for effective cleaning

020 X-Fine dia. 5 mm, 6 pcs/pack  
021 Fine dia. 8 mm, 6 pcs/pack  
022 Medium dia. 14 mm, 6 pcs/pack  
023 Large dia. 20 mm, 6 pcs/pack


### AUTOCLAVABLE SUCTION TIP SLEEVES



- For protection from sharp edges on suction tubes
- Fully autoclavable up to 134°C
- Colour: Light Blue

290 Suction Tip Sleeves  
6 pcs/pack

### AUTOCLAVABLE SUCTION TUBE ADAPTORS



Saliva Ejector 6.5 mm

11 mm

306 Saliva Ejector Adaptor reduces 11 mm to 6.5 mm

Standard HVE Hose 11 mm

Suction Tube 11 mm

16 mm

11 mm

305 Suction Tube Adaptor reduces 16 mm to 11 mm

European HVE Hose 16 mm

## OBLONG ARTICULATING MOUNTING PLATES



- Fit Denar type articulators
- Fully disposable
- Colours: Blue or White

490 Mounting Plates, 100 pcs/box

## AUTOCLAVABLE ACRYLIC MIXING CUPS



- Fully autoclavable up to 134°C

280 Mixing Cup Small dia. 25 mm  
281 Mixing Cup Medium dia. 35 mm  
282 Mixing Cup Large dia. 50 mm  
283 Mixing Cup X-Large dia. 70 mm

## DISPOSABLE EVACUATION TRAPS



dia. 1 7/8" fits Belmont, Den-Tal-Ez Pelton & Crane, Proma etc

041




dia. 2 1/8" fits Adec models with plastic canisters

042

041 Disposable Evacuator Traps dia. 48 mm Green, 20 pcs/pack  
042 Disposable Evacuator Traps dia. 54 mm Blue, 20 pcs/pack

## AUTOCLAVABLE ALUMINIUM BUR BLOCKS



- Ideal for storage of dental burs
- Specially designed for ultrasonic cleaning and autoclaving of burs
- Sturdy construction, fully autoclavable up to 134°C
- Colours: Silver and Blue

**Small Bur Block**  
038-FG 15 FG W58 mm x D25 mm x H31 mm  
039-MIX 10 FG + 5 RA W58 mm x D25 mm x H31 mm  
037-HP 6 HP W54 mm x D28 mm x H62 mm  
**Large Bur Block**  
2832-FG 32 FG W80 mm x D42 mm x H33 mm  
2832-MIX 16 FG + 16 RA W80 mm x D42 mm x H33 mm  
2880-FG 80 FG W100 mm x D75 mm x H33 mm  
2880-MIX 40 FG + 40 RA W100 mm x D75 mm x H33 mm


## DISPOSABLE EVACUATION TRAPS



- Fits Adec type unit

044 Evacuator Traps  
20 pcs/pack

## AUTOCLAVABLE PLASTIC BUR DISPENSERS



- Specially designed for cleaning and autoclaving of burs
- Locks open or closed for processing
- Medical grade plastic, autoclavable up to 134°C
- Size: 75 mm x 35 mm x 10 mm

050-FG, Bur Dispenser 12 FG  
050-MIX Bur Dispenser 6 FG + 6 RA

COLOURS	Yellow
Blue	White
Light Blue	Light Beige
Green	Lavender

## AUTOCLAVABLE MAGNETIC BUR STANDS



COLOURS	Yellow
Blue	White
Light Blue	Light Beige
Green	Lavender

- Stands autoclavable up to 134°C
- Fits both FG and RA burs
- Size: Small Square 42 mm x 42 mm x 12 mm  
Large Rectangular 80 mm x 42 mm x 12 mm

033 Bur Stand Square, 8 Holes 034 Bur Stand Rectangular, 14 Holes  
033A Bur Stand Cover Square 034A Bur Stand Cover Rectangular

## AUTOCLAVABLE BRASS BUR BRUSH



- Telescopic handle with brass wire bristles
- Length: 12 cm

1110 Brass Bur Brush

## DISPOSABLE CURVED TIP SYRINGES



078  
Curved Tip  
Syringes  
50 pcs/box

- Made of polypropylene plastic
- Full 12 ml capacity
- Self sealing medical grade plunger
- Designed for perio irrigation, impression materials, post surgical home use, fluoride treatments etc

## AUTOCLAVABLE FRACTURE DETECTOR

The instrument for locating cracks in teeth



- Fully autoclavable to 134°C
- Simple tool for practice
- Quick and precise to locate cracks
- Can be used for crown and bridge seating

025 Autoclavable  
Fracture Detector Blue

## GLASS CEMENT MIXING SLAB



- Polished glass mixing slab with bevelled edges
- Size: 65 mm x 135 mm x 15 mm

045 Glass Slab

## COTTON ROLL DISPENSER



1140 Roll Dispenser

- Provides an aseptic barrier against operator aerosols and accidental splatter
- Top flips back for easy loading, sliding drawer dispenses up to 4 rolls with a single pull
- Foam rubber strips inserted into base prevent slipping on table top or counter

COLOURS	Blue	Light Blue	Green
Yellow	White	Light Beige	Lavender

### AUTOCLAVABLE IMPLANT SCALER SET



- Fully autoclavable up to 134°C
  - Replaceable carbon reinforced resin tips, can be resharpended
  - Extra light hollow handle
- Implant Scaler H6/H7 Hollow Handle*  
1128 Handle + 24 Tips (12 each H6/H7)  
1128A 24 Tips (12 each H6/H7)

### DISPOSABLE SELF-SEALING STERILIZATION POUCHES



- With internal and external steam indicators
- Triple adhesive seal to ensure sterile storage
- Blue tinted plastic window helps detect puncture holes that may be caused by jagged instruments

**Standard Pack, 200 pcs/box**  
6001 2 1/4" x 5" (55 mm x 130 mm)  
6003 3 1/2" x 5 3/4" (90 mm x 145 mm)  
6005 2 3/4" x 9" (70 mm x 230 mm)  
6007 3 1/2" x 9" (90 mm x 230 mm)  
6009 5 1/4" x 10" (135 mm x 260 mm)  
6011 7 1/2" x 13" (190 mm x 330 mm)

**Standard Pack, 100 pcs/box**  
6015 10" x 15" (255 mm x 380 mm)



### AUTOCLAVABLE NAPKIN HOLDERS

090 Napkin Holder Adult (L47 cm)  
090C Napkin Holder Child (L41 cm)  
090XS Napkin Holder Extra Short (L37 cm)

- Medical grade silicone
  - Fully autoclavable up to 134°C
  - Stainless steel alligator clips
  - Non-slip design
  - Bright, attractive colours
- Colours available  
Adult: Light Blue, Blue, Green, Red, Yellow, White, Pink and Lavender  
Child: Light Blue, Lavender, Blue and Pink  
XS (for neck contour bibs): Blue only



### DENTAL BIB HOLDER AND CLIP

- Cold disinfection only
- Cute design

091 Clip  
092 Bib Holder, White or Red



### PEN / TOOTHBRUSH HOLDER

- Cute design

2005  
Pen / Toothbrush Holder



### JUMPING TEETH

- Cute design
- Economical giveaway for kids

8001 Jumping Teeth  
24 pcs/box

### DENTAL MANDRELS STAINLESS STEEL



**Dental Mandrels, 6 pcs/box**  
M01 HP  
M02 RA  
M03 Pop On

M04 Super-Snap® Type  
M05 OneGloss® Type  
M06 Snap-On Prophy  
M08 FG/RA Adaptors

### AUTOCLAVABLE X-RAY FILM HOLDER PLASTIC



- One end for anterior x-rays, the other end for posterior
- Both ends securely hold film without bending it
- Improves quality and accuracy of intraoral films
- Helps to standardize the bisecting angle technique
- Colours: Blue and Light Beige

1188 X-Ray Film Holder



### DISPOSABLE BITE WING TABS

- Attach to adult and children's x-ray films or sensor sleeves
- Holds securely when patient closes mouth for accurate distortion-free bite wing radiographs
- Convenient pop-out dispenser

2010 Bite Wing Tabs, 500 pcs/box

### DISPOSABLE BITE WING LOOPS



- Convenient pop-open design ensures quick and easy use
- Manufactured from quality ten point paper stock
- Quick, simple and economical


**Bite Wing Loops, 500 pcs/box**  
A. 2011-1 Film Child  
B. 2011-2 Film Adult  
C. 2012-1 Sensor Size 1  
D. 2012-2 Sensor Size 2



## DELUXE ACRYLIC COMPOSITE ORGANIZERS

- Smart, flexible design
- Sturdy, durable construction
- Clear material makes contents easy to see
- For storage of composites, adhesive, micro applicators etc.

**Composite Organizer**  
 M060 Mini *W19.5 cm x H9.5 cm x D14.5 cm*  
 S060 Standard *W29 cm x H9.5 cm x D14.5 cm*



## DISPOSABLE MICRO APPLICATOR DISPENSERS

**Applicator Dispenser**  
 18-900 *Standard Handle*  
 18-9900 *Mini Handle*

- Small, sturdy, durable construction
- Provides aseptic barrier against operator aerosols and accidental splatter

## ACRYLIC ORTHO CHAIN DISPENSER



- Long lasting dispenser that will fit the decor of any office
- Chain is dispensed easily with an interior roller
- Holds up to 15 spools


059  
 Chain Dispenser



## ACRYLIC COMPOSITE DISPENSER

- Smart, durable construction
- Clear material makes contents easy to see
- For storage of 20 composite syringes, etch and bonding material etc.

H060 Composite Dispenser



## ACRYLIC COMPOSITE SYRINGE HOLDERS

551 Syringe Holder *Small*  
 552 Syringe Holder *Medium*

- Smart, durable construction
- Size: 551 Small (11 holes) *W10 cm x D12.5 cm x H12 cm*  
 552 Medium (23 holes) *W20 cm x D12.5 cm x H12 cm*



## LIGHT PROTECTION HAND SHIELD

- Hand held shield for curing light
- Protection for both dentist and assistant
- Cold disinfection only
- Size: *12.5 cm x 26 cm x 0.3 cm*

063 Light Protection Hand Shield



## LIGHT PROTECTION COVERS

- Extends working time of composite materials
- Sturdy, durable construction
- Cold disinfection only
- Size: Small *W10 cm x D7.5 cm x H3 cm*  
 Medium *W15 cm x D10 cm x H3 cm*

062-0 Light Protection Cover *Small*  
 062-1 Light Protection Cover *Medium*



## ACRYLIC IMPRESSION GUN HOLDERS

- For both loaded or empty guns
  - Store all guns in upright position
  - Include storage compartments
  - Sturdy, durable construction
  - Sizes:  
072 W145 mm x D140 mm x H150 mm  
074 W265 mm x D140 mm x H150 mm  
076 W395 mm x D140 mm x H150 mm
- 072 Gun Holder for 2 guns  
074 Gun Holder for 4 guns  
076 Gun Holder for 6 guns



## ACRYLIC GRINDING BOX

- Ideal for grinding and sand blasting
  - Retains flying debris
  - Sturdy durable construction
  - Size: W30 cm x D20 cm x H18 cm
- 065 Grinding Box

## ACRYLIC GLOVE AND TISSUE DISPENSERS



- Smart and attractive acrylic design
  - Sturdy, durable construction
  - Pre-drilled holes for wall mounting
  - Internal measurements:-  
064 Glove Dispenser Standard W25 cm x D14 cm x H9 cm  
066 Glove Dispenser Small W25 cm x D14 cm x H7 cm  
068 Glove/Tissue Dispenser Large W24 cm x D12.5 cm x H10 cm  
069 Glove/Tissue Dispenser Extra Large W24 cm x D12.5 cm x H11 cm
- 064 Glove Dispenser *Standard*  
066 Glove Dispenser *Small*  
068 Glove/Tissue Dispenser *Large*  
069 Glove/Tissue Dispenser *Extra Large*



## ACRYLIC CUP DISPENSER STAINLESS STEEL


- Wall mounted deluxe cup dispenser
  - Holds approx. one hundred cups
  - Can be adjusted to hold cups from 120 mls to 150 mls (5 oz to 7 oz)
  - Sturdy design
  - Size:  
H30 cm x W8 cm x D8 cm
- 1150 Cup Dispenser



## ACRYLIC CUP DISPENSER PLASTIC

- Magnetic, stick on or screw-in
  - Holds approx. one hundred cups
  - Can be adjusted to hold cups from 120 mls to 150 mls (5 oz to 7 oz)
  - Size: Dia. 80 mm x L360 mm
- 1250 Cup Dispenser

## ACRYLIC BIB AND C-FOLD TOWEL DISPENSERS




- Sturdy, durable construction
  - Smart, attractive acrylic design
  - Sizes:  
057 W15 cm x D11 cm x H12 cm  
058 W23.5 cm x D10 cm x H13.5 cm
- 057 Bib Dispenser  
058 C-Fold Towel Dispenser  
057CV Bib Dispenser With Cover  
058CV C-Fold Towel Dispenser With Cover



## ACRYLIC MASK DISPENSER

095 Mask Dispenser

- Sturdy, durable construction
- Clear material makes contents easy to see
- Pre-drilled holes for wall mounting
- Size: W94 mm x D106 mm x H108 mm



## ACRYLIC POUCH DISPENSERS

- Can be used on table top or wall mounted in either vertical or horizontal position
- Sizes:-  
061-0 (7.5 cm x 14.5 cm x 21 cm)  
061-1 (8.5 cm x 27.5 cm x 21 cm)  
061-2 (10 cm x 18.5 cm x 21 cm)  
061-3 (10 cm x 27.5 cm x 21 cm)  
061-4 (15 cm x 29.5 cm x 21 cm)

061-0 Pouch Dispenser for 5.5 cm x 10 cm Pouches  
061-1 Pouch Dispenser for 7 cm x 23 cm Pouches  
061-2 Pouch Dispenser for 9 cm x 13.5 cm Pouches  
061-3 Pouch Dispenser for 9 cm x 23 cm Pouches  
061-4 Pouch Dispenser for 13.5 cm x 26 cm Pouches

## ACRYLIC BARRIER FILM DISPENSER



- Sturdy durable construction
- Clear material makes contents easy to see
- Size: dia. W14 cm x D23 cm x H19 cm

1155 Barrier Film Dispenser

# Orthodontic Brackets



## ORTHODONTIC MACHINING BRACKETS

- Wire cutting technology
- Colour I.D. for easy identification
- Generous tie wing area for easy ligation

- Torque-in-base design techniques to provide level slot line-up and helps to avoid occlusal interference
- Smooth surface and rounded corners for enhanced patient comfort

MACHINING BRACKETS STANDARD ROTH® TYPE (20 pcs/set)				
ROTH® TYPE STANDARD	Slot Size	w/o Hook	w Hook 3	w Hooks 3,4,5
	18	BR011-20	BR011-23	BR011-25
	22	BR012-20	BR012-23	BR012-25
MACHINING BRACKETS MINI ROTH® TYPE (20 pcs/set)				
ROTH® TYPE MINI	Slot Size	w/o Hook	w Hook 3	w Hooks 3,4,5
	18	BR011-10	BR011-13	BR011-15
	22	BR012-10	BR012-13	BR012-15

MACHINING BRACKETS STANDARD MBT® TYPE (20 pcs/set)				
MBT® TYPE STANDARD	Slot Size	w/o Hook	w Hook 3	w Hooks 3,4,5
	18	BR021-20	BR021-23	BR021-25
	22	BR022-20	BR022-23	BR022-25
MACHINING BRACKETS MINI MBT® TYPE (20 pcs/set)				
MBT® TYPE MINI	Slot Size	w/o Hook	w Hook 3	w Hooks 3,4,5
	18	BR021-10	BR021-13	BR021-15
	22	BR022-10	BR022-13	BR022-15



## ICE CRYSTAL CERAMIC BRACKETS

- Ceramic brackets made from polycrystalline alumina ensures clear and transparent brackets
- Low contour design for maximum patient comfort
- Mechanical base design for better bonding strength and convenience when debonding
- Generous tie wing space for easy ligating
- Enhanced base design for easy and safe debonding



## SAPPHIRE ULTRA CLEAR BRACKETS

- Lowest friction aesthetic sapphire ceramic brackets
- Superior monocrystalline clear appearance
- Low contour design for maximum patient comfort

ICE CRYSTAL CERAMIC BRACKETS ROTH® TYPE (20 pcs/set)			
ROTH® TYPE	Slot Size	w Hook 3	w Hooks 3,4,5
	18	BR211-23	BR211-25
	22	BR212-23	BR212-25
ICE CRYSTAL CERAMIC BRACKETS MBT® TYPE (20 pcs/set)			
MBT® TYPE	Slot Size	w Hook 3	w Hooks 3,4,5
	18	BR221-23	BR221-25
	22	BR222-23	BR222-25

SAPPHIRE ULTRA CLEAR BRACKETS ROTH® TYPE (20 pcs/set)			
ROTH® TYPE	Slot Size	w Hook 3	w Hooks 3,4,5
	18	BR211-33	BR211-35
	22	BR212-33	BR212-35
SAPPHIRE ULTRA CLEAR BRACKETS MBT® TYPE (20 pcs/set)			
MBT® TYPE	Slot Size	w Hook 3	w Hooks 3,4,5
	18	BR221-33	BR221-35
	22	BR222-33	BR222-35

## NOBLE CLASSIC MINI CERAMIC BRACKETS WITH METAL SLOT



- Poly-crystalline aesthetic ceramic
- Wire slot enhanced with metal, low friction
- Twin wing
- Contoured designed for maximum comfort
- Triple polished smooth surface
- D80 Mesh base ensures good bonding

## ICE METAL MINI BRACKETS LOW PROFILE NICKEL FREE



- Hard **nickel free** stainless steel
- One piece construction mesh base (MIM)
- Superb appearance for the maximum patient comfort
- Low profile wing design for easy ligating
- D80 mesh base ensures good bonding

ROTH® TYPE MINI CERAMIC METAL SLOT	w Hook 3	018	022
MBT® TYPE MINI CERAMIC METAL SLOT	w Hook 3	018	022

ROTH® TYPE MINI NICKEL FREE	w Hook 3	022
	w Hooks 3, 4, 5	
MBT® TYPE MINI NICKEL FREE	w Hook 3	022
	w Hooks 3, 4, 5	

## CLEVER M SELF-LIGATING BRACKETS ACTIVE LOW NICKEL



- MIM advanced technology, hard 17-4 stainless steel
- Newest aesthetic design metal brackets with Interactive Self-Ligating Clip
- Contoured appearance for maximum comfort
- Low friction and constant gentle force
- Smooth surface for maximum comfort
- Laser marking for easy identification
- Mini profile design
- D80 Mesh base ensures good bonding

## THOMAS M SELF-LIGATING BRACKETS PASSIVE LOW NICKEL



- Stainless steel 17-4 moulded construction for exceptional strength and durability
- Passive smooth clip design, opens and closes easily
- Small size with low profile, smooth round corners for patients' maximum comfort
- Vertical and horizontal auxiliary slots, can be used for traction, and provide greater versatility
- Laser welding mesh base for strong bonding
- Advanced passive self-ligation technology with light force orthodontic system, can effectively reduce treatment time

ROTH® TYPE ACTIVE	w Hook 3	022
	w Hooks 3, 4, 5	

MBT® TYPE PASSIVE	w Hook 3	022
	w Hooks 3, 4, 5	

## COSY-II BONDABLE BUCCAL TUBES TRUMPET MONOBLOCK



- Teeth number marks are easy for identification
- Trumpet design at mesial entrance for easy wire engagement
- Comfort ball hook enhances patient comfort and avoids gingival impingement
- Low profile, round corners with smooth surface enhances patient comfort
- Torque-in-base design techniques to provide level slot line-up, and help avoid occlusal interference
- Compound contour base anatomically conforms to the tooth and provides proper fit

COSY-II BONDABLE NON-CONVERTIBLE ROTH® TYPE (4 pcs/kit)			
ROTH® TYPE 1 <sup>st</sup> MOLAR	U/L single	018	BT611-4310
		022	BT612-4310
ROTH® TYPE 2 <sup>nd</sup> MOLAR	U/L single	018	BT711-6310
		022	BT712-6310

COSY-II BONDABLE NON-CONVERTIBLE MBT® TYPE (4 pcs/kit)			
MBT® TYPE 1 <sup>st</sup> MOLAR	U/L single	018	BT621-4310
		022	BT622-4310
MBT® TYPE 2 <sup>nd</sup> MOLAR	U/L single	018	BT721-6310
		022	BT722-6310



## BONDABLE BUCCAL TUBES CONVERTIBLE / NON-CONVERTIBLE

- Comfort ball hook enhances patient comfort and avoids gingival irritation
- Round corners and smooth surface for enhanced patient comfort
- Torque-in-base design techniques to provide level slot line-up, and help avoid occlusal interference
- Compound contour base anatomically conforms to the tooth and provides proper fit

BONDABLE CONVERTIBLE ROTH® TYPE (4 pcs/kit)			
ROTH® TYPE 1 <sup>st</sup> MOLAR	0.018"	0.022"	Designation
	BT611-1110	BT612-1110	upper 1 tube / lower 1 tube
	BT611-1120	BT612-1120	upper 2 tube / lower 1 tube
	BT611-1140	BT612-1140	upper / lower double rectangular
	BT611-1130	BT612-1130	upper 3 / lower double rectangular
	BT611-1160	BT612-1160	upper 3 tube / lower 2 tube

BONDABLE CONVERTIBLE MBT® TYPE (4 pcs/kit)			
MBT® TYPE 1 <sup>st</sup> MOLAR	0.018"	0.022"	Designation
	BT621-1110	BT622-1110	upper 1 tube / lower 1 tube
	BT621-1120	BT622-1120	upper 2 tube / lower 1 tube
	BT621-1140	BT622-1140	upper / lower double rectangular
	BT621-1130	BT622-1130	upper 3 / lower double rectangular
	BT621-1160	BT622-1160	upper 3 tube / lower 2 tube

BONDABLE NON-CONVERTIBLE ROTH® TYPE (4 pcs/kit)			
ROTH® TYPE 1 <sup>st</sup> MOLAR	0.018"	0.022"	Designation
	BT611-2010	BT612-2010	upper 1 tube / lower 1 tube
	BT611-2020	BT612-2020	upper 2 tube / lower 1 tube
ROTH® TYPE 2 <sup>nd</sup> MOLAR	0.018"	0.022"	Designation
BT711-5010	BT712-5010	upper 1 tube / lower 1 tube	

BONDABLE NON-CONVERTIBLE MBT® TYPE (4 pcs/kit)			
MBT® TYPE 1 <sup>st</sup> MOLAR	0.018"	0.022"	Designation
	BT621-2010	BT622-2010	upper 1 tube / lower 1 tube
	BT621-2020	BT622-2020	upper 2 tube / lower 1 tube
MBT® TYPE 2 <sup>nd</sup> MOLAR	0.018"	0.022"	Designation
BT721-5010	BT722-5010	upper 1 tube / lower 1 tube	

## MOLAR BANDS ROUGH INNER SURFACE



1 <sup>st</sup> MOLAR BANDS (ROUGH INNER SURFACE) (4 pcs/kit)																					
No	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
Peri-meter mm	30	30.5	31	31.5	32	32.5	33	33.5	34	34.5	35	35.5	36	36.5	37	37.5	38	38.5	39	39.5	40

2 <sup>nd</sup> MOLAR BANDS (ROUGH INNER SURFACE) (4 pcs/kit)																	
No	32+	33	33+	34	34+	35	35+	36	36+	37	37+	38	38+	39	39+	40	40+
Peri-meter mm	32.5	33	33.5	34	34.5	35	35.5	36	36.5	37	37.5	38	38.5	39	39.5	40	40.5



## WELDABLE BUCCAL TUBES CONVERTIBLE / NON-CONVERTIBLE

- Comfort ball hook enhances patient comfort and avoids gingival irritation
- Round corners and smooth surface for enhanced patient comfort
- Torque-in-base design techniques to provide level slot line-up, and help avoid occlusal interference
- Compound contour base anatomically conforms to the tooth and provides proper fit

WELDABLE CONVERTIBLE ROTH® TYPE (4 pcs/kit)		
ROTH® TYPE 1 <sup>st</sup> MOLAR	0.022"	Designation
	BT612-0110	upper 1 tube / lower 1 tube
	BT612-0120	upper 2 tube / lower 1 tube
	BT612-0140	upper / lower double rectangular
	BT612-0130	upper 3 / lower double rectangular
	BT612-0160	upper 3 tube / lower 2 tube




WELDABLE CONVERTIBLE MBT® TYPE (4 pcs/kit)		
MBT® TYPE 1 <sup>st</sup> MOLAR	0.022"	Designation
	BT622-0110	upper 1 tube / lower 1 tube
	BT622-0120	upper 2 tube / lower 1 tube
	BT622-0140	upper / lower double rectangular
	BT622-0130	upper 3 / lower double rectangular
	BT622-0160	upper 3 tube / lower 2 tube




WELDABLE NON-CONVERTIBLE ROTH® TYPE (4 pcs/kit)		
ROTH® TYPE 1 <sup>st</sup> MOLAR	0.022"	Designation
	BT612-0010	upper 1 tube / lower 1 tube
	BT612-0020	upper 2 tube / lower 1 tube
ROTH® TYPE 2 <sup>nd</sup> MOLAR	0.022"	Designation
BT712-0010	upper 1 tube / lower 1 tube	

WELDABLE NON-CONVERTIBLE MBT® TYPE (4 pcs/kit)		
MBT® TYPE 1 <sup>st</sup> MOLAR	0.022"	Designation
	BT622-0010	upper 1 tube / lower 1 tube
	BT622-0020	upper 2 tube / lower 1 tube
MBT® TYPE 2 <sup>nd</sup> MOLAR	0.022"	Designation
BT722-0010	upper 1 tube / lower 1 tube	









## NITI SUPER ELASTIC ARCHWIRES ROUND AND RECTANGULAR

SUPER ELASTIC ARCHWIRES ROUND			
 SQUARE	.012"	UPPER	Super Elastic Niti Wires Round Upper / Lower 10 pcs/pack
	.014"		
	.016"		
 OVOID	.018"		
	.020"		
	.012"		
.014"			
.016"			
 TAPERED	.018"		
	.020"		



SUPER ELASTIC ARCHWIRES RECTANGULAR					
 SQUARE	.016" x .016"	UPPER	Super Elastic Niti Wires Rectangular Upper / Lower 10 pcs/pack		
	.016" x .022"				
	.017" x .022"				
	.017" x .025"				
	.018" x .022"				
 OVOID	.018" x .025"				
	.019" x .025"				
	.021" x .025"				
	.016" x .016"			LOWER	Nitinol Coated Ivory Super Elastic Niti Wire Rectangular Upper / Lower 1 piece/pack
	.016" x .022"				
.017" x .022"					
.017" x .025"					
 TAPERED	.018" x .022"				
	.018" x .025"				
	.019" x .025"				
	.021" x .025"				



## NITI HEAT ACTIVATED ARCHWIRES ROUND AND RECTANGULAR

HEAT ACTIVATED ARCHWIRES ROUND			
 SQUARE	.012"	UPPER	Heat Activated Niti Wires Round Upper / Lower 5 pcs/pack
	.014"		
	.016"		
 OVOID	.018"		
	.020"		
	.012"		
.014"			
.016"			
 TAPERED	.018"		
	.020"		


HEAT ACTIVATED ARCHWIRES RECTANGULAR				
 SQUARE	.016" x .016"	UPPER	Heat Activated Niti Wires Rectangular Upper / Lower 5 pcs/pack	
	.016" x .022"			
	.017" x .022"			
	.017" x .025"			
	.018" x .022"			
 OVOID	.018" x .025"			
	.019" x .025"			
	.021" x .025"			
	.016" x .016"			LOWER
	.016" x .022"			
.017" x .022"				
.017" x .025"				
 TAPERED	.018" x .022"			
	.018" x .025"			
	.019" x .025"			
	.021" x .025"			


## NITI REVERSE CURVE ARCHWIRES ROUND AND RECTANGULAR


REVERSE CURVE ARCHWIRES ROUND			
 UPPER	.012"	Reverse Curve Niti Wires Round Upper / Lower 2 pcs/pack	
	.014"		
	.016"		
	.018"		
	.020"		
 LOWER	.012"		
	.014"		
	.016"		
	.018"		
	.020"		

REVERSE CURVE ARCHWIRES RECTANGULAR			
 UPPER	.016" x .016"	Reverse Curve Niti Wires Rectangular Upper / Lower 2 pcs/pack	
	.016" x .022"		
	.017" x .022"		
	.017" x .025"		
	.018" x .022"		
	.018" x .025"		
	.019" x .025"		
	.021" x .025"		
	 LOWER		.016" x .016"
			.016" x .022"
.017" x .022"			
.017" x .025"			
.018" x .022"			
.018" x .025"			
.019" x .025"			
.021" x .025"			




## NITI WIRES AND SPRINGS




SUPER ELASTIC WIRES 5 M			
	ROUND	.012"	Super Elastic Niti Wires Round 5M
		.014"	
		.016"	
		.018"	
		.020"	

OPEN SPRINGS			
 180 mm	.010"	.012"	Open Springs 2 pcs/pack




CLOSE SPRINGS				
	6 mm	.010"	.012"	Close Springs 10 pcs/pack
	9 mm			
	12 mm			
	15 mm			




## NITI SUPER ELASTIC ARCHWIRES ROUND AND RECTANGULAR LOW FRICTION

NITI SUPER ELASTIC ARCHWIRES ROUND LOW FRICTION			
 SQUARE	.012"	UPPER	Super Elastic Niti Wires <b>Round</b> <b>Low Friction</b> Upper / Lower 10 pcs/pack
	.014"		
	.016"		
 OVOID	.018"		
	.020"		
	.012"		
.014"			
 TAPERED	.016"		
	.018"		
	.020"		



NITI SUPER ELASTIC ARCHWIRES RECTANGULAR LOW FRICTION			
 SQUARE	.016" x .016"	UPPER	Super Elastic Niti Wires <b>Rectangular</b> <b>Low Friction</b> Upper / Lower 10 pcs/pack
	.016" x .022"		
	.017" x .022"		
	.017" x .025"		
	.018" x .022"		
.018" x .025"			
 OVOID	.019" x .025"		
	.021" x .025"		
	.016" x .016"	LOWER	
	.016" x .022"		
.017" x .022"			
 TAPERED	.017" x .025"		
	.018" x .022"		
	.018" x .025"		
	.019" x .025"		
	.021" x .025"		




## NITI HEAT ACTIVATED ARCHWIRES ROUND AND RECTANGULAR LOW FRICTION

NITI HEAT ACTIVATED ARCHWIRES ROUND LOW FRICTION			
 SQUARE	.012"	UPPER	Heat Activated Niti Wires <b>Round</b> <b>Low Friction</b> Upper / Lower 5 pcs/pack
	.014"		
	.016"		
 OVOID	.018"		
	.020"		
	.012"		
.014"			
 TAPERED	.016"		
	.018"		
	.020"		


NITI HEAT ACTIVATED ARCHWIRES RECTANGULAR LOW FRICTION			
 SQUARE	.016" x .016"	UPPER	Heat Activated Niti Wires <b>Rectangular</b> <b>Low Friction</b> Upper / Lower 5 pcs/pack
	.016" x .022"		
	.017" x .022"		
	.017" x .025"		
	.018" x .022"		
.018" x .025"			
 OVOID	.019" x .025"		
	.021" x .025"		
	.016" x .016"	LOWER	
	.016" x .022"		
.017" x .022"			
 TAPERED	.017" x .025"		
	.018" x .022"		
	.018" x .025"		
	.019" x .025"		
	.021" x .025"		


## STAINLESS STEEL ARCHWIRES ROUND AND RECTANGULAR

STAINLESS STEEL ARCHWIRES ROUND			
 OVOID	.012"	UPPER	Stainless Steel Wires <b>Round</b> [ only available in Ovoid and Natural ] Upper / Lower 10 pcs/pack
	.014"		
	.016"		
	.018"		
	.020"		
 NATURAL	.012"		
	.014"		
	.016"		
	.018"		
	.020"		

STAINLESS STEEL ARCHWIRES RECTANGULAR			
 OVOID	.016" x .016"	UPPER	Stainless Steel Wires <b>Rectangular</b> [ only available in Ovoid and Natural ] Upper / Lower 10 pcs/pack
	.016" x .022"		
	.017" x .022"		
	.017" x .025"		
	.018" x .022"		
	.018" x .025"		
 NATURAL	.019" x .025"	LOWER	
	.021" x .025"		
	.016" x .016"		
	.016" x .022"		
	.017" x .022"		
	.017" x .025"		
 NATURAL	.018" x .022"		
	.018" x .025"		
	.019" x .025"		
	.021" x .025"		

## CuNITI SUPER ELASTIC ARCHWIRES with STOPPERS

CuNITI SUPER ELASTIC ARCHWIRES ROUND W STOPPERS		
	.013"	CuNiti Super Elastic Wires <b>Round</b> 1 piece/pack
	.014"	
	.016"	

CuNITI SUPER ELASTIC ARCHWIRES RECTANGULAR W STOPPERS		
	.014" x .025"	CuNiti Super Elastic Wires <b>Rectangular</b> 1 piece/pack
	.016" x .025"	
	.018" x .025"	

## INTRAORAL ACCESSORIES



### LINGUAL BUTTONS

- Semi-Flat
- Round base

IA04-001 (WELDABLE)  
IA04-051 (BONDABLE)  
10 pcs/pack



### SHORT CRIMPABLE HOOKS

- Easily place onto any archwire inside or outside the mouth
- Secure anchor for any auxiliary attachment
- Bi-Directional hook for easy elastic hanging

IA05-001  
10 pcs/pack



### LONG CRIMPABLE HOOKS

- Easily place onto any archwire inside or outside the mouth
- Secure anchor for any auxiliary attachment
- Bi-Directional hook for easy elastic hanging

IA05-101  
10 pcs/pack



### LINGUAL BUTTONS

- Semi-Flat
- Elliptic base

IA04-101 (WELDABLE)  
IA04-151 (BONDABLE)  
10 pcs/pack




### LONG CURVED CRIMPABLE HOOKS

- Excellent for use with TAD mini screw attachments
- Securely attach to closed coil springs or any auxiliary attachments
- Can be placed onto archwires inside or outside of the mouth and fits archwires up to .021"x.025"

IA05-201 (LEFT)  
IA05-202 (RIGHT)  
10 pcs/pack

### SPIRAL CRIMPABLE HOOKS



IA05-301 (LEFT)  
IA05-302 (RIGHT)  
10 pcs/pack



### HI02-001 ARCH TURRET STAINLESS STEEL

- Torque 0° 5 slots
- Slot sizes: 0.016"; 0.017"; 0.018"; 0.019"; 0.021";



### HI02-101 ARCH TURRET STAINLESS STEEL

- Torque 5°, 7°, 6 slots
- Slot sizes: 5° : 0.017"; 0.018"; 0.019"  
7° : 0.017"; 0.018"; 0.019"

## HAND INSTRUMENTS



### BRACKET POSITIONING GAUGE

HI01-101 (2.0 mm - 2.5 mm)  
HI01-102 (3.0 mm - 3.5 mm)  
HI01-103 (4.0 mm - 4.5 mm)  
HI01-104 (5.0 mm - 5.5 mm)  
HI01-105 (6.0 mm - 6.5 mm)



### HI01-001 BRACKET POSITIONING CROSS GAUGE



### HI06-001 FORCE GAUGE

- Precisely measures forces released by elastics, archwires, and coil springs. Double-ended calibrated shaft
- Max.16 ounces, laser marked



### RETAINER BOXES

- Size: W75 mm x D61 mm x H28 mm
- Colours: Assorted- Blue, Lime, Green, Yellow, Light Blue, White and Lavender


600 Retainer Boxes  
Assorted Colours  
10 pcs/pack



### ACRYLIC ORTHO CHAIN DISPENSER

- Long lasting dispenser that will fit the decor of any office
- Chain is dispensed easily with an interior roller
- Holds up to 15 spools


## LIGATURE TIES LATEX FREE



Ligature Ties  
26 x 40  
(1,040 pcs)

COLOURS		
LAKE BLUE	PURPLE	YELLOW
BLUE	PINK	WHITE
GREEN	RED	SILVER GREY
PEARL GREEN	ORANGE	CLEAR

## ELASTOLINK CHAIN LATEX FREE



Elastolink Chain  
5 M/spool

COLOURS		
LAKE BLUE	PURPLE	YELLOW
BLUE	PINK	WHITE
GREEN	RED	SILVER GREY
PEARL GREEN	ORANGE	CLEAR



### 'S' STICK SEPARATORS LATEX FREE

ES-S01 Separators Blue  
70 x 15 (1,050 pcs)

## ELASTICS LATEX

50 x 100  
(5,000 pcs)



Type	1/8" (3.2 mm)	3/16" (4.6 mm)	1/4" (6.4 mm)	5/16" (7.9 mm)	3/8" (9.5 mm)
LIGHT 2.5oz/71 gm	L01	L02	L03	L04	L05
MEDIUM 3.5oz/99 gm	M01	M02	M03	M04	M05
HEAVY 6.5oz/184 gm	H01	H02	H03	H04	H05

## ORTHODONTIC PLIERS



**HI07-018 Light Wire Cutters Standard**  
Hardness: HRC 60-65. Cutting edge used for ligature wires or elastic rings' diameter up to 0.021" (0.30 mm) (can not be used for cutting arch wires)




**HI08-018 Distal End Cutters Standard**  
Hardness: HRC 60-65. Cutting edge used to cut the long distal end of arch wires and safely holding the cut part. Arch wires up to 0.021" x 0.025" (0.53 mm x 0.64 mm)




with flush end

**HI09-018 Distal End Cutters Standard**  
Hardness: HRC 60-65. Cutting edge used to cut the distal end of arch wires and safely holding the cut part, the distal end of arch wire will be flush with tube slot's distal end. Arch wires up to 0.021" x 0.025" (0.53 mm x 0.64 mm)



**HI10-018 Light Wire Pliers**  
Hardness of tip: HRC 60-65. Cutting edge used to make different fine loops for arch wires, and to make spring loops. Wires, diameter up to 0.030" (0.76 mm)



**HI11-018 Light Wire Pliers**  
Hardness of tip: HRC 60-65. Cutting edge used to make different fine loops for arch wires, and to make spring loops. Tip with cutting edge for wire cutting



**HI12-018 Heavy Cutters**  
Hardness of tip: HRC 60-65. Cutting edge used to cut Australia wires, NiTi wires' stainless steel wires. Wires' diameter up to 0.030" (0.76 mm)



**HI14-018 V Stop Pliers**  
Hardness of tip: HRC 60-65. To make V stop along arch wires to prevent any movement



**HI16-018 Multi-Function Pliers**  
Hardness of tip: HRC 52-55  
Used to bend, flatten and cut wires



**HI17-018 Torque Pliers**  
Hardness of tip: HRC 60-65. Cutting edge used in pair to make torques on rectangle wires, thickness of tips: 0.050" (1.27 mm)



**HI18-018 Tweed Loop Pliers**  
Hardness of tip: HRC 60-65. Cutting edge used to make small loops, distal stops and Ω loops on arch wires



**HI19-018 Band Remover**  
Hardness of tip: HRC 52-55. Used to remove bands and other bondable attachment



HI19-02  
Stopper



anterior

**HI20-018 Bracket Removing Pliers**  
Hardness of tip: HRC 52-55. Used to remove anterior brackets



posterior

**HI21-018 Bracket Removing Pliers Standard**  
Hardness of tip: HRC 52-55  
Used to remove posterior brackets




**20-077 Adam Pliers Flat / Flat**  
Tungsten Carbide Tips  
Hardness of tip: HRC 60-65.  
General applicator



**20-074 Adam Pliers Flat / Round**  
Tungsten Carbide Tips  
Hardness of tip: HRC 60-65.  
General applicator



**HI13-018 NiTi Cinch Back Pliers Standard**  
Hardness of tip: HRC 60-65. Cutting edge used to cinch back distal end of arch wires to avoid any injury of soft tissue. NiTi wires' diameter up to 0.021" (0.53 mm)




**HI22-018 Weingart Pliers**  
Hardness of tip: HRC 52-55.  
Used to cinch back distal end of arch wires to avoid any injury of soft tissue




**HI23-018 Crimpable Hook Pliers Standard**  
Hardness of tip: HRC 52-55.  
To fix the crimpable hooks along arch wires




**HI24-018 Loop Bending Pliers**  
Hardness of tip: HRC 52-55. To make different small loops for multi-loop edgewise arch wires, diameter of loops: 2.00, 3.20, 4.20 mm




**HI25-018 HOW's Pliers**  
Hardness of tip: HRC 52-55. To hold arch wires when guiding into buccal tubes, to hold ligature wires, or to form anterior band



**HI26-018 Kim Pliers**  
Hardness of tip: HRC 52-55  
To bend multi-loop edgewise wires, with cutter on tips



**HI27-018 Nance Loop Bending Pliers**  
Hardness of tip: HRC 52-55. To make vertical loops of different height along arch wires, 3.0, 4.0, 5.0, 6.0mm. Wires diameter up to: 0.022" (0.56 mm)



**HI28-018 Band Removing Pliers**  
Hardness of tip: HRC 52-55. Used to remove bond on the crown without any injury on enamel




HI28-02  
Stopper



**HI29-018 Three Jaw Pliers**  
Hardness of tip: HRC 52-55. To make V loops and bend arch wires. Wires' diameter up to 0.030" (0.76 mm)



**HI30-018 Omega Loop Bending Pliers**  
Hardness of tip: HRC 52-55. To make fine Ω loops. Three different loops: diameter 0.063", 0.075", 0.10" (1.6 mm, 2.0 mm, 2.6 mm)



**HI31-018 Elastic Pliers**  
Hardness of tip: HRC 52-55  
For expanding elastic separating rings before separating molar gap




**ONE YEAR WARRANTY** (except scaler tips)

MaxPiezo 7+ (EMS Type) Portable LED Scaler Set  
MaxPiezo DS7+ (SATELEC Type) Portable LED Scaler Set

## REFINE MaxPiezo 7+ / DS7+ PORTABLE LED SCALERS

- Automatic frequency tracking, real-time feedback
- Extra large, capacity 600 ml water tank for water, chlorhexidine, hydrogen peroxide, sodium hypochlorite etc...
- Power and water controlled by touch panel, easy to clean for infection control
- Detachable handpiece
- Frequency: 28 kHz ± 3 kHz
- Complete with 8 tips
- Size: W180 mm x D180 mm x H150 mm
- Weight: 1.6 Kg (unit only)
- Colour: White



1. EMS is a registered trademark of Electro Medical System S.A. Switzerland  
2. SATELEC is a registered trademark of Satelec/Acteon Group France



**ONE YEAR WARRANTY** (except scaler tips)

MaxPiezo 6+ (EMS Type) LED Scaler Set  
MaxPiezo DS6+ (SATELEC Type) LED Scaler Set

## REFINE MaxPiezo 6+ / DS6+ LED SCALERS

- Automatic frequency tracking, real-time feedback
- Power and water controlled by touch panel, easy to clean for infection control
- Detachable handpiece
- Frequency: 28 kHz ± 3 kHz
- Complete with 6 tips
- Size: W155 mm x D168 mm x H80 mm
- Weight: 0.6 Kg (unit only)
- Colour: White



1. EMS is a registered trademark of Electro Medical System S.A. Switzerland  
2. SATELEC is a registered trademark of Satelec/Acteon Group France



**ONE YEAR WARRANTY** (except scaler tips)

N3000L (EMS Type) Built-In LED Scaler Set  
V3000L (SATELEC Type) Built-In LED Scaler Set

## REFINE N3000L / V3000L BUILT-IN LED SCALERS

- Automatic frequency tracking, real-time feedback
- Detachable handpiece
- Frequency: 28 KHz ± 3 KHz
- Complete with 6 tips



1. EMS is a registered trademark of Electro Medical System S.A. Switzerland  
2. SATELEC is a registered trademark of Satelec/Acteon Group France

## REFINE SCALER HANDPIECES



HP-3H Non-LED Scaler Handpiece (EMS<sup>1</sup> Type)



HS-7H Non-LED Scaler Handpiece (SATELEC<sup>2</sup>, NSK<sup>3</sup> Type)



HP-5L LED Scaler Handpiece (EMS<sup>1</sup> Type)



HS-7L LED Scaler Handpiece (SATELEC<sup>2</sup>, NSK<sup>3</sup> Type)

1. EMS is a registered trademark of Electro Medical System S.A. Switzerland  
 2. SATELEC is a registered trademark of Satelec/Acteon Group France  
 3. NSK is a registered trademark of Nakanishi Inc. Japan

## REFINE SCALER TORQUE WRENCHES



**TW-1L**  
For REFINE,  
EMS<sup>1</sup>, SATELEC<sup>2</sup>,  
NSK<sup>3</sup> Type scaler



**TS-1L**  
for SIRONA<sup>4</sup>  
Type scaler



**TK-1L**  
for KAVO<sup>5</sup>  
SonicFlex  
Type scaler

1. EMS is a registered trademark of Electro Medical System S.A. Switzerland  
 2. SATELEC is a registered trademark of Satelec/Acteon Group France  
 3. NSK is a registered trademark of Nakanishi Inc. Japan  
 4. SIRONA is a registered trademark of Sirona Dental System  
 5. KAVO is a registered trademark of Kavo Dental GmbH Germany



3H Non-LED (EMS<sup>1</sup> Type)  
Handpiece Hose



7H Non-LED (SATELEC<sup>2</sup> Type)  
Handpiece Hose



5L LED (EMS<sup>1</sup> Type)  
Handpiece Hose



7L LED (SATELEC<sup>2</sup> Type)  
Handpiece Hose



## REFINE SCALER HANDPIECE HOSES

1. EMS is a registered trademark of Electro Medical System S.A. Switzerland  
 2. SATELEC is a registered trademark of Satelec/Acteon Group France

## REFINE ENDO WRENCH



Multi function wrench  
for both endo chuck & bur

TW-A



## REFINE SOLENOID VALVE



## REFINE FOOT SWITCH



## REFINE POWER SUPPLY

## REFINE SCALER TIP HOLDERS

NO. 1



NO. 2



## E04 ULTRASONIC FILES NiTi alloy

- For use with ultrasonic scaler
- Excellent tenacity
- Length: 33 mm

E04 Ultrasonic Files  
Single Size  
6 pcs/pack




E04 ULTRASONIC FILES NiTi alloy	
#15 ○	#30 ●
#20 ●	#35 ●
#25 ●	#40 ●


## REFINE SCALER TIPS STAINLESS STEEL

FOR EMS<sup>1</sup>, SATELEC<sup>2</sup>, NSK<sup>3</sup> and SIRONA<sup>4</sup> REFINE Type Scalers


### SCALING




**G1 (EMS<sup>1</sup> Type)**  
**GD1 (SATELEC<sup>2</sup>, NSK<sup>3</sup> Type)**  
**GS1 (SIRONA<sup>4</sup> Type)**  
 Supragingival deposits  
 in all quadrants




**G4 (EMS<sup>1</sup> Type)**  
**GD4 (SATELEC<sup>2</sup>, NSK<sup>3</sup> Type)**  
**GS4 (SIRONA<sup>4</sup> Type)**  
 Supragingival and  
 interproximal deposits




**G2 (EMS<sup>1</sup> Type)**  
**GD2 (SATELEC<sup>2</sup>, NSK<sup>3</sup> Type)**  
**GS2 (SIRONA<sup>4</sup> Type)**  
 Heavy supragingival  
 deposits on anterior  
 teeth



**G5 (EMS<sup>1</sup> Type)**  
**GD5 (SATELEC<sup>2</sup>, NSK<sup>3</sup> Type)**  
**GS5 (SIRONA<sup>4</sup> Type)**  
 Supragingival deposits  
 in all quadrants




**G3 (EMS<sup>1</sup> Type)**  
**GD3 (SATELEC<sup>2</sup>, NSK<sup>3</sup> Type)**  
**GS3 (SIRONA<sup>4</sup> Type)**  
 Supragingival deposits  
 in all quadrants, including the  
 interproximal and sulcus areas




**G6 (EMS<sup>1</sup> Type)**  
**GD6 (SATELEC<sup>2</sup>, NSK<sup>3</sup> Type)**  
**GS6 (SIRONA<sup>4</sup> Type)**  
 Voluminous supragingival  
 deposits


### ENDO




**E1 (EMS<sup>1</sup> Type)**  
**ED1 (SATELEC<sup>2</sup>, NSK<sup>3</sup> Type)**  
**ES1 (SIRONA<sup>4</sup> Type)**  
 120°file holder for anterior  
 root canal irrigation




**E4 (EMS<sup>1</sup> Type)**  
**ED4 (SATELEC<sup>2</sup>, NSK<sup>3</sup> Type)**  
**ES4 (SIRONA<sup>4</sup> Type)**  
 Removing root filling obstructions  
 during root canal retreatment




**E2 (EMS<sup>1</sup> Type)**  
**ED2 (SATELEC<sup>2</sup>, NSK<sup>3</sup> Type)**  
**ES2 (SIRONA<sup>4</sup> Type)**  
 95°file holder for molar teeth  
 root canal irrigation



**E4D (EMS<sup>1</sup> Type)**  
**ED4D (SATELEC<sup>2</sup>, NSK<sup>3</sup> Type)**  
**ES4D (SIRONA<sup>4</sup> Type)**  
 Removing calcification and  
 filling material during  
 root canal retreatment




**E3 (EMS<sup>1</sup> Type)**  
**ED3 (SATELEC<sup>2</sup>, NSK<sup>3</sup> Type)**  
**ES3 (SIRONA<sup>4</sup> Type)**  
 Gutta-percha melting




**E3D (EMS<sup>1</sup> Type)**  
**ED3D (SATELEC<sup>2</sup>, NSK<sup>3</sup> Type)**  
**ES3D (SIRONA<sup>4</sup> Type)**  
 Diamond coated instrument for  
 root canal enlarging




**E5 (EMS<sup>1</sup> Type)**  
**ED5 (SATELEC<sup>2</sup>, NSK<sup>3</sup> Type)**  
**ES5 (SIRONA<sup>4</sup> Type)**  
 Removing root filling obstructions  
 during root canal retreatment



**P4D (EMS<sup>1</sup> Type)**  
**PD4D (SATELEC<sup>2</sup>, NSK<sup>3</sup> Type)**  
**PS4D (SIRONA<sup>4</sup> Type)**  
 Diamond coated instrument for  
 root canal enlarging



**E14 (EMS<sup>1</sup> Type)**  
**ED14 (SATELEC<sup>2</sup>, NSK<sup>3</sup> Type)**  
**ES14 (SIRONA<sup>4</sup> Type)**  
 Removing root filling obstructions  
 during root canal retreatments



**E15 (EMS<sup>1</sup> Type)**  
**ED15 (SATELEC<sup>2</sup>, NSK<sup>3</sup> Type)**  
**ES15 (SIRONA<sup>4</sup> Type)**  
 Removing root filling obstructions  
 during root canal retreatments



**E5D (EMS<sup>1</sup> Type)**  
**ED5D (SATELEC<sup>2</sup>, NSK<sup>3</sup> Type)**  
**ES5D (SIRONA<sup>4</sup> Type)**  
 Removing calcification and  
 filling material during  
 canal retreatment


1. EMS is a registered trademark of Electro Medical System S.A. Switzerland  
 2. SATELEC is a registered trademark of Satelec/Acteon Group France  
 3. NSK is a registered trademark of Nakanishi Inc. Japan  
 4. SIRONA is a registered trademark of Sirona Dental System

## REFINE SCALER TIPS STAINLESS STEEL


### FOR EMS<sup>1</sup>, SATELEC<sup>2</sup>, NSK<sup>3</sup> and SIRONA<sup>4</sup> Type Scalers




**P1 (EMS<sup>1</sup> Type)**  
**PD1 (SATELEC<sup>2</sup>, NSK<sup>3</sup> Type)**  
**PS1 (SIRONA<sup>4</sup> Type)**  
 Shallow pocket subgingival deposits



**P3 (EMS<sup>1</sup> Type)**  
**PD3 (SATELEC<sup>2</sup>, NSK<sup>3</sup> Type)**  
**PS3 (SIRONA<sup>4</sup> Type)**  
 Deep pocket subgingival deposits



**P4 (EMS<sup>1</sup> Type)**  
**PD4 (SATELEC<sup>2</sup>, NSK<sup>3</sup> Type)**  
**PS4 (SIRONA<sup>4</sup> Type)**  
 Shallow pocket subgingival deposits



**P3D (EMS<sup>1</sup> Type)**  
**PD3D (SATELEC<sup>2</sup>, NSK<sup>3</sup> Type)**  
**PS3D (SIRONA<sup>4</sup> Type)**  
 Diamond coated instrument for root planning after periodontal flap surgery


### PERIO



**P2L (EMS<sup>1</sup> Type)**  
**PD2L (SATELEC<sup>2</sup>, NSK<sup>3</sup> Type)**  
 Left angle instrument to remove calculus from very narrow inter-root space and furcation



**P2R (EMS<sup>1</sup> Type)**  
**PD2R (SATELEC<sup>2</sup>, NSK<sup>3</sup> Type)**  
 Right angle instrument to remove calculus from very narrow inter-root space and furcation



**P2LD (EMS<sup>1</sup> Type)**  
**PD2LD (SATELEC<sup>2</sup>, NSK<sup>3</sup> Type)**  
 Removal of calculus from narrow inter-root spaces (diamond-coated)



**P2RD (EMS<sup>1</sup> Type)**  
**PD2RD (SATELEC<sup>2</sup>, NSK<sup>3</sup> Type)**  
 Diamond coated instrument for root planning

### FOR AMDENT<sup>6</sup> Type Scalers



**A1 (AMDENT<sup>6</sup> Type)**  
 Supragingival scaling (#37)



**A2 (AMDENT<sup>6</sup> Type)**  
 Supragingival scaling (#39)

### FOR KAVO SONICflex<sup>5</sup> Type Scalers



**GK1 (KAVO<sup>5</sup> Type)**  
 Supragingival deposits in all quadrants




**GK2 (KAVO<sup>5</sup> Type)**  
 Supragingival heavy calculus and plaque



**GK3 (KAVO<sup>5</sup> Type)**  
 Supragingival interdental calculus and calculus at the neck of the tooth




**GK4 (KAVO<sup>5</sup> Type)**  
 Supragingival deposits, including the interproximal and sulcus areas



**GK5 (KAVO<sup>5</sup> Type)**  
 Subgingival calculus




**GK6 (KAVO<sup>5</sup> Type)**  
 Left angle instrument for removing subgingival calculus



**GK7 (KAVO<sup>5</sup> Type)**  
 Right angle instrument for removing subgingival calculus

### FOR KAVO SONOsoft<sup>5</sup> Type Scalers



**GC1 (KAVO<sup>5</sup> Type)**  
 Supragingival deposits in all quadrants



**GC2 (KAVO<sup>5</sup> Type)**  
 Supragingival interdental calculus and calculus at the neck of the tooth



**PC1 (KAVO<sup>5</sup> Type)**  
 Subgingival calculus, including the interproximal and sulcus areas

### KAVO PiezoLED<sup>5</sup> Type Scalers



**201 (KAVO<sup>5</sup> Type)**  
 Supragingival and interproximal deposits



**202 (KAVO<sup>5</sup> Type)**  
 Supragingival deposits in all quadrants, including the interproximal and sulcus areas



**203 (KAVO<sup>5</sup> Type)**  
 Deep pocket subgingival deposits

1. EMS is a registered trademark of Electro Medical System S.A. Switzerland  
 2. SATELEC is a registered trademark of Satelec/Acteon Group France  
 3. NSK is a registered trademark of Nakanishi Inc. Japan  
 4. SIRONA is a registered trademark of Sirona Dental System  
 5. KAVO, SONICflex, SONOsoft, PiezoLED are registered trademarks of Kavo Dental GmbH Germany  
 6. AMDENT is a registered trademark of Amdent AB Sweden





No.3 (20 mm)  
No.4 (22 mm)  
No.5 (24 mm)

## MOUTH MIRRORS STAINLESS STEEL

**Mouth Mirrors Stem**  
10 pcs/box  
Plain No.3, 4, 5  
Magnifying No.3, 4, 5  
Front Surface No.3, 4, 5



**SOLID HANDLE**

47-050-003 Mirror Handle **Round Stem**

47-050-004 Mirror Handle **Octagonal Stem**

47-050-008 Mirror Handle **with Rule**

**HOLLOW HANDLE**

HH-80-006 Hollow Mirror Handle **XL Stem**

HH-80-005 Hollow Mirror Handle **XL Titanium Coated Stem**

47-051-002 Hollow Mirror Handle **7.5 mm Extra Light**

**MIRROR  
HANDLES  
STAINLESS  
STEEL**

## SCALERS AND CURETTES **HOLLOW HANDLE EXTRA LIGHT**



49-155-001 Hollow Handle Curette **XL Gracey 1-2**

49-155-003 Hollow Handle Curette **XL Gracey 3-4**


49-155-005 Hollow Handle Curette **XL Gracey 5-6**

49-155-007 Hollow Handle Curette **XL Gracey 7-8**

49-155-009 Hollow Handle Curette **XL Gracey 9-10**

49-155-011 Hollow Handle Curette **XL Gracey 11-12**

49-155-013 Hollow Handle Curette **XL Gracey 13-14**



49-175-002 Hollow Handle Curette **XL Columbia 2L-2R**

49-175-004 Hollow Handle Curette **XL Columbia 4L-4R**

49-160-013 Hollow Handle Curette **XL McCall 13S-14S**


49-006-007 Hollow Handle Scaler **XL Sickle H6-H7**

49-005-015 Hollow Handle Scaler **XL Sickle H5-U15**

49-005-033 Hollow Handle Scaler **XL Sickle H5-33**

49-015-030 Hollow Handle Scaler **XL Sickle U15-30**

49-015-033 Hollow Handle Scaler **XL Sickle U15-33**



49-190-006 Hollow Handle Probe **XL Fig 6**

49-190-008 Hollow Handle Probe **XL Fig 8**

49-190-009 Hollow Handle Probe **XL Fig 9**




49-190-016 Hollow Handle Probe **XL Fig 16**

49-190-023 Hollow Handle Probe **XL Fig 23**

**PROBES  
HOLLOW  
HANDLE  
EXTRA LIGHT**

## PERIO PROBES **LASER MARKED HOLLOW HANDLE EXTRA LIGHT**



**DOUBLE ENDED**

49-180-001 Hollow Handle Perio Probe **XL 23/Williams**

49-180-002 Hollow Handle Perio Probe **XL 23/WHO**

49-180-003 Hollow Handle Perio Probe **XL 23/11**

49-180-004 Hollow Handle Perio Probe **XL 23/12**

10-9-8-7-5-3-2-1  
11-5-5-5-5-3-5  
11-8-6-3  
12-9-6-3

**SINGLE ENDED**

49-181-001 Hollow Handle Perio Probe **XL Williams**


49-181-002 Hollow Handle Perio Probe **XL WHO**

49-181-003 Hollow Handle Perio Probe **XL 11**


49-181-004 Hollow Handle Perio Probe **XL 12**

1-2-3-5-7-8-9-10  
3-5-5-5-5-11-5  
3-6-9-11  
3-6-9-12


## AMALGAM INSTRUMENTS HOLLOW HANDLE EXTRA LIGHT



## FILLING INSTRUMENTS HOLLOW HANDLE EXTRA LIGHT



## EXCAVATORS HOLLOW HANDLE EXTRA LIGHT




## SCISSORS

STAINLESS STEEL

TUNGSTEN CARBIDE



NEEDLE HOLDERS  
TUNGSTEN CARBIDE TIPS




NEEDLE HOLDERS  
STAINLESS STEEL

HEMOSTATIC FORCEPS  
STAINLESS STEEL




# Hand Instruments


## TWEEZERS STAINLESS STEEL




## CEMENT SPATULAS STAINLESS STEEL



## AMALGAM INSTRUMENTS STAINLESS STEEL




## PROBES STAINLESS STEEL




## EXCAVATORS STAINLESS STEEL




## FILLING INSTRUMENTS STAINLESS STEEL



## FILLING INSTRUMENTS STAINLESS STEEL



## MATRIX RETAINERS STAINLESS STEEL



## CHEEK RETRACTORS STAINLESS STEEL

20-121-140 Sternberg 14 cm  
20-121-160 Sternberg 16 cm  
20-122-145 Roux-Magnus 14.5 cm



## IMPRESSION TRAYS STAINLESS STEEL PERF-LOCK

UPPER	NO.	WIDTH (mm)	DEPTH (mm)
Extra-Large	60-120-001	76	64
Large	60-120-002	71	60
Medium	60-120-003	71	58
Small	60-120-004	65	54
Extra-Small	60-120-005	65	52

LOWER	NO.	WIDTH (mm)	DEPTH (mm)
Extra-Large	60-130-001	75	57
Large	60-130-002	75	54
Medium	60-130-003	72	54
Small	60-130-004	62	50
Extra-Small	60-130-005	55	45

PERF-PARTRIAL	NO.	WIDTH (mm)	DEPTH (mm)
Upper Anterior	61-180-004	66	38
Lower Anterior	61-180-003	64	36
Quadrant	61-180-002	32	54
Quadrant	61-180-001	32	54



## IMPRESSION TRAYS STAINLESS STEEL PERF-LOCK EDENTULOUS


UPPER	No.	WIDTH (mm)	DEPTH (mm)
Extra-Large	60-160-001	76	70
Large	60-160-002	73	66
Medium	60-160-003	70	66
Small	60-160-004	68	64

LOWER	No.	WIDTH (mm)	DEPTH (mm)
Extra-Large	61-170-001	82	62
Large	61-170-002	70	62
Medium	61-170-003	68	60
Small	61-170-004	70	56

## MIXING SPATULAS AND PLASTER KNIVES



## WAXING UP INSTRUMENTS STAINLESS



## SUNDRY INSTRUMENTS STAINLESS STEEL

Crown Remover





56-300-001



## AMALGAM WELL STAINLESS STEEL

53-030-001  
Amalgam well


# Hand Instruments




## SEMKEN TISSUE FORCEPS




## INJECTION SYRINGES



## EXTRACTION FORCEPS ANATOMICAL HANDLE



ANATOMICAL HANDLE




EXTRACTION FORCEPS  
PEDO  
ANATOMICAL HANDLE

## SURGICAL INSTRUMENTS STAINLESS STEEL



STRAIGHT

CURVED



SPLINTER  
FORCEPS  
STAINLESS  
STEEL

## AUTOCLAVABLE TRAYS PLASTIC




STANDARD COLOURS
Blue
Light Blue
Green
Yellow
White
Light Beige
Lavender

- Fully autoclavable up to 134°C

650 Standard Flat 'B' Tray 34 cm x 24 cm x 2 cm  
 651 Standard Divided 'B' Tray 34 cm x 24 cm x 2 cm  
 652 Mini Tray Flat 'F' Tray 24 cm x 16.5 cm x 2 cm

## AUTOCLAVABLE CABINET TRAYS PLASTIC



- Autoclavable up to 134°C
- Indispensable drawer organizer
- Scientifically designed for maximum convenience and efficiency
- Easy to clean, sterilize and disinfect
- Resists breaking, chipping and chemical attack
- Colours: Light Blue, Blue, White, Yellow, Lavender, Green and Light Beige
- Size: 9.5 cm x 19.5 cm x 2 cm

653 Cabinet Tray

STANDARD COLOURS	Yellow	Light Blue	Light Beige	Blue	White	Green	Lavender
Blue	White	Green	Lavender				



## INSTRUMENT TRAY RACKS

- Sturdy, chromium plated steel construction
- Size: Standard (8 trays) W37 cm x D24 cm x H37 cm  
Mini (6 trays) W28 cm x D17 cm x H29 cm

658 Standard Rack 659 Mini Rack

658 Tray Rack Standard 'B'  
 659 Tray Rack Mini 'F'

## AUTOCLAVABLE INSTRUMENT MAT



655 Instrument Mat

- Fully autoclavable
- Medical grade silicone material
- Size: L172 mm x W50 mm x H10 mm

STANDARD COLOURS	Yellow	Light Blue	Light Beige
Blue	White	Green	Lavender

## AUTOCLAVABLE NEEDLE STATION STAINLESS STEEL



- Holds any dental needle sheath
- Solid stainless construction, fully autoclavable
- Silicone ring base resists sliding
- Size: Dia. D55 mm x H38 mm

031 Needle Station

## AUTOCLAVABLE STERILIZATION CASSETTES PLASTIC



- Sturdy, durable construction
- Fully autoclavable

3881 Small (5 instruments)  
 3882 Medium (8 instruments)

STANDARD COLOURS	Blue	White	Light Blue	Green	Lavender
------------------	------	-------	------------	-------	----------

## AUTOCLAVABLE SILICONE INSTRUMENT CODE RINGS



- Fully autoclavable to 134°C
- Sizes: Small - for small instruments  
Large - for pliers and forceps etc
- Available in single colors: Blue, Light Blue, Green, Yellow, Pink, Red, Orange, Lavender, Black and White

Silicone Instrument Rings  
 657S Small, 50 pcs/box  
 657L Large, 25 pcs/box

## AUTOCLAVABLE BUR/ENDO ORGANIZERS



4810 (78 holes)  
 for both hand/engine  
 endo instruments  
 21 mm - 30 mm

4811 (48 holes)  
 for both hand/engine  
 endo instruments  
 21 mm - 30 mm

4812 (36 holes)  
 for both hand/engine  
 endo instruments  
 21 mm - 30 mm

4813  
 12 endo vials

4818  
 78 FG

4819  
 42 FG + 36 CA

- Smart, flexible design for customize setups
- Fully autoclavable up to 134°C
- Colour: Blue

Endo Stand, W115 mm x D52 mm x H40 mm  
 4810 - 78 Holes  
 4811 - 48 Holes + 4 Endo Vials  
 4812 - 36 Holes + Endo Clean  
 4813 - 12 Endo Vials  
 Bur Stand, W115 mm x D52 mm x H40 mm  
 4818 - 78 FG Holes  
 4819 - 42 FG + 36 CA Holes

Stainless Steel Cover, W115 mm x D52 mm x H15 mm  
 4830 - Sterilization Cover

Bur/Endo Box, (box only)  
 4821 Organizer Box for 1 Stand W120 mm x D55 mm x H56 mm  
 4823 Organizer Box for 3 Stands W160 mm x D120 mm x H56 mm  
 4825 Organizer Box for 5 Stands W265 mm x D120 mm x H56 mm

### PERFORATED BUR/ENDO ORGANIZER BOXES



4823

4825

4821

### STAINLESS STEEL COVER



4830

**Remarks:** \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_





**Premium Plus International Ltd.**

1001, Yuen Long Trading Centre  
No.33, Wang Yip Street West  
Yuen Long, N.T.

**Hong Kong**

Tel: 852-2761-3211

Fax: 852-2761-3210

Email: [mailbox@premiumplushk.com](mailto:mailbox@premiumplushk.com)

Web: [www.premiumplushk.com](http://www.premiumplushk.com)



**Premium Plus GmbH**

Bgm.-Finsterwalder Ring 10  
D-82515 Wolfratshausen

**Germany**

Email: [info@premiumplusde.com](mailto:info@premiumplusde.com)

**Premium Plus UK Ltd.**

2, Knighton Heath Industrial Estate  
847 Ringwood Road, Bournemouth  
Dorset BH11 8NE

**U.K.**

Tel: 44-1202-611011

Fax: 44-1202-611010

Email: [info@premiumplusuk.com](mailto:info@premiumplusuk.com)

Web: [www.premiumplusuk.com](http://www.premiumplusuk.com)

**Premium Plus (Dongguan) Limited**

No. 1 Industrial Area, Tutang,  
Changping, Dongguan,  
Guangdong 523581

**China**

Tel: 86-769-83397277

Fax: 86-769-83818081

Email: [info@premiumpluschina.com](mailto:info@premiumpluschina.com)

**Premium Plus Dental Supplies Inc.**

155 Liberty Avenue  
Brooklyn NY 11212

**U.S.A.**

Tel: 1-718-366-8686

Fax: 1-718-366-1288

Email: [info@premiumplususa.com](mailto:info@premiumplususa.com)

Web: [www.premiumplususa.com](http://www.premiumplususa.com)

**Premium Plus Japan Co., Ltd.**

2nd Floor, 1-4-2, Jonan, Ikedashi  
Osaka, 563-0025

**JAPAN**

Tel: 81-72-735-7101

Fax: 81-72-735-7102

Email: [info@premiumplus.jp](mailto:info@premiumplus.jp)

Web: [www.premiumplus.jp](http://www.premiumplus.jp)