


Premium Plus International

dental products 2022



innovators, designers, manufacturers,
partners and much, much more...

premium
plus™

P.1-2 INTRODUCTION

Who are we
Why choose us
What do we make
Where are we

P.3 Personal Protection

Disposable Face Masks
Disposable Face Shields
Face Shields
Mask Extenders
Safety Goggles

P.4-10 SMALL EQUIPMENT

Apex Locators
Curing Lights
Curing Light Meter
High and Low Speed Handpieces

P.11-17 PROPHY

Prophy Angles
Prophy Brushes
Prophy Cups
Prophy Handpieces
Prophy Points

P.18 SYRINGE TIPS

Disposable Syringe Tips
Disposable Syringe Tip Shields
Metal Syringe Tip

P.19-21 Green Line

Disposable Micro Applicator
Disposable Micro Applicator Tips
Disposable Syringe Tips
Interdental Brushes
PSP Barrier Envelopes

P.22 DISPOSABLE

Disposable Bibs
Disposable Face Masks
Disposable Napkin Holders
Disposable Tray Cover Papers

P.23-24 BARRIERS

Barrier Film
Bite Block Sleeves
Chair Sleeves
Curing Light Guide Sleeves
Curing Light Sleeves
Flowable Composite Sleeves
Handpiece & Pen Sleeves
Headrest Covers Non Woven
Headrest Sleeves

P.23-24 BARRIERS

Impression Gun Sleeves
Instrument Tray Sleeves
Intraoral Camera Sleeves
Keyboard Sleeves
Mouth Mirror Sleeves
PC Mouse Sleeves
PSP Barrier Envelopes
Suction Tube Sleeves
Syringe Sleeves
T-Style Light Handle Sleeves
Tubing Sleeves
X-Ray Head Sleeves
X-Ray Sensor Sleeves

P.25 APPLICATORS

Brush Applicators
Brush Tips and Handle
Composite Gun and Tips
Micro Applicators
Micro Applicator Dispensers
Micro Applicator Tips

P.26 EVACUATION

HVE Suction Tubes
Oral Evacuators
Saliva Ejectors
Surgical Aspirators

P.27-28 IMPRESSION

Alginate and Plaster Spatulas
Autoclavable Impression Trays
Autoclavable Mixing Bowls
Disposable 3-in-1 Trays
Disposable Impression Syringes
Disposable Impression Trays
Disposable Mixing Tips
Intraoral Syringe Tips

P.29 INSTRUMENTS

Aluminium Endo Holder
Endo Dispensers and Holders
Endo Measuring Blocks
Endo Rulers
Endo Sponges
Rubber Dam Frame

P.30-34 SUNDRIES

Acrylic Mixing Cups
Aluminium Bur Blocks
Articulating Mounting Plates
Aspirator Brushes
Autoclavable Bur/Endo Organizers
Autoclavable Cabinet Trays
Autoclavable Cheek Retractors
Autoclavable Instrument Trays

P.30-34 SUNDRIES

Autoclavable Mouth Lippers
Autoclavable Napkin Holders
Autoclavable Prophylaxis Rings
Brass Bur Brush
Cotton Roll Dispenser
Dental Appliance Boxes
Dental Bib Holder and Clip
Dental Mandrels
Disposable Bite Loops
Disposable Cheek Retractors
Disposable Composite Tips
Disposable Curved Tip Syringes
Disposable Dappen Dishes
Disposable Evacuation Traps
Disposable Eye Shields
Disposable Fluoride Trays
Disposable Mixing Pads
Disposable Mixing Wells
Disposable Mouth Lippers
Disposable Mouth Mirrors
Disposable Plastic Wedges
Disposable Prophylaxis Rings
Disposable Sticky Applicators
Fracture Detector
Glass Dappen Dishes
Glass Mixing Slab
Implant Scaler
Jumping Teeth
Magnetic Bur Stands
Mouth Props
Needle Station
Patient Chair Cushions
Patient Tongue Cleaners
Pen/Toothbrush Holder
Plastic Bur Dispensers
Protective Eyewears
Retainer/Mouth Boxes/Denture Baths
Silicone Instrument Mat
Silicone Instrument Rings
Suction Tip Sleeves and Adaptors
X-Ray Film Holders

P.35-36 ACRYLIC PRODUCTS

Barrier Film Dispenser
Bib and C-Fold Towel Dispensers
Chain Dispenser
Composite Dispensers and Organizers
Composite Syringe Holders
Cup Dispensers
Glove and Tissue Dispensers
Grinding Box
Impression Gun Holders
Light Protection Covers
Light Protection Hand Shield
Mask Dispenser
Micro Applicator Dispensers
Pouch Dispensers

Who are we ?



At **Premium Plus**, manufacturing dental products is what we do. Only for the past decade we've not just made dental products... we've made better dental products.

Through innovation, research, investment in manufacturing, technical know-how and over 40 years experience in the dental trade, we've learned how to control and improve every step of the manufacturing process... and that's a huge advantage to our customers.


With a global view of the dental industry and knowledge of the different geographical markets, we have learned to understand what you expect and are committed to far exceeding your requirements.



Why choose us ?

Whether its under our own **Premium Plus** brand, private labels or OEM projects, our emphasis is on quality for not only the products, but also its presentation. We not only manufacture in-house – We innovate and develop products, research raw material selection, manufacture moulds, construct production lines, create and design artwork and packaging for our partners, warehouse goods and turn our customers' ideas and requests into reality... providing the complete package.

Our continual investment in purpose-built machinery, technical training and market research along with our flexible approach keeps us at the forefront of dental manufacturers.



What do we make ?

Over the past decade our range has increased to cover small equipment and the majority of plastic and paper disposable dental products. This gives us, and more importantly you, the ability to offer a huge range that includes LED curing lights, apex locators, handpieces, disposable prophy cups and angles, micro applicators, syringe tips, barrier sleeves, saliva ejectors and suction tubes etc.



disposable oral evacuators

disposable saliva ejectors



LED curing lights



disposable syringe tips



disposable trays for implants



disposable 3-in-1 trays



disposable prophy angles




disposable prophy cups



disposable prophy brushes



disposable micro applicators



disposable PSP barrier envelopes

Where are we ?

Premium Plus now supplies to over 90 countries and regions throughout the world and has regional offices in the USA, Japan, Poland and the United Kingdom. Our corporate head office in Hong Kong oversees the production of our range in our purpose-built manufacturing site. It is also home to our dedicated regulatory affairs team that are constantly working on not only maintaining but exceeding standards set by ISO, CE and the FDA. We are proud of the fact that we hold **ISO13485: 2016** certification for the **Quality Management System of Medical Devices** and have all of the class IIa products manufactured by us certified with the **CE₀₁₂₃** mark by **TUV SUD**. We know that with our experience, our integrated range of production processes and our ever-expanding range we have the ability to cover all our customers' requirements.

DISPOSABLE IIR FACE MASKS FLUID RESISTANT LATEX FREE

PERFORMANCE TYPE IIR

EN14683

BFE ≥99% PFE ≥99%

LEVEL 3 BARRIER

ASTM F2100



- Non-irritating
- Fluid resistant, excellent filtration (BFE 99%)
- Permits air to pass through mask to prevent glasses from fogging
- Colours: Blue, Green, Pink, Lavender or Black (ear-loop only)

Level 3 (IIR) Face Masks,
50 pcs/box
2288 Ear-Loop
2288T Tie-On

MASK EXTENDERS



098
Mask Extenders
5 pcs/pack

- Relieve ear fatigue and pain
- Adjustable buckle strip for earloop masks
- Autoclavable soft silicone
- Colours: Light Blue, Light Green, Lavender, Pink, White or Assorted

PERFORMANCE TYPE IIR

EN14683

BFE ≥99% PFE ≥99%

LEVEL 3 BARRIER

ASTM F2100



DISPOSABLE IIR FACE MASKS WITH SHIELDS EAR-LOOP FLUID RESISTANT LATEX FREE

- Non-irritating
- Fluid resistant, excellent filtration (BFE 99%)
- Anti-fog shields
- Permits air to pass through mask to prevent shield from fogging
- Colours: Blue, Green, Pink or Lavender

FS2288 Level 3 (IIR) Face Masks with Shields
Ear-Loop, 25 pcs/box

DISPOSABLE FACE SHIELDS WITH HEADBAND



- For personal protection against droplets and splashes
- Light weight
- Colour: Blue
- Optically clear plastic
- Vented for air flow
- No pressure on temples

1148
Disposable Face Shields
20 pcs/box

PPE standards applied: EN 166:2001, ANSI / ISEA Z87 1-2015




SAVE 50% MATERIALS USED
MORE ENVIRONMENTALLY
FRIENDLY

1168 Clip-on Face Shields, 12 pcs/pack
1168 Clip-on Face Shields, 5 x 12 pcs/box
1170 Glasses Frame (without lenses) with 100 Sleeves, Black or Silver
1171 Disposable Glasses Frame Silicone Sleeves, 100 pcs/pack
[only used for thin metal frames arms]

PPE standards applied: EN 166:2001, ANSI / ISEA Z87 1-2015

DISPOSABLE CLIP-ON FACE SHIELDS

- For personal protection against droplets and splashes
- Lightweight face shields that attach directly to your prescription glasses or glasses frame
- No more bulky and uncomfortable headgear



FACE SHIELDS

- For personal protection against droplets and splashes
- For personal protection against droplets and splashes
- Light weight
- Adjustable design to suit every user
- High quality replaceable shields
- Colours: Blue, White or Lavender

1145 Face Shield complete with 3 shields
1146 Replacement Shields, 5 pcs/pack

PPE standards applied:
EN 166:2001, ANSI / ISEA Z87 1-2015

REFINE SAFETY GOGGLES



EP-01
Safety Goggles

- For personal protection against droplets and splashes
- Light weight
- Soft TPU for all day comfort
- Soft nosepiece
- Compatible with wearing glasses
- Polycarbonate lenses for clear vision
- Flexible head strap

PPE standards applied: EN 166:2001, ANSI / ISEA Z87 1-2015

C01-M / C02-M LED CURING LIGHTS MINI (with USB Charger)




- 360° Rotatable tips
- Automatic compensation circuit for constant power output
- Fully covered membrane panel for easy cleaning and infection control
- Light weight: approx. 85 gms
- Dual modes: FULL / SOFT
- 12 Operation times
- Fully charged battery provides over 100 x 10 second curing cycles

C01-M / C02-M SPECIFICATIONS		
Model	C01-M (F.O. Light Guide)	C02-M (Direct Light Source)
LIGHT TIP	8 mm (dia.)	11 mm (dia.)
LIGHT WAVELENGTH	440~480 nm	440~480 nm
MODE	FULL	1,200 mW/cm ² for 5, 10, 15, 20, 25, 30 seconds±10%
	SOFT	1 st five seconds from 0-900 mW/cm ² next 5, 10, 15, 20, 25 seconds 900 mW/cm ² ±10%
WEIGHT	86 gms	84 gms
LENGTH	230 mm	220 mm

C01-M / C02-M CURING LIGHT DEPTH OF CURE							
MODE	MODEL	C01-M (F.O. Light Guide) Curing Time	Depth Measurements	ISO 4049 Depth Measurements (-50%)	C02-M (Direct Light Source) Curing Time	Depth Measurements	ISO 4049 Depth Measurements (-50%)
FULL		10 seconds	4.3 mm	2.15 mm	10 seconds	5.1 mm	2.55 mm
SOFT		15 seconds	4.6 mm	2.30 mm	15 seconds	4.8 mm	2.40 mm

REMARKS: Above data is based on test results using an A2 shade composite. The light tip was placed above and pointed perpendicular to the composite. Tests were done as per ISO 4049 Standard. The samples were measured and the data was divided by two (50% of the total composite depth cured). ISO 4049 Standard requires that the results for class II restorative materials should be at least 1.5 mm for non-opaque shades and 1.0 mm for opaque shades. All resins and composites have different curing times. It is recommended that you should test your resins and composites performance to determine the curing times required.



C10 LED & HALOGEN CURING LIGHT METER



- Simple, compact design for easy use
- Measures intensity up to 2500 mW/cm²
- Size: 110 mm x 68 mm x 33 mm

C10 Light Meter



C01-M Curing Light (F.O. Light Guide)
 C02-M Curing Light (Direct Light Head)
 C01-M-1 F.O. Light Guide 8 mm
 C02-M-1 Direct Light Head 11 mm
 C01/2-M-1 Replacement Battery
 228 Light Sleeves Mini 500 pcs/box

24 MONTH WARRANTY

C01 LED CURING LIGHTS with fibre optic light guides

C01 SPECIFICATIONS			
Model	C01-C (classic)	C01-S (super power)	C01-D (dual range)
F.O. LIGHT GUIDE	8 mm (dia.)	11 mm (dia.)	
LIGHT WAVELENGTH	440~480 nm	440~480 nm	390~430 nm 440~480 nm
APPLICATION	General	Ortho / General	General / Ortho
MODE	TACK	N/A	900 mW/cm ² for 1 second±10%
	TURBO	1,500 mW/cm ² for 5 seconds ±10%	2,000 mW/cm ² for 4 seconds±10%
	3' + 3'	N/A	2,000 mW/cm ² for 3 + 3 seconds±10% (with 1 second blink in between)
	6' + 6'	N/A	1,500 mW/cm ² for 6 + 6 seconds±10% (with 1 second blink in between)
	FULL	1,200 mW/cm ² for 5, 10, 15, 20 seconds±10%	
	RAMP	1 st five seconds from 0-1,200 mW/cm ²	next 5, 10, 15, 20 seconds 1,200 mW/cm ² ±10%
	SOFT CURE	1 st five seconds from 0-900 mW/cm ²	next 5, 10, 15, 20, 25 seconds 900 mW/cm ² ±10%
WEIGHT	205 gms	200 gms	200 gms
LENGTH	260 mm	265 mm	265 mm

CURING LIGHT MODE APPLICATIONS	
TACK	Tack cure, a very short initial curing for temporary processing before final cure. For example, veneer bonding or temporary cement for in-position and easy removal of excess material.
TURBO	Super fast curing for resin and composite under 2 mm thickness. Also recommended for ceramic brackets, veneer bonding etc.
3' + 3'	3 Seconds Mesial + 3 Seconds Distal for metal brackets resin bonding. (with 1 second blink in between)
6' + 6'	6 Seconds Mesial + 6 Seconds Distal for metal brackets with bonding resin requiring more curing time, and buccal tubes bonding. (with 1 second blink in between)
FULL	For general applications.
RAMP	With first 0-5 seconds ramping for general application.
SOFT CURE	With first 0-5 seconds ramping for more gentle processing, designed for the patients sensitive to the heat.

AUTOCLAVABLE FIBRE OPTIC LIGHT GUIDE CAN BE ROTATED 360°

AUTOMATIC COMPENSATION CIRCUIT FOR CONSTANT POWER OUTPUT

FULLY COVERED MEMBRANE PANEL FOR EASY CLEANING AND INFECTION CONTROL

SIZE: L265 x Dia.27 (mm)

REPLACEABLE BATTERY OPENING

A FULL CHARGE OF BATTERY CAN BE USED MORE THAN 100 TIMES WHEN SET AT 10 SECOND TIME MODE

ON/OFF SWITCH

DISPLAY SCREEN

TIME SWITCH

MODE INDICATORS

MODE SWITCH


C01 Complete Set with Built-In Light Intensity Indicator

063 Protective Hand Light Shield

C01-C Curing Light (Classic) with 8 mm Light Guide
 C01-S Curing Light (Super Power) with 11 mm Light Guide
 C01-D Curing Light (Dual Range) with 11 mm Light Guide
 C01-C-1 8 mm Fibre Optic Light Guide
 C01-S/D-1 11 mm Fibre Optic Light Guide
 C01/2-1 Replacement Battery
 123S Light Sleeves Small (for C01-C) 500 pcs/box
 123L Light Sleeves Large (for C01-S/C01-D) 500 pcs/box



24 MONTH WARRANTY



C01-C / C01-S / C01-D CURING LIGHT DEPTH OF CURE									
MODEL	C01-C (classic) Curing Time	Depth Measurements	ISO 4049 Depth Measurements (-50%)	C01-S (super power) Curing Time	Depth Measurements	ISO 4049 Depth Measurements (-50%)	C01-D (dual range) Curing Time	Depth Measurements	ISO 4049 Depth Measurements (-50%)
TURBO	5 seconds	4.2 mm	2.10 mm	4 seconds	4.6 mm	2.30 mm	4 seconds	4.3 mm	2.15 mm
FULL	10 seconds	4.3 mm	2.15 mm	10 seconds	5.1 mm	2.55 mm	10 seconds	4.4 mm	2.20 mm
RAMP	10 seconds	4.1 mm	2.05 mm	10 seconds	5.0 mm	2.50 mm	10 seconds	4.2 mm	2.10 mm
SOFT CURE	15 seconds	4.6 mm	2.30 mm	15 seconds	4.8 mm	2.40 mm	15 seconds	4.6 mm	2.30 mm

REMARKS: Above data is based on test results using an A2 shade composite. The light guide tip was placed above and pointed perpendicular to the composite. The light guide tip was placed above and pointed perpendicular to the composite. Tests were done as per ISO 4049 Standard. The samples were measured and the data was divided by two (50% of the total composite depth cured). ISO 4049 Standard requires that the results for class II restorative materials should be at least 1.5 mm for non-opaque shades and 1.0 mm for opaque shades. All resins and composites have different curing times. It is strongly recommended that you should test your composite and resin performance to determine the curing times required.


Data is only for reference, always follow the instructions from your resin and composite manufacturers.

- There are many variables (distance and angle between the light tip and composite, area, shade, thickness and type etc.) affecting composite curing times. It is recommended to place and cure a maximum of 2 mm for each composite placement.
- Using a barrier sleeve on the curing light will reduce 5-10% of light intensity.
- The light guide tip should be placed approximately 1-2 mm above the composite or resin and pointed perpendicular to the tooth surface.
- High intensity curing lights produce more heat. This is a direct indication of the power of the light being emitted. When using the device for more than 10 seconds, do not keep the light tip in one position. It is always recommended to use a dental syringe to blow a cold air flow on the area during long periods of curing. This will maintain a lower temperature for the composite and surrounding area, particularly when using rubber dam so as to avoid any rubber dam breakage due to heat built-up.
- Do not expose soft tissue to the light being emitted for more than 2 seconds or burns may occur.
- Avoid the light guide tip coming into contact with any soft tissue in the mouth.

10 colours available:




C02 LED CURING LIGHTS 90° right angle direct light source head




C02 SPECIFICATIONS			
Model	C02-C (classic)	C02-S (super power)	C02-D (dual range)
DIRECT LIGHT HEAD	10 mm (dia.)		
LIGHT WAVELENGTH	440~480 nm	440~480 nm	390~430 nm 440~480 nm
APPLICATION	General	Ortho / General	General / Ortho
MODE	TACK	N/A	900 mW/cm ² for 1 second±10%
	TURBO	1,500 mW/cm ² for 5 seconds ±10%	2,500 mW/cm ² for 3 seconds±10%
	3' + 3'	N/A	2,000 mW/cm ² for 3 + 3 seconds±10% (with 1 second blink in between)
	6' + 6'	N/A	1,500 mW/cm ² for 6 + 6 seconds±10% (with 1 second blink in between)
	FULL	1,200 mW/cm ² for 5, 10, 15, 20 seconds±10%	
	RAMP	1 st five seconds from 0-1,200 mW/cm ²	next 5, 10, 15, 20 seconds 1,200 mW/cm ² ±10%
SOFT CURE	1 st five seconds from 0-900 mW/cm ²	next 5, 10, 15, 20, 25 seconds 900 mW/cm ² ±10%	
WEIGHT	185 gms	185 gms	185 gms
LENGTH	255 mm	255 mm	255 mm

CURING LIGHT MODE APPLICATIONS	
TACK	Tack cure, a very short initial curing for temporary processing before final cure. For example, veneer bonding or temporary cement for in-position and easy removal of excess material.
TURBO	Super fast curing for resin and composite under 2 mm thickness. Also recommended for ceramic brackets, veneer bonding etc.
3' + 3'	3 Seconds Mesial + 3 Seconds Distal for metal brackets resin bonding. (with 1 second blink in between)
6' + 6'	6 Seconds Mesial + 6 Seconds Distal for metal brackets with bonding resin requiring more curing time, and buccal tubes bonding. (with 1 second blink in between)
FULL	For general applications.
RAMP	With first 0-5 seconds ramping for general application.
SOFT CURE	With first 0-5 seconds ramping for more gentle processing, designed for the patients sensitive to the heat.



C02-C (Classic) Curing Light with 10 mm Light Head
 C02-S (Super Power) Curing Light with 10 mm Light Head
 C02-D (Dual Range) Curing Light with 10 mm Light Head
 C02-C-1 (Classic) Direct Light Source Head
 C02-S-1 (Super Power) Direct Light Source Head
 C02-D-1 (Dual Range) Direct Light Source Head
 C01/2-1 Replacement Battery
 C02-3 Heat Protect Rings, 6 pcs/pack
 C02-4 Tack Heat Protective Rings, 6 pcs/pack
 125 Light Sleeves (for C02-C/S/D) 500 pcs/box

24 MONTH WARRANTY



C02-C / C02-S / C02-D CURING LIGHT DEPTH OF CURE									
MODEL	C02-C (classic) Curing Time	Depth Measurements	ISO 4049 Depth Measurements (-50%)	C02-S (super power) Curing Time	Depth Measurements	ISO 4049 Depth Measurements (-50%)	C02-D (dual range) Curing Time	Depth Measurements	ISO 4049 Depth Measurements (-50%)
TURBO	5 seconds	4.4 mm	2.20 mm	3 seconds	4.6 mm	2.30 mm	3 seconds	4.5 mm	2.25 mm
FULL	10 seconds	5.1 mm	2.55 mm	10 seconds	5.3 mm	2.65 mm	10 seconds	5.2 mm	2.60 mm
RAMP	10 seconds	4.8 mm	2.40 mm	10 seconds	5.3 mm	2.65 mm	10 seconds	4.9 mm	2.45 mm
SOFT CURE	15 seconds	5.0 mm	2.50 mm	15 seconds	5.4 mm	2.70 mm	15 seconds	5.1 mm	2.55 mm

REMARKS: Above data is based on test results using an A2 shade composite. The light tip was placed above and pointed perpendicular to the composite. Tests were done as per ISO 4049 Standard. The samples were measured and the data was divided by two (50% of the total composite depth cured). ISO 4049 Standard requires that the results for class II restorative materials should be at least 1.5 mm for non-opaque shades and 1.0 mm for opaque shades. All resins and composites have different curing times. It is strongly recommended that you should test your composite and resin performance to determine the curing times required.

Data is only for reference, always follow the instructions from your resin and composite manufacturers.

- There are many variables (distance and angle between the light tip and composite, area, shade, thickness and type etc.) affecting composite curing times. It is recommended to place and cure a maximum of 2 mm for each composite placement.
- Using a barrier sleeve on the curing light will reduce 5-10% of light intensity.
- The light tip should be placed approximately 1-2 mm above the composite or resin and pointed perpendicular to the tooth surface.
- High intensity curing lights produce more heat. This is a direct indication of the power of the light being emitted. When using the device for more than 10 seconds, do not keep the light tip in one position. It is always recommended to use a dental syringe to blow a cold air flow on the area during long periods of curing. This will maintain a lower temperature for the composite and surrounding area, particularly when using rubber dam so as to avoid any rubber dam breakage due to heat built-up.
- Do not expose soft tissue to the light being emitted for more than 2 seconds or burns may occur.
- Avoid the light tip coming into contact with any soft tissue in the mouth. As an option, you may cover the entire curing light with a disposable sleeve and then place the Heat Protective Ring (included inside the packing box) onto the light tip covered with sleeve as shown below. Ensure the Heat Protective Ring is securely positioned on the curing light before starting the cure.

10 colours available:





L10-1
Straight
Handpiece
1:1 Standard

L20-1L
Latch
CA Handpiece
1:1 Standard

L20-1P
Push Button
CA Handpiece
1:1 Standard

L20-1PM
Push Button Mini Head
CA Handpiece
1:1 Standard

LOW SPEED HANDPIECES AND AIR MOTOR

precisionplus™



L20-4P
Push Button
CA Handpiece
4:1 Reduction

L20-4PM
Push Button
Mini Head
CA Handpiece
4:1 Reduction



L30
Air Motor (2 or 4 holes)
20,000 r.p.m.

12 MONTH WARRANTY

PROPHY LOW SPEED HANDPIECES



L20-1PY (SS)
Screw-In / Snap-On
CA Handpiece
1:1 Standard

L20-4PY (SS)
Screw-In / Snap-On
CA Handpiece
4:1 Reduction



L20-4PY (R)
60° Turn Reciprocating
Push Button
Standard Head
CA Handpiece
4:1 Reduction



L20-4PY (RM)
60° Turn Reciprocating
Push Button
Mini Head
CA Handpiece
4:1 Reduction

precisionplus™

12 MONTH WARRANTY




L38-M4 Prophy Angle Air Motor Complete with 3 Nose Cones
(4 holes) 4:1 Reduction 3,000 r.p.m.



L38-QM Prophy Angle Air Motor complete with 3 Nose Cones
[QM Quick Connect Coupling not included] 4:1 Reduction 3,000 r.p.m.

PROPHY ANGLE HANDPIECE AND MOTORS

precisionplus™



L10-4PY Prophy Angle
Straight Handpiece 4:1 Reduction

12 MONTH WARRANTY

A08 DENTAL APEX LOCATOR

precisionplus™



10 colours available:

- Bright LCD screen with the trajectory of the file being clearly indicated by different colours
- Foldable design with adjustable view angle, saves space and protects the screen
- The unit is designed with a smooth seamless body and membrane buttons to enhance cold sterilization and help avoid cross contamination
- Advanced multiple frequency technology with automatic calibrating ensures accurate measurement
- Measurement results are not affected by tooth type, size of root canal or the internal environment of root canal
- Different beep tones indicate the position of the file in the root canal
- Screen size: 110 mm x 70 mm Overall size: L122 mm x W118 mm x H25 mm
- Accessories autoclavable up to 134°C
- Weight: 300 gms (including rechargeable lithium battery)

A08 Apex Locator
A08-1 Battery
A08-2 Spare Probe

A08-3 Spare Hook
A08-4 Spare Cord
A08-5 Spare Clip

24 MONTH WARRANTY

HIGH SPEED TURBINE WRENCH STANDARD CONNECTIONS



T80 WRENCH


Torque Power: 18 Watts
Speed: 400,000 r.p.m. ±10%

- Ceramic Bearings
- Stainless steel head
- Ergonomic design
- Quiet operation noise level, less than 65 dB
- Autoclavable up to 134°C
- Midwest 4 holes or Borden 2 holes

Turbine Wrench Standard Connection
T80W Turbine (2 or 4 holes)
T80W-A Turbine Cartridge

12 MONTH WARRANTY

HIGH SPEED TURBINE PUSH BUTTON STANDARD CONNECTIONS



T90P HIGH TORQUE PUSH BUTTON

Torque Power: 22 Watts
Speed: 350,000 r.p.m. ±10%

- Anti-retraction
- Ceramic Bearings
- Stainless steel head
- Ergonomic design
- Quiet operation noise level, less than 65 dB
- Autoclavable up to 134°C
- Midwest 4 holes or Borden 2 holes

Turbine Push Button Standard Connection
T90P High Torque Push Button (2 or 4 holes)
T90P-A High Torque PB Cartridge

12 MONTH WARRANTY

HIGH SPEED TURBINE QUICK DISCONNECT



T90+QK HIGH TORQUE PUSH BUTTON

Torque Power: 22 Watts
Speed: 350,000 r.p.m. ±10%

NON LED QUICK COUPLING



QK KaVo¹ Type
with Spray Adjustment

- Anti-retraction
- Ceramic Bearings
- Stainless steel head
- Quiet operation noise level, less than 65 dB
- Ergonomic design
- Autoclavable up to 134°C
- KaVo¹ type fitting

NON LED High Torque Turbine Quick Disconnect
T90+QK KaVo¹ Type (without Coupling)
T90+-A High Torque PB Cartridge

NON LED Quick Coupling
QK KaVo¹ Type (non LED) w Spray Adjustment

24 MONTH WARRANTY

1. KAVO is a registered trademark of Kavo Dental GmbH Germany

LED HIGH SPEED TURBINE QUICK DISCONNECT



T90+LK HIGH TORQUE PUSH BUTTON

Torque Power: 22 Watts
Speed: 350,000 r.p.m. ±10%

LED QUICK COUPLING



LK KaVo¹ Type
with Spray Adjustment


- LED Fibre Optic
- Anti-retraction
- Ceramic Bearings
- Stainless steel head
- Quiet operation noise level, less than 65 dB
- Ergonomic design
- Autoclavable up to 134°C
- KaVo¹ type fitting

LED High Torque Turbine Quick Disconnect
T90+LK KaVo¹ Type (without Coupling)
T90+-A High Torque PB Cartridge

LED Quick Coupling
LK KaVo¹ Type (LED) w Spray Adjustment

24 MONTH WARRANTY

1. KAVO is a registered trademark of Kavo Dental GmbH Germany



DISPOSABLE 60° TURN RECIPROCATING LATEX FREE PROPHY ANGLES

NO SPLATTER - The centrifugal force generated by a normal 360° turn prophy angle results in splattering and mess. Our **NEW 60° turn reciprocating** prophy angles eliminates potential splattering of prophy paste, blood and saliva, keeping the prophy paste in the cup where it belongs.

NO HEAT - Our **NEW 60° turn reciprocating** prophy angles build up much less heat during use, assuring a smooth, cool treatment for your patient.

EFFECTIVE - Due to less heat generated during use, if necessary more pressure may be applied to remove stubborn stains.

FASTER - Thanks to our innovatively designed **NEW 60° turn reciprocating** prophy angles, the prophy cup or brush can be used continuously on the tooth surface, saving valuable chair time.

SAFER - The rotational motion of traditional prophy cups and brushes can cause laceration and bleeding of the gums. Our **NEW 60° turn reciprocating** prophy angles provide a much safer and pain-free prophylaxis treatment.

WARNING: USE OF THESE 60° TURN RECIPROCATING PROPHY ANGLES IS EXTREMELY ADDICTIVE
Once you've tried them, you will never want to use any other prophy angles again

DISPOSABLE 60° TURN RECIPROCATING PROPHY ANGLES STANDARD CUPS LATEX FREE TURBINE BLADES

STANDARD CUPS (8.5 mm)

REGULAR STANDARD CUPS TURBINE BLADES

1,500 - 3,000 r.p.m.

Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	CT3399-60WT	CT3399-60LB	CT3399-60LV	CT3399-60OR
Regular Assorted Colours Order No.	CT3399-60RA Assorted: White, L.Blue, Lavender, Orange			

STANDARD CUPS (8.5 mm)

SOFT STANDARD CUPS TURBINE BLADES

1,500 - 3,000 r.p.m.

Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	CT3399-60PK	CT3399-60LM	CT3399-60GR	CT3399-60DB
Soft Assorted Colours Order No.	CT3399-60SA Assorted: Pink, Lime, Green, D.Blue			

STANDARD CUPS (8.5 mm)

X-SOFT STANDARD CUPS TURBINE BLADES

1,500 - 3,000 r.p.m.

X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	CT3399-60RD	CT3399-60YL
X-Soft Assorted Colours Order No.	CT3399-60XA Assorted: Red, Yellow	

DISPOSABLE 60° TURN RECIPROCATING PROPHY ANGLES SHORT CUPS LATEX FREE TURBINE BLADES

SHORT CUPS (7.0 mm)

REGULAR SHORT CUPS TURBINE BLADES

1,500 - 3,000 r.p.m.

Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	CTS3399-60WT	CTS3399-60LB	CTS3399-60LV	CTS3399-60OR
Regular Assorted Colours Order No.	CTS3399-60RA Assorted: White, L.Blue, Lavender, Orange			

SHORT CUPS (7.0 mm)

SOFT SHORT CUPS TURBINE BLADES

1,500 - 3,000 r.p.m.

Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	CTS3399-60PK	CTS3399-60LM	CTS3399-60GR	CTS3399-60DB
Soft Assorted Colours Order No.	CTS3399-60SA Assorted: Pink, Lime, Green, D.Blue			

SHORT CUPS (7.0 mm)

X-SOFT SHORT CUPS TURBINE BLADES

1,500 - 3,000 r.p.m.

X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	CTS3399-60RD	CTS3399-60YL
X-Soft Assorted Colours Order No.	CTS3399-60XA Assorted: Red, Yellow	

DISPOSABLE 60° TURN RECIPROCATING PROPHY ANGLES LATEX FREE BRUSHES REGULAR / SOFT / X-SOFT

REGULAR WHITE (nylon brush dia 0.2 mm)

1,500 - 3,000 r.p.m.

Regular Type	Flat	Tapered	Mini Cup
Regular White Colour Order No.	C3399-60-R261F	C3399-60-R261T	C3399-60-R269
Regular Assorted Colours Order No.	C3399-60-R261FA C3399-60-R261TA C3399-60-R269A Assorted: White, L.Blue, Lavender, Orange		

SOFT CLEAR (nylon brush dia 0.15 mm)

1,500 - 3,000 r.p.m.


Soft Type	Flat	Tapered	Mini Cup
Soft Clear Colour Order No.	C3399-60-S261F	C3399-60-S261T	C3399-60-S269
Soft Assorted Colours Order No.	C3399-60-S261FA C3399-60-S261TA C3399-60-S269A Assorted: Pink, Lime, Green, D.Blue		

X-SOFT YELLOW (nylon brush dia. 0.10 mm)

1,500 - 3,000 r.p.m.

X-Soft Type	Flat	Tapered	Mini Cup
X-Soft Yellow Colour Order No.	C3399-60-X261F	C3399-60-X261T	C3399-60-X269
X-Soft Assorted Colours Order No.	C3399-60-X261FA C3399-60-X261TA C3399-60-X269A Assorted: Yellow, Red		

PROPHY ANGLE STRAIGHT HANDPIECE 4:1




L10-4PY Prophy Angle Handpiece Short 4:1 Reduction

STANDARD STRAIGHT HANDPIECE 1:1




L10-1 Standard Straight Handpiece 1:1

DISPOSABLE PROPHY ANGLES STANDARD CUPS LATEX FREE RIBBED AND WEBBED



Type	90° Right Angle		105° Contra Angle	
Cup Colours	Regular Lavender	Soft Green	Regular Lavender	Soft Green
Order No.	3399-090RL	3399-090SG	3399-105RL	3399-105SG

DISPOSABLE PROPHY ANGLES BRUSHES LATEX FREE FLAT AND TAPERED



Type	90° Right Angle		105° Contra Angle	
Brushes	Regular Flat	Regular Tapered	Regular Flat	Regular Tapered
Order No.	3399-090F	3399-090T	3399-105F	3399-105T

DISPOSABLE CLEAR PROPHY ANGLES STANDARD CUPS LATEX FREE RIBBED AND WEBBED




Type	90° Right Angle		90° Right Angle	
Cup Colours	Regular Light Blue	Regular Lavender	Soft Green	Soft Pink
Order No.	C3399-090RB	C3399-090RL	C3399-090SG	C3399-090SP



Type	105° Contra Angle		105° Contra Angle	
Cup Colours	Regular Light Blue	Regular Lavender	Soft Green	Soft Pink
Order No.	C3399-105RB	C3399-105RL	C3399-105SG	C3399-105SP

DISPOSABLE CLEAR PROPHY ANGLES SHORT CUPS LATEX FREE RIBBED AND WEBBED




Type	90° Right Angle		90° Right Angle	
Cup Colours	Regular Light Blue	Regular Lavender	Soft Green	Soft Pink
Order No.	CS3399-090RB	CS3399-090RL	CS3399-090SG	CS3399-090SP




Type	105° Contra Angle		105° Contra Angle	
Cup Colours	Regular Light Blue	Regular Lavender	Soft Green	Soft Pink
Order No.	CS3399-105RB	CS3399-105RL	CS3399-105SG	CS3399-105SP

DISPOSABLE CLEAR PROPHY ANGLES STANDARD CUPS LATEX FREE TURBINE BLADES



Type	90° Right Angle		90° Right Angle	
Cup Colours	Regular Light Blue	Regular Lavender	Soft Green	Soft Pink
Order No.	CT3399-090RB	CT3399-090RL	CT3399-090SG	CT3399-090SP



Type	105° Contra Angle		105° Contra Angle	
Cup Colours	Regular Light Blue	Regular Lavender	Soft Green	Soft Pink
Order No.	CT3399-105RB	CT3399-105RL	CT3399-105SG	CT3399-105SP



60° TURN RECIPROCATING PROPHY HANDPIECE 4:1 REDUCTION PUSH BUTTON STANDARD HEAD

NO SPLATTER - The centrifugal force generated by a normal 360° turn prophy handpiece results in splattering and mess. Our **NEW 60° turn reciprocating** prophy handpiece eliminates potential splattering of prophy paste, blood and saliva, keeping the prophy paste in the cup where it belongs.

NO HEAT - Our **NEW 60° turn reciprocating** prophy handpieces build up much less heat during use, assuring a smooth, cool treatment for your patient.

EFFECTIVE - Due to less heat generated during use, if necessary more pressure may be applied to remove stubborn stains.

FASTER - Thanks to our innovatively designed **NEW 60° turn reciprocating** prophy handpiece, the prophy cup or brush can be used continuously on the tooth surface, saving valuable chair time.

SAFER - The rotational motion of traditional prophy cups and brushes can cause laceration and bleeding of the gums. Our **NEW 60° turn reciprocating** prophy handpieces provide a much safer and pain-free prophylaxis treatment.

L20-4PY (R)
60° Turn Reciprocating
Push Button
Standard Head
Prophy Handpiece
4:1 Reduction

METAL SHANKS CUPS (8.5 mm)



REGULAR STANDARD CUPS
RIBBED AND WEBBED

Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	M222WT	M222LB	M222LV	M222OR
Regular Assorted Colours Order No.	M222RA Assorted: White, L.Blue, Lavender, Orange			

METAL SHANKS CUPS (8.5 mm)



SOFT STANDARD CUPS
RIBBED AND WEBBED

Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	M222PK	M222LM	M222GR	M222DB
Soft Assorted Colours Order No.	M222SA Assorted: Pink, Lime, Green, D.Blue			

METAL SHANKS CUPS (8.5 mm)



X-SOFT STANDARD CUPS
RIBBED AND WEBBED

X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	M222RD	M222YL
X-Soft Assorted Colours Order No.	M222XA Assorted: Red, Yellow	

DISPOSABLE LATEX FREE STANDARD PROPHY CUPS CA METAL SHANKS RIBBED AND WEBBED

PLASTIC SHANKS CUPS (8.5 mm)



REGULAR STANDARD CUPS
RIBBED AND WEBBED

Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	P222WT	P222LB	P222LV	P222OR
Regular Assorted Colours Order No.	P222RA Assorted: White, L.Blue, Lavender, Orange			

PLASTIC SHANKS CUPS (8.5 mm)



SOFT STANDARD CUPS
RIBBED AND WEBBED

Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	P222PK	P222LM	P222GR	P222DB
Soft Assorted Colours Order No.	P222SA Assorted: Pink, Lime, Green, D.Blue			

PLASTIC SHANKS CUPS (8.5 mm)



X-SOFT STANDARD CUPS
RIBBED AND WEBBED

X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	P222RD	P222YL
X-Soft Assorted Colours Order No.	P222XA Assorted: Red, Yellow	

DISPOSABLE LATEX FREE STANDARD PROPHY CUPS CA PLASTIC SHANKS RIBBED AND WEBBED

METAL SHANKS SHORT CUPS (7.0 mm)



REGULAR SHORT CUPS
RIBBED AND WEBBED

Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	S222WT	S222LB	S222LV	S222OR
Regular Assorted Colours Order No.	S222RA Assorted: White, L.Blue, Lavender, Orange			

METAL SHANKS SHORT CUPS (7.0 mm)



SOFT SHORT CUPS
RIBBED AND WEBBED

Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	S222PK	S222LM	S222GR	S222DB
Soft Assorted Colours Order No.	S222SA Assorted: Pink, Lime, Green, D.Blue			

METAL SHANKS SHORT CUPS (7.0 mm)



X-SOFT SHORT CUPS
RIBBED AND WEBBED

X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	S222RD	S222YL
X-Soft Assorted Colours Order No.	S222XA Assorted: Red, Yellow	

DISPOSABLE LATEX FREE SHORT PROPHY CUPS CA METAL SHANKS RIBBED AND WEBBED

DISPOSABLE LATEX FREE PROPHY BRUSHES CA REGULAR / SOFT / X-SOFT



Regular Type	Flat	Tapered	Micro	Junior Cup	Mini Cup
Regular White Colour Order No.	R-262F	R-262T	R-265	R-M268N	R-269N
Regular Assorted Colours Order No.	Assorted: White, L.Blue, Lavender, Orange				



Soft Type	Flat	Tapered	Micro	Junior Cup	Mini Cup
Soft Clear Colour Order No.	S-262F	S-262T	S-265	S-268	S-269
Soft Assorted Colours Order No.	Assorted: Pink, Lime, Green, D.Blue				



X-Soft Type	Flat	Tapered	Micro	Junior Cup	Mini Cup
X-Soft Yellow Colour Order No.	X-262F	X-262T	X-265	X-268	X-269
X-Soft Assorted Colours Order No.	Assorted: Red, Yellow				



Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	225WT	225LB	225LV	225OR
Regular Assorted Colours Order No.	225RA Assorted: White, L.Blue, Lavender, Orange			



Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	225PK	225LM	225GR	225DB
Soft Assorted Colours Order No.	225SA Assorted: Pink, Lime, Green, D.Blue			



X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	225RD	225YL
X-Soft Assorted Colours Order No.	225XA Assorted: Red, Yellow	

DISPOSABLE LATEX FREE PROPHY POINTS CA METAL SHANKS [for lingual and interproximal]



Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	T-M222 WT	T-M222 LB	T-M222 LV	T-M222 OR
Regular Assorted Colours Order No.	T-M222RA Assorted: White, L.Blue, Lavender, Orange			



Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	T-M222 PK	T-M222 LM	T-M222 GR	T-M222 DB
Soft Assorted Colours Order No.	T-M222SA Assorted: Pink, Lime, Green, D.Blue			



X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	T-M222 RD	T-M222 YL
X-Soft Assorted Colours Order No.	T-M222XA Assorted: Red, Yellow	

DISPOSABLE LATEX FREE STANDARD PROPHY CUPS CA METAL SHANKS TURBINE BLADES



Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	T-P222 WT	T-P222 LB	T-P222 LV	T-P222 OR
Regular Assorted Colours Order No.	T-P222RA Assorted: White, L.Blue, Lavender, Orange			



Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	T-P222 PK	T-P222 LM	T-P222 GR	T-P222 DB
Soft Assorted Colours Order No.	T-P222SA Assorted: Pink, Lime, Green, D.Blue			



X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	T-P222 RD	T-P222 YL
X-Soft Assorted Colours Order No.	T-P222XA Assorted: Red, Yellow	

DISPOSABLE LATEX FREE STANDARD PROPHY CUPS CA PLASTIC SHANKS TURBINE BLADES



Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	T-S222 WT	T-S222 LB	T-S222 LV	T-S222 OR
Regular Assorted Colours Order No.	T-S222RA Assorted: White, L.Blue, Lavender, Orange			



Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	T-S222 PK	T-S222 LM	T-S222 GR	T-S222 DB
Soft Assorted Colours Order No.	T-S222SA Assorted: Pink, Lime, Green, D.Blue			



X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	T-S222 RD	T-S222 YL
X-Soft Assorted Colours Order No.	T-S222XA Assorted: Red, Yellow	

DISPOSABLE LATEX FREE SHORT PROPHY CUPS CA METAL SHANKS TURBINE BLADES

PROPHY HANDPIECE SNAP-ON / SCREW-IN 4:1 REDUCTION

SNAP-ON SCREW-IN



L20-4PY (SS)
Prophy Handpiece CA
Snap On / Screw In
4:1 Reduction

DISPOSABLE LATEX FREE PROPHY POINTS SCREW-IN [for lingual and interproximal]



REGULAR POINTS SCREW-IN

Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	226WT	226LB	226LV	226OR
Regular Assorted Colours Order No.	226RA Assorted: White, L.Blue, Lavender, Orange			



SOFT POINTS SCREW-IN

Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	226PK	226LM	226GR	226DB
Soft Assorted Colours Order No.	226SA Assorted: Pink, Lime, Green, D.Blue			



X-SOFT POINTS SCREW-IN

X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	226RD	226YL
X-Soft Assorted Colours Order No.	226XA Assorted: Red, Yellow	

DISPOSABLE LATEX FREE PROPHY CUPS SNAP-ON RIBBED AND WEBBED



REGULAR CUPS RIBBED AND WEBBED

Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	220WT	220LB	220LV	220OR
Regular Assorted Colours Order No.	220RA Assorted: White, L.Blue, Lavender, Orange			



SOFT CUPS RIBBED AND WEBBED

Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	220PK	220LM	220GR	220DB
Soft Assorted Colours Order No.	220SA Assorted: Pink, Lime, Green, D.Blue			



X-SOFT CUPS RIBBED AND WEBBED

X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	220RD	220YL
X-Soft Assorted Colours Order No.	220XA Assorted: Red, Yellow	

DISPOSABLE LATEX FREE STANDARD PROPHY CUPS SCREW-IN RIBBED AND WEBBED



REGULAR STANDARD CUPS RIBBED AND WEBBED

Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	221WT	221LB	221LV	221OR
Regular Assorted Colours Order No.	221RA Assorted: White, L.Blue, Lavender, Orange			



SOFT STANDARD CUPS RIBBED AND WEBBED

Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	221PK	221LM	221GR	221DB
Soft Assorted Colours Order No.	221SA Assorted: Pink, Lime, Green, D.Blue			



X-SOFT STANDARD CUPS RIBBED AND WEBBED

X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	221RD	221YL
X-Soft Assorted Colours Order No.	221XA Assorted: Red, Yellow	

DISPOSABLE LATEX FREE SHORT PROPHY CUPS SCREW-IN RIBBED AND WEBBED



REGULAR SHORT CUPS RIBBED AND WEBBED

Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	S221WT	S221LB	S221LV	S221OR
Regular Assorted Colours Order No.	S221RA Assorted: White, L.Blue, Lavender, Orange			



SOFT SHORT CUPS RIBBED AND WEBBED

Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	S221PK	S221LM	S221GR	S221DB
Soft Assorted Colours Order No.	S221SA Assorted: Pink, Lime, Green, D.Blue			



X-SOFT SHORT CUPS RIBBED AND WEBBED

X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	S221RD	S221YL
X-Soft Assorted Colours Order No.	S221XA Assorted: Red, Yellow	



REGULAR WHITE
(nylon brush dia. 0.20 mm)

Regular Type	Flat	Tapered
Regular White Brushes Order No.	R-260F	R-260T



SOFT CLEAR
(nylon brush dia. 0.15 mm)

Soft Type	Flat	Tapered
Soft Clear Brushes Order No.	S-260F	S-260T



X-SOFT YELLOW
(nylon brush dia. 0.10 mm)

X-Soft Type	Flat	Tapered
X-Soft Yellow Brushes Order No.	X-260F	X-260T

DISPOSABLE LATEX FREE PROPHY BRUSHES SNAP-ON REGULAR / SOFT / X-SOFT

DISPOSABLE LATEX FREE PROPHY BRUSHES SCREW-IN REGULAR / SOFT / X-SOFT



REGULAR WHITE (nylon brush dia. 0.20 mm)

Regular Type	Flat	Tapered	Micro	Junior Cup	Mini Cup
Regular White Colour Order No.	R-261F	R-261T	R-265-1	R-268-1	R-269-1
Regular Assorted Colours Order No.	R-261FA R-261TA R-265-1A R-268-1A R-269-1A Assorted: White, L.Blue, Lavender, Orange				



SOFT CLEAR (nylon brush dia. 0.15 mm)

Soft Type	Flat	Tapered	Micro	Junior Cup	Mini Cup
Soft Clear Colour Order No.	S-261F	S-261T	S-265-1	S-268-1	S-269-1
Soft Assorted Colours Order No.	S-261FA S-261TA S-265-1A S-268-1A S-269-1A Assorted: Pink, Lime, Green, D.Blue				



X-SOFT YELLOW (nylon brush dia. 0.10 mm)

X-Soft Type	Flat	Tapered	Micro	Junior Cup	Mini Cup
X-Soft Yellow Colour Order No.	X-261F	X-261T	X-265-1	X-268-1	X-269-1
X-Soft Assorted Colours Order No.	X-261FA X-261TA X-265-1A X-268-1A X-269-1A Assorted: Red, Yellow				



SNAP-ON CUPS (8.0 mm)



REGULAR CUPS
TURBINE BLADES

Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	T-220WT	T-220LB	T-220LV	T-220OR
Regular Assorted Colours Order No.	T-220RA Assorted: White, L.Blue, Lavender, Orange			



SNAP-ON CUPS (8.0 mm)



SOFT CUPS
TURBINE BLADES

Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	T-220PK	T-220LM	T-220GR	T-220DB
Soft Assorted Colours Order No.	T-220SA Assorted: Pink, Lime, Green, D.Blue			



SNAP-ON CUPS (8.0 mm)



X-SOFT CUPS
TURBINE BLADES

X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	T-220RD	T-220YL
X-Soft Assorted Colours Order No.	T-220XA Assorted: Red, Yellow	

DISPOSABLE LATEX FREE PROPHY CUPS SNAP-ON TURBINE BLADES



SCREW-IN CUPS (8.5 mm)



REGULAR STANDARD CUPS
TURBINE BLADES

Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	T-221WT	T-221LB	T-221LV	T-221OR
Regular Assorted Colours Order No.	T-221RA Assorted: White, L.Blue, Lavender, Orange			



SCREW-IN CUPS (8.5 mm)



SOFT STANDARD CUPS
TURBINE BLADES

Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	T-221PK	T-221LM	T-221GR	T-221DB
Soft Assorted Colours Order No.	T-221SA Assorted: Pink, Lime, Green, D.Blue			



SCREW-IN CUPS (8.5 mm)



X-SOFT STANDARD CUPS
TURBINE BLADES

X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	T-221RD	T-221YL
X-Soft Assorted Colours Order No.	T-221XA Assorted: Red, Yellow	

DISPOSABLE LATEX FREE STANDARD PROPHY CUPS SCREW-IN TURBINE BLADES



SCREW-IN SHORT CUPS (7.0 mm)



REGULAR SHORT CUPS
TURBINE BLADES

Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	T-S221WT	T-S221LB	T-S221LV	T-S221OR
Regular Assorted Colours Order No.	T-S221RA Assorted: White, L.Blue, Lavender, Orange			



SCREW-IN SHORT CUPS (7.0 mm)



SOFT SHORT CUPS
TURBINE BLADES

Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	T-S221PK	T-S221LM	T-S221GR	T-S221DB
Soft Assorted Colours Order No.	T-S221SA Assorted: Pink, Lime, Green, D.Blue			



SCREW-IN SHORT CUPS (7.0 mm)



X-SOFT SHORT CUPS
TURBINE BLADES

X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	T-S221RD	T-S221YL
X-Soft Assorted Colours Order No.	T-S221XA Assorted: Red, Yellow	

DISPOSABLE LATEX FREE SHORT PROPHY CUPS SCREW-IN TURBINE BLADES

Prophy Brushes and Cups



60° TURN RECIPROCATING PROPHY HANDPIECE 4:1 REDUCTION PUSH BUTTON MINI HEAD

NO SPLATTER - The centrifugal force generated by a normal 360° turn prophy handpiece results in splattering and mess. Our **NEW 60° turn reciprocating** prophy handpiece eliminates potential splattering of prophy paste, blood and saliva, keeping the prophy paste in the cup where it belongs.

NO HEAT - Our **NEW 60° turn reciprocating** prophy handpieces build up much less heat during use, assuring a smooth, cool treatment for your patient.

EFFECTIVE - Due to less heat generated during use, if necessary more pressure may be applied to remove stubborn stains.

FASTER - Thanks to our innovatively designed **NEW 60° turn reciprocating** prophy handpiece, the prophy cup or brush can be used continuously on the tooth surface, saving valuable chair time.

SAFER - The rotational motion of traditional prophy cups and brushes can cause laceration and bleeding of the gums. Our **NEW 60° turn reciprocating** prophy handpieces provide a much safer and pain-free prophylaxis treatment.

RIBBED AND WEBBED
1,500 - 3,000 r.p.m.

SHORT CUPS (7.0 mm)

REGULAR SHORT CUPS
SHORT SHANKS
RIBBED AND WEBBED

Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	SS-222WT	SS-222LB	SS-222LV	SS-222OR
Regular Assorted Colours Order No.	SS-222RA Assorted: White, L.Blue, Lavender, Orange			

RIBBED AND WEBBED
1,500 - 3,000 r.p.m.

SHORT CUPS (7.0 mm)

SOFT SHORT CUPS
SHORT SHANKS
RIBBED AND WEBBED

Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	SS-222PK	SS-222LM	SS-222GR	SS-222DB
Soft Assorted Colours Order No.	SS-222SA Assorted: Pink, Lime, Green, D.Blue			

RIBBED AND WEBBED
1,500 - 3,000 r.p.m.

SHORT CUPS (7.0 mm)

X-SOFT SHORT CUPS
SHORT SHANKS
RIBBED AND WEBBED

X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	SS-222RD	SS-222YL
X-Soft Assorted Colours Order No.	SS-222XA Assorted: Red, Yellow	

DISPOSABLE LATEX FREE
SHORT
PROPHY CUPS CA
METAL SHANKS
SHORT SHANKS
RIBBED AND WEBBED
[fit mini head CA handpieces only]

TURBINE BLADES
1,500 - 3,000 r.p.m.

SHORT CUPS (7.0 mm)

REGULAR SHORT CUPS
SHORT SHANKS
TURBINE BLADES

Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	SS-T222WT	SS-T222LB	SS-T222LV	SS-T222OR
Regular Assorted Colours Order No.	SS-T222RA Assorted: White, L.Blue, Lavender, Orange			

TURBINE BLADES
1,500 - 3,000 r.p.m.

SHORT CUPS (7.0 mm)

SOFT SHORT CUPS
SHORT SHANKS
TURBINE BLADES

Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	SS-T222PK	SS-T222LM	SS-T222GR	SS-T222DB
Soft Assorted Colours Order No.	SS-T222SA Assorted: Pink, Lime, Green, D.Blue			

TURBINE BLADES
1,500 - 3,000 r.p.m.

SHORT CUPS (7.0 mm)

X-SOFT SHORT CUPS
SHORT SHANKS
TURBINE BLADES

X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	SS-T222RD	SS-T222YL
X-Soft Assorted Colours Order No.	SS-T222XA Assorted: Red, Yellow	

DISPOSABLE LATEX FREE
SHORT
PROPHY CUPS CA
METAL SHANKS
SHORT SHANKS
TURBINE BLADES
[fit mini head CA handpieces only]

DISPOSABLE LATEX FREE SHORT SHANKS PROPHY BRUSHES CA [fit mini head CA handpieces only]

REGULAR WHITE (nylon brush dia. 0.20 mm)

1,500 - 3,000 r.p.m.

SOFT CLEAR (nylon brush dia. 0.15 mm)

1,500 - 3,000 r.p.m.

X-SOFT YELLOW (nylon brush dia. 0.10 mm)

1,500 - 3,000 r.p.m.

Regular Type	Flat	Tapered	Micro	Junior Cup	Mini Cup
Regular White Colour Order No.	SS-R262F	SS-R262T	SS-R265	SS-R268	SS-R269
Regular Assorted Colours Order No.	SS-R262FA SS-R262TA SS-R265A SS-R268A SS-R269A Assorted: White, L.Blue, Lavender, Orange				

Soft Type	Flat	Tapered	Micro	Junior Cup	Mini Cup
Soft Clear Colour Order No.	SS-S262F	SS-S262T	SS-S265	SS-S268	SS-S269
Soft Assorted Colours Order No.	SS-S262FA SS-S262TA SS-S265A SS-S268A SS-S269A Assorted: Pink, Lime, Green, D.Blue				

X-Soft Type	Flat	Tapered	Micro	Junior Cup	Mini Cup
X-Soft Yellow Colour Order No.	SS-X262F	SS-X262T	SS-X265	SS-X268	SS-X269
X-Soft Assorted Colours Order No.	SS-X261FA SS-X261TA SS-X265A SS-X268A SS-X269A Assorted: Red, Yellow				

DISPOSABLE AIR-WATER SYRINGE TIPS (SHORT) EURO TYPE



- Directly interchangeable with American type short syringe tips (syringes which have been converted to use American type short syringe tips)
- No sharp edges for patient comfort
- Unique quad air core design to ensure dry air

Syringe Tips Short, 250 pcs/pack
 028S Syringe Tips **Short** Spectrum
 028W Syringe Tips **Short** White



DISPOSABLE AIR-WATER 'e' SYRINGE TIPS WITHOUT CORE

ready to use, no converter required for American style syringe

- Directly interchangeable with most American style metal air/water syringe tips
- Economically designed without separate water core, in exchange for a split second of moisture when switching from water to air
- Rounded tip end provides maximum patient comfort

Disposable 'e' Air-Water Syringe Tips without core 250 pcs/pack
 0120S Opaque **Colour** Tips
 012W Opaque **White** Tips




DISPOSABLE AIR-WATER 'a' SYRINGE TIPS

ready to use, no converter required for American style syringe

- Directly interchangeable with most American style metal air/water syringe tips
- No sharp edges for patient comfort
- Unique quad air core design to ensure dry air

Disposable 'a' Air-Water Syringe Tips, 250 pcs/pack
 018S **Spectrum** Colour Tips/White Cores
 018C **Clear** Tips/Colour Cores
 018OS **Opaque Colour** Tips/White Cores
 018W **Opaque White** Tips/White Cores



DISPOSABLE TRANSPARENT SYRINGE TIPS

- Directly interchangeable with *Sani-Tip® type transparent tip (syringes which have been converted to use transparent tip)
- Costs less than sterilizing metal tips
- Size: Long(American type) Short(European type)
- Colour: Clear Transparent

Clear Tips, 250 pcs/pack
 011A Standard Length (*American Type Syringe*)
 011E Short Length (*European Type Syringe*)

*Sani-Tip is a registered trademark of Dentsply Sirona

AUTOCLAVABLE METAL AIR-WATER SYRINGE TIP



- Designed with a positive locking feature between the inner and outer tubes for increased safety
- Improved spray pattern
- Special alloy is employed to ensure a virtually flawless finish
- Compatible with all American type standard quick change connectors



DISPOSABLE SYRINGE TIP SHIELDS

- Prevent spattering of saliva and blood
- Increases the effectiveness of all rinsing and dry
- Fully disposable

015 Syringe Tip Shields
 100 pcs/pack

Welcome to our new Greenline. We believe that making simple changes to the everyday products you purchase is the most sustainable way of helping our planet. This range is our first step towards offering you the choice of more environmentally friendly products.

After scratching our heads, we've created our interdental brushes with a few 'ez' changes:

ez interdental brushes




1 We've reduced the size of our handle to a smaller, more practical design – this uses up to 50% less plastic than other leading brands.

2 Instead of using new plastic, we use 100% recycled plastic for the brush handle and 50% recycled plastic for the protective cap. This recycled plastic is made using waste materials left over from other manufacturing processes.

3 Unlike other brands who package interdental brushes in unnecessary plastic boxes and bags, we use cardboard which is 100% plastic free and recyclable.

- ISO size (colour coded)
- Latex and BPA free
- Vegan Friendly
- Reusable
- Non-slip, easy grip design
- Plastic free packaging
- Small, compact and easy to carry

Interdental Brushes, ISO Size No. 1-6
Single or Assorted Sizes, 6 pcs/box



ISO Size	1	2	3	4	5	6
Reference Size	XX Fine	X Fine	Fine	Fine -Medium	Medium	Large
Made to fit hole (mm)	<0.8	0.9-1.0	1.1-1.2	1.3-1.5	1.6-1.8	>1.9
Wire Dia. (mm)	0.35	0.35	0.45	0.45	0.45	0.45
Order Code	G800-1	G800-2	G800-3	G800-4	G800-5	G800-6
Assorted Order Code	G800-A No.1-6 Assorted					

Greenline micro applicators *mini*

In an effort to limit the environmental effect of one of the most used items in a dental surgery we came up with some great ideas:



ACTUAL SIZE

G-9901 REGULAR

G-9902 FINE

G-9903 ULTRA FINE

- Bends into any angle for the placement of small controlled amounts of material in limited access areas.
- Economical - saving up to 80% by the reduction of excess materials used.

*Greenline Micro Applicators Mini, 100 pcs/box
(Length 85 mm, dia 2.5 mm handles)
G-9901 Regular (dia. 2.5 mm) Light Green
G-9902 Fine (dia. 2 mm) Medium Green
G-9903 Ultra Fine (dia. 1.5 mm) Dark Green*

- 1 We designed the handle to be a smaller, more practical size that reduces the amount of plastic used by 50% in comparison to competitive products.
- 2 Instead of new plastic we used 100% recycled plastic. This recycled plastic is made using waste materials from other manufacturing processes.
- 3 No more plastic tube that can use up to a third of the amount of plastic as the actual product - just a 100% recyclable box.



ACTUAL SIZE

ø2.5 mm

ø2 mm

G-921 REGULAR

G-922 FINE

7502W

EASY TO ATTACH TO HANDLE

- For the placement of small controlled amounts of material in limited access areas.
- Economical - saving up to 80% by the reduction of excess materials used.

*Greenline Micro Applicator Tips, 100 pcs/box
G-921 Regular (dia. 2.5 mm) Light Green
G-922 Fine (dia. 2 mm) Medium Green
Autoclavable (134°C) Micro Applicator Tip Handle
7502W White
7502T Transparent*

Greenline micro applicator *tips*

Our most environmentally friendly micro applicators:

- 1 By using these tips you save 80% of the plastic that you would normally use in the handle of a standard micro applicator.
- 2 An autoclavable handle saves plastic. Clean it, autoclave it and re-use it. Simple and effective.
- 3 100% recyclable environmentally friendlier box. Recycling means preventing pollution by reducing the need to collect new raw materials.
- 4 Instead of new plastic we use 100% recycled plastic. The recycled plastic is made using waste materials from other manufacturing processes.



Greenline air-water syringe tips

Reusable products are a great idea. Unfortunately, trying to clean the inside of a reusable metal air-water syringe tip is just not possible and if you can't clean it, you can't guarantee the effectiveness of its sterilisation. This means that using a disposable plastic air-water syringe tip is the only way to guarantee avoiding cross contamination.

We have however created a more environmentally friendly air-water syringe tip by making a few changes.


- Directly interchangeable with American type syringe tips (long), or European syringes which have been converted to use American type syringe tips (short)
- Made from 100% recycled plastic
- No sharp edges for patient comfort
- Unique quad air core design to ensure dry air
- Environmentally friendly packaging

Greenline Air-water Syringe Tips
 G-018 (American Type) 200 pcs/box
 G-028 (European Type) 200 pcs/box

1 Our Greenline air-water syringe tips are made from 100% recycled plastic, using waste materials left over from other manufacturing processes.


2 The Greenline air-water syringe tips are packaged in 100% recyclable cardboard boxes, instead of large plastic bags, helping to reduce the amount of plastic packaging in your surgery.

Greenline PSP barrier envelopes



In order to improve on an already great product and in an effort to reduce waste we've: Removed the plastic packaging and supplied them in a 100% recyclable paper box, which conserves natural resources like timber and water.





- Fits most phosphor plates
- Fully disposable envelopes
- Soft vinyl with rounded corners
- Environmentally friendly packaging



Mid-Opening Barrier Envelopes
 G188-0 No.0 (I.D. 39 mm x 28 mm) 100 pcs/box
 G188-1 No.1 (I.D. 44 mm x 29 mm) 100 pcs/box
 G188-2 No.2 (I.D. 45 mm x 35 mm) 100 pcs/box

Top-Opening Barrier Envelopes
 G198-0 No.0 (I.D. 27 mm x 40 mm) 100 pcs/box
 G198-1 No.1 (I.D. 28 mm x 45 mm) 100 pcs/box
 G198-2 No.2 (I.D. 34 mm x 46 mm) 100 pcs/box

DISPOSABLE FACE MASKS EAR-LOOP WITH PATTERNS FLUID RESISTANT LATEX FREE






LEVEL 3 BARRIER
INTERNATIONAL STANDARD ASTM F2100
BFE ≥99% PFE ≥99%
PERFORMANCE TYPE IIR
EUROPEAN STANDARD EN14683

- Non-irritating
- Fluid resistant, excellent filtration (**BFE 99%**)
- Permits air to pass through mask to prevent glasses from fogging

2288 Face Masks Level 3 (IIR) Ear-Loop with Patterns, 50 pcs/box


DISPOSABLE FLUID RESISTANT DENTAL BIBS WITH PATTERNS



- Two plies of facial quality tissue plus a layer of polyethylene film offers quick absorbency with a positive moisture barrier
- Paper specially treated for wet strength
- Embossed horizontal pattern, water contact disperses horizontally instead of "channelling" down onto patients clothing
- Size: 45 cm x 33 cm

Dental Bibs with Patterns
2200FA Flowers Assorted, 500 pcs/box
2200LA Leaves Assorted, 500 pcs/box
2200KID Kids, 500 pcs/box
2200MON Monsters, 500 pcs/box


DISPOSABLE PROPHY BIBS
EXTRA LONG FLUID RESISTANT



- Ideal for prophy, where extra protection is needed
- Features contour neck design for additional coverage to shoulders and arm
- Two plies of facial quality tissue plus a layer of polyethylene
- Size: **W45 cm x L60 cm**
- Colour: **Light Blue**

2208 Prophy Bibs Extra Long
100 pcs/box + Bib Holder XS


DISPOSABLE FLUID RESISTANT
TRAY COVER PAPERS



- Two plies of facial quality tissue plus a layer of polyethylene film for positive barrier
- Colours: **Blue** and white barrier

2220 "B" Plastic Tray Covers 200'S (29.5 cm X 21.5 cm)
2221 "F" Plastic Tray Covers 200'S (20 cm X 13.5 cm)
2223 "ISO" Aluminium Tray Covers 200'S (27.5 cm X 17.5 cm)


DISPOSABLE NAPKIN HOLDERS



- Fully disposable
- Easy to apply
- Width: **15 mm**
- Length: **44 cm** extendable up to **60 cm**

2209 Napkin Holders, 100 pcs/pack

DISPOSABLE PSP BARRIER ENVELOPES



- Quick tear-off opening
- Fits most phosphor plates
- Fully disposable envelope
- Soft vinyl with rounded corners

DURR® TYPE

Mid-Opening Barrier Envelopes

188-0 No. 0 (I.D. 39 mm x 28 mm) 100 pcs/box
 188-1 No. 1 (I.D. 44 mm x 29 mm) 100 pcs/box
 188-2 No. 2 (I.D. 45 mm x 35 mm) 100 pcs/box
 188-3 No. 3 (I.D. 58 mm x 30 mm) 50 pcs/box
 188-4 No. 4 (I.D. 82 mm x 60 mm) 50 pcs/box

Mid-Opening Barrier Envelopes Bulk

188-2B No. 2 (I.D. 45 mm x 35 mm) 500 pcs/box

Top-Opening Barrier Envelopes

198-0 No. 0 (I.D. 27 mm x 40 mm) 100 pcs/box
 198-1 No. 1 (I.D. 28 mm x 45 mm) 100 pcs/box
 198-2 No. 2 (I.D. 34 mm x 46 mm) 100 pcs/box

Top-Opening Barrier Envelopes Bulk

198-2B No. 2 (I.D. 33 mm x 44 mm) 500 pcs/box

DIGORA® (SOREDEX®) TYPE

Top-Opening Barrier Envelopes

189-0 No. 0 (I.D. 24 mm x 34 mm) 100 pcs/box
 189-1 No. 1 (I.D. 26 mm x 43 mm) 100 pcs/box
 189-2 No. 2 (I.D. 33 mm x 44 mm) 100 pcs/box

KODAK® (CARESTEAM®) TYPE

Top-Opening Barrier Envelopes

199-0 No. 0 (I.D. 31 mm x 36 mm) 100 pcs/box
 199-1 No. 1 (I.D. 33 mm x 44 mm) 100 pcs/box
 199-2 No. 2 (I.D. 40 mm x 46 mm) 100 pcs/box

Protective Paper Covers

178-0 (2 x 23 mm x 36 mm) 100 pcs/box
 178-1 (2 x 24 mm x 42 mm) 100 pcs/box
 178-2 (2 x 30 mm x 42 mm) 100 pcs/box



133

DISPOSABLE PC MOUSE SLEEVES

- Universal size, fit both cordless and corded optical mouse
- Soft, pliable polypropylene sleeves
- Easy to apply

PC Mouse Sleeves, 500 pcs/box
 133S Small 15 cm x 8.5 cm (5 7/8" x 3 3/8")
 133L Large 16 cm x 10 cm (6 1/4" x 3 7/8")

DISPOSABLE X-RAY SLEEVES



108

109

- Full x-ray unit coverage
- Disposable yet economical
- Soft, pliable polypropylene sleeve holds position until removed
- Easy application and removal

X-Ray Head Sleeves, 250 pcs/box

108 Small 20 cm x 17 cm x 60 cm (7 7/8" x 6 5/8" x 23 1/2")
 109 Large 35 cm x 24 cm x 75 cm (13 3/4" x 9 1/2" x 29 1/2")

DISPOSABLE HVE TUBE SLEEVES




135S
135L

- Pre-cut opening for easier application
- Made from very supple poly plastic for greater tactile feel

HVE Suction Tube Sleeves, 500 pcs/box

135S Small (11 mm) 35 cm x 5.5 cm (13 3/4" x 2 1/8")
 135L Large (16 mm) 38 cm x 6 cm (15" x 2 3/8")

A. STANDARD



140XS
140S
140M

B. TWO STEP



183-1
183-2

DISPOSABLE X-RAY SENSOR SLEEVES

- Made from very supple poly plastic for greater tactile feel
- Easy to apply

A. Standard Sensor Sleeve, 500 pcs/box
 140XS X-Small 3.8 cm x 20 cm (1 1/2" x 7 7/8")
 140S Small 4.2 cm x 20 cm (1 5/8" x 7 7/8")
 140M Medium 4.5 cm x 20 cm (1 3/4" x 7 7/8")
 140L Large 6 cm x 20 cm (2 3/8" x 10")

B. Two Step Sensor Sleeves, 500 pcs/box
 183-1 Size 1 3.4 cm / 5 cm x 22 cm (1 3/8" / 2" x 8 5/8")
 183-2 Size 2 3.8 cm / 5.5 cm x 22 cm (1 1/2" / 2 1/8" x 8 5/8")



101

101

102

DISPOSABLE INSTRUMENT TRAY SLEEVES

- Instruments are placed directly on top of the disposable plastic sleeve
- Reduces the possibility of cross-contamination, saves cleaning time

Tray Sleeves, 500 pcs/box
 101 Standard 'B' 27 cm x 36 cm (10 1/2" x 14")
 102 Mini 'F' 19 cm x 28 cm (7 1/2" x 11")



DISPOSABLE ISO INSTRUMENT TRAY SLEEVES

- Instruments are placed directly on top of the disposable plastic sleeve
- Reduces the possibility of cross-contamination, saves cleaning time

ISO Metal Tray Sleeves

500 pcs/box
 190 Standard 31 cm x 21 cm (12 1/8" x 8 1/4")



107

DISPOSABLE IMPRESSION GUN SLEEVES

- Plastic sleeve covers cartridge and handle
- Provides protection in case of cartridge failure or leakage

Impression Gun Sleeves, 500 pcs/pack
 107 Universal 25 cm x 21 cm (9 7/8" x 8 1/4")




140XS

DISPOSABLE MOUTH MIRROR SLEEVES

- Made from very supple poly plastic for greater tactile feel

Mouth Mirror Sleeves, 500 pcs/box
 140XS X-Small 3.8 cm x 20 cm (1 1/2" x 7 7/8")

DISPOSABLE BARRIER FILM



8008-B 8008-C

- With 5 mm non-adhesive edge for simple placement and removal
- Residue free
- Colours: Blue and Clear

Barrier Film, 1200 pcs/roll

8008-C Clear 10 cm x 15 cm (3 7/8" x 6")
 8008-B Blue 10 cm x 15 cm (3 7/8" x 6")

DISPOSABLE T-SHAPED LIGHT HANDLE SLEEVES



106
T-SHAPE

- Universal sizing to fit most T-style handles
- Full handle coverage
- Disposable yet economical
- Soft, pliable polypropylene sleeve holds position until removed
- Easy application and removal

T-Shape Light Handle Sleeves, 500 pcs/box

106 Universal 10 cm / 15 cm x 8.5 cm (3 7/8" / 6" x 3 3/8")



145

DISPOSABLE BITE BLOCK SLEEVES

- Polypropylene sleeves
- Easy application and removal
- Fully disposable

Bite Block Sleeves, 5 x 200 pcs/box

145 Universal 5.5 cm x 2.6 cm (2 1/8" x 1")

DISPOSABLE PEN TYPE CURING LIGHT SLEEVES



- Made from very supple poly plastic for greater tactile feel
- Easy to apply

Pen Type Curing Light Sleeves, 500 pcs/box
228 Mini Short 25 cm x 3.8 cm (9" x 1 1/2")
123SS Small Short 28 cm x 5 cm (11" x 2")

123S Small 32 cm x 5 cm (12 1/2" x 2")
123L Large 35 cm x 6.5 cm (13 3/4" x 2 1/2")
123XL X-Large 38 cm x 7.2 cm (15" x 2 7/8")

DISPOSABLE INTRAORAL CAMERA /RA CURING LIGHT SLEEVES



- Made from very supple poly plastic for greater tactile feel
- Easy to apply

Intraoral Camera Sleeves / RA Curing Light Sleeves
228 Mini Short 25 cm x 3.8 cm (9" x 1 1/2"), 500 pcs/box
125 Standard 33 cm x 4.8 cm (13" x 1 7/8"), 500 pcs/box
125SS Short 28 cm x 4.8 cm (11" x 1 7/8"), 500 pcs/box



DISPOSABLE PEEL-N-STICK LIGHT GUIDE SLEEVES

- Soft, pliable polypropylene Sleeve
- Easy to apply
- Economical

Light Guide Sleeves, 200 pcs/box
128 Small (fits 7-8 mm Guide)
129 Large (fits 9-11 mm Guide)



DISPOSABLE PEN AND HANDPIECE SLEEVES

- Disposable sleeves for pen, handpieces, and ultrasonic scalars
- Made from soft pliable poly plastic for greater tactile feel

500 pcs/box
103 Pen Shape Instrument 22 cm x 2.4 cm (8 1/2" x 1")
104S Low Speed/Universal Slim 22 cm x 3.8 cm (8 1/2" x 1 1/2")
104L Low Speed/Universal Large 22 cm x 4 cm (8 1/2" x 1 5/8")



DISPOSABLE SYRINGE/HVE SLEEVES

- Pre-cut opening for easier application
- Ideal for 3-way syringes, saliva ejectors and HVE valves etc

Syringe Sleeves, 500 pcs/box
100C Clear 26 cm x 6.5 cm (10" x 2 1/2")
100B Blue 26 cm x 6.5 cm (10" x 2 1/2")

DISPOSABLE TUBING SLEEVES



- For straight and coiled tubing
- Affordable protection from infectious aerosols and minimizes time spent on cleaning
- 1200 ft/roll (366 m)

118 Tubing Sleeve for straight tubing 2" (5 cm)
119 Tubing Sleeve for coiled tubing 4" (10 cm)

DISPOSABLE FLOWABLE COMPOSITE SYRINGE SLEEVES



- Made from very supple poly plastic for greater tactile feel
- Easy to apply
- Economical

Flowable Composite Syringe Sleeves, 200 pcs/box
138 Universal 1.8 cm / 5 cm x 12.5 cm (5/8" / 2" x 5")



DISPOSABLE KEYBOARD SLEEVES

- Soft, pliable polypropylene sleeves

Keyboard Sleeves
158S Standard 16 cm x 55 cm (6 1/4" x 21 5/8"), 500 pcs/box
158M Medium 21.5 cm x 55 cm (8 1/2" x 21 5/8"), 250 pcs/box
108 Large 20 cm + 17 cm x 60 cm (7 7/8" + 6 5/8" x 23 1/2"), 250 pcs/box




DISPOSABLE CHAIR SLEEVES

- Eliminates the need for separate headrest covers
- Avoids the need for harmful disinfectants and aerosols prolonging the life of dental chairs/stools

Chair Sleeves, 125 pcs/box
111 Full Chair 180 cm x 72 cm (70" x 28 3/8")
112 Half Chair 100 cm x 65 cm (39 3/8" x 25 5/8")

DISPOSABLE HEADREST SLEEVES



- Polypropylene sleeves
- Reduces possible cross-contamination from infectious aerosol spray

Headrest Sleeves, 500 pcs/box
113 Small 26 cm x 27.5 cm (10 1/4" x 10 7/8")
114 Large 35 cm x 28 cm (13 3/4" x 11")



DISPOSABLE NON-WOVEN HEADREST COVERS

- Breathable non-woven fabric
- Soft, provides maximum patient comfort
- Colours: Light Blue and White

Non Woven Headrest Covers, 500 pcs/box
780S Small 26 cm x 28 cm (10 1/4" x 11")
780M Medium 35 cm x 28 cm (13 3/4" x 11")

Applicators and Tips




DISPOSABLE MICRO APPLICATORS

- Economical and disposable, saving up to 80% by the reduction of excess materials used
- Bends into any angle for the placement of small controlled amounts of material in limited access areas



Available in :-
Standard Handles
 100 mm x dia. 3 mm
Mini Handles
 85 mm x dia. 2.5 mm



Available in :-
Standard Handles
 100 mm x dia. 3 mm
Mini Handles
 85 mm x dia. 2.5 mm



Available in :-
Standard Handles
 100 mm x dia. 3 mm


Standard Handles
 Tube Pack 400 pieces
 901T Regular
 902T Fine
 903T Ultrafine
 904T Cylinder

Mini Handles
 Tube Pack 400 pieces
 9901T Regular
 9902T Fine
 9903T Ultrafine

Standard Handles
 Bulk Pack 1000 pieces
 901B Regular
 902B Fine
 903B Ultrafine
 904B Cylinder

Mini Handles
 Bulk Pack 1000 pieces
 9901B Regular
 9902B Fine
 9903B Ultrafine

Standard Handles
 Dispenser Pack 400 pieces
 901D Regular
 902D Fine
 903D Ultrafine
 904D Cylinder




DISPOSABLE MICRO APPLICATOR DISPENSERS

- Small, sturdy, durable construction
- Provides aseptic barrier against operator aerosols and accidental splatter

Applicator Dispenser
 18-900 **Standard Handle**
 18-9900 **Mini Handle**

DISPOSABLE MICRO APPLICATOR TIPS




- Economical and disposable saving up to 80% by the reduction of excess materials used
- For the placement of small controlled amounts of material in limited access areas
- Can be used for etchants, bonding agents, cavity lining, fluoride varnishing, hemostatic and sealers etc

Micro Applicator Tips, 100 pcs/box
 (Blue, Green, Lavender or Yellow)
 921 Regular
 922 Fine

Brush Handle Autoclavable
 7502W - White
 7502T - Transparent

DISPOSABLE BRUSH APPLICATORS



- Perfect for placement of etchants, bonding agents, sealers and cements etc
- Bendable to reach around corners and the tightest spaces
- Handles are colour coded to eliminate confusion between materials
- Fully disposable, eliminates cross-contamination
- Available colours: Blue, Green, Lavender or Yellow

Brush Applicator Single Colour, 100 pcs/pack
 7500 Regular
 7510 Mini

DISPOSABLE BRUSH TIPS AND HANDLES



- Perfect for placement of etchants, bonding agents, sealers and cements etc
- Available Colours: Blue, Green, Lavender, or Yellow


7501 Brush Tips **Regular** Standard Pack, 100 pcs/box
 7511 Brush Tips **Mini** Standard Pack, 100 pcs/box
 7502W Brush Handle Autoclavable **White**
 7502T Brush Handle Autoclavable **Clear**

COMPOSITE TIP GUN AND TIPS



- For placement of composites, cements etc
- Reusable handle with disposable tubes and plugs
- Autoclavable plastic handle with pistol grip design
- Tips are photoblock

AUTOCLAVABLE HVE SUCTION TUBES



- Fits European style 16 mm or standard 11 mm suction tubing
 - Autoclavable up to 135°C
 - Special high quality plastic, sturdy enough for cheek retraction
 - Fully polished internal surface and edge for quiet operation and patient comfort
 - Assortment of bright, lively colours: Lavender, Light Blue, Light Green, Lime and Yellow
- Assorted Colours, 10 pcs/pack**
 335L Adult Large (16 mm)
 335M Adult Medium (16 mm)
 336P Pedo (16 mm)
 005L Adult Large (16 mm)
 005M Adult Medium (16 mm)
 006P Pedo (16 mm)
 205M Adult Medium (11 mm)
 206P Pedo (11 mm)

VENTED

SOFT TOUCH™



ST001 WHITE

DISPOSABLE SOFT TOUCH™ ORAL EVACUATORS (SLOTTED) VENTED LATEX FREE

- Soft Touch edge eliminates the risk of lacerating the oral tissue
- Rigid tubing strong enough for cheek retraction
- Vented to reduce unnecessary tissue grab and increase evacuation cover area
- Thoroughly polished surface for smooth and quiet operation
- Fits standard 11 mm suction holders
- Colour: White (Blue or Lavender Tip)

PATENT PENDING

ST001 Soft Touch™ Oral Evacuators, 50 pcs/pack

VENTED (SLOTTED) ORAL EVACUATORS



- Vented to reduce unnecessary tissue grab
- Soft rigid tubing, allows for greater patient comfort yet is strong enough for cheek retraction
- Thoroughly polished surface for smooth and quiet operation
- Fits standard 11 mm suction holders

Oral Evacuators, 100 pcs/pack
 001 White Disposable
 001A Single Colour Disposable
 003 Grey Autoclavable


VENTED/NON-VENTED ORAL EVACUATORS 'S' CUT



- Vented and non-vented combination
- Soft rigid tubing, allows for greater patient comfort yet is strong enough for cheek retraction
- Thoroughly polished surface for smooth and quiet operation
- Fits standard 11 mm suction holders

Oral Evacuators, 100 pcs/pack
 301 White Disposable
 301A Single Colour Disposable
 303 Grey Autoclavable

VENTED ORAL EVACUATORS 'S' CUT



- Vented to reduce unnecessary tissue grab
- Soft rigid tubing, allows for greater patient comfort yet is strong enough for cheek retraction
- Thoroughly polished surface for smooth and quiet operation
- Fits standard 11 mm suction holders
- Colours: Blue, Yellow, Green, Lime or Lavender

Oral Evacuators, 100 pcs/pack
 308 White Disposable
 308A Single Colour Disposable
 309 Grey Autoclavable

AUTOCLAVABLE SURGICAL ASPIRATOR TIPS

- Non-sterile
- Autoclavable up to 135°C
- Colours: Standard - White Large orifice - Blue



Autoclavable Surgical Tips, 25 pcs/pack
 007 Mini Tip
 009 Standard Non-Vented
 009V Standard Vented
 008 Large Orifice Non-Vented
 008V Large Orifice Vented

DISPOSABLE SALIVA EJECTORS



- Tip design that eliminates suction of tissue
- Smooth fixed tip for patient safety and comfort

Saliva Ejectors, 100 pcs/pack
 300W White Tip/Clear Tube
 300B Blue Tip/Clear Tube
 300A White Tip/Coloured Tube

DISPOSABLE SALIVA EJECTORS SOFT TIPS



- Flexible soft cushion tip for patient comfort
- Latex free tip design that eliminates suction of soft tissue
- Fixed tip for safety

Saliva Ejectors, 100 pcs/pack
 333W White Tip/Coloured Tube
 333B Blue Tip/Clear Tube
 333A White Tip/Coloured Tube




DISPOSABLE 3-IN-1 TRAYS

- Wide occlusal plane
- Strong and tough plastic
- Specially designed walls for extra retention
- Tear resistant quality mesh
- Colour: Blue

3-in-1 Trays
 18-N401 Posterior 36 pcs/box
 18-N402 Quadrant 36 pcs/box
 18-N403 Sideless 36 pcs/box
 18-N404 Anterior 36 pcs/pack
 18-N405 Full Arch 24 pcs/pack

DISPOSABLE 3-IN-1 TRAYS

Stainless Steel Wire Connectors



- RIGID Stainless Steel Wire Connector
- Stainless steel wire connector for greater access in limited space
- Wide occlusal plane
- Special rigid walls designed for extra retention
- Tear resistant quality mesh
- Colour: Blue

3-in-1 Trays
 18-S401 Posterior 36 pcs/box
 18-S402 Quadrant 36 pcs/box
 18-S403 Sideless 36 pcs/box

DISPOSABLE 3-IN-1 TRAYS

Rigid Bar Connectors




- RIGID Aluminium Bar Connector
- Strong and tough
- Wide occlusal plane
- Special rigid walls designed for extra retention
- Tear resistant quality mesh
- Colour: Blue


3-in-1 Trays
 18-T401 Posterior 36 pcs/box
 18-T402 Quadrant 36 pcs/box
 18-T403 Sideless 36 pcs/box

DISPOSABLE IMPRESSION TRAYS WITH RIM LOCK

for Implant Post and General Impressions



base sections can be easily removed with a bur for implant posts



rim lock reduced thickness trim lines on labial wall for easy removal in inclination post cases

Disposable Impression Trays, 10 pcs/pack
 410 Intro (one size each)
 411-420 Single pack

- Save time and expensive laboratory charges
- Strong, rigid distortion free plastic construction
- Anatomically designed
- Smooth edges for patient comfort
- Excellent retention for impression materials whilst providing sufficient pressure for all details
- Colour: Blue



AUTOCLAVABLE IMPRESSION TRAYS WITH RIM LOCK

- Autoclavable up to 135°C (275°F)
- Strong, rigid distortion free plastic construction
- Anatomically designed
- Smooth edges for patient comfort
- Excellent retention for impression materials whilst providing sufficient pressure for all details
- Colour: Light Grey

Autoclavable Impression Trays
 430 Intro 10 pcs/pack (one size each)
 431-440 Single 5 pcs/pack



DISPOSABLE 3-IN-1 TRAYS

- Wide occlusal plane
- Strong and tough plastic
- Specially designed walls for extra retention
- Tear resistant quality mesh
- Individually wrapped

3-in-1 Trays

N401 Posterior 36 pcs/box
 N402 Quadrant 36 pcs/box
 N403 Sideless 36 pcs/box
 N404 Anterior 36 pcs/pack
 N405 Full Arch 24 pcs/pack

DISPOSABLE 3-IN-1 TRAYS

Wire Connectors Stainless Steel



- Stainless steel wire connector for greater access in limited space
- Wide occlusal plane
- Special rigid walls designed for extra retention
- Tear resistant quality mesh

3-in-1 Trays

S401 Posterior 36 pcs/box
 S402 Quadrant 36 pcs/box
 S403 Sideless 36 pcs/box

DISPOSABLE 3-IN-1 TRAYS

Rigid Bar Connectors




- Rigid Aluminium Bar Connector
- Strong and tough
- Wide occlusal plane
- Special rigid walls designed for extra retention
- Tear resistant quality mesh

3-in-1 Trays

T401 Posterior 36 pcs/box
 T402 Quadrant 36 pcs/box
 T403 Sideless 36 pcs/box

DISPOSABLE PLASTIC COATED METAL Pliable IMPRESSION TRAYS



- Fully disposable
- Pliable metal

451-460 Pliable Metal Tray

AUTOCLAVABLE ALGinate/PLASTER CURVED SPATULA

- For alginate or plaster
- Light weight, 58 gms only
- Fully autoclavable up to 135°C
- Stainless steel blade
- Handle colour: Blue

MS455 Curved Spatula

DISPOSABLE IMPRESSION SYRINGES

- Pre-bent for easy application
- Fully disposable
- 5 ml capacity

089 Impression Syringes 50 pcs/pack

FLEXIBLE MIXING BOWLS



- For effective mixing of alginates and plasters
- Sizes: 7460 Small dia. 11 cm x H8 cm
7470 Medium dia. 12.5 cm x H10 cm
- Available in Transparent Red, Green, Blue, Purple, Yellow and solid Blue

7460 Flexible Mixing Bowl Small
 7470 Flexible Mixing Bowl Medium


AUTOCLAVABLE ALGinate AND PLASTER SPATULAS



- Fully autoclavable up to 134°C
- Fit all flexible mixing bowls
- Will not scratch the surface of bowls and clean up easily
- Flexible control in mixing with easy grip handle
- Size: Length: 190 mm

450R Regular Blue
 450XS Extra Stiff Grey

DISPOSABLE INTRAORAL SYRINGE TIPS



- Universal fitting
- Available in Yellow or White

Intraoral Tips, 100 pcs/pack
 085 Yellow for Yellow mixing tips
 086 Clear for Green and Pink mixing tips

DISPOSABLE DYNAMIC MIXING TIPS

- Fits Pentamix[®] type mixing machine
- Economically priced
- Colours: Red or White

7055 Dynamic Mixing Tips, 50 pcs/pack

¹ Pentamix is a registered trademark of 3M or 3M ESPE AG.


AUTOCLAVABLE ENDO STAND



- Fully autoclavable plastic stand with disposable sponge insert
- Cleans and holds endo instruments
- Silicone ring base resists sliding
- Fully autoclavable up to 134°C
- Colour: Blue

4833 Endo Stand Large (dia. 45 mm)
812 Endo Sponges, 50 pcs/pack

DISPOSABLE ENDO SPONGES



811
for round endo stand
size: 2" dia. x 3/8" thick
(50mm dia. x 8mm)

812
for oval and round endo stand
size: 2 1/2" dia. x 3/8" thick
(65mm dia. x 8mm)

- Fully disposable blue polyurethane sponges
- Colour: Blue

Endo Sponges, 50 pcs/pack
811 for Round Endo Stand
812 for Oval and Large Endo Stand


DISPOSABLE PEEL-AND-STICK ENDO CUSHIONS



802 Peel-and-Stick
Endo Cushions
50 pcs/box

- Hold, arrange and clean small instruments during endodontic or implant procedures
- Ideal for files, screws, drills and pluggers etc
- Peel-and-stick back secures to patient napkins, tray covers or anywhere handy
- Pre-arranging saves set up time
- Dense blue polyurethane foam 75 mm long

AUTOCLAVABLE ENDO CUSHION HOLDER



- Holder for used files with disposable foam cushions
- Built in endo file ruler
- Stainless steel construction
- Fully autoclavable up to 134°C
- Sizes: 50 mm x 50 mm x 25 mm

1180 Endo Cushion Holder
(without Cushion)
801 Square Endo Cushions
50 pcs/pack

AUTOCLAVABLE TRIANGULAR ENDO RING AND CUSHIONS



- Fully disposable cushions
- Colour: Blue
- Autoclavable up to 134°C

088 Autoclavable Cushion Holders
with Large and Small Rings
803 Disposable Cushions 50 pcs/box



AUTOCLAVABLE ENDO HOLDER ALUMINIUM

2800 Endo Holder


- Specially designed for ultrasonic cleaning and autoclaving of root canal instruments
- Holds 18 hand or engine endo instrument up to size 100 all lengths
- Millimeter rule for rubber stop measurements for both right or left hand use
- Adjustable cover with silicone stop usage counter for sterilization
- Size: W65 mm x D30 mm x H55 mm
- Fully autoclavable up to 134°C
- Colour: Blue

AUTOCLAVABLE RUBBER DAM FRAME



4866
Rubber
Dam Frame

- Proven design
- Better grip
- Enhance patient comfort




AUTOCLAVABLE ENDO DISPENSER PLASTIC

055 Endo Dispenser

- Specially designed for cleaning and autoclaving of endo instruments
- Holds 18 hand or engine instruments up to size 100 (all lengths)
- Medical grade plastic, autoclavable up to 134°C
- Size: 85 mm x 52 mm x 14 mm
- Colour: Blue

AUTOCLAVABLE ENDO MEASURING BLOCK PLASTIC



4888
Endo Measuring Block

- Accurate settings
- Durable construction
- Fully autoclavable up to 134°C
- Size: L100 mm x D24 mm x H32 mm
- Colour: Blue

AUTOCLAVABLE ENDO RULER STAINLESS STEEL



- Durable Stainless steel
- Accurate setting
- Size: L80 mm x D12.5 mm x H3 mm


4899 Endo Ruler

AUTOCLAVABLE ENDO RULER ALUMINIUM



- For both hand and engine instruments
- Gauges from 1 to 30 mm for right or left hand

4898 Endo Ruler



10005
Endo Measuring Block

AUTOCLAVABLE ENDO MEASURING BLOCK ALUMINIUM

- Durable aluminium construction
- Accurate settings
- Fully autoclavable up to 134°C
- Size: L150 mm x W24 mm x H33 mm
- Colour: Blue

AUTOCLAVABLE TRAYS PLASTIC




- Fully autoclavable up to 134°C

650 Standard Flat 'B' Tray 34 cm x 24 cm x 2 cm
 651 Standard Divided 'B' Tray 34 cm x 24 cm x 2 cm
 652 Mini Tray Flat 'F' Tray 24 cm x 16.5 cm x 2 cm

STANDARD COLOURS
Blue
Light Blue
Green
Yellow
White
Light Beige
Lavender

AUTOCLAVABLE CABINET TRAYS PLASTIC



- Autoclavable up to 134°C
- Indispensable drawer organizer
- Scientifically designed for maximum convenience and efficiency
- Easy to clean, sterilize and disinfect
- Resists breaking, chipping and chemical attack
- Colours: Light Blue, Blue, White, Yellow, Lavender, Green and Light Beige
- Size: 9.5 cm x 19.5 cm x 2 cm

653 Cabinet Tray

STANDARD COLOURS	Yellow	Light Blue	Light Beige	Blue	White	Green	Lavender
Blue	White	Green	Lavender				



658 Standard Rack 659 Mini Rack

658 Tray Rack Standard 'B'
 659 Tray Rack Mini 'F'

INSTRUMENT TRAY RACKS

- Sturdy, chromium plated steel construction
- Size: Standard (8 trays) W37 cm x D24 cm x H37 cm
 Mini (6 trays) W28 cm x D17 cm x H29 cm

AUTOCLAVABLE INSTRUMENT MAT



- Fully autoclavable
- Medical grade silicone material
- Size: L172 mm x W50 mm x H10 mm

STANDARD COLOURS	Yellow	Light Blue	Light Beige
Blue	White	Green	Lavender


AUTOCLAVABLE NEEDLE STATION STAINLESS STEEL



- Holds any dental needle sheath
- Solid stainless construction, fully autoclavable
- Silicone ring base resists sliding
- Size: Dia. D55 mm x H38 mm

031 Needle Station

AUTOCLAVABLE STERILIZATION CASSETTES PLASTIC



- Sturdy, durable construction
- Fully autoclavable

3881 Small (5 instruments)
 3882 Medium (8 instruments)

STANDARD COLOURS	Blue	White	Light Blue	Green	Lavender
------------------	------	-------	------------	-------	----------

AUTOCLAVABLE SILICONE INSTRUMENT CODE RINGS



- Fully autoclavable to 134°C
- Sizes: Small - for small instruments
 Large - for pliers and forceps etc
- Available in single colors: Blue, Light Blue, Green, Yellow, Pink, Red, Orange, Lavender, Black and White

Silicone Instrument Rings
 657S Small, 50 pcs/box
 657L Large, 25 pcs/box

AUTOCLAVABLE BUR/ENDO ORGANIZERS



4810 (78 holes)
 for both hand/engine endo instruments
 21 mm - 30 mm

4811 (48 holes)
 for both hand/engine endo instruments
 21 mm - 30 mm

4812 (36 holes)
 for both hand/engine endo instruments
 21 mm - 30 mm

4813
 12 endo vials

4818
 78 FG

4819
 42 FG + 36 CA



PERFORATED BUR/ENDO ORGANIZER BOXES

- Smart, flexible design for customize setups
- Fully autoclavable up to 134°C
- Colour: Blue

Endo Stand, W115 mm x D52 mm x H40 mm
 4810 - 78 Holes
 4811 - 48 Holes + 4 Endo Vials
 4812 - 36 Holes + Endo Clean
 4813 - 12 Endo Vials
 Bur Stand, W115 mm x D52 mm x H40 mm
 4818 - 78 FG Holes
 4819 - 42 FG + 36 CA Holes

Stainless Steel Cover, W115 mm x D52 mm x H15 mm
 4830 - Sterilization Cover

Bur/Endo Box, (box only)
 4821 Organizer Box for 1 Stand W120 mm x D55 mm x H56 mm
 4823 Organizer Box for 3 Stands W160 mm x D120 mm x H56 mm
 4825 Organizer Box for 5 Stands W265 mm x D120 mm x H56 mm



STAINLESS STEEL COVER

4830



DISPOSABLE EYE SHIELDS

- Fully disposable
- Available in tinted or clear lenses
- Light weight, optical plastic shields
- Assorted colours frames: **Blue, Light Blue, Lavender, Green and White**

Disposable Eye Shield
20 Lenses + 20 Frames (Assorted Colours)
720T, Tinted
720C, Clear

Replacement Eye Shield
Lens Only, 50 pcs/pack
720T-R, Tinted
720C-R, Clear



505

PROTECTIVE EYEWEARS

- Adjustable side arms
- Contemporary design
- Economical price



500

500 Protective Eyewear **Clear**
505 Protective Eyewear **Tinted**




DISPOSABLE MOUTH MIRRORS

- Fully disposable plastic mouth mirrors for patients
- Anti-fog front surface fiber glass mirror
- Mirror size: dia. **23 mm**

800 Disposable Mouth Mirrors, 10 pcs/box

DISPOSABLE PLASTIC WEDGES



actual size

- Anatomically designed
- Soft tough plastic

Disposable Plastic Wedges, 100 pcs/box
7650X **X-Large (17 mm) White**
7650L **Large (15 mm) Yellow**
7650M **Medium (13 mm) Green**
7650S **Small (11 mm) Blue**



TONGUE CLEANERS

- Effective in the treatment of halitosis
- Multiple cleaning ridges
- Assorted colours

9000 Tongue Cleaners
20 pcs/pack

RETAINER BOXES/MOUTH GUARD BOXES/DENTURE BATHS



600 Retainer Box W75 mm x D61 mm x H28 mm
601 Mouth Guard Box W77 mm x D64 mm x H40 mm
602 All Guard Box W77 mm x D77 mm x H38 mm
603 Denture Bath W83 mm x D68 mm x H52 mm

Colours: **Blue, Green, Yellow, White, Lime, Light Blue, Lavender** or Assorted

600 Retainer Boxes, 10 pcs/box
601 Mouth Guard Boxes, 10 pcs/box

602 All Guard Boxes, 10 pcs/box
603 Denture Baths, 10 pcs/box

DENTAL APPLIANCE BOXES WITH MIRRORS




PASTEL COLOURS	
Vanilla	Sea Green
Sand	Lilac
Light Mauve	Baby Blue
White	Assorted

- Sturdy strong plastic
- Size: **W80 mm x D72 mm x H32 mm**

612
Dental Appliance Boxes
10 pcs/box

INTRAORAL CHEEK RETRACTORS



- Holds mouth open, depresses tongue movement
- Anatomically designed
- Available in both autoclavable and disposable

Autoclavable Intraoral Cheek Retractor
2060A Adult Large, 1 pc/box
2060B Adult Small, 1 pc/box
2060C Child, 1 pc/box



Disposable Intraoral Cheek Retractors
2080A Adult Large, 10 pcs/box
2080B Adult Small, 10 pcs/box
2080C Child, 10 pcs/box

INTRAORAL MOUTH LIPPERS



- Expands the mouth vertically rather than horizontally
- Ideal for working in the gingival areas for anteriors, ortho procedures, bleaching etc...
- Available in both autoclavable and disposable

Autoclavable Intraoral Lippers, 2 pcs/box
2040A Adult, Semi-Transparent
2040C Child, Semi-Transparent

Disposable Intraoral Lippers, 10 pcs/box
2030A Adult, White
2030C Child, White



AUTOCLAVABLE CHEEK RETRACTORS



- Autoclavable up to **134°C**

Mouth Retractor
2050C Child, 2 pcs/pack
2050A Adult, 2 pcs/pack
2055C Child, 2 pairs/pack
2055A Adult, 2 pairs/pack



AUTOCLAVABLE MOUTH PROPS LATEX FREE

2090A Mouth Prop Adult
2090P Mouth Prop Pedo

- Made of chemical resistant, autoclavable silicone
- Durable
- Available in: - Adult (30 mm x 40 mm x 20 mm)
Pedo (30 mm x 25 mm x 20 mm)



DISPOSABLE MULTI-PURPOSE DAPPEN DISHES

004B Standard Pack Assorted Colours 1000 pcs/box

- Fully disposable to save time and prevent cross-contamination
- For prophylaxis, etchant and bonding etc
- Colours: Lavender, White, Green, Light blue and Yellow

DISPOSABLE PROPHY RINGS



- Fully disposable to save time and prevent cross-contamination
- For prophylaxis, endo etc
- Colour: Blue

010 Disposable Propyl Rings, 200 pcs/box

AUTOCLAVABLE PROPHY RINGS



- Autoclavable up to 134°C
- For prophylaxis, endo etc
- Colour: Grey

A010 Autoclavable Propyl Rings, 10 pcs/pack


GLASS DAPPEN DISH



- Colours: Light Blue, Green, Clear or Brown

040 Glass Dappen Dish

DOUBLE ARCH



LARGE blue

MEDIUM yellow

SMALL green

SINGLE ARCH



LARGE blue

MEDIUM yellow


SMALL green

DISPOSABLE FLUORIDE TRAYS

- Extremely comfortable, anatomically designed trays
- Economical price
- Colour coded

Double Arch, 50 pcs/box
7800S Small (green)
7800M Medium (yellow)
7800L Large (blue)

Single Arch, 100 pcs/box
7801S Small (green)
7801M Medium (yellow)
7801L Large (blue)



DISPOSABLE MIXING PADS

Mixing Pads, 3X100 pcs/pack
3220 Small 35 mm X 45 mm
3225 Large 70 mm X 80 mm



DISPOSABLE MIXING WELLS

- Fully disposable
- Sturdy design

906 Mixing Wells 4 Wells 200 pcs/box
907 Mixing Wells 2 Wells 500 pcs/box



DISPOSABLE FLOWABLE COMPOSITE AND ETCHING GEL TIPS

- Universal Luer-lock design
- Bendable
- Gauges: 080 G.20 (0.9 mm)
081 G.23 (0.6 mm)
082 G.18 (1.2 mm)

100 pcs/pack
080 Flow-Tips Black G.20
081 Etch-Tips Light Blue G.23
082 Liner-Tips Pink G.18

DISPOSABLE VALUSTICK STICKY APPLICATORS



- Bends into any angle for placement of Veneers, Inlays, Onlays, Crowns etc
- Strong adhesive for multiple applications
- Fully disposable
- Length: 110 mm
- Colour: White

909 Valustick Sticky Applicators 50 pcs/box

AUTOCLAVABLE ASPIRATOR BRUSHES



- Autoclavable up to 134°C
- Rigid bristles for effective cleaning

020 X-Fine dia. 5 mm, 6 pcs/pack
021 Fine dia. 8 mm, 6 pcs/pack
022 Medium dia. 14 mm, 6 pcs/pack
023 Large dia. 20 mm, 6 pcs/pack


AUTOCLAVABLE SUCTION TIP SLEEVES



- For protection from sharp edges on suction tubes
- Fully autoclavable up to 134°C
- Colour: Light Blue

290 Suction Tip Sleeves 6 pcs/pack

AUTOCLAVABLE SUCTION TUBE ADAPTORS



Saliva Ejector 6.5 mm
306 Saliva Ejector Adaptor reduces 11 mm to 6.5 mm
Standard HVE Hose 11 mm

Suction Tube 11 mm
305 Suction Tube Adaptor reduces 16 mm to 11 mm
European HVE Hose 16 mm

OBLONG ARTICULATING MOUNTING PLATES



- Fit Denar type articulators
- Fully disposable
- Colours: Blue or White

490 Mounting Plates, 100 pcs/box


AUTOCLAVABLE ACRYLIC MIXING CUPS



- Fully autoclavable up to 134°C


280 Mixing Cup Small dia. 25 mm
281 Mixing Cup Medium dia. 35 mm
282 Mixing Cup Large dia. 50 mm
283 Mixing Cup X-Large dia. 70 mm

DISPOSABLE EVACUATION TRAPS



041 Disposable Evacuator Traps dia. 48 mm Green, 20 pcs/pack
042 Disposable Evacuator Traps dia. 54 mm Blue, 20 pcs/pack

AUTOCLAVABLE ALUMINIUM BUR BLOCKS



- Ideal for storage of dental burs
- Specially designed for ultrasonic cleaning and autoclaving of burs
- Sturdy construction, fully autoclavable up to 134°C
- Colours: Silver and Blue

Small Bur Block
038-FG 15 FG W58 mm x D25 mm x H31 mm
039-MIX 10 FG + 5 RA W58 mm x D25 mm x H31 mm
037-HP 6 HP W54 mm x D28 mm x H62 mm
Large Bur Block
2832-FG 32 FG W80 mm x D42 mm x H33 mm
2832-MIX 16 FG + 16 RA W80 mm x D42 mm x H33 mm
2880-FG 80 FG W100 mm x D75 mm x H33 mm
2880-MIX 40 FG + 40 RA W100 mm x D75 mm x H33 mm


DISPOSABLE EVACUATION TRAPS



- Fits Adec type unit

044 Evacuator Traps
20 pcs/pack

AUTOCLAVABLE PLASTIC BUR DISPENSERS



- Specially designed for cleaning and autoclaving of burs
- Locks open or closed for processing
- Medical grade plastic, autoclavable up to 134°C
- Size: 75 mm x 35 mm x 10 mm

050-FG, Bur Dispenser 12 FG
050-MIX Bur Dispenser 6 FG + 6 RA

COLOURS	Yellow
Blue	White
Light Blue	Light Beige
Green	Lavender

AUTOCLAVABLE MAGNETIC BUR STANDS



COLOURS	Yellow
Blue	White
Light Blue	Light Beige
Green	Lavender

- Stands autoclavable up to 134°C
- Fits both FG and RA burs
- Size: Small Square 42 mm x 42 mm x 12 mm
Large Rectangular 80 mm x 42 mm x 12 mm

033 Bur Stand Square, 8 Holes 034 Bur Stand Rectangular, 14 Holes
033A Bur Stand Cover Square 034A Bur Stand Cover Rectangular

AUTOCLAVABLE BRASS BUR BRUSH



- Telescopic handle with brass wire bristles
- Length: 12 cm

1110 Brass Bur Brush

DISPOSABLE CURVED TIP SYRINGES



078
Curved Tip
Syringes
50 pcs/box

- Made of polypropylene plastic
- Full 12 ml capacity
- Self sealing medical grade plunger
- Designed for perio irrigation, impression materials, post surgical home use, fluoride treatments etc

AUTOCLAVABLE FRACTURE DETECTOR

The Instrument for locating cracks in teeth



- Fully autoclavable to 134°C
- Simple tool for practice
- Quick and precise to locate cracks
- Can be used for crown and bridge seating

025 Autoclavable
Fracture Detector Blue

GLASS CEMENT MIXING SLAB



- Polished glass mixing slab with bevelled edges
- Size: 65 mm x 135 mm x 15 mm

045 Glass Slab

COTTON ROLL DISPENSER



- Provides an aseptic barrier against operator aerosols and accidental splatter
- Top flips back for easy loading, sliding drawer dispenses up to 4 rolls with a single pull
- Foam rubber strips inserted into base prevent slipping on table top or counter

1140 Roll Dispenser

COLOURS	Blue	Light Blue	Green
Yellow	White	Light Beige	Lavender



AUTOCLAVABLE IMPLANT SCALER SET

- Fully autoclavable up to 134°C
- Replaceable carbon reinforced resin tips, can be resharpened
- Extra light hollow handle

Implant Scaler H6/H7 Hollow Handle
 1128 Handle + 24 Tips (12 each H6/H7)
 1128A 24 Tips (12 each H6/H7)



AUTOCLAVABLE NAPKIN HOLDERS

090 Napkin Holder Adult (L47 cm)
 090C Napkin Holder Child (L41 cm)
 090XS Napkin Holder Extra Short (L37 cm)

- Medical grade silicone
 - Fully autoclavable up to 134°C
 - **Stainless steel** alligator clips
 - Non-slip design
 - Bright, attractive colours
- Colours available
 Adult: **Light Blue, Blue, Green, Red, Yellow, White, Pink** and **Lavender**
 Child: **Light Blue, Lavender, Blue** and **Pink**
 XS (for neck contour bibs): **Blue** only



DENTAL BIB HOLDER AND CLIP

- Cold disinfection only
- Cute design

091 Clip
 092 Bib Holder, **White** or **Red**



PEN / TOOTHBRUSH HOLDER

- Cute design


2005
 Pen / Toothbrush Holder



JUMPING TEETH

- Cute design
- Economical giveaway for kids

8001 Jumping Teeth
 24 pcs/box



DENTAL MANDRELS STAINLESS STEEL

- HP
- RA
- Pop On Disc
- Snap-On Prophy Cup
- Super-Snap® Type
- OneGloss® Type
- FG/RA Adaptor
- M01
- M02
- M03
- M06
- M04
- M05
- M08

Dental Mandrels, 6 pcs/box
 M01 HP
 M02 RA
 M03 Pop On


M04 Super-Snap® Type
M05 OneGloss® Type
M06 Snap-On Prophy
M08 FG/RA Adaptors



AUTOCLAVABLE X-RAY FILM HOLDER PLASTIC

- One end for anterior x-rays, the other end for posterior
- Both ends securely hold film without bending it
- Improves quality and accuracy of intraoral films
- Helps to standardize the bisecting angle technique
- Colours: **Blue** and **Light Beige**

1188 X-Ray Film Holder



DISPOSABLE BITE WING LOOPS

- Convenient pop-open design ensures quick and easy use
- Manufactured from quality ten point paper stock
- Quick, simple and economical


Bite Wing Loops, 500 pcs/box
 A. 2011-1 Film **Child**
 B. 2011-2 Film **Adult**
 C. 2012-1 Sensor Size 1
 D. 2012-2 Sensor Size 2



DELUXE ACRYLIC COMPOSITE ORGANIZERS

- Smart, flexible design
- Sturdy, durable construction
- Clear material makes contents easy to see
- For storage of composites, adhesive, micro applicators etc.

Composite Organizer
 M060 Mini *W19.5 cm x H9.5 cm x D14.5 cm*
 S060 Standard *W29 cm x H9.5 cm x D14.5 cm*



DISPOSABLE MICRO APPLICATOR DISPENSERS

Applicator Dispenser
 18-900 *Standard Handle*
 18-9900 *Mini Handle*

- Small, sturdy, durable construction
- Provides aseptic barrier against operatory aerosols and accidental splatter

ACRYLIC ORTHO CHAIN DISPENSER




- Long lasting dispenser that will fit the decor of any office
- Chain is dispensed easily with an interior roller
- Holds up to 15 spools



ACRYLIC COMPOSITE DISPENSER

- Smart, durable construction
- Clear material makes contents easy to see
- For storage of 20 composite syringes, etch and bonding material etc.

H060 Composite Dispenser



ACRYLIC COMPOSITE SYRINGE HOLDERS

551 Syringe Holder Small
552 Syringe Holder Medium

- Smart, durable construction
- Size: 551 Small (11 holes) *W10 cm x D12.5 cm x H12 cm*
 552 Medium (23 holes) *W20 cm x D12.5 cm x H12 cm*



LIGHT PROTECTION HAND SHIELD

- Hand held shield for curing light
- Protection for both dentist and assistant
- Cold disinfection only
- Size: *12.5 cm x 26 cm x 0.3 cm*

063 Light Protection Hand Shield



LIGHT PROTECTION COVERS

- Extends working time of composite materials
- Sturdy, durable construction
- Cold disinfection only
- Size: Small *W10 cm x D7.5 cm x H3 cm*
 Medium *W15 cm x D10 cm x H3 cm*

062-0 Light Protection Cover Small
062-1 Light Protection Cover Medium



ACRYLIC IMPRESSION GUN HOLDERS

- For both loaded or empty guns
 - Store all guns in upright position
 - Include storage compartments
 - Sturdy, durable construction
 - Sizes:
 - 072 W145 mm x D140 mm x H150 mm
 - 074 W265 mm x D140 mm x H150 mm
 - 076 W395 mm x D140 mm x H150 mm
- 072 Gun Holder for 2 guns
074 Gun Holder for 4 guns
076 Gun Holder for 6 guns



ACRYLIC GRINDING BOX

- Ideal for grinding and sand blasting
 - Retains flying debris
 - Sturdy durable construction
 - Size: W30 cm x D20 cm x H18 cm
- 065 Grinding Box

ACRYLIC GLOVE AND TISSUE DISPENSERS



- Smart and attractive acrylic design
 - Sturdy, durable construction
 - Pre-drilled holes for wall mounting
 - Internal measurements:-
 - 064 Glove Dispenser **Standard** W25 cm x D14 cm x H9 cm
 - 066 Glove Dispenser **Small** W25 cm x D14 cm x H7 cm
 - 068 Glove/Tissue Dispenser **Large** W24 cm x D12.5 cm x H10 cm
 - 069 Glove/Tissue Dispenser **Extra Large** W24 cm x D12.5 cm x H11 cm
- 064 Glove Dispenser **Standard**
066 Glove Dispenser **Small**
068 Glove/Tissue Dispenser **Large**
069 Glove/Tissue Dispenser **Extra Large**




PLASTIC CUP DISPENSER

- Holds approx. 50 disposable cups from 150 ml to 200 ml (5 oz to 7 oz)
- Self adhesive
- Made of ABS plastic (acrylonitrile butadiene styrene)
- Size: dia. 100 mm x L280 mm

1280
Plastic Cup
Dispenser

ACRYLIC BIB AND C-FOLD TOWEL DISPENSERS




- Sturdy, durable construction
 - Smart, attractive acrylic design
 - Sizes:
 - 057 W15 cm x D11 cm x H12 cm
 - 058 W23.5 cm x D10 cm x H13.5 cm
- 057 Bib Dispenser
058 C-Fold Towel Dispenser
057CV Bib Dispenser With Cover
058CV C-Fold Towel Dispenser With Cover



ACRYLIC MASK DISPENSER

095 Mask Dispenser

- Sturdy, durable construction
- Clear material makes contents easy to see
- Pre-drilled holes for wall mounting
- Size: W94 mm x D106 mm x H108 mm



ACRYLIC POUCH DISPENSERS

- Can be used on table top or wall mounted in either vertical or horizontal position
- Sizes:-
 - 061-0 (7.5 cm x 14.5 cm x 21 cm)
 - 061-1 (8.5 cm x 27.5 cm x 21 cm)
 - 061-2 (10 cm x 18.5 cm x 21 cm)
 - 061-3 (10 cm x 27.5 cm x 21 cm)
 - 061-4 (15 cm x 29.5 cm x 21 cm)

- 061-0 Pouch Dispenser for 5.5 cm x 10 cm Pouches
061-1 Pouch Dispenser for 7 cm x 23 cm Pouches
061-2 Pouch Dispenser for 9 cm x 13.5 cm Pouches
061-3 Pouch Dispenser for 9 cm x 23 cm Pouches
061-4 Pouch Dispenser for 13.5 cm x 26 cm Pouches

ACRYLIC BARRIER FILM DISPENSER



- Sturdy durable construction
- Clear material makes contents easy to see
- Size: dia. W14 cm x D23 cm x H19 cm

1155 Barrier Film Dispenser

premium
plus™



Premium Plus International Ltd.

1001, Yuen Long Trading Centre
No.33, Wang Yip Street West
Yuen Long, N.T.

Hong Kong

Tel: 852-2761-3211
Email: mailbox@premiumplushk.com
Web: www.premiumplushk.com



Premium Plus Poland sp.z o.o.

ul. Bukowska 27
62-081 Wysogotowo

Poland

Tel: 48-61-880-1094
Email: info@premiumpluspl.com

Premium Plus UK Ltd.

2, Knighton Heath Industrial Estate
847 Ringwood Road, Bournemouth
Dorset BH11 8NE

U.K.

Tel: 44-1202-611011
Email: info@premiumplusuk.com
Web: www.premiumplusuk.com

Premium Plus (Dongguan) Limited

No. 1 Industrial Area, Tutang,
Changping, Dongguan,
Guangdong 523581

China

Tel: 86-769-83397277
Email: info@premiumpluschina.com

Premium Plus Dental Supplies Inc.

155 Liberty Avenue
Brooklyn NY 11212

U.S.A.

Tel: 1-718-366-8686
Email: info@premiumplususa.com
Web: www.premiumplususa.com

Premium Plus Japan Co., Ltd.

2nd Floor, 1-4-2, Jonan, Ikedashi
Osaka, 563-0025

JAPAN

Tel: 81-72-735-7101
Email: info@premiumplus.jp
Web: www.premiumplus.jp