

# premium plus™ 2020 Greenline

**Premium Plus International Ltd.**  
1001, Yuen Long Trading Centre  
No.33, Wang Yip Street West  
Yuen Long, N.T.  
**Hong Kong**  
Web: [www.premiumplushk.com](http://www.premiumplushk.com)

**Premium Plus UK Ltd.**  
2, Knighton Heath Industrial Estate  
847 Ringwood Road, Bournemouth  
Dorset BH11 8NE  
**U.K.**  
Web: [www.premiumplusuk.com](http://www.premiumplusuk.com)

Welcome to our new Greenline. We believe that making simple changes to the everyday products you purchase is the most sustainable way of helping our planet. This range is our first step towards offering you the choice of more environmentally friendly products.

## ez interdental brushes


After scratching our heads, we've created our interdental brushes with a few 'ez' changes:

- 1 We've reduced the size of our handle to a smaller, more practical design – this uses up to 50% less plastic than other leading brands.
- 2 Instead of using new plastic, we use 100% recycled plastic for the brush handle and 50% recycled plastic for the protective cap. This recycled plastic is made using waste materials left over from other manufacturing processes.
- 3 Unlike other brands who package interdental brushes in unnecessary plastic boxes and bags, we use cardboard which is 100% plastic free and recyclable.




- ISO size (colour coded)
- Latex and BPA free
- Vegan Friendly
- Reusable
- Non-slip, easy grip design
- Plastic free packaging
- Small, compact and easy to carry

*Interdental Brushes, ISO Size No. 1-6  
Single or Assorted Sizes, 6 pcs/box*



ISO Size	1	2	3	4	5	6
Reference Size	XX Fine	X Fine	Fine	Fine -Medium	Medium	Large
Made to fit hole (mm)	<0.8	0.9-1.0	1.1-1.2	1.3-1.5	1.6-1.8	>1.9
Wire Dia. (mm)	0.35	0.35	0.45	0.45	0.45	0.45
Order Code	G800-1	G800-2	G800-3	G800-4	G800-5	G800-6
Assorted Order Code	G800-A No.1-6 Assorted					

# Greenline biodegradable PSP barrier envelopes



- Fits most phosphor plates
- Fully disposable envelopes
- Soft vinyl with rounded corners
- Environmentally friendly packaging

In order to improve on an already great product and in an effort to reduce waste we've:

- 1 Made our product from oxy-biodegradable plastic, which means the material will breakdown more easily than conventional plastics.
- 2 Removed the plastic packaging and supplied them in a 100% recyclable paper box, which conserves natural resources like timber and water.

#### Mid-Opening Barrier Envelopes

G188-0 No.0 (I.D. 39 mm x 28 mm) 100 pcs/box  
 G188-1 No.1 (I.D. 44 mm x 29 mm) 100 pcs/box  
 G188-2 No.2 (I.D. 45 mm x 35 mm) 100 pcs/box

#### Top-Opening Barrier Envelopes

G198-0 No.0 (I.D. 27 mm x 40 mm) 100 pcs/box  
 G198-1 No.1 (I.D. 28 mm x 45 mm) 100 pcs/box  
 G198-2 No.2 (I.D. 34 mm x 46 mm) 100 pcs/box

# Greenline soft touch™ reusable oral evacuators vented (slotted) latex free



We've kept all the features of our market leading SOFT TOUCH™ oral evacuator BUT added a few greener touches.

- 1 Our newest version is reusable - clean it, autoclave it and use it at least 10 times. This reduces the production and disposal of a lot of plastic.
- 2 To avoid disposing of the tube too early we have built in an indicator - when the white soft edge turns yellow its time to dispose of it.
- 3 Instead of a plastic bag we've used a 100% recyclable paper box, reducing the amount of waste sent to landfill and incinerators.



- Autoclavable up to at least 10 times
- Dispose of the evacuator when soft white tip turns light yellow
- Soft Touch™ edge eliminates the risk of lacerating the oral tissue
- Rigid tubing, strong enough for cheek retraction
- Vented to reduce unnecessary tissue grab
- Thoroughly polished surface for smooth and quiet operation

G-ST003 Greenline Oral Evacuators, 20 pcs/pack  
 G-022 Aspirator Cleaning Brush 14 mm

## Greenline micro applicators mini

In an effort to limit the environmental effect of one of the most used items in a dental surgery we came up with some great ideas:

- 1 We designed the handle to be a smaller, more practical size that reduces the amount of plastic used by 50% in comparison to competitive products.
- 2 Instead of new plastic we used 100% recycled plastic. This recycled plastic is made using waste materials from other manufacturing processes.
- 3 No more plastic tube that can use up to a third of the amount of plastic as the actual product - just a 100% recyclable box.



- Bends into any angle for the placement of small controlled amounts of material in limited access areas.
- Economical - saving up to 80% by the reduction of excess materials used.

**Greenline Micro Applicators Mini, 100 pcs/box**  
 ( Length 85 mm, dia 2.5 mm handles )  
 G-9901 Regular ( dia. 2.5 mm ) Light Green  
 G-9902 Fine ( dia. 2 mm ) Medium Green  
 G-9903 Ultra Fine ( dia. 1.5 mm ) Dark Green

## Greenline micro applicator tips

Our most environmentally friendly micro applicators:

- 1 By using these tips you save 80% of the plastic that you would normally use in the handle of a standard micro applicator.
- 2 An autoclavable handle saves plastic. Clean it, autoclave it and re-use it. Simple and effective.
- 3 100% recyclable environmentally friendlier box. Recycling means preventing pollution by reducing the need to collect new raw materials.
- 4 Instead of new plastic we use 100% recycled plastic. The recycled plastic is made using waste materials from other manufacturing processes.



- For the placement of small controlled amounts of material in limited access areas.
- Economical - saving up to 80% by the reduction of excess materials used.

**Greenline Micro Applicator Tips, 100 pcs/box**  
 G-921 Regular ( dia. 2.5 mm ) Light Green  
 G-922 Fine ( dia. 2 mm ) Medium Green  
**Autoclavable Micro Applicator Tip Handle**  
 7502W White  
 7502T Transparent



## Greenline air-water syringe tips

Reusable products are a great idea. Unfortunately, trying to clean the inside of a reusable metal air-water syringe tip is just not possible and if you can't clean it, you can't guarantee the effectiveness of its sterilisation. This means that using a disposable plastic air-water syringe tip is the only way to guarantee avoiding cross contamination.

We have however created a more environmentally friendly air-water syringe tip by making a few changes.

- Directly interchangeable with American type syringe tips (long), or European syringes which have been converted to use American type syringe tips (short)
- Made from 100% recycled plastic
- No sharp edges for patient comfort
- Unique quad air core design to ensure dry air
- Environmentally friendly packaging

**Greenline Air-water Syringe Tips**  
**G-018 (American Type) 200 pcs/box**  
**G-028 (European Type) 200 pcs/box**

- 1 Our Greenline air-water syringe tips are made from 100% recycled plastic, using waste materials left over from other manufacturing processes.
- 2 The Greenline air-water syringe tips are packaged in 100% recyclable cardboard boxes, instead of large plastic bags, helping to reduce the amount of plastic packaging in your surgery.




## Greenline bamboo cups

Looking to cut down on plastic use? Switching to bamboo cups is a simple and easy change that can have a big impact on reducing plastic waste in your surgery. Here are some reasons to make the switch to using bamboo:

- Leakproof construction, non-slip design
- Dispenses easily from cup dispenser due to their rolled lip
- Conveniently designed 170ml (6oz) size

**Bamboo Cups 170 ml ( 6 oz ), 1000 pcs/box**  
**6800BB Blank**  
**6800BH Healthy Teeth ( 4 assorted designs )**

### Healthy Teeth Designs



- 1 Bamboo is the world's quickest growing plant and provides a sustainable source of material. When bamboo is cut down to make paper, it is able to quickly recover to be harvested again.
- 2 Bamboo absorbs a large amount of CO<sup>2</sup> emissions as it grows, which helps to offset human carbon emissions.
- 3 Our Greenline outer carton is 100% recyclable and the inner sleeves, which are needed to keep the product sanitary, are biodegradable.