

High Quality Needles | Painless Steel®

Illustrations scale 1:1 - REF: Transcodent item numbers

transcodent®

A Sulzer Brand

Painless Steel® Needles


- Siliconized cannula for smooth gliding
- 3-edge lancet grinding
- Plastic hub, for use with cylinder ampoule syringes
- G 30 for the intraligamentary anaesthesia
- Precisely threaded hub for secure fit and easy handling
- Heat-sealed container
- EO sterilized



Painless Steel® Injection

Sizes
G 25 / 0.5 mm
G 27 / 0.4 mm
G 30 / 0.3 mm

Packaging
100 pieces per box



Paste Needles

Sizes
G 20 / 0.9 mm
G 21 / 0.8 mm
G 23 / 0.6 mm

Packaging
100 pieces per box



Luminject® Thin Wall Needles

- G 30, G 31 and G 33 intraligamentary and infiltration anaesthesia.
- Especially recommended for the use in pediatrics and for sensitive patients
- Hub colour coding according to DIN EN ISO 6009
- Thin-Wall design for increased flow




Increased flow. Enlarged lumen.

G 31/G 33: extra thin

Luminject® Thin Wall Injection

Sizes
G 27 / 0.4 mm
G 30 / 0.3 mm
G 31 / 0.28 mm
G 33 / 0.26 mm


Packaging
100 pieces



Injection Luer

Sizes
G 25 / 0.5 mm
G 27 / 0.4 mm
G 30 / 0.3 mm

Packaging
100 pieces per box



Paste Needles

- Threaded hub (metric)
- Suitable for cylinder ampoule syringes
- Rounded tip for comfortable usage in the gingival sulcus
- Bendable without significant decrease of the flow rate
- EO sterilized

Intraligamentalia®

- Cylinder ampoule syringe for the intraligamentary and intraseptal anaesthesia
- High quality stainless steel. autoclavable and sterilizable
- Single tooth anaesthesia for ~30 minutes. Multiple quadrants within one session
- Comfortable and precise dosage: 0.06 ml per click
- Pressure limitation

Intraligamentalia®

Guarantee
1 year

Material
Stainless Steel

Packaging
1 piece per box



Handpiece for Luer Lock Suction Needles

Guarantee
2 years

Material
Stainless Steel


Packaging
1 piece per box



For use with suction needles – sterilizable. Included adapter, usable with suction unit adapter.

Suction Needles Luer Lock

- Thin-Wall tubes
- Annealed steel
- Easy to bend
- For local suction during endodontic and surgical treatments
- For usage with Luer Lock handpiece




Annealed steel. Easy to bend!

Suction Needles Luer Lock

Sizes
G 14 / 2.1 mm
G 16 / 1.6 mm
G 18 / 1.2 mm
G 20 / 0.9 mm

Packaging
100 pieces




Apical Side Vent

Endo Evacuation Needles

Sizes
G 27 / 0.4 mm
G 30 / 0.3 mm

Packaging
25 pieces per box



Laser Marks

Endo Evacuation Needles

- Apical side vent
- Four laser marks
- 32 mm length
- Easy to bend
- For usage with Luer Lock handpiece



transcodent.com

You do precision work.
We create your tools.

High Quality Needles | Micro-Application

Illustrations scale 1:1 - REF: Transcodent item numbers

Endo Irrigation Needles

- Soft and effective rinsing of root canal & dental pocket
- Luer Lock hub
- Container-packaging
- Thin-Wall – extra thin tube for more flexibility
- EO sterilized



Double Side Vent

Irrigation Luer Lock Endo Thin Wall

Sizes
G 23 / 0.6 mm
G 25 / 0.5 mm
G 27 / 0.4 mm
G 30 / 0.3 mm
G 31 / 0.28 mm
PU 25 pieces



Double Side Vent

Irrigation Luer Lock Endo Thin Wall


Sizes
G 27 / 0.4 mm
Packaging
25 pieces per box



Single Side Vent

Irrigation Luer Lock Endo Thin Wall

Sizes
G 25 / 0.5 mm
G 30 / 0.3 mm
Packaging
25 pieces per box



Double Vent

- Two lateral openings
- Closed end / rounded tip
- Side openings prevent injuries of the apex

Single Vent

- One lateral opening
- Side opening and rounded tip for apex protection
- Closed end | rounded tip

Paste Needles

- Luer Lock hub
- EO sterilized
- Rounded tip for usage in dental pockets
- Bendable without significant decrease of the flow rate


Paste Needles Luer Lock

Sizes
G 20 / 0.9 mm
G 23 / 0.6 mm
Packaging
100 pieces per box



Irrigation blunt Luer Lock Endo

Sizes
G 25 / 0.5 mm
G 27 / 0.4 mm
G 30 / 0.3 mm
Packaging
100 pieces per box



Elephant-Tip™

Dimensions
OD: 1.3 mm
ID: 0.6 mm
Material
Polypropylene
Packaging
100 pieces per box




Elephant-Tip™

- Luer Lock Application tip for gel & paste
- Two versions: straight | bent 45°
- Etching gel
- Bleaching gel
- Fissure sealant

Sterile single-use Syringes

Sizes Luer
2 ml, 5 ml, 10 ml
Size Luer Lock
5 ml
Packaging
100 pieces per box



Luer Lock:
Recommended for Endo Irrigation Needles



Luer:
Recommended for Luer Injection Needles