


REFINE MaxPiezo 7+ / DS7+ PORTABLE LED SCALERS



- Automatic frequency tracking, real-time feedback
- Extra large, capacity 600 ml water tank for water, chlorhexidine, hydrogen peroxide, sodium hypochlorite etc...
- Power and water controlled by touch panel, easy to clean for infection control
- Detachable handpiece
- Frequency: 28 kHz \pm 3 kHz
- Complete with 8 tips
- Size: W180 mm x D180 mm x H150 mm
- Weight: 1.6 Kg (unit only)
- Colour: White

ONE YEAR WARRANTY
(except scaler tips)

MaxPiezo 7 (EMS Type) Portable **Non-LED** Scaler Set
 MaxPiezo 7+ (EMS Type) Portable **LED** Scaler Set
 MaxPiezo DS7 (SATELEC Type) Portable **Non-LED** Scaler Set
 MaxPiezo DS7+ (SATELEC Type) Portable **LED** Scaler Set



1. EMS is a registered trademark of Electro Medical System S.A. Switzerland
 2. SATELEC is a registered trademark of Satelec/Acteon Group France


REFINE MaxPiezo 6+ / DS6+ LED SCALERS



- Automatic frequency tracking, real-time feedback
- Power and water controlled by touch panel, easy to clean for infection control
- Detachable handpiece
- Frequency: 28 kHz \pm 3 kHz
- Complete with 6 tips
- Size: W155 mm x D168 mm x H80 mm
- Weight: 0.6 Kg (unit only)
- Colour: White

ONE YEAR WARRANTY
(except scaler tips)

MaxPiezo 6 (EMS Type) **Non-LED** Scaler Set
 MaxPiezo 6+ (EMS Type) **LED** Scaler Set
 MaxPiezo DS6 (SATELEC Type) **Non-LED** Scaler Set
 MaxPiezo DS6+ (SATELEC Type) **LED** Scaler Set



1. EMS is a registered trademark of Electro Medical System S.A. Switzerland
 2. SATELEC is a registered trademark of Satelec/Acteon Group France

REFINE N3000L / V3000L BUILT-IN LED SCALERS



- Automatic frequency tracking, real-time feedback
- Detachable handpiece
- Frequency: 28 KHz \pm 3 KHz
- Complete with 6 tips

ONE YEAR WARRANTY
(except scaler tips)

N3000 (EMS Type) Built-In **Non-LED** Scaler Set
 N3000L (EMS Type) Built-In **LED** Scaler Set
 V3000 (SATELEC Type) Built-In **Non-LED** Scaler Set
 V3000L (SATELEC Type) Built-In **LED** Scaler Set



1. EMS is a registered trademark of Electro Medical System S.A. Switzerland
 2. SATELEC is a registered trademark of Satelec/Acteon Group France

REFINE AUTOCLAVABLE SCALER HANDPIECE KITS



- Complete with:-
- 1x Handpiece 3H / 5L / 7H / 7L
- 3x Perio scaler tip P1 / PD1
- 1x Torque wrench TW-1L
- 1x Scaler tip holder No.1
- 1x Sterilization box



3H Non-LED Handpiece Kit (EMS Type)
7H Non-LED Handpiece Kit (SATELEC Type)
5L LED Handpiece Kit (EMS Type)
7L LED Handpiece Kit (SATELEC Type)
 (for handpiece connections, please refer to Page 37)

ONE YEAR WARRANTY
(except scaler tips)

1. EMS is a registered trademark of Electro Medical System S.A. Switzerland
 2. SATELEC is a registered trademark of Satelec/Acteon Group France

Scaler Handpieces and Accessories

REFINE SCALER HANDPIECES



HP-3H Non-LED Scaler Handpiece (EMS¹ Type)



HS-7H Non-LED Scaler Handpiece (SATELEC² Type)



HP-5L LED Scaler Handpiece (EMS¹ Type)



HS-7L LED Scaler Handpiece (SATELEC² Type)

1. EMS is a registered trademark of Electro Medical System S.A. Switzerland
2. SATELEC is a registered trademark of Satelec/Acteon Group France
3. NSK is a registered trademark of Nakanishi Inc. Japan

REFINE SCALER TORQUE WRENCHES



TW-1L
For EMS¹,
KAVO PiezoLED⁵
SATELEC², NSK³
type scalers



TS-1L
For SIRONA⁴
type scaler



TK-1L
For KAVO⁵
SONICflex⁵
SONOsoft⁵
type scalers

1. EMS is a registered trademark of Electro Medical System S.A. Switzerland
2. SATELEC is a registered trademark of Satelec/Acteon Group France
3. NSK is a registered trademark of Nakanishi Inc. Japan
4. SIRONA is a registered trademark of Sirona Dental System
5. KAVO, SONICflex, SONOsoft, PiezoLED are registered trademarks of Kavo Dental GmbH Germany

REFINE SCALER HANDPIECE HOSES



3H Non-LED (EMS¹ Type)
Handpiece Hose



7H Non-LED (SATELEC² Type)
Handpiece Hose




5L LED (EMS¹ Type)
Handpiece Hose



7L LED (SATELEC² Type)
Handpiece Hose

1. EMS is a registered trademark of Electro Medical System S.A. Switzerland
2. SATELEC is a registered trademark of Satelec/Acteon Group France

REFINE ENDO WRENCH



Multi function wrench
for both endo chuck & bur

TW-A

REFINE FOOT SWITCH



REFINE SCALER TIP HOLDER



NO. 2

E04 ULTRASONIC FILES NiTi alloy

- For use with ultrasonic scaler
- Excellent tenacity
- Length: 33 mm

E04 Ultrasonic Files
Single or
Assorted Sizes
6 pcs/pack

E04 ULTRASONIC FILES NiTi alloy	
#15 ○	#30 ●
#20 ●	#35 ●
#25 ●	#40 ●


REFINE SCALER TIPS STAINLESS STEEL

FOR EMS¹, SATELEC², NSK³ and SIRONA⁴ REFINE Type Scalers


SCALING




G1 (EMS¹ Type)
GD1 (SATELEC², NSK³ Type)
GS1 (SIRONA⁴ Type)
 Supragingival deposits in all quadrants




G2 (EMS¹ Type)
GD2 (SATELEC², NSK³ Type)
GS2 (SIRONA⁴ Type)
 Heavy supragingival deposits on anterior teeth




G3 (EMS¹ Type)
GD3 (SATELEC², NSK³ Type)
GS3 (SIRONA⁴ Type)
 Supragingival deposits in all quadrants, including the interproximal and sulcus areas



G4 (EMS¹ Type)
GD4 (SATELEC², NSK³ Type)
GS4 (SIRONA⁴ Type)
 Supragingival and interproximal deposits




G5 (EMS¹ Type)
GD5 (SATELEC², NSK³ Type)
GS5 (SIRONA⁴ Type)
 Supragingival deposits in all quadrants




G6 (EMS¹ Type)
GD6 (SATELEC², NSK³ Type)
GS6 (SIRONA⁴ Type)
 Voluminous supragingival deposits


ENDO



E1 (EMS¹ Type)
ED1 (SATELEC², NSK³ Type)
ES1 (SIRONA⁴ Type)
 120° file holder for anterior root canal irrigation



E2 (EMS¹ Type)
ED2 (SATELEC², NSK³ Type)
ES2 (SIRONA⁴ Type)
 95° file holder for molar teeth root canal irrigation




E3 (EMS¹ Type)
ED3 (SATELEC², NSK³ Type)
ES3 (SIRONA⁴ Type)
 Gutta-percha melting




E3D (EMS¹ Type)
ED3D (SATELEC², NSK³ Type)
ES3D (SIRONA⁴ Type)
 Diamond coated instrument for root canal enlarging




P4D (EMS¹ Type)
PD4D (SATELEC², NSK³ Type)
PS4D (SIRONA⁴ Type)
 Diamond coated instrument for root canal enlarging




E14 (EMS¹ Type)
ED14 (SATELEC², NSK³ Type)
ES14 (SIRONA⁴ Type)
 Removing root filling obstructions during root canal retreatments




E15 (EMS¹ Type)
ED15 (SATELEC², NSK³ Type)
ES15 (SIRONA⁴ Type)
 Removing root filling obstructions during root canal retreatments




E4 (EMS¹ Type)
ED4 (SATELEC², NSK³ Type)
ES4 (SIRONA⁴ Type)
 Removing root filling obstructions during root canal retreatment



E4D (EMS¹ Type)
ED4D (SATELEC², NSK³ Type)
ES4D (SIRONA⁴ Type)
 Removing calcification and filling material during root canal retreatment



E5 (EMS¹ Type)
ED5 (SATELEC², NSK³ Type)
ES5 (SIRONA⁴ Type)
 Removing root filling obstructions during root canal retreatment



E5D (EMS¹ Type)
ED5D (SATELEC², NSK³ Type)
ES5D (SIRONA⁴ Type)
 Removing calcification and filling material during canal retreatment


1. EMS is a registered trademark of Electro Medical System S.A. Switzerland
 2. SATELEC is a registered trademark of Satelec/Acteon Group France
 3. NSK is a registered trademark of Nakanishi Inc. Japan
 4. SIRONA is a registered trademark of Sirona Dental System

REFINE SCALER TIPS STAINLESS STEEL

FOR EMS¹, SATELEC², NSK³ and SIRONA⁴ Type Scalers




P1 (EMS¹ Type)
PD1 (SATELEC², NSK³ Type)
PS1 (SIRONA⁴ Type)
 Shallow pocket subgingival deposits



P3 (EMS¹ Type)
PD3 (SATELEC², NSK³ Type)
PS3 (SIRONA⁴ Type)
 Deep pocket subgingival deposits



P4 (EMS¹ Type)
PD4 (SATELEC², NSK³ Type)
PS4 (SIRONA⁴ Type)
 Shallow pocket subgingival deposits



P3D (EMS¹ Type)
PD3D (SATELEC², NSK³ Type)
PS3D (SIRONA⁴ Type)
 Diamond coated instrument for root planning after periodontal flap surgery

PERIO



P2L (EMS¹ Type)
PD2L (SATELEC², NSK³ Type)
 Left angle instrument to remove calculus from very narrow inter-root space and furcation



P2R (EMS¹ Type)
PD2R (SATELEC², NSK³ Type)
 Right angle instrument to remove calculus from very narrow inter-root space and furcation



P2LD (EMS¹ Type)
PD2LD (SATELEC², NSK³ Type)
 Removal of calculus from narrow inter-root spaces (diamond-coated)



P2RD (EMS¹ Type)
PD2RD (SATELEC², NSK³ Type)
 Diamond coated instrument for root planning

FOR AMDENT⁶ Type Scalers



A1 (AMDENT⁶ Type)
 Supragingival scaling (#37)



A2 (AMDENT⁶ Type)
 Supragingival scaling (#39)

FOR KAVO SONICflex⁵ Type Scalers




GK1 (KAVO⁵ Type)
 Supragingival deposits in all quadrants



GK2 (KAVO⁵ Type)
 Supragingival heavy calculus and plaque



GK3 (KAVO⁵ Type)
 Supragingival interdental calculus and calculus at the neck of the tooth




GK4 (KAVO⁵ Type)
 Supragingival deposits, including the interproximal and sulcus areas



GK5 (KAVO⁵ Type)
 Subgingival calculus




GK6 (KAVO⁵ Type)
 Left angle instrument for removing subgingival calculus




GK7 (KAVO⁵ Type)
 Right angle instrument for removing subgingival calculus

FOR KAVO SONOsoft⁵ Type Scalers



GC1 (KAVO⁵ Type)
 Supragingival deposits in all quadrants



GC2 (KAVO⁵ Type)
 Supragingival interdental calculus and calculus at the neck of the tooth



PC1 (KAVO⁵ Type)
 Subgingival calculus, including the interproximal and sulcus areas

KAVO PiezoLED⁵ Type Scalers



201 (KAVO⁵ Type)
 Supragingival and interproximal deposits



202 (KAVO⁵ Type)
 Supragingival deposits in all quadrants, including the interproximal and sulcus areas



203 (KAVO⁵ Type)
 Deep pocket subgingival deposits

1. EMS is a registered trademark of Electro Medical System S.A. Switzerland
 2. SATELEC is a registered trademark of Satelec/Acteon Group France
 3. NSK is a registered trademark of Nakanishi Inc. Japan
 4. SIRONA is a registered trademark of Sirona Dental System
 5. KAVO, SONICflex, SONOsoft, PiezoLED are registered trademarks of Kavo Dental GmbH Germany
 6. AMDENT is a registered trademark of Amdent AB Sweden



Premium Plus International Ltd.
1001, Yuen Long Trading Centre
No.33, Wang Yip Street West
Yuen Long, N.T.
Hong Kong
Tel: 852-2761-3211
Email: mailbox@premiumplushk.com
Web: www.premiumplushk.com



Premium Plus Poland sp. z o.o.
ul. Bukowska 27
62-081 Wysogotowo
Poland
Tel: 48-61-880-1094
Email: info@premiumpluspl.com
Web: www.premiumpluspl.com

Premium Plus UK Ltd.
2, Knighton Heath Industrial Estate
847 Ringwood Road, Bournemouth
Dorset BH11 8NE
U.K.
Tel: 44-1202-611011
Email: info@premiumplusuk.com
Web: www.premiumplusuk.com

Premium Plus (Dongguan) Limited
No. 1 Industrial Area, Tutang,
Changping, Dongguan,
Guangdong 523581
China
Tel: 86-769-83397277
Email: info@premiumpluschina.com

Premium Plus Dental Supplies Inc.
155 Liberty Avenue
Brooklyn NY 11212
U.S.A.
Tel: 1-718-366-8686
Email: info@premiumplususa.com
Web: www.premiumplususa.com

Premium Plus Japan Co., Ltd.
5F, 1-1-31, Nishimidorigaoka,
Toyonakashi, Osaka, 560-0005
JAPAN
Tel: 81-6-6845-0066
Email: info@premiumplus.jp
Web: www.premiumplus.jp