

PREMIUM DENTAL PRODUCTS

PRODUCT CATALOG

-
The complete assortment of
Cavex premium products

EN
edition

YOUR
IMPRESSION
IS OUR
SPECIALTY

CAVEX

PREMIUM ALGINATE SYSTEM

THE BEST ALGINATES IN THE WORLD

Proven by facts, stated by dentists



YOUR IMPRESSION IS OUR SPECIALTY

ABOUT US

Since 1908, Cavex has been developing and manufacturing premium solutions for the dental industry. Right from the start, all the company's operations have been performed in-house in Haarlem, The Netherlands. Every new Cavex product is a result of intensive collaboration between R&D, Regulatory, Production, Marketing, Sales and Logistics. This enables Cavex to develop, produce and distribute high-quality products all over the world.

MANUFACTURER FOR DENTAL PROFESSIONALS SINCE 1908

With Haarlem as its permanent home base, Cavex has been inventing and manufacturing innovative solutions within the dental industry since 1908. We have been known as a manufacturer of top products with guaranteed quality from the offset. Find out how Cavex' story began when two ambitious Haarlem pharmacists laid the foundations for one of the largest dental product manufacturers in the world.

HAARLEM BASED PHARMACISTS KEUR & SNELTJES

Back in 1908, there was one single solution to dental problems: extraction. The Haarlem pharmacists Keur and Sneltsjes were convinced that extraction could be made a lot more pleasant for patients under local anesthesia. Hence the first product for dentistry was born: Selectocaine anesthetic fluid.

'CAVEX SILVER AMALGAM'

Over time and due to the increasing availability of dental materials, the goal of treatment shifted from 'extraction' to 'conservation and restoration'. The established dental product company Keur & Sneltsjes developed a revolutionary filler based on amalgam. As tooth decay (cavity) now belonged to the past (ex), the name Cavex was born. Later, the company split into two brands: 'Keur & Sneltsjes' and 'Cavex Holland'.

YOUR IMPRESSION IS OUR SPECIALTY

Today, we have grown to become a leading manufacturer in the world of dental products. We help dental professionals to optimize daily performance by supplying them reliable, top quality products. In addition to supporting dental professionals globally, Cavex offers consumers effective and safe home-solutions in the field of halitosis and whitening. Stating that 'Your Impression Is Our Specialty', our dental products also enable consumers to always make a perfect impression. Ever since the introduction of the CA37 alginate, the golden combination of human determination and professional perfectionism has transformed us into a benchmark in several areas in the field of dentistry.

THE POWER OF COLLABORATION

We are in close contact with dentists, universities and opinion leaders to develop new products, but also to perfect existing products. In this way Cavex remains a leader in the dental industry and we continue to be able to support dental professionals worldwide in their daily work.

CAVEX' SECRET

Our experienced team of well-trained people has always been the key to our success. As the world of dentistry is changing rapidly, our professionals adapt just as quickly – if not faster. We maintain close partnerships with dealers all over the world, providing them the best products and support with promotions and product trainings. To keep thriving in the field of innovative product development and optimize our product line, we are closely connected to dentists, universities, and opinion leaders. This is how we maintain our leader status in dentistry and how we continue to be able to support dental professionals worldwide in their daily work.

We value the opinion of the dentist. That's why we let The Dental Advisor regularly test our products. After clinical evaluations by a large team of dentists, the tested products receive a score of up to five pluses. Several Cavex products, including Cavex Cream Alginate, Cavex ColorChange and the Cavex Bite&White ABC Masterkit, received very high ratings. Cavex Cream even achieved the title 'Top Award Winner' for several years in a row.



YOUR IMPRESSION IS OUR SPECIALTY

CAVEX

CAVEX CREAM ALGINATE



Cavex Cream Alginate is the latest addition to the Cavex alginate range. Due to a combination of 5 µm accuracy, high tear strength and excellent elasticity, the properties of this dental impression material challenge those of silicone quality. Add extremely rapid water absorption, a user-friendly creamy consistency during mixing, snap-set, and five-year shelf life. Cavex Cream Alginate is an unrivaled alginate.

This alginate impression comes with a beyond pleasant bubble-gum scent to sink one's teeth into, hence avoiding patient bite resistance. Cavex Cream Alginate is available in Normal (90 sec.) and Fast Set (60 sec.).

Features

- Detail reproduction of at least 5 µm
- Easy mixing
- Highly elastic
- High tear strength

Art. Nr.	Description	Content
AA600	Pack Cavex Cream Alginate NS	20 x 500 gr
AA605	Pack Cavex Cream Alginate FS	20 x 500 gr

CAVEX IMPRESSIONAL

Cavex Impressional enables you to take alginate impressions down to the last detail. The alginate has a unique combination of very high elasticity, tear strength, and detail reproduction (25 µm). Consequently, this alginate product is very suitable for working around metal frames and/or deep undercuts.

Features

- Maximum elasticity and tear strength
- Clearly legible blue color
- Ideal for the alginate hydrocolloid technique
- Optimal gypsum compatibility

Cavex Impressional is available in Normal (90 sec.) and Fast Set (60 sec.). Due to the combination of fast setting and a pleasant spearmint aroma, patients enjoy more pleasant procedures.

Art. Nr.	Description	Content
AA096	Pack Impressional NS	20 x 500 gr
AA097	Pack Impressional FS	20 x 500 gr



CAVEX COLORCHANGE



User-friendly impression material with virtual intelligence: Cavex ColorChange. This product belongs to the latest generation of Cavex alginates. Color changes indicate different procedural stages.

Features

- Indicates three color changes
- Short setting time
- Highly elastic
- Maximum compatibility with plaster

Does the alginate's color change from violet to pink? Mixing time is over. Does the alginate change from pink to white as the impression is taken? The alginate has hardened, and the impression tray can be removed from the patient's mouth. After thirty minutes, the completely cured alginate turns pink again.

Art. Nr.	Description	Content
AA323	Pack ColorChange FS	20 x 500 gr

CAVEX ORTHOTRACE

Cavex Orthotrace is an orthodontic alginate impression material with a natural fruit scent. The product has been developed based on the needs of orthodontists specifically. Additionally, Cavex Orthotrace contains many features that provide patients with a more pleasant procedure.

avoided. Moreover, the combination of pleasant scent and color invites one to sink one's teeth into the impression material.

Features

- Extra fast setting
- Pleasant fruit aroma
- Highly elastic
- Clearly legible pink color

Due to reduced setting time, Cavex Orthotrace alginate does not get into the back of the mouth. Hence, by using Cavex Orthotrace a gag reflex in patients is

Art. Nr.	Description	Content
AA380	Pack Orthotrace EFS	20 x 500 gr



CAVEX CA37

Cavex CA37 is an all-purpose alginate that is used by tens of thousands of dental clinics worldwide daily. It is the very first alginate impression material to become a worldwide benchmark.

This strong alginate is available in Normal (90 sec.) and Fast Set (60 sec.). Due to the short setting time and a pleasant peppermint flavor, patients will experience a pleasant procedure. Furthermore, Cavex CA37 is very suitable for prosthetic work.

Features

- Optimal tear strength and deformation resistant
- Peppermint flavor
- Natural pink color
- Compatible with all plaster types



Art. Nr.	Description	Content
AA075	Pack CA37 NS	20 x 500 gr
AA055	Pack CA37 NS	24 x 453 gr
AA025	Pack CA37 FS	20 x 500 gr
AA035	Pack CA37 FS	24 x 453 gr

CAVEX ALGINATE MIXER 3

The Cavex Alginate Mixer 3 allows dental professionals to maximize impression results by standardizing a homogeneous alginate mixture. Thanks to optimized mixing technology and smart design, this device is a showpiece in every dental clinic.

Hand-mixed alginates cause product variations. Differentiations in mixing time, intensity, and usage of various techniques affect both curing time and alginate quality. The Cavex Alginate Mixer 3 eliminates these variables and consistently delivers evenly mixed alginates with perfectly predictable

properties. In addition to consistent quality, the Cavex Alginate Mixer 3 also offers significant time savings. Preparation, processing, and cleaning become a matter of just seconds.

4 menu's

- Extra Fast Set alginate
- Fast Set alginate
- Normal Set alginate
- Gypsum

Art. Nr.	Description
	Cavex Alginate Mixer 3, 4 x mixing cup, spatula, mixingrod, water dosing bottle, accessory holder
AT540	incl. 10 x 500 gr Cream Alginate NS
AT541	incl. 10 x 500 gr Cavex CA37 NS
AT542	incl. 10 x 500 gr Cavex Impressional NS
AT543	incl. 10 x 500 gr Cavex ColorChange
AT544	incl. 10 x 500 gr Cavex Orthotrace



CAVEX ALGINATE MIXER II

Improve impression results with perfectly mixed alginate at the touch of a button: Cavex Alginate Mixer II. Mixing alginate into a smooth paste requires a lot of experience. The user's mixing skills influence the texture, consistency, and hence the quality of the alginate mix. Less accurately mixed alginates result in less accurate dental impressions.

The Cavex Alginate Mixer II enhances dental impression performance by standardizing

the alginate mixing procedure and eliminating even the slightest occasional imperfections. Consequently, dental professionals are assured of a smooth and homogeneous mixture, free of lumps or air bubbles. Time after time.

Features

- Hygienic
- User friendly
- Easy to clean
- More accurate impression

Art. Nr.	Description
	Cavex Alginate Mixer II, 4 x mixing cup, spatula, mixingrod, water dosing bottle, accessory holder
AT520	Incl. 10 x 500 gr Cavex Cream NS
AT521	incl. 10 x 500 gr Cavex CA37 NS
AT522	incl. 10 x 500 gr Cavex Impressional NS
AT523	incl. 10 x 500 gr Cavex ColorChange
AT524	incl. 10 x 500 gr Cavex Orthotrace



Accessories for Cavex Alginate Mixer

AT530	Mixing bowl for Cavex Alginate Mixer II	1 pce
AT550	Mixing bowl for Cavex Alginate Mixer 3	1 pce
AT511	Mixing spatula for Cavex Alginate Mixer	1 pce
AT513	Mixing rod for Cavex Alginate Mixer	1 pce

CAVEX ALGINATE ACCESSORIES

Various external factors can influence the quality of alginate impressions. This includes the water temperature, mixing ratio, mixing technique and storage conditions. Cavex has developed various alginate accessories to optimize and standardize these external factors. With our alginate accessories, you can always rely on a perfect alginate impression.

The following accessories are available:

- The alginate container for controlled storage and exact dosing of the alginate powder
- The water dosing bottle for a controlled, constant water temperature and exact dosing of the water
- The flexible mixing bowl with mixing spatula (spatula also available in metal)



Art. Nr.	Description	Content
AT021	Cavex Water dosing bottle	1 pce
AT041	Cavex mixing bowl blue	1 pce
AT042	Cavex mixing spatula, plastic	1 pce
AT060	Cavex mixing spatula, metaal	1 pce
AT040/043	Cavex plastic container (box + lid)	1 pce
AT076	Cavex powder scoop + water measure	50 sets
AT240	Cavex Alginate container blue	1 pce
AT260	Cavex Alginate container pink	1 pce



CAVEX IMPRESAFE



Hygiene and sterilization are essential in dental clinics and laboratories. That is why Cavex developed a fast, safe, and highly effective dental impression disinfectant: Cavex ImpreSafe.

In just three minutes, this revolutionary disinfectant liquid kills a minimum of 99.99% of all bacteria, fungi, and viruses, including the coronavirus. The complete system for dental clinics contains:

- A disinfection container
- One liter of disinfection concentrate
- A digital timer
- 50 zip lock bags

Features

- Quick and efficient
- Concentrated
- Aldehyde free
- Complete system
- Reliable and secure
- Pleasant smell



Art. Nr.	Description	Content
HA053	Cavex ImpreSafe Starterkit Cavex ImpreSafe 1000 ml, Container, Timer, 50 plastic bags Refill Cavex ImpreSafe	1 pce 1 x 1.000 ml
HA050	Cavex ImpreSafe Refill pack	6 x 1.000 ml
HT020	Cavex ImpreSafe Container	1 pce
HT025	Cavex ImpreSafe plastic bags	60 x 50 pcs

CAVEX GREENCLEAN

Quickly and easily remove alginate and gypsum residues with Cavex GreenClean. Who would be better to produce an alginate dissolver than the makers of "The Best Alginates in the World"?

Cavex GreenClean is a smart solution for cleaning dental impression trays, instruments, and other accessories. Not by scratching or sanding, but simply by dissolving alginate and gypsum residues. Simply submerge for

15-30 minutes, rinse under running water and the sparkling clean impression tray is ready for use.

Features

- PH neutral
- Skin-friendly
- Biodegradable
- 1 kg > 500 clean impression trays
- Refreshing smell

Art. Nr.	Description	Content
HA001	Cavex GreenClean jar	1 kg



CAVEX ALGINATE ADHESIVE

Cavex Alginate Adhesive optimizes both the mechanical and chemical connection between the impression tray and the alginate. The product prevents detachment of the impression from the tray when it is taken out of the patients mouth.

Thanks to a user-friendly brush in the cap and a striking blue color, Cavex Alginate Adhesive is very easy to apply.

Features

- Specifically developed for use with alginates
- Compatible with both metal and plastic impression trays
- Very easy to apply
- No special ventilation measures are necessary
- A minimum of 50 applications per bottle
- Simple removal of any residue with isopropyl alcohol or ethanol (95%)

Art. Nr.	Description	Content
AB010	Cavex alginate adhesive	2 x 14 ml

CAVEX SCANSPRAY

Cavex ScanSpray was developed to optimize the scannability of alginate impressions. The addition of titanium dioxide enables even the smallest details to be clearly scannable. This powder-based spray is suitable for any type of dental impression material.

Features

- Highly elaborated and detailed digital impression models
- Compatible with all types of dental impression materials
- Powder-based instead of oil-based

Art. Nr.	Description	Content
AT001	Cavex ScanSpray	1 x 200 ml



CAVEX VACUFORMER SYSTEM

The Cavex VacuFormer System is a complete vacuum thermoforming system consisting of the Cavex VacuFormer, laminates, a unique gypsum composite and accessories.

The Cavex VacuFormer is a no-nonsense all round thermal vacuum former. This robust machine has an electric high-performance vacuum motor that is suitable for any conventional pressure forming process. The Cavex VacuFormer is easy to operate and protected against overheating.

The Cavex VacuFormer is a no-nonsense all round thermal vacuum former to produce bleaching trays, splints, and mouthguards.

The Starter kit equips the dental professional with all the equipment and materials needed to create ready-to-use pressure forming appliances.



Art. Nr.	Description	Content
VF001	Cavex VacuFormer Starterkit	1 pce
	12 Cavex Mouth Protector Transparent 3,8 mm	
	25 Cavex Splint+ Transparent 1,0 mm	
	25 Cavex Bleach Heavy Transparent 1,5 mm	
	18 sets Instant Stone with Activator and scissors	
VT036	Cavex Mouth Protector Transparent 3,8 mm	12 plates
VT033	Cavex Mouth Protector X-pro Transparent 5 mm	12 plates
VT028	Cavex Splint+ Transparent 1,0 mm	25 plates
VT026	Cavex Splint X-heavy Transparent 1,5 mm	25 plates
VT022	Cavex Bleach Transparent 1,0 mm	25 plates
VT030	Cavex Bleach Heavy Transparent 1,5 mm	25 plates
VT027	Cavex Bruxism 2,0 mm	25 plates



CAVEX INSTANT STONE

Many dental clinics lack a separate area for pouring dental impressions and only a few clinics use a gypsum (vacuum) mixer. Cavex Instant Stone is a gypsum composite to cast perfect molds for bleaching trays, splints and mouthguards in just 10 minutes. Treatments that usually require multiple dental visits can now be performed in just one session. Cavex Instant Stone increases clinic productivity and improves the patient experience worldwide.

Features

- Cavex Instant Stone is a unique gypsum composite to make a perfect gypsum model in just 10 minutes. Equipped with a highly sophisticated mixing system, Cavex Instant Stone is very clean and user-friendly



Art. Nr.	Description	Content
VT010	Cavex Instant Stone with Activator (4 x 18 sets)	72 sets

CAVEX DISPOTRAYS

Hygienic disposable trays for use in standard trays. The trays are designed for disposal after use and can therefore be used once. The disposable tray is the ideal solution for the storage or preparation of instruments in the dental clinic.

Features

- Size: 18 x 28 cm
- Made of (environmentally friendly) polystyrene
- Perfect fit for standard trays



Art. Nr.	Description	Content
HA061	Cavex Dispotrays (18 x 28 cm), box 8 x 50	400 pcs



CAVEX OUTLINE

Cavex Outline is an impression paste developed specifically for taking accurate denture impressions. As the impression paste is eugenol-free, the patient will not experience any stinging or burning sensation upon tissue contact during the impression procedure.

The impression paste is non-elastic and flows under light pressure, resulting in a very accurate impression of the soft tissue.

Features

- Optimal consistency
- Eugenol free
- Non-elastic
- Blue color
- Non-sticky impression paste
- Natural pleasant smell
- Three-year shelf life

An appropriate denture perfectly fits the denture bearing soft tissue. Cavex Outline was developed to enhance elaborated denture impressions.

Art. Nr.	Description	Content
AP030	Cavex Outline	
	white paste	1 x 140 gr
	blue paste	1 x 65 gr
	mixing pad	1 pce



CAVEX SILICONA & SILICONA BITE

Cavex SiliconA is a carefully selected system made of superior silicone material. Even under wet surface conditions, perfect detail reproduction is guaranteed. As non-sticky dental impression system and hence particularly easy to work, Cavex SiliconA is unique in its kind.

Cavex SiliconA bite ensures perfect bite registration time after time.

Features

- User friendly
- 100% accuracy
- Non-sticky
- Excellent long-term dimensional stability
- Highly detailed reproduction
- Optimal time to cure
- Thixotropic, flows under light pressure
- Three-year shelf life

Combine the product with Cavex SiliconA bite. This user-friendly bite registration material offers 100% accuracy, combined with short processing and setting times.

Art. Nr.	Description	Content
AS40	Cavex SiliconA Soft Putty (yellow)	500 ml
AS52	Cavex SiliconA Heavy Body (green)	2 x 50 ml
AS53	Cavex SiliconA Injection Type (blue)	2 x 50 ml

Cavex SiliconA Accessoires

AS70	Cavex SiliconA Dispensing gun	1 pce
AS94	Cavex SiliconA Injection - tips, yellow	48 pcs
AS95	Cavex SiliconA Heavy Body - tips, green	48 pcs
AS96	Cavex SiliconA Injection - intra-oral tips, yellow	96 pcs
AS97	Cavex SiliconA Heavy Body - intra-oral tips, white	96 pcs
AS54	Cavex SiliconA Bite	2 x 50 ml



CAVEX BITE&WHITE IN-OFFICE STARTER PACK



Maximum whitening result with minimal sensitivity? That is the Cavex Bite&White In-Office System. The whitening gel is based on revolutionary Hydrogen Peroxide Superior (HPS) technology, which achieves unprecedented whitening results with 25% hydrogen peroxide. Using the extremely quick and effective HPS gel, the maximum whitening result is achieved in just two or three sessions!

Cavex' exclusive In-Office whitening system is the only system that applies whitening gel with a brush, facilitating maximum control during application. The Whitening Barrier/ Spacer protects the surrounding soft tissue and all areas that are not to be whitened.

Features

- The most effective whitening gel on the market
- Maximum whitening result after only 2 to 3 treatments
- The thermogelling HPS-gel is fixated during treatment thus does not run
- Ready-to-use product, no mixing required
- Unique application technique with a brush pen
- HPS-gel can be stored outside of the refrigerator



Art. Nr.	Description	Content
BW050	Cavex Bite&White in-Office HP25 starter pack	
	Whitening pen	2 x 2 ml
	Barrier / Spacer	1 x 1 ml
	Application tip	5 pcs

CAVEX BITE&WHITE IN-OFFICE PEN BULK PACK



A bulk package of 12 whitening pens for in-office use. The whitening gel is based on 25% Hydrogen Peroxide Superior (HPS) technology for unprecedented whitening results.

The Cavex In-Office Pen contains a brush top which makes it easy applicable on the teeth and gives maximum control during application.

Features

- 12 Cavex In-Office Whitening Pens containing 2 ml gel per pen.
- Applicable with a brush instead of a syringe.
- Can be dosed by turning the back of the pen for easy handling.
- Just two or three in-office treatments is sufficient for an optimal whitening result.

Art. Nr.	Description	Content
BW051	Cavex Bite&White in-Office Pen Bulk Pack	
	Whitening pen	12 x 2 ml



CAVEX BITE&WHITE IN-OFFICE BARRIER BULK PACK

Bulk package of 3 Barriers/Spacers to prevent soft tissue contact and to shield areas that should not come in contact with the whitening gel.

The barrier is easy applicable with the included application tips and is cured within 20 seconds with a curing lamp. Easily removable afterwards.

Features

- 3 Cavex In-Office Barriers of 1 ml and 15 application tips.
- Easily removable after treatment.
- Curing time of 20 seconds with a curing lamp.
- Distinctive blue color.

Art. Nr.	Description	Content
BW052	Cavex Bite&White in-Office Barrier Bulk Pack	
	Barrier / Spacer	3 x 1 gr
	Application tip	15 pcs



CAVEX BITE&WHITE ABC MASTERKIT



Cavex offers the only three-step whitening solution for safe, effective, and above all, pain-free teeth whitening: The Cavex Bite&White ABC Masterkit.

As the syringes are packed in a separate tray there is no need to refrigerate the entire kit. Furthermore, the 3 x 3 ml whitening gel is sufficient for a 15-day treatment. Does it take less than 15 days for the patient to be satisfied? The resealable syringes can be used for a touch-up later.

Activate with StainLess

Clean the teeth thoroughly with Cavex Bite&White StainLess. Other than cleaning the teeth, this paste raises oral acidity levels and optimizes the whitening gel's performance.

Whitening with Bite&White whitening gel

Keep the whitening tray with whitening gel in the mouth for 30 to 60 minutes to achieve optimal whitening results quickly and safely.

Art. Nr.	Description	Content
BW090	Cavex Bite&White ABC Masterkit Cavex Bite&White StainLess, tube Cavex Bite&White refill, syringe Cavex Bite&White ExSense, tube tray cassette, shade guide and instruction manual	1 x 35 ml / 50 gr 3 x 35 ml 1 x 42 ml / 50 gr
BW010	Cavex Bite&White syringe	1 x 3 ml
BW015	Cavex Bite&White bulk Refill Pack	10 x 3 ml

Cavex whitening gel has been proven to be effective.

Conditioning with ExSense

After whitening, Cavex Bite&White ExSense provides essential aftercare. The whitening tray with this conditioning gel is kept in the mouth for 5-10 minutes. ExSense re-mineralizes the teeth, fills the tubules and micro-cracks and so prevents sensitivity.

Features

- A complete whitening treatment under dentist-supervision
- Maximal results, minimal sensitivity
- Professional packaging with clear instructions
- Additional video instructions available on YouTube
- 3 x 3 ml whitening gel for a 15-day treatment
- Teeth remain whiter and brighter for at least 2.5 to 3 years



CAVEX BITE&WHITE 3 SYRINGES



The compact, budget-friendly Cavex Bite&White 3 Syringes package contains three syringes with Cavex Bite&White whitening gel. The amount is sufficient for a 15-day whitening treatment. Cavex Bite&White 3 Syringes is the perfect solution to prolong whitening treatment or simply for a touch-up as well.

The thixotropic gel ensures that the teeth are covered with whitening gel evenly. Furthermore, the gel is prevented from coming into contact with the gums. Cavex Bite&White 3 Syringes is available with 16% carbamide peroxide (CP16= 6% hydrogen

peroxide) for the European market. For other regions, the product is available with 22% carbamide peroxide (CP22= 9% hydrogen peroxide).

Features

- 3 x 3 ml Cavex Bite&White whitening gel
- Thixotropic gel flows under pressure
- Whitening sessions of just 30 to 60 minutes per day
- Three-year shelf life (refrigerated)
- Convenient compact packaging
- Instruction manual included

Art. Nr.	Description	Content
BW007	Cavex Bite&White 3 Syringes CP16	3 x 3 ml
BW008	Cavex Bite&White 3 Syringes CP22	3 x 3 ml



CAVEX BITE&WHITE IMPACT KIT



The Cavex Bite&White Impact Kit is a complete do-it-yourself dental whitening set for use at home. No need for syringes, as the whitening gel is applied using a tube. With the gel being stable at room temperature, the Cavex Bite&White Impact Kit is a one-of-a-kind home whitening set.

The Cavex Bite&White Impact Kit consists of a thermoforming whitening tray and a tube with enough whitening gel for a complete whitening treatment and additional touch-up. Clear instructions are included.

Features

- Stable at room temperature
- No refrigeration needed
- Thermally deformable whitening tray
- Sufficient whitening gel for potential touch up
- Clear instructions included

Art. Nr.	Description	Content
BW040	Cavex Bite&White Impact Kit 1 tube of Cavex Bite&White PVP 1 Boil&Bite tray, 1 tray and instruction manual	1 x 40 gr



CAVEX BITE&WHITE WHITENING PEN



The Cavex Bite&White Whitening Pen is an incredibly user-friendly must-have whitening tool that can be used anywhere and at any time. Whether the patient wants to whiten the top and bottom of the teeth, or just touch up hard-to-reach spots, the Cavex Bite&White Whitening Pen does the trick.

How does it work? Remove the cap and twist the bottom of the pen until a small drop of

whitening gel appears. Apply directly to the teeth.

Features

- No refrigerated storage required
- Effective and safe
- 6% hydrogen peroxide
- Leak-proof whitening pen

Art. Nr.	Description	Content
BW080	Cavex Bite&White Whitening Pen	1 x 2 ml



CAVEX BITE&WHITE PATIENT KIT



The Cavex Bite&White Patient Kit was the first whitening product that Cavex developed. This professional home whitening set is user-friendly and ensures impressive whitening results in the blink of an eye.

As the thixotropic gel must cover the patient's teeth perfectly, the use of the Cavex Bite&White Patient Kit requires a tailor-made whitening tray.

Features

- 2 x 3 ml syringe Cavex Bite&White whitening gel
- Whitening tray cassette
- Shade guide
- User manual

Art. Nr.	Description	Content
BW000	Cavex Bite&White Patient Kit	
	Cavex Bite&White refill, syringe	2 x 3 ml
	tray cassette, shade guide and instruction manual	



CAVEX BITE&WHITE READY 2 USE



Cavex Bite&White Ready 2 Use is a professional, unique, quick, safe, and above all, easy to use whitening kit for at home. This second generation Ready 2 Use teeth whitening kit has been upgraded with Hydrogen Peroxide Superior (HPS) technology. Extremely quick and effective HPS gel boosts performance so a fantastic whitening result is reached after bleaching for just 8 days, 15 minutes per day. On top of that, this home whitening kit is immediately ready for usage and pH neutral.

Due to temperature stability, the new Cavex Bite&White Ready 2 Use does not need to be stored refrigerated. Also, the product is very user-friendly. Unlike other systems that require whitening trays to be worn overnight, treatment with Cavex Ready 2 Use takes just 15 minutes per day.

Features

- Consistent whitening results with 6% Hydrogen Peroxide Superior (HPS)
- 8 flexible whitening trays for 8 days
- One size fits all tray that simultaneously whitens the upper and lower teeth
- HPS technology: complete whitening treatment in 8 days
- Only 15 minutes per day
- Stable at room temperature
- PH neutral

Art. Nr.	Description	Content
BW031	Cavex Bite&White Ready 2 Use, 8 pre-filled trays	8 x 3 ml
BW036	Cavex Bite&White Ready 2 Use Bulk	50 pre-filled trays
		50 x 3 ml



CAVEX BITE&WHITE BARRIER/SPACER

A blue colored flowable resin especially designed for creating reservoirs in the bleach tray as well as a barrier for In-Office treatments to prevent gum contact.

Features

- Usable as barrier and spacer.
- Distinctive blue color for accurate application.
- Curing time of just 20 seconds.

With the included dispensing tips the barrier/spacer is easy to apply. After the treatment the material is easy to remove.

Art. Nr.	Description	Content
VT029	Cavex Bite&White Whitening Barrier/Spacer	2 x 1 gr
	Application tip	10 pcs





CAVEX BITE&WHITE STAINLESS

Cavex Bite&White StainLess is one of the best extreme cleaning, non-abrasive pastes, and has been specifically developed to achieve optimal whitening results. The paste is used as a pre-treatment activator in the ABC Masterkit, before the whitening gel, but is also available separately.

StainLess paste cleans the teeth thoroughly, allowing the whitening gel to come into direct contact with the tooth enamel.

Also, by temporarily raising oral pH values, the StainLess paste optimizes hydrogen peroxide's performance.

Features

- Thorough cleansing and removal of plaque
- Raises oral pH and optimized whitening conditions
- Fresh, pleasant taste

Art. Nr.	Description	Content
BW070	Cavex Bite&White StainLess, tube	10 x 50 gr / 35 ml

CAVEX ORAL PRE RINSE

Cavex developed Oral Pre Rinse by combining scientific findings with years of experience in the field of hydrogen peroxide and mouthwashes. Cavex Oral Pre Rinse is a pre-procedural mouth rinse with 1.5% stabilized hydrogen peroxide. 30 to 60 seconds of flushing before treatment significantly reduces the risk of any viruses spreading, including the coronavirus.

Cavex Oral Pre Rinse contains a bio-adhesive that creates an ideal viscosity and foam condition. Erythritol has an antibacterial effect and strengthens cleansing.

NEW

The addition of glycerin prevents dry mouth. Cavex Oral Pre Rinse has a nice, slightly sweet mint flavor, offering patients a pleasant rinsing experience.

Features

- Pre-procedural oral rinse
- Contains 1.5% stabilized hydrogen peroxide
- Kills viruses, including COVID-19
- Ready to use
- Effective in 30-60 seconds
- Patient-friendly, pleasant taste
- 500 ml for about 50 rinses
- Bottle with integrated dosing cup

Art. Nr.	Description	Content
MW001	Cavex Oral Pre Rinse	6 x 500 ml



CAVEX EXSENSE

Cavex ExSense is the perfect solution for sensitive teeth. Sensitive teeth often occur as a result of dental erosion or exposed dentin. Cavex ExSense is based on 2% hydroxyapatite that seals the tubules and microcracks in the enamel using a hydro-dispersing sol. In doing so, it completely closes off sensitive areas. The addition of 0.1% potassium nitrate also restores the micro-hardness of the elements. In addition to complete closure, crystallization takes place which supports and accelerates the remineralization process.

Cavex ExSense contains natural ingredients, hence patients can use the product without restriction at home. This desensitizer

NEW



is a Medical Device Class IIa that meets all European requirements for safety, health, and environment in Europe. Say goodbye to sensitive teeth!

Features

- Instant and long-lasting relief from sensitive teeth
- Helps to prevent sensitive teeth
- Accelerates and reinforces the process of remineralization
- Polishes teeth surface, luster effect
- Slows down the process of plaque formation
- Suited for (unlimited) home use

Art. Nr.	Description	Content
BW160	Cavex Bite&White ExSense, tube	1 x 50 gr / 42 ml



CAVEX HEALTHYPOPS



Lollipops with a dental function! Fight dry mouth and stimulate fresh breath with these unique lollipops. Healthypops are the ideal approach to halitosis and xerostomia in dental clinics. Additionally, the lollipops promote patients' oral care. Furthermore, the lollipops are very suitable as a tasty, yet responsible snack!

Cavex Healthypops are sugar-, gluten-, dairy- soy- and nuts-free and do not contain added aromas or flavors. Consequently, they can safely be used by diabetic patients.

Art. Nr.	Description	Content
BW900	Cavex Healthypops Cherry	60 x Lollipop (454 gr)
BW901	Cavex Healthypops Mint	60 x Lollipop (454 gr)
BW902	Cavex Healthypops Orange	60 x Lollipop (454 gr)

NEW

The lollipops are available in three flavors:

- Cherry** – mouth moistening in case of xerostomia
- Orange** – fresh breath in case of halitosis
- Mint** – fresh breath in case of halitosis

Features

- Sugar-free lollipops with a dental function
- Available in cherry, orange, and mint flavor
- Relieves dry mouth
- Promotes fresh breath
- No sugar, gluten, dairy, soy, nuts, or added fragrances and flavors
- Suitable for diabetics

CAVEX RUSH BRUSH



The Cavex Rush Brush allows patients to brush their teeth easily and quickly before dental treatment or on the go. The toothbrush is immediately ready for use and no toothpaste is needed. The brush is invisibly coated with fresh, fluoride-free toothpaste.

The small brush head reaches every area in the mouth and the soft, rounded bristles are safe for the gums. Also, the toothbrush is packaged in a individual wrapper. This makes the brush easy to take anywhere while maintaining optimal hygiene.

Art. Nr.	Description	Content
HA090	Cavex Rush Brush 100 pieces per box (red, green, yellow and blue 25 of each)	10 x 100 pcs

Features

- Impregnated with toothpaste
- Hygienic packaging
- On-the-go solution
- Fresh mint flavor
- Safe, soft brush
- Small brush head for optimal precision

CAVEX FRESH4SURE

Prevent bad breath with Cavex Fresh4Sure. Cavex Fresh4Sure is a unique 3-step oral care system that restores the balance between 'good' and 'bad' bacteria in the mouth in a natural way. Stabilizing oral bacteria is essential for a healthy mouth and fresh breath. The 3-step oral care system is therefore effective against halitosis and a dry mouth.

The Fresh4Sure kit consists of a tongue cleaner, mouthwash, and tongue gel.

Features

- Contains a tongue cleaner, mouthwash and tongue gel
- Long-lasting and durable effect
- Mouthwash also effective against dry mouth
- Does not contain SLS
- Does not contain alcohol

Art. Nr.	Description	Content
MW030	Cavex Fresh4Sure Kit	1 pce
	1 x tongue cleaner	
	1 x bottle Fresh4Sure Mouthwash 500 ml	
	1 x tube Fresh4Sure Tongue Gel 45 gr	
	1 x instruction manual	
MW010	Fresh4Sure Mouthwash	6 x 500 ml
MW020	Fresh4Sure Tongue Gel Tube	1 x 45 gr



CAVEX BOIL&BITE

Cavex Boil&Bite is a double-sided thermoplastic tray for the treatment of tooth sensitivity. This dental tray is an upper and lower tray in one and can be used for different purposes. The revolutionary tray is well tolerated hence suitable for everyone!

The Cavex Boil&Bite can be easily used by both dentists and patients. Consequently, Cavex Boil&Bite is can be used for treatment in dental clinics and use at home.

Features

- Perfect match with Cavex ExSense
- For use in dental clinic as well as at home
- One box contains eight mouthguards

Art. Nr.	Description	Content
BW085	Cavex Boil&Bite	8 pcs



QUADRANT UNIVERSAL LC



Quadrant Universal Living Colors (QULC) is a universal hybrid composite, enriched with submicron particles. Thanks to this unique and balanced composition, QULC is easy to apply and very easy to process. The composite filling is easy to model, does not stick to instruments, and reaches a mesmerizing high gloss when polished. The product contains a unique fluorescence pigment, hence mimicking organic tooth issues. As a result, QULC grants teeth a natural and lively appearance.

Features

- Prism-shaped filler particles for optimal light scattering
- Natural fluorescence
- Suitable for mono layering and multi-layering technique
- Perfect chameleon effect
- Durable and wear-resistant
- Perfect aesthetics
- Radiopaque
- BPA free (Bisphenol A)

Art. Nr.	Description	Content
FU02	Quadrant Universal LC starterkit syringes Content: 5 x 4 g (A1, A2, A3, A3.5, A4) + Uni-SE-Bond + shade guide & accessoires	1 pce
FU20	Quadrant Universal LC syringe A1	4 gr
FU21	Quadrant Universal LC syringe A2	4 gr
FU22	Quadrant Universal LC syringe A3	4 gr
FU23	Quadrant Universal LC syringe A3.5	4 gr
FU28	Quadrant Universal LC syringe A4	4 gr
FU24	Quadrant Universal LC syringe B2	4 gr
FU29	Quadrant Universal LC syringe B3	4 gr
FU25	Quadrant Universal LC syringe C2	4 gr
FU30	Quadrant Universal LC syringe C3	4 gr
FU26	Quadrant Universal LC syringe OA2	4 gr
FU27	Quadrant Universal LC syringe Incisal	4 gr
FU52	Quadrant Universal LC starterkit capsules Content: 5 x 5 x 0,25 g (A1, A2, A3, A3.5, A4) + Uni-SE-Bond + shade guide & accessoires	1 pce
FU60	Quadrant Universal LC capsule A1	20 x 0,25 gr
FU61	Quadrant Universal LC capsule A2	20 x 0,25 gr
FU62	Quadrant Universal LC capsule A3	20 x 0,25 gr
FU63	Quadrant Universal LC capsule A3.5	20 x 0,25 gr
FU68	Quadrant Universal LC capsule A4	20 x 0,25 gr
FU64	Quadrant Universal LC capsule B2	20 x 0,25 gr
FU69	Quadrant Universal LC capsule B3	20 x 0,25 gr
FU65	Quadrant Universal LC capsule C2	20 x 0,25 gr
FU70	Quadrant Universal LC capsule C3	20 x 0,25 gr
FU66	Quadrant Universal LC capsule OA2	20 x 0,25 gr
FU67	Quadrant Universal LC capsule Incisal	20 x 0,25 gr

Accessoires for Quadrant Universal LC

FG73	Shade guide Quadrant Universal LC 11 colors	1 pce
FG57	Quadrant mixingdish for single use	25 pcs



QUADRANT ANTERIOR SHINE



Give your patients a healthy, natural look with our light-curing, micro-hybrid composite Quadrant Anterior Shine. A glossy, easy-to-polish surface combined with natural fluorescence makes it very hard to distinguish a restoration from a natural tooth. This extremely durable micro-hybrid composite guarantees a wear and fracture-resistant surface. Additionally, low liquid absorption results in a high color-stability.

Features

- Barium glass technology
- High gloss polishable
- Natural fluorescence
- Perfect chameleon effect
- Black light effect
- High elasticity
- Superior processing properties
- Thixotropic
- Color stable
- High fracture toughness
- Does not stick to instruments
- BPA free (Bisphenol A)

Quadrant Anterior Shine is a perfect composite filling for:

- Front restorations
- Treatment of enamel hypoplasia
- Discolorations
- Tooth corrections
- Diastema closure



Art. Nr.	Description	Content
FA60	Quadrant Anterior Shine capsule A1	10 x 0,25 gr
FA61	Quadrant Anterior Shine capsule A2	20 x 0,25 gr
FA62T	Quadrant Anterior Shine capsule A3	20 x 0,25 gr
FA63	Quadrant Anterior Shine capsule A3.5	20 x 0,25 gr
FA64	Quadrant Anterior Shine capsule A4	10 x 0,25 gr
FA65	Quadrant Anterior Shine capsule B2	10 x 0,25 gr
FA66	Quadrant Anterior Shine capsule B3	20 x 0,25 gr
FA67	Quadrant Anterior Shine capsule C2	10 x 0,25 gr
FA68	Quadrant Anterior Shine capsule C4	10 x 0,25 gr
FA69	Quadrant Anterior Shine capsule D3	10 x 0,25 gr
FA70	Quadrant Anterior Shine capsule Incisal	20 x 0,25 gr
FA71	Quadrant Anterior Shine capsule OA2	10 x 0,25 gr
FA72	Quadrant Anterior Shine capsule OA3.5	10 x 0,25 gr
FA73	Quadrant Anterior Shine capsule Dark Shade	10 x 0,25 gr
FA74	Quadrant Anterior Shine capsule Bleach	10 x 0,25 gr

Accessoires for Quadrant Anterior Shine

FG57	Quadrant mixingdish for single use	25 pcs
------	------------------------------------	--------

QUADRANT QUAST

The Quadrant Quast, or 'Quadrant Anterior Shine Tool', is a brush specifically selected for the precise application and distribution of Quadrant Anterior Shine. This revolutionary brush realizes superior composite facings. The Cavex Quadrant Quast was developed in collaboration with Dutch dentist and prominent opinion leader Hein de Kloet.

Features

- Specifically developed for the application and distribution of Quadrant Anterior Shine
- Developed in collaboration with Hein de Kloet
- For superior composite facings

Art. Nr.	Description	Content
FF30	Quadrant Anterior Shine Tool (QUAST)	12 pcs

QUADRANT FLOW

Quadrant Flow is a liquid, light-curing composite. This versatile product is used for various treatments. The syringe is easy to use. Thanks to the small syringes and included tips, the composite can be applied directly.

Quadrant Flow is compatible with all modern composite systems. The most unique feature of this product is its thixotropic character. It flows easily, but does not run.

Features

- Easy to apply
- Available in four shades
- Barium-glass composite
- Minimally invasive composite fillings
- Does not stick to instruments
- Compatible with other composite materials
- Thixotropic
- Three-year shelf life
- BPA-free (Bisphenol A)
- Light-curing
- Radiopaque

Art. Nr.	Description	Content
FD09	Quadrant Flow A2 syringe + 5 tips	1,8 gr / 1 ml
FD10	Quadrant Flow A3 syringe + 5 tips	1,8 gr / 1 ml
FD11	Quadrant Flow A3.5 syringe + 5 tips	1,8 gr / 1 ml
FD12	Quadrant Flow OA2 syringe + 5 tips	1,8 gr / 1 ml

QUADRANT TOTAL ETCH

Quadrant Total Etch is a 35% phosphoric acid gel for the total etching technique.

Features

- Ultimate concentration for total etching technique
- Optimized glaze preparation
- No excessive dentin demineralization

- Distinctive blue color
- Perfect visual overview of the etching process
- Quick and easy residue detection
- Excellent base for adhesive bonding agents
- Complete hybridization of exposed collagen fibrils

Art. Nr.	Description	Content
FE10	Quadrant Total Etch + 10 tips	2,5 ml
FG55	Tips Quadrant Total Etch	10 pcs

QUADRANT UNI-SE-BOND

Quadrant Uni-SE-Bond is a seventh-generation self-etch bonding system. The system can be combined with light-curing restorative materials. A unique integrated desensitizer prevents post-operative sensitivity.

Features

- Strong and durable bonding due to 4-META and mild etch
- Integrated desensitizer prevents sensitivity
- Simple user-friendly procedure
- Deep surface penetration
- Can be stored at room temperature

Art. Nr.	Description	Content
FB20	Quadrant Uni-SE-Bond	4 ml



QUADRANT UNI-1-BOND

Quadrant Uni-1-Bond is a one-bottle primer and sealer for direct restorations. The light-curing composite adhesive has excellent bonding properties due to the 4-META monomer.

Features

- 4-META monomer quality optimization
- Excellent adhesion
- Optimal dentin penetration abilities
- Prevents bacterial invasion
- Ethanol water base
- Constant composition
- Low technique-sensitivity
- Easy to use
- Humidification capacity

Art.-Nr.	Description	Content
FB30	Quadrant Uni-1-Bond	4 ml



QUADRANT FINISENSE

Quadrant FiniSense is a quick and effective desensitizer. The liquid provides pain relief caused by exposed dentin. Also, it relieves post-operative discomfort after crown preparation and inlays within just 30 seconds. Quadrant FiniSense is easy to use and offers a long-lasting result. One bottle contains sufficient liquid for about 150 applications.

Features

- Effective within 30 seconds
- Durable 12-month effect
- Surface covering with film
- Easy to use
- About 150 applications per bottle

Art. Nr.	Description	Content
FS10	Quadrant FiniSense	5 ml



QUADRANT Q-SYSTEMISER CAPSULES



Functional storage system for composite capsules and accessories.

Features

- 16 compartments, each with a recess for a used capsule.
- Made from orange acrylic; exclusion of all light influences.

- Supplied incl. stickers with color code.
- Built-in brush dispenser with unique rotation system. With every turn, one brush is released; practical and hygienic.
- Can be cleaned in the thermo disinfector. Dimensions: 33 x 17 x 4 (lxwxh).

Art. Nr.	Description	Content
FG40	Quadrant Q-Systemiser Capsules, incl. 100 Quadrant Q-Brushes XL	1 pce

QUADRANT Q-SYSTEMISER SYRINGES

Functional storage system for composite syringes and accessoires.

Features

- Made from transparent acrylic.
- 14 slots for composite syringes.

- Built-in brush dispenser with unique rotation system. With every turn, one brush will be released; practical and hygienic.
- Can be cleaned in the thermo disinfector.

Art. Nr.	Description	Content
FG41	Quadrant Q-Systemiser Syringes, incl. 100 Quadrant Q-Brushes XL	1 pce



QUADRANT Q-BRUSHES XL


Quadrant Q-Brushes XL are blue disposable applicators for the application of primer, bonding, and desensitizer, among other things. These bendable brushes can be used for multiple purposes. Once bent, the brush will maintain a fixed position. The extra-large application brushes significantly aid when applying larger amounts of liquid.

Features

- For the application of etch, bonding and desensitizer
- Brushes can be bent in different positions
- Brush maintains desired position
- XL size



Art. Nr.	Description	Content
FF21	Quadrant Q-brushes XL, borstels for single use	100 pcs



QUADRANT ACCESSORIES

Art. Nr.	Description	Content
FG73	Shade guide Quadrant Universal LC 11 colors	1 pce
FG55	Tips Total Etch	10 pcs
FG57	Mixingtrays	25 pcs



CAVEX SET UP WAX



Cavex Set Up Wax is a top-shelf modeling wax. The wax is available in Soft, Regular, Hard, and Extra Hard. Processing is identical but melting points and carving wax hardness vary. In addition to these variants, the wax is now available in Regular Flex. This dental wax is characterized by extra heat transition flexibility, offering technicians a Cavex Set-Up Wax that perfectly matches the processing of the wax.

Set Up wax plates are large and usually sufficient for two prostheses. In addition to setting up the bite wall, Cavex dental wax adheres well to the spoon and is therefore

also used as a boxing wax. Due to its ideal consistency, this wax is one of the top material choices for bite registration in dental clinics. A prosthesis must both be beautiful and fit perfectly. To that end, Cavex Set Up Wax holds flow when molded by the dental technician and can be processed as desired.

Features

- For complete and partial dentures
- Low casting shrinkage
- Dimensionally stable
- Elastic, holds the flow
- Modeling and boxing wax

Art. Nr.	Description	Content
WA004	Cavex Set Up Soft 500 gr, 23 sheets	500 gr
WA005	Cavex Set Up Regular 500 gr, 23 sheets	500 gr
WA002	Cavex Set Up Regular Flex 500 gr, 23 sheets	500 gr
WA006	Cavex Set Up Hard 500 gr, 23 sheets	500 gr
WA008	Cavex Set Up Extra Hard 500 gr, 23 sheets	500 gr
WA014	Cavex Set Up Soft 2500 gr, 115 sheets	2500 gr
WA015	Cavex Set Up Regular 2500 gr, 115 sheets	2500 gr
WA012	Cavex Set Up Regular Flex 2500 gr, 500 sheets	2500 gr
WA016	Cavex Set Up Hard 2500 gr, 115 sheets	2500 gr
WA018	Cavex Set Up Extra Hard 2500 gr, 115 sheets	2500 gr

CAVEX YELLOW WAX

Cavex Yellow Wax is a softer dental wax and even easier to mold. Perfect for eliminating undercuts and adjusting the edges of impression trays.

Features

- For adjustment to impression trays
- Blocking out undercuts
- Flexible, easy to mold
- Cost-efficient
- Five-year shelf life

Art.-Nr.	Description	Content
WB006	Cavex Yellow Wax	17 sheets



CAVEX BASEPLATES

Cavex Baseplates are the most reliable basis for a perfectly fitting prosthesis. Shellac is one of the natural components of the product. Years of technological perfection within Cavex have led to the possibility to form an extremely stable three-dimensional model from a two-dimensional plate.

After a short heat-through in the flame, the pink plate with a standard thickness of 1.4 mm can easily be molded onto the gypsum model. Furthermore, the user-friendly Cavex Baseplates are easy to cut.

The Mica-free version further improves the performance of the Cavex Dental Baseplates. Aluminum hydroxide ensures the release of cooling crystal water once

temperatures exceed 200°C, preventing burning of the plate. Additionally, the result is smoother, non-porous, and shiny, exactly like the prosthetic plastic will turn out.

Features

- Laboratory material for prosthesis manufacturing
- Simple technique
- Strong with high dimensional stability
- Standardizes thickness
- No deformation at mouth temperature
- Packed per 100 pieces
- Mica-free baseplates available

Art. Nr.	Description	Content
BP004	Cavex Baseplates uppers	100 pcs
BP054	Cavex Baseplates lowers	100 pcs
BP024	Cavex Baseplates uppers, double thickness 2.70 - 2.80 mm	50 pcs
BP005	Cavex Baseplates mica-free uppers	100 pcs
BP055	Cavex Baseplates mica-free lowers	100 pcs



CAVEX LC DENTAL TRAY

Cavex Dental Trays have been used for many years in dental laboratories all over the world to make well-fitting impression trays. We differentiate between Cavex LC Dental Tray and Cavex Individual Tray.

The Cavex Light Curing Dental Tray is a light-curing synthetic resin spoon plate for creating individual or functional impression trays. The material is also suitable for making base plates for wax-ups and bite templates.

The Cavex LC Dental Tray offers various benefits, among which its extra-wide shape. Consequently, the hassle with sticking up

pieces belongs to the past. Furthermore, the processing is simplified as the material can be cured with halogen and UV. Once properly cured, the material is dimensionally stable.

Features

- Light-curing tray material
- Curing time with halogen or UV: 3 to 5 minutes
- Extended width
- Three-year shelf life
- Fifty sheets per unit

Art. Nr.	Description	Content
BP010	Cavex LC Dental Tray	50 pcs

CAVEX INDIVIDUAL TRAYS

The Cavex Individual Tray is a gray-colored tray with a standard thickness of 1.9 mm. It guarantees an optimal fit even when the spoon is stored for a longer period. Fast processing and a favorable price/quality ratio are reasons why the Cavex Individual Tray has been popular in dental laboratories all over the world for many years.

Features

- Simple technique
- Rapid processing
- Dimensionally stable
- Cost-efficient
- Standard thickness of 1.9 mm
- Metal handles included

Art. Nr.	Description	Content
BL006	Cavex Individual Trays uppers	100 pcs
BL056	Cavex individual Trays lowers	100 pcs

CAVEX GLASS IONOMER CEMENT

Cavex Glass Ionomer Cement is a dental cement for the crown- bridge- and in-lay- cementation. Prepare the cement by gradually adding small amounts of powder to the fluid and mix to enhance absorption. Formulas of both powder and liquid are developed and optimized to attain a cement with aesthetic properties.

Features

- Biocompatible
- Optimal dental tissue adhesion
- User-friendly processing- and curing-time
- Increased compressive and tensile strength
- Radiopaque
- No over-acidification
- 5.5-year shelf life

Art.-Nr.	Description	Content
CE020	Cavex Glass Ionomer Cement Powder	35 gr
	Liquid	15 ml
	Measuring spoon	1 pce



CAVEX TEMPORARY CEMENT



A temporary dental cement must be reliable and provide effective protection. It should also be easy to apply and remove without leaving any residue. Cavex developed the ideal solution to meet these needs: Cavex Temporary Cement. The properties of the cement are perfect for cementing temporary crowns and bridges. Cavex Temporary Cement is user-friendly, odorless, aesthetically pleasing, and also eugenol-free.

Features

- Easy to use
- Eugenol free
- Distinctive color
- Perfect sealing
- Low solubility
- Short curing time
- Easy to remove

Art.-Nr.	Description	Content
CE090	Cavex Temporary Cement Base	35 gr
	Catalyst	16 gr
	Mixing pad	1 pce



CAVEX NON GAMMA-2



3 spill

2 spill

1 spill

Cavex NON-Gamma-2 is a dispersed, high silver gamma-2-free amalgam. The amalgam has been proven effective and has already proven its reliability and wear resistance in dental clinics.

The self-activating and resealable capsule offers optimum safety both during and after mixing. Cavex amalgams avoid the corrosion susceptible gamma-2 phase and hence can

be used to produce durable restorations with a permanent shine. In short: a safe product for a new future.

Features

- High silver content
- Exceptionally durable
- Easy to process
- Relatively short curing time

Art. Nr.	Description	Content
LK021	jar with capsules, 1-spill (400 mg alloy)	50 caps.
LK022	jar with capsules, 2-spill (600 mg alloy)	50 caps.
LK023	jar with capsules, 3-spill (800 mg alloy)	50 caps.
LK031	jar with capsules, 1-spill (400 mg alloy)	300 caps.
LK032	jar with capsules, 2-spill (600 mg alloy)	300 caps.
LK033	jar with capsules, 3-spill (800 mg alloy)	300 caps.

CAVEX AVALLOY



3 spill

2 spill

1 spill

Cavex has developed a future-proof amalgam with improved mechanical properties: Cavex Avalloy. An innovative high-copper alloy that contains lathe-cut particles.

Cavex Avalloy has already proven its reliability, usability, and wear resistance in the dental clinic. And now, thanks to the new self-activating resealable capsule, Cavex amalgams offer optimum safety. Both during and after mixing. By avoiding the gamma 2 phase, Cavex amalgams significantly decrease corrosion susceptibility. As a result, Cavex amalgams can be used to produce durable restorations with a permanent shine. In short: a safe product for a new future.

Features

- Safe, self-activating, resealable capsule
- High-quality restoration work
- Color stable
- Easy to process
- Low silver and high copper alloy
- Gamma 2 phase free
- Complies with the latest ISO and ADA requirements

Art. Nr.	Description	Content
LM021	jar with capsules, 1-spill (400 mg alloy)	50 caps.
LM022	jar with capsules, 2-spill (600 mg alloy)	50 caps.
LM023	jar with capsules, 3-spill (800 mg alloy)	50 caps.
LM031	jar with capsules, 1-spill (400 mg alloy)	300 caps.
LM032	jar with capsules, 2-spill (600 mg alloy)	300 caps.
LM033	jar with capsules, 3-spill (800 mg alloy)	300 caps.



**YOUR
IMPRESSION
IS OUR
SPECIALTY**

Fustweg 5 | 2031 CJ Haarlem | The Netherlands | cavex.nl

The Cavex logo is presented within a white square. The word "CAVEX" is written in a bold, black, sans-serif font, with a distinctive dot above the letter 'A'.

CAVEX

ZY200E.14