



CREATED BY SDI
POLA
ADVANCED TOOTH
WHITENING SYSTEM

WHITER
BRIGHTER
YOU

Australia

POLA OFFICE	POLA OFFICE+	POLA RAPID	POLA DAY	POLA DAY CP	POLA NIGHT	POLA LUMINATE	POLA LIGHT	SOOTHE	GINGIVAL BARRIER
<ul style="list-style-type: none"> Minimal chair time No bleaching light required. Can be used with or without a curing light Economical Flexible barrier Mix as you go Only thin layer required Contains potassium nitrate for the reduction of sensitivity Sticky gel formula for easy application Gel is freshly activated for direct application 	<ul style="list-style-type: none"> Fast and simple to use: from 3 x 8 minute applications Can be used with or without a light Auto-mix syringe – direct application Contains potassium nitrate to inhibit sensitivity Sticky gel formula maximises adhesion and whitening results 37.5% hydrogen peroxide 1 and 3 patient kits Optional retractors available in kits 	<ul style="list-style-type: none"> Fast and simple to use: from 3 x 8 minute applications 38% hydrogen peroxide Non-stick gel for faster application and removal Enhanced blue gel for superior application visibility Auto-mix syringe for direct application Contains potassium nitrate and fluoride to reduce sensitivity Compatible with or without a whitening light Available in 1 and 3 patient kits Optional retractors available in kits 	<ul style="list-style-type: none"> High Viscosity - prevents saliva washout Fluoride releasing Reduces dehydration of the enamel Built-in desensitiser for the reduction of the gels drying effect on the teeth High water content to reduce sensitivity Whitening trays required Spearmint flavour 	<ul style="list-style-type: none"> High Viscosity - prevents saliva washout Fluoride releasing Reduces dehydration of the enamel Reduction of gels drying effect on the teeth High water content to reduce sensitivity Whitening trays required Spearmint flavour 	<ul style="list-style-type: none"> High Viscosity - prevents saliva washout Fluoride releasing Reduces dehydration of the enamel Built-in desensitiser for the reduction of the gels drying effect on the teeth High water content to reduce sensitivity Whitening trays required Reduction of post-operative sensitivity Spearmint flavour 	<ul style="list-style-type: none"> Quick and easy to use Brush on whitening system Touch up and maintain Pola whitening Fluoride releasing Over 60 applications in the pen Spearmint flavour 	<ul style="list-style-type: none"> Fast acting, delivering whiter teeth in 5 days Returns a patient's smile to their whitest natural shade possible Extra hydration to minimise sensitivity Remineralization to strengthen your teeth Fully rechargeable LED light for future tooth whitening treatments 	<ul style="list-style-type: none"> Desensitises the teeth Will not interfere with the whitening process Can be easily applied Can be applied before and after treatment Re-mineralises the teeth Reduces sensitivity 	<ul style="list-style-type: none"> Applied to the gums to protect them while whitening Blue in color Easy to apply and remove Flexible

ORDER CODES

<p>3 Patient Kit - 7700011 1 Patient Kit - 7700031</p>	<p>3 Patient Kit 7700412 with retractors 7700411 without retractors</p> <p>1 Patient Kit 7700427 with retractors 7700426 without retractors</p>	<p>3 Patient Kit 7700518 with retractors 7700519 without retractors</p> <p>1 Patient Kit 7700516 with retractors 7700517 without retractors</p>	<p>10 x 1.3g Syringe Kits 7710020 - 6% HP 7700151 - 7.5% HP 7700154 - 9.5 HP</p> <p>4 x 1.3g Mini Kits 7710021- 6% HP 7700150 - 7.5% HP 7700153 9% - 9.5%</p> <p>50 x 3g Dispenser Kits 7700075 - 6% HP 7700081 - 7.5% HP 7700082 - 9.5% HP</p>	<p>10 x 1.3g Syringe Kits 7700328 - 35% CP</p> <p>4 x1.3g Mini Kits 7700326 - 35% CP</p> <p>50 x 3g Dispenser Kits 7700084 - 35% CP</p>	<p>10 x 1.3g Syringe Kits 7700157 - 10% CP 7700160 - 16% CP 7710023 - 18% CP 7700163 - 22% CP</p> <p>4 x 1.3g Mini Kits 7700156 - 10% CP 7700159 - 16% CP 7710022 - 18% CP 7700162 - 22 % CP</p> <p>50 x 3g Dispenser Kits 7700076 - 10% CP 7700077 - 16% CP 7700078 - 18% CP 7700083 - 22% CP</p>	<p>Pen 7700092 1 x 5.5ml</p>	<p>Pola Light 6% Hydrogen Peroxide Kit 7700982</p>	<p>4 x 1.2ml Syringe 8150002</p> <p>50 x 1.2ml Syringes 8150033</p>	<p>2x 1g Gingival Barrier Syringes - 7750006</p> <p>Accessories</p> <p>Tray Cases - 7750008</p> <p>Space Resin 2 syringe refill - 7700012</p> <p>Tray Material 2 refill - 7750005 x 20</p> <p>Cosmetic Pola pouches - 7700999 x 10</p>
--	---	---	---	---	--	----------------------------------	--	---	---

