



CEMENT PRODUCTS INDICATIONS



RIVA LUTING
(Self-Curing GI Cement)



RIVA LUTING PLUS
(Resin Modified GI Cement)



RIVA CEM
(Resin Modified GI Cement)



RIVA CEM AUTOMIX
(Resin Modified GI Cement)



SET
(Self Etching, Self Adhesive Resin Cement)



SET PP
(Self etching, Self adhesive Resin Cement)

	RIVA LUTING	RIVA LUTING PLUS	RIVA CEM	RIVA CEM AUTOMIX	SET	SET PP
Metal/Metal Based						
Crowns	✓	✓	✓	✓	✓	✓
Bridges		✓	✓	✓	✓	✓
Inlays/Onlays	✓	✓	✓	✓	✓	✓
Orthodontic appliances	✓ (need conditioner/etch)			✓	✓ (need etch)	✓ (need etch)
Posts and screws	✓	✓	✓	✓	✓	✓
Resin Based						
Posts (fibre reinforced)		✓	✓	✓	✓	✓
High Strength Ceramics (e.g. Zirconia, Lithium disilicate)						
Crowns		✓	✓	✓	✓	✓
Bridges		✓	✓	✓	✓	✓
Inlays/Onlays		✓	✓	✓	✓	✓
Veneers						✓ (need etch)
Low Strength Ceramics (e.g. Porcelain, Glass ceramic, Feldspar)						
Crowns					✓	✓
Bridges					✓	✓
Inlays/Onlays		✓	✓	✓	✓	✓
Veneers						✓ (need etch)
Other Indications						
Base/liner	✓					