

premium  
plus™

# the sterilisation range

our full range of sterilisation  
rolls and pouches



# our sterilisation range

Autoclave pouches and rolls are essential for effective sterilisation across various fields, including dental, medical, veterinary, hospital, podiatry, tattoo, and body piercing.

In this catalogue, you will find detailed information about our sterilisation pouches, rolls, and cassettes to help you get to know our products better and provide you with options to suit your customer's every need and meet their strictest hygiene standards.

## contents

<b>Page 3</b>	About our sterilisation range
<b>Page 4</b>	Sterilisation pouches
<b>Page 5</b>	Sterilisation rolls
<b>Page 6-8</b>	Sterilisation pouches - specifications
<b>Page 9</b>	Sterilisation rolls - specifications
<b>Page 10-11</b>	Sterilisation cassettes - specifications
<b>Page 12-13</b>	Packaging
<b>Page 14</b>	Order information
<b>Page 15</b>	About Premium Plus


## the importance of sterilisation

In the world of modern healthcare, maintaining a sterile environment is critical to ensuring patient safety. Proper sterilisation of instruments is a fundamental aspect of infection control, safeguarding both the patient and the practitioner from harmful microorganisms that could be transmitted through contaminated instruments. Sterilisation is not just a recommended practice but a professional obligation that protects the integrity of treatments and builds trust with patients.

Effective sterilisation helps to prevent cross-contamination, reduces the risk of infections, and supports compliance with health regulations and standards. By ensuring that all instruments are thoroughly sterilised, healthcare providers demonstrate a commitment to high standards of care and operational excellence.

## expanding your product line

At Premium Plus, we offer a comprehensive Sterilisation Range designed to meet the highest standards of hygiene and safety. Our selection of sterilisation pouches, rolls, and cassettes ensures reliable performance, integrating seamlessly into daily workflow and delivering efficient performance across a variety of instruments and procedures. By choosing our products you are offering trusted, high-performance sterilisation solutions that meet the needs of your customers and their practices.




# sterilisation pouches

Sterilisation pouches provide an easy sealing system that eliminates the need for expensive heat-sealing equipment. This ensures a uniform seal that withstands the sterilisation process and keeps instruments sterile.

- Manufactured from 70 gsm medical grade paper and 60 µm blue tinted film to help detect puncture holes
- Easy to read steam indicators change from blue to dark grey to ensure sterile conditions have been met
- Triple adhesive seal to protect contents
- Internal and external indicators
- Shelf life: 3 years before use, 6 months after processing
- Each pouch is printed with LOT no. for traceability
- Large range of sizes for all applications
- Complies with ISO 11607-1 and EN 868-5


<b>6000</b>	57 mm x 90 mm	(200 pcs/box)	<b>6007</b>	90 mm x 230 mm	(200 pcs/box)
<b>6001</b>	57 mm x 130 mm	(200 pcs/box)	<b>6009</b>	135 mm x 260 mm	(200 pcs/box)
<b>6003</b>	90 mm x 145 mm	(200 pcs/box)	<b>6011</b>	190 mm x 330 mm	(200 pcs/box)
<b>6005</b>	70 mm x 230 mm	(200 pcs/box)	<b>6015</b>	250 mm x 380 mm	(100 pcs/box)



# sterilisation rolls

Sterilisation rolls are essential in healthcare and laboratory settings for their cost-effective and customizable sterilisation solutions. These rolls allow users to create pouches of varying lengths, accommodating instruments of different sizes while minimizing waste. Designed to withstand autoclaving temperatures, sterilisation rolls ensure a secure, uniform seal that maintains sterility until instruments are needed. Their transparent panels enable easy identification of contents without compromising sterility.


- Designed for use with a suitable heat-sealing equipment or autoclave tape to seal each end of the roll after being cut to size.
- Manufactured from 70 gsm medical grade paper and 60 µm blue tinted film to help detect puncture holes
- Triple adhesive seal to protect contents
- Easy to read external steam indicators change from blue to dark grey to ensure sterile conditions have been met
- Three universal sizes for all applications
- Complies with ISO 11607-1 and EN 868-5



6052 50 mm (200 m roll)

6053 75 mm (200 m roll)

6054 100 mm (200 m roll)



## specifications sterilisation pouches



**6000** 57 mm x 90 mm (200 pcs/box)  
*diamond and carbide burs, scaler tips, orthodontic appliances etc.*




**6001** 57 mm x 130 mm (200 pcs/box)  
*mirror heads, matrix retainers, orthodontic appliances, rubber dam clamps etc.*




**6003** 90 mm x 145 mm (200 pcs/box)  
*slow speed handpieces, scissors, small instruments etc.*






**6005** 70 mm x 230 mm (200 pcs/box)  
probes, currettes, scalers, tweezers and other hand instruments etc.




**6007** 90 mm x 230 mm (200 pcs/box)  
forceps, elevators, scissors and other hand instruments etc.





**6009** 135 mm x 260 mm (200 pcs/box)  
ruber dam forceps, rubber dam frames and other large hand instruments etc.





**6011** 190 mm x 330 mm (200 pcs/box)  
*large instruments and sterilisation cassettes*



**6015** 250 mm x 380 mm (100 pcs/box)  
*large instruments and sterilisation cassettes*





## specifications

# sterilisation cassettes


- Lightweight plastic for easy handling and transport
- More economical than expensive metal cassettes
- Rust and corrosion resistant
- Protects delicate instruments during sterilisation and reduces the risk of user injuries caused by heat



unlocked



locked



**3881 small (fits 5 instruments)**  
190 mm x 90 mm x 25 mm (1 pc/box)



**3882 medium (fits 8 instruments)**  
190 mm x 135 mm x 25 mm (1 pc/box)

**6011**

190 mm x 330 mm (200 pcs/box)

**3881**190 mm x 90 mm x 25 mm  
(1 pc/box)**3882**190 mm x 135 mm x 25 mm  
(1 pc/box)**6015**

250 mm x 380 mm (100 pcs/box)

**Available in the following colours:**

white

light blue

blue

lavender

green

# packaging



**6000**  
10 x boxes per carton



**6001**  
10 x boxes per carton



**6003**  
10 x boxes per carton



**6005**  
10 x boxes per carton



**6007**  
10 x boxes per carton



**6009**  
10 x boxes per carton



**6011**  
10 x boxes per carton



**6015**  
10 x boxes per carton



**3881**  
50 x sterilisation cassettes per carton



**3882**  
50 x sterilisation cassettes per carton



**6052**  
10 x rolls per carton



**6053**  
10 x rolls per carton



**6054**  
10 x rolls per carton

# order information

Order Code	Product Description	Packaging	Carton QTY	Carton Dimensions L x W x H	Carton Weight (kg)	Cartons per EU Pallet
6000	6000 Pouches 200s 57 mm x 90 mm	200 pcs/box	10 boxes/carton	33 cm x 15 cm x 15 cm	2.35 kg	208 cartons
6001	6001 Pouches 200s 57 mm x 130 mm	200 pcs/box	10 boxes/carton	35 cm x 19 cm x 15 cm	3.05 kg	156 cartons
6003	6003 Pouches 200s 90 mm x 145 mm	200 pcs/box	10 boxes/carton	35 cm x 21 cm x 21 cm	4.95 kg	90 cartons
6005	6005 Pouches 200s 70 mm x 230 mm	200 pcs/box	10 boxes/carton	35 cm x 29 cm x 16 cm	5.80 kg	96 cartons
6007	6007 Pouches 200s 90 mm x 230 mm	200 pcs/box	10 boxes/carton	35 cm x 29 cm x 21 cm	7.30 kg	72 cartons
6009	6009 Pouches 200s 135 mm x 260 mm	200 pcs/box	10 boxes/carton	33 cm x 32 cm x 30 cm	11.70 kg	36 cartons
6011	6011 Pouches 200s 190 mm x 330 mm	200 pcs/box	10 boxes/carton	42 cm x 38 cm x 35 cm	20.60 kg	20 cartons
6015	6015 Pouches 100s 250 mm x 380 mm	100 pcs/box	10 boxes/carton	44 cm x 34 cm x 28 cm	15.75 kg	28 cartons
6052	6052 Sterilisation Roll 200 m x 50 mm	200 m/box	10 boxes/carton	44 cm x 24 cm x 26 cm	13.75 kg	42 cartons
6053	6053 Sterilisation Roll 200 m x 75 mm	200 m/box	10 boxes/carton	44 cm x 22 cm x 40 cm	20.65 kg	30 cartons
6054	6054 Sterilisation Roll 200 m x 100 mm	200 m/box	10 boxes/carton	43 cm x 22 cm x 52 cm	27.65 kg	18 cartons
3881	3881 Plastic Sterilisation Cassette Small	1 pc/box	50 boxes/carton	48 cm x 30 cm x 21 cm	4.60 kg	36 cartons
3882	3882 Plastic Sterilisation Cassette Medium	1 pc/box	50 boxes/carton	71 cm x 32 cm x 21 cm	6.15 kg	27 cartons

# about Premium Plus

With over 50 years of experience in the dental industry, we specialise in the manufacturing and distribution of high-quality dental products that combine precision, affordability, and innovation. Our extensive knowledge ensures that every product we offer meets the highest standards of quality, supporting dental professionals worldwide.

Since 2011, we've been serving the European market, initially from the UK and, since 2020, from our distribution centre located in Poland. We are proud to be part of the global Premium Plus partnership, with strong connections across Hong Kong, the USA, Australia, Japan, and Europe. Our 10,000-square-meter facility allows us to deliver premium products at competitive prices.

## Mission Statement

Remember your first visit to the dentist? For many, it meant a stark room, unsettling sounds, and a wave of anxiety.

At Premium Plus, we're on a mission to change that. We believe a bright, healthy smile can uplift your entire day. That's why we've created a vibrant range of products to make dental care a positive, joyful experience. Our high-quality designs bring warmth and comfort to practices, helping patients relax and smile with confidence.

# why buy from us?

We understand that in today's fast-paced environment, you need more than just a product – you need a partner who supports your business every step of the way. That's why we offer:

## Regulatory Compliance

Full adherence to regulatory requirements, so you can be assured of the quality and safety of our products.

## Marketing Support

Tailored marketing assets such as flyers, catalogues, product photos and more, to help you stand out and grow your business with our support.

## Product Training

Comprehensive training to ensure your team is fully equipped to use and promote our products with confidence.

## Logistics Support

Reliable and efficient logistics solutions, ensuring timely delivery of your products every time.

## Manufacturing Experience

Decades of expertise in producing dental products that consistently meet industry standards and exceed expectations.





[www.premiumpluspoland.com](http://www.premiumpluspoland.com)

**European Office**  
**Premium Plus Poland sp. z o. o.**  
Bukowska 27  
62-081 Wysogotowo  
**Poland**

+48 61 880 1094  
info@premiumpluspl.com  
www.premiumpluspoland.com

**UK Office**  
**Premium Plus UK Ltd.**  
Unit 2, Knighton Heath Ind. Estate  
847 Ringwood Road,  
BH11 8NE Bournemouth, Dorset  
**United Kingdom**

+44 1202 611011  
info@premiumplusuk.com  
www.premiumplusuk.com

**Corporate Head Office**  
**Premium Plus International Ltd.**  
Flat 1001, Yuen Long Trading Centre,  
No.33 Wang Yip Street West,  
Yuen Long, N. T.  
**Hong Kong**

+852 2761 3211  
mailbox@premiumplushk.com  
www.premiumplushk.com

**Manufacturing Facility**  
**Premium Plus (Dongguan) Limited**  
Flat 101, No. 122, Tianqiao Road  
Changping, Dongguan,  
Guandong 523581  
**China**

+86 769 83397277  
info@premiumpluschina.com

**USA Office**  
**Premium Plus Dental Supplies Inc.**  
155 Liberty Avenue  
Brooklyn  
NY 11212  
**United States of America**

+1 718 366 8686  
info@premiumplusdentalusa.com  
www.premiumplususa.com

**Japan Office**  
**Premium Plus Japan Co., Ltd.**  
5th Floor, 1-1-31, Nishimidorigaoka,  
Toyonakashi,  
Osaka 560-0005  
**Japan**

+81 6 6845 0066  
premiumplus.info@premiumplus.jp  
www.premiumplus.jp