



intraoral photographic mirrors

Our titanium intraoral mirrors offer superior image clarity with a hard, flat surface that eliminates distortion, ensuring that every photograph is sharp and accurate. Crafted with 3mm thickness, these mirrors are designed for strength and precision. The front titanium coated surface provides maximum durability and reflection, while remaining lightweight and comfortable for both clinician and patient during use.

◆ KEY FEATURES FEATURES KEY

Double-Sided Design

Offering reflection on both sides, these mirrors provide flexibility in positioning and ensure better access to all areas in the oral cavity.

Durable Titanium Surface

The hard, flat surface prevents distortion in images, providing the clarity needed for accurate diagnostics and patient documentation.

Lightweight and Comfortable

Despite their durability, these mirrors are lightweight, reducing hand fatigue for clinicians and discomfort for patients.

Scratch-Resistant

The titanium surface is less prone to scratches or damage, ensuring long-lasting use and clear reflections over time.

Autoclavable

Designed for repeated use, our mirrors are autoclavable up to 134°C, making them safe for sterilisation in busy dental practices.

Versatile Shapes

Available in five different shapes, these mirrors cover all areas of the oral cavity, including occlusal, buccal, and lingual surfaces.

3mm Thickness

Built with strength, while being thin enough for easy maneuverability in the mouth.

◆ WHY CHOOSE OUR INTRAORAL MIRRORS? INTRAORAL MIRRORS? WHY CHOOSE OUR

These reversible mirrors are designed for both clinical precision and patient comfort, providing an essential tool for capturing high-quality images across a range of dental procedures and photographic applications.

◆
**OUR
INTRAORAL
MIRRORS**
MIRROR
INTRAORAL
OUR TITANIUM

PM001



Buccal

Designed for use in the buccal area of the mouth to view hard to see areas of the mouth.

PM002



Lingual

Designed for use in the lingual area of the mouth and to examine the tongue, cheeks and other soft tissue for abnormalities and lesions.

PM003



Child Occlusal

Designed to inspect the occlusal surfaces of children's molars and premolars, particularly for early detection of cavities or dental issues.


PM004



Adult Occlusal

Designed to inspect the occlusal surfaces of the adult teeth.

PM005




XL Occlusal

Designed to inspect the occlusal surfaces of adults with larger oral cavities.



INTRAORAL MIRROR CARE

- Clean with soapy water
- Use soft cloth or microfibre to avoid scratching
- Autoclavable to 134°C
- Avoid using paper towels, strong cleaning agents and ultrasonic cleaning as these can damage the mirror surfaces



PACKAGING

PM001	Buccal mirror	1 pc/pack
PM002	Lingual mirror	1 pc/pack
PM003	Child Occlusal mirror	1 pc/pack
PM004	Adult Occlusal mirror	1 pc/pack
PM005	XL Occlusal mirror	1 pc/pack



www.premiumpluspoland.com

Minimum order value of €500 applies. Pricing is exclusive to dealers.
Minimum order quantity 1 piece.