


premium  
pluspl™

# paper points & gutta percha


Reliable canal drying  
and accurate obturation.




radiopacity - clear  
visibility on x-ray



ISO consistent sizing  
for predictable fit



pre-measured depth markings  
for accurate placement



accurate tip sizing  
and consistent taper

[www.premiumpluspoland.com](http://www.premiumpluspoland.com)

# Paper Points

High absorbency paper points designed for predictable canal drying and efficient workflow.

## Use

For drying and moisture control of root canals.

## Key Features

- ISO colour coded sizes for easy identification
- High absorbency for efficient canal drying
- Pre-measured depth markings
- Non-sterile
- Supplied in a convenient slide box



## PP Points 2%

- Standard ISO paper points with a 0.02 taper.
- Suitable for use with conventional hand instrumentation.

| PP POINTS 2% |       |       |       |       |       |
|--------------|-------|-------|-------|-------|-------|
| #15 ○        | #20 ● | #25 ● | #30 ● | #35 ● | #40 ● |
| #45 ○        | #50 ● | #55 ● | #60 ● | #70 ● | #80 ● |

## PP Points 4%

- Standard ISO paper points with a 0.04 taper.
- Suitable for use with conventional hand instrumentation.
- Compatible with **Premium Plus E03, E18, E25, E28** rotary file instruments.

| PP POINTS 4% |       |       |       |       |       |
|--------------|-------|-------|-------|-------|-------|
| #15 ○        | #20 ● | #25 ● | #30 ● | #35 ● | #40 ● |
| #45 ○        | #50 ● | #55 ● | #60 ● | #70 ● | #80 ● |

## PP Points 6%

- Standard ISO paper points with a 0.06 taper.
- Suitable for use with conventional hand instrumentation.
- Compatible with **Premium Plus E03, E18, E25, E28** rotary file instruments.

| PP POINTS 6% |       |       |       |       |       |
|--------------|-------|-------|-------|-------|-------|
| #15 ○        | #20 ● | #25 ● | #30 ● | #35 ● | #40 ● |
| #45 ○        | #50 ● | #55 ● | #60 ● | #70 ● | #80 ● |

## PP ONE Points


- Matched taper paper points designed for reciprocating file systems. Available in commonly used sizes R25, R40 and R50.
- Compatible with **Premium Plus E05 & E06 ONE FILE** instruments.

| PP ONE POINTS |        |        |
|---------------|--------|--------|
| #R25 ●        | #R40 ● | #R50 ● |

## PP SUPER Points

- Matched taper paper points designed for variable taper file systems. Available in commonly used sizes F1, F2, F3, F4 and F5.
- Compatible with **Premium Plus E01 & E02 SUPER FILE** instruments.

| PP SUPER POINTS |      |      |      |      |
|-----------------|------|------|------|------|
| F1 ●            | F2 ● | F3 ● | F4 ● | F5 ● |



# Gutta Percha


Reliable obturation points designed for accuracy, consistency and compatibility across traditional and modern endodontic systems.

## Use

For obturation of root canals using cold obturation techniques such as single cone and lateral compaction or as a master cone for warm vertical condensation.

## Key Features

- Hand rolled construction
- Cadmium free and latex free
- Radiopaque for clear visibility on radiographs
- ISO colour coded sizes for easy identification
- Accurate tip sizing and consistent taper
- Suitable for warm and cold obturation techniques
- Low melting point for improved handling
- Pre-measured depth markings
- Supplied in a convenient slide box



## GP Points 2%

- Standard ISO gutta percha points with a 0.02 taper.
- Suitable for use with conventional hand instrumentation.

| GP 2% |       |       |       |       |       |
|-------|-------|-------|-------|-------|-------|
| #15 ○ | #20 ● | #25 ● | #30 ● | #35 ● | #40 ● |
| #45 ○ | #50 ● | #55 ● | #60 ● | #70 ● | #80 ● |

## GP Points 4%

- Standard ISO gutta percha points with a 0.04 taper.
- Suitable for use with conventional hand instrumentation.
- Compatible with **Premium Plus E03, E18, E25, E28** rotary file instruments.

| GP 4% |       |       |       |       |       |
|-------|-------|-------|-------|-------|-------|
| #15 ○ | #20 ● | #25 ● | #30 ● | #35 ● | #40 ● |
| #45 ○ | #50 ● | #55 ● | #60 ● | #70 ● | #80 ● |

## GP Points 6%

- Standard ISO gutta percha points with a 0.06 taper.
- Suitable for use with conventional hand instrumentation.
- Compatible with **Premium Plus E03, E18, E25, E28** rotary file instruments.

| GP 6% |       |       |       |       |       |
|-------|-------|-------|-------|-------|-------|
| #15 ○ | #20 ● | #25 ● | #30 ● | #35 ● | #40 ● |
| #45 ○ | #50 ● | #55 ● | #60 ● | #70 ● | #80 ● |

## GP ONE Points


- Matched taper gutta percha points designed for reciprocating file systems. Available in commonly used sizes R25, R40 and R50.
- Compatible with **Premium Plus E05 & E06 ONE FILE** instruments.

| GP ONE |        |        |
|--------|--------|--------|
| #R25 ● | #R40 ● | #R50 ● |

## GP SUPER Points

- Matched taper gutta percha points designed for variable taper file systems. Available in commonly used sizes F1, F2, F3, F4 and F5.
- Compatible with **Premium Plus E01 & E02 SUPER FILE** instruments.

| GP SUPER |      |      |      |      |
|----------|------|------|------|------|
| F1 ●     | F2 ● | F3 ● | F4 ● | F5 ● |



## packaging



**PP 2%**  
200 pcs/pack  
10 packs/bundle



**PP 4%**  
100 pcs/pack  
10 packs/bundle



**PP 6%**  
100 pcs/pack  
10 packs/bundle



**PP ONE**  
100 pcs/pack  
10 packs/bundle



**PP SUPER**  
100 pcs/pack  
10 packs/bundle



**GP 2%**  
120 pcs/pack  
10 packs/bundle



**GP 4%**  
60 pcs/pack  
10 packs/bundle



**GP 6%**  
60 pcs/pack  
10 packs/bundle



**GP ONE**  
60 pcs/pack  
10 packs/bundle



**GP SUPER**  
60 pcs/pack  
10 packs/bundle

| product range   | description               | sizes                                                                      | packaging    |
|-----------------|---------------------------|----------------------------------------------------------------------------|--------------|
| <b>PP 2%</b>    | 2% Taper                  | 15, 20, 25, 30, 35, 40, 45, 50, 55, 60, 70, 80,<br>Assorted 15-40 or 45-80 | 200 pcs/pack |
| <b>GP 2%</b>    |                           |                                                                            | 120 pcs/pack |
| <b>PP 4%</b>    | 4% Taper                  | 15, 20, 25, 30, 35, 40, 45<br>Assorted 15-40 or 45-80                      | 100 pcs/pack |
| <b>GP 4%</b>    |                           |                                                                            | 60 pcs/pack  |
| <b>PP 6%</b>    | 6% Taper                  | 15, 20, 25, 30, 35, 40, 45<br>Assorted 15-40 or 45-80                      | 100 pcs/pack |
| <b>GP 6%</b>    |                           |                                                                            | 60 pcs/pack  |
| <b>PP ONE</b>   | matched<br>variable taper | R25, R40, R50, Assorted R25-R50                                            | 100 pcs/pack |
| <b>GP ONE</b>   |                           |                                                                            | 60 pcs/pack  |
| <b>PP SUPER</b> | matched<br>variable taper | F1, F2, F3, F4, F5, Assorted F1-F3                                         | 100 pcs/pack |
| <b>GP SUPER</b> |                           |                                                                            | 60 pcs/pack  |